



**Technical Publications**

**DICOM Conformance Statement**

**Hipax DICOM Print Broker 1.0**

REV 0.1 140405a

April 2005

**Steinhart** Medizinsysteme GmbH  
**Medical Imaging, Archiving, and Communication**

© 1995-2005 Steinhart Medizinsysteme GmbH, Grubstraße 6-8, D-78279 Vörstetten





### Document Information

<b>Author</b>	Thomas Wittmer	Tel: +49 7666 9007-0 <a href="mailto:wittmer@hipax.de">wittmer@hipax.de</a>
<b>Editing &amp; Layout</b>	Thomas Wittmer	
<b>Project Manager</b>	Thomas Wittmer	
<b>Product Manager</b>	Thomas Wittmer	
<b>Service-related issues</b>	Steinhart Medizinsysteme GmbH Vörstetten	Tel: +49 7666 9007-0 <a href="mailto:support@hipax.de">support@hipax.de</a>

Issued by:  
Steinhart Medizinsysteme GmbH  
Grubstraße 6-8  
D-79279 Vörstetten  
Germany

Steinhart shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Steinhart reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Steinhart software by dedicated personnel. The information presented herein is sensitive and is classified Company Confidential. Without written authority from the proprietor, further distribution outside the company is not allowed.

Tel: +49 7666 9007 - 0  
Fax: +49 7666 9007 - 11  
Email: [info@hipax.de](mailto:info@hipax.de)

**Copyright © July, 05**  
Steinhart Medizinsysteme GmbH  
All rights reserved

### Revision Record

<b>Revision</b>	<b>Responsible</b>	<b>Date</b>	<b>Remark</b>
1.0	Thomas Wittmer	2005-04-14	Initial revision



### Conformance Statement Overview

This product, Hipax Print Broker, implements the necessary DICOM services to facilitate the Print Jobs (SCP/SCU) for Management in the healthcare departments. It enables the capabilities to capture images at any networked DICOM modality, administrate them and then send them anywhere they were needed in the medical facility.

**Table 1.1:** Network Services Supported

<b>SOP Classes</b>	<b>User of Service</b>	<b>Provider of Service</b>
Verification SOP Class	No	Yes
Basic Grayscale Print Management Meta SOP class	Yes	Yes
Basic Color Print Management Meta SOP class	Yes	Yes
Basic Annotation Box SOP Class	Yes	Yes
Print Job SOP Class	Yes	Yes
Presentation LUT SOP Class	Yes	Yes
Print Queue Management SOP Class	Yes	Yes



## Table of Contents

<b>1. USE OF THE DICOM PRINT BROKER .....</b>	<b>7</b>
<b>2. INTRODUCTION TO THIS DOCUMENT.....</b>	<b>7</b>
<b>2.1 PURPOSE AND INTENDED AUDIENCE OF THIS DOCUMENT .....</b>	<b>7</b>
<b>2.2 GENERAL REMARKS.....</b>	<b>7</b>
<b>2.2.1 Integration and Validation Activities .....</b>	<b>7</b>
<b>2.2.2 Future Evolution.....</b>	<b>7</b>
<b>2.2.3 Acronyms and Abbreviations .....</b>	<b>8</b>
<b>2.2.4 Related Documents .....</b>	<b>8</b>
<b>3. NETWORKING.....</b>	<b>9</b>
<b>3.1 IMPLEMENTATION MODEL .....</b>	<b>9</b>
<b>3.1.1 Application Data Flow Diagram .....</b>	<b>9</b>
<b>4. DICOM CONFORMANCE.....</b>	<b>10</b>
<b>4.1 HIPAX PRINT BROKER .....</b>	<b>10</b>
<b>4.1.1 Implementation Identifying Information.....</b>	<b>10</b>
<b>4.1.2 Supported SOP Classes .....</b>	<b>10</b>
<b>4.1.3 Supported Transfer Syntaxes .....</b>	<b>10</b>
<b>4.1.4 Number of Associatons .....</b>	<b>10</b>
<b>4.1.5 Asynchronous Nature .....</b>	<b>10</b>
<b>4.1.6 Accepted Presentation Contexts .....</b>	<b>11</b>
<b>4.1.7 SOP Specific Conformance.....</b>	<b>12</b>
<b>4.1.7.1 Specific Conformance for Verification SOP Class.....</b>	<b>12</b>
<b>4.1.7.2 Specific Conformance to Grayscale Print Management Meta SOP Class:.....</b>	<b>12</b>
<b>4.1.7.2.1 Specific Conformance for Basic Film Session SOP Class.....</b>	<b>12</b>
<b>4.1.7.2.1.1 Film Session SOP Class Operations for N-CREATE.....</b>	<b>13</b>
<b>4.1.7.2.2 Specific Conformance for Basic Film Box SOP Class .....</b>	<b>15</b>
<b>4.1.7.2.2.1 Basic Film Box SOP Class Operations for N-CREATE .....</b>	<b>15</b>
<b>4.1.7.2.3 Specific Conformance for Basic Grayscale Image Box SOP Class.....</b>	<b>17</b>
<b>4.1.7.2.3.1 Basic Grayscale Image Box SOP Class Operations for N-SET.....</b>	<b>17</b>
<b>4.1.7.2.4 Specific Conformance for Printer SOP Class .....</b>	<b>18</b>
<b>4.1.7.2.4.1 Printer SOP Class Operations for N-EVENT-REPORT.....</b>	<b>18</b>
<b>4.1.7.2.4.2 Printer SOP Class Operations for N-GET.....</b>	<b>19</b>
<b>4.1.7.3 Specific Conformance to Color Print Management Meta SOP Class:.....</b>	<b>19</b>
<b>4.1.7.3.1 Specific Conformance for Basic Film Session SOP Class.....</b>	<b>20</b>
<b>4.1.7.3.2 Specific Conformance for Basic Film Box SOP Class .....</b>	<b>20</b>
<b>4.1.7.3.3 Specific Conformance for Basic Color Image Box SOP Class .....</b>	<b>20</b>
<b>4.1.7.3.3.1 Basic Color Image Box SOP Class Operations for N-SET .....</b>	<b>20</b>
<b>4.1.7.3.4 Specific Conformance for Basic Color Printer SOP Class.....</b>	<b>21</b>
<b>4.1.7.4 Specific Conformance to Basic Annotation Box SOP Class .....</b>	<b>21</b>
<b>4.1.7.5 Specific Conformance to Print Job SOP Class.....</b>	<b>22</b>
<b>4.1.7.5.1 Print Job SOP Class N-EVENT-REPORT .....</b>	<b>22</b>
<b>4.1.7.5.2 Print Job SOP Class N-GET Request .....</b>	<b>23</b>
<b>4.1.7.6 Specific Conformance to Presentation LUT SOP Class.....</b>	<b>23</b>
<b>4.1.7.6.1 Presentation LUT SOP Class N-CREATE Request.....</b>	<b>23</b>
<b>4.1.7.6.2 Presentation LUT SOP Class N-DELETE Request .....</b>	<b>24</b>
<b>4.1.7.7 Specific Conformance to Print Queue Management SOP Class .....</b>	<b>24</b>
<b>4.1.7.7.1 Print Queue SOP Class N-EVENT-REPORT .....</b>	<b>24</b>
<b>4.1.7.7.2 Print Queue SOP Class N-GET Request .....</b>	<b>25</b>
<b>4.1.7.7.3 Print Queue SOP Class N-ACTION Request .....</b>	<b>26</b>
<b>4.2 NETWORK INTERFACES .....</b>	<b>26</b>
<b>5. MEDIA INTERCHANGE .....</b>	<b>26</b>





## 1. Use of the DICOM Print Broker

The Hipax DICOM Print Broker can be used for DICOM print job management. The program is able to administrate DICOM print jobs and to emit them bit by bit corresponding to the capability of the used printer. Print jobs can also be buffered and redirected. Thus, the Print Broker prevents that print jobs are cancelled and lost because of a capacity overload of the printer.

Furthermore, the Hipax DICOM Print Broker allows the user to process the same print job on different DICOM printers. To make this, the physical printers are related to one or several groups, called logical printers, using mapping and constraints, e.g. depending on the film size ID, medium type, or destination.

Using the Hipax DICOM Print Broker together with the Print Server, images can also be printed on paper via Windows printer.

## 2. Introduction to this Document

### 2.1 Purpose and Intended Audience of this Document

This document is a DICOM Conformance Statement for the DICOM Services of the Hipax Print Broker Software. The user of this document is involved with system integration and/or software design. We assume that the reader is familiar with the terminology and concepts that are used in the DICOM 3.0 standard and the IHE Technical Framework.

Readers not familiar with DICOM 3.0 terminology should first read the appropriate parts of the DICOM standard itself, prior to reading this conformance statement. Although the use of this conformance statement in conjunction with the DICOM 3.0 standard is intended to facilitate communication with Steinhart Medizinsysteme GmbH software, it is not sufficient to guarantee, by itself, the inter-operation of the connection.

### 2.2 General Remarks

#### 2.2.1 Integration and Validation Activities

The integration of any device into a system of interconnected devices goes beyond the scope of the DICOM 3.0 standard and this conformance statement when *interoperability* is desired. The responsibility for analyzing the applications requirements and developing a solution that integrates the Hipax software with other vendors. systems is the user.s responsibility and should not be underestimated. In some circumstances it might be necessary to perform a validation to make sure that functional interoperability between the Hipax software and non-Hipax devices works as expected. The user should ensure that any non-Hipax provider accepts responsibility for any validation required for their connection with the Hipax software.

#### 2.2.2 Future Evolution

As the DICOM 3.0 standard evolves to meet the user.s growing requirements and to incorporate new features and technologies, Steinhart Medizinsysteme GmbH will follow the evolution of the standard. This evolution of the standard may require changes to devices that have implemented DICOM 3.0. The user should ensure that any non-Hipax provider, who connects with Hipax devices, also plans for future evolution of the DICOM standard. A refusal to do so may result in the loss of functionality and/or connectivity between the different products.



### 2.2.3 Acronyms and Abbreviations

Definitions, terms and abbreviations used in this document are defined within the different parts of the DICOM standard. Abbreviations and terms are as follows:

AE	DICOM Application Entity
AET	Application Entity Title
ASCE	Association Control Service Element
DICOM	Digital Imaging and Communications in Medicine
FSC	File-Set Creator
FSU	File-Set Updater
FSR	File-Set Reader
GSDF	Grayscale Standard Display Function
GSPP	Grayscale Softcopy Presentation State
IE	Information Entity
IOD	(DICOM) Information Object Definition
ISO	International Standard Organization
MPPS	Modality Performed Procedure Step
MSPS	Modality Scheduled Procedure Step
PDU	DICOM Protocol Data Unit
SCU	DICOM Service Class User (DICOM client)
SCP	DICOM Service Class Provider (DICOM server)
SOP	DICOM Service-Object Pair
UID	Unique Identifier
VR	Value Representation

### 2.2.4 Related Documents

- ACR-NEMA Digital Imaging and Communications in Medicine (DICOM) V3.0. 2003.
- IHE Radiology Technical Framework Revision 5.5 . Final Text, November 20, 2003



### 3. Networking

#### 3.1 Implementation Model

##### 3.1.1 Application Data Flow Diagram

Print Composer SCU sends images and print management information.

----- DICOM Standard

Hipax Print Broker receives images, print management information and sends them to printers.

----- DICOM Standard

Printer receives images and print management information.



## 4. DICOM Conformance

### 4.1 Hipax Print Broker

#### 4.1.1 Implementation Identifying Information

The information listed in Table 4.1 is used to identify the Hipax Print Broker.

**Table 4.1:** Hipax Print Broker identifying information.

<b>Implementation Version Name</b>	HIPAXDICOM3V5000
<b>Implementation Class UID</b>	1.2.276.0.5420030101.5.0.0.0
<b>DICOM Application Context Name</b>	1.2.840.10008.3.1.1.1

#### 4.1.2 Supported SOP Classes

Hipax Print Broker provides standard conformance to the DICOM 3.0 Meta SOP Classes and DICOM 3.0 Optional SOP Classes as SCP listed in Table 4..

**Table 4.2:** Supported SOP Classes.

<b>SOP Class UID</b>	<b>SOP Classes Name</b>	<b>SCU</b>	<b>SCP</b>
1.2.840.10008.5.1.1.9	Basic Grayscale Print Management Meta SOP class	Yes	Yes
1.2.840.10008.5.1.1.18	Basic Color Print Management Meta SOP class	Yes	Yes
1.2.840.10008.1.1	Verification SOP Class	No	Yes
1.2.840.10008.5.1.1.15	Basic Annotation Box SOP Class	Yes	Yes
1.2.840.10008.5.1.1.14	Print Job SOP Class	Yes	Yes
1.2.840.10008.5.1.1.23	Presentation LUT SOP Class	Yes	Yes
1.2.840.10008.5.1.1.26	Print Queue Management SOP Class	Yes	Yes

#### 4.1.3 Supported Transfer Syntaxes

Hipax Print Broker provides support for each of the SOP Classes listed in Table 4. with the following Transfer Syntaxes:

**Table 4.3:** Supported Transfer Syntaxes.

<b>Name</b>	<b>UID</b>
DICOM Implicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2
DICOM Explicit VR Little Endian Transfer Syntax	1.2.840.10008.1.2.1
DICOM Explicit VR Big Endian Transfer Syntax	1.2.840.10008.1.2.2

#### 4.1.4 Number of Associatons

Unlimited on Plus Version, 4 SCP / 4 SCU on Standard Version

#### 4.1.5 Asynchronous Nature

Maximum Number of outstanding asynchronous transactions: 0



### 4.1.6 Accepted Presentation Contexts

**Table 4.4:** Presentation Contexts Accepted/Proposed by Hipax Print Broker.

<i>Abstract Syntax</i>		<i>Transfer Syntax</i>		<i>Role</i>	<i>Extended Negotiation</i>
<i>Name</i>	<i>UID</i>	<i>Name</i>	<i>UID</i>		
Basic Grayscale Print Management Meta SOP class	1.2.840.10008.5.1.1.9	DICOM Implicit VR Little Endian	1.2.840.10008.1.2	SCP SCU	None
		DICOM Explicit VR Little Endian	1.2.840.10008.1.2.1		
		DICOM Explicit VR Big Endian	1.2.840.10008.1.2.2		
Basic Color Print Management Meta SOP class	1.2.840.10008.5.1.1.18	DICOM Implicit VR Little Endian	1.2.840.10008.1.2	SCP SCU	None
		DICOM Explicit VR Little Endian	1.2.840.10008.1.2.1		
		DICOM Explicit VR Big Endian	1.2.840.10008.1.2.2		
Verification SOP Class	1.2.840.10008.1.1	DICOM Implicit VR Little Endian	1.2.840.10008.1.2	SCP	None
		DICOM Explicit VR Little Endian	1.2.840.10008.1.2.1		
		DICOM Explicit VR Big Endian	1.2.840.10008.1.2.2		
Basic Annotation Box SOP Class	1.2.840.10008.5.1.1.15	DICOM Implicit VR Little Endian	1.2.840.10008.1.2	SCP SCU	None
		DICOM Explicit VR Little Endian	1.2.840.10008.1.2.1		
		DICOM Explicit VR Big Endian	1.2.840.10008.1.2.2		
Print Job SOP Class	1.2.840.10008.5.1.1.14	DICOM Implicit VR Little Endian	1.2.840.10008.1.2	SCP SCU	None
		DICOM Explicit VR Little Endian	1.2.840.10008.1.2.1		
		DICOM Explicit VR Big Endian	1.2.840.10008.1.2.2		
Presentation LUT SOP Class	1.2.840.10008.5.1.1.23	DICOM Implicit VR Little Endian	1.2.840.10008.1.2	SCP SCU	None
		DICOM Explicit VR Big Endian	1.2.840.10008.1.2.1		
		DICOM Explicit VR Big Endian	1.2.840.10008.1.2.2		



**Table 4.4 (continue):** Presentation Contexts Accepted/Proposed by Hipax Print Broker.

<b>Abstract Syntax</b>		<b>Transfer Syntax</b>		<b>Role</b>	<b>Extended Negotiation</b>
<b>Name</b>	<b>UID</b>	<b>Name</b>	<b>UID</b>		
Print Queue Management SOP Class	1.2.840.10008.5.1.1.26	DICOM Implicit VR Little Endian	1.2.840.10008.1.2	SCP SCU	None
		DICOM Explicit VR Little Endian	1.2.840.10008.1.2.1		
		DICOM Explicit VR Big Endian	1.2.840.10008.1.2.2		

### 4.1.7 SOP Specific Conformance

#### 4.1.7.1 Specific Conformance for Verification SOP Class

The Hipax Print Broker provides standard conformance to the DICOM Verification Service Class as a SCP. The following status code is returned in response to a C-ECHO.

**Table 4.5:** C-ECHO Response Status.

<b>Service Status</b>	<b>Further Meaning</b>	<b>Error Code</b>	<b>Reason</b>
Success	Success	\$0000	

#### 4.1.7.2 Specific Conformance to Grayscale Print Management Meta SOP Class:

The Hipax Print Brokersupports the following mandatory SOP classes as defined by the Basic Grayscale Print Management Meta SOP Class:

**Table 4.6:** SOP Classes for Basic Grayscale Print Management Meta SOP Class.

<b>SOP Class Name</b>	<b>SOP Class UID</b>	<b>SCU</b>	<b>SCP</b>
Basic Film Session	1.2.840.10008.5.1.1.1	Yes	Yes
Basic Film Box	1.2.840.10008.5.1.1.2	Yes	Yes
Basic Grayscale Image Box	1.2.840.10008.5.1.1.4	Yes	Yes
Printer	1.2.840.10008.5.1.1.16	Yes	Yes

#### 4.1.7.2.1 Specific Conformance for Basic Film Session SOP Class

The Hipax Print Broker accepts all DICOM Print related attributes. Any DICOM Print related attributes that are received by Hipax Print Broker are sent to the associated printer without change. Attributes not related to DICOM Print are not sent to the associated printers.

The Hipax Print Broker provides support for the following DIMSE Services:

- N-CREATE
- N-SET
- N-ACTION
- N-DELETE

Hipax Print Broker only supports one Basic Film Session per Association. However, a sequential Film Session on the same Association is allowed after deleting the previous Film Session.



### 4.1.7.2.1.1 Film Session SOP Class Operations for N-CREATE

The Hipax Print Broker provides the following support for the Film Session attributes sent by the N-CREATE DIMSE service:

**Table 4.7:** Film Session SOP Class Operations for N-CREATE.

<b>Attribute</b>	<b>Tag</b>	<b>Valid Range</b>	<b>Default Value If not sent by SCU or invalid value received</b>
Number of Copies	(2000,0010)		
Print Priority	(2000,0020)		
Medium Type	(2000,0030)		
Film Destination	(2000,0040)		
Film Session Label	(2000,0050)		
Memory Allocation	(2000,0060)		
Owner ID	(2100,0160)		
Proposed Study Seq.	(2130,00A0)		
>Patient.s Name	(0010,0010)		
>Patient ID	(0010,0020)		
>Patient.s Birth Date	(0010,0030)		
>Patient.s Birth Time	(0010,0032)		
>Patient.s Sex	(0010,0040)		
>Other Patient ID	(0010,1000)		
>Other Patient Names	(0010,1001)		
>Patient.s Age	(0010,1010)		
>Patient.s Size	(0010,1020)		
>Patient Weight	(0010,1030)		
>Ethnic Group	(0010,2160)		
>Occupation	(0010,2180)		
>Add. Patient.s History	(0010,21B0)		
>Patient Comments	(0010,4000)		
>Study ID	(0020,0010)		
>Series Number	(0020,0011)		
>Study Instance UID	(0020,000D)		
>Study Date	(0008,0020)		
>Study Time	(0008,0030)		
>Accession Number	(0008,0050)		
>Ref. Physician.s Name	(0008,0090)		
>Study Description	(0008,1030)		
>Name of Physician Reading Study	(0008,1060)		
>Admitting Diagnosis Description	(0008,1080)		



The following status codes are returned in response to a N-CREATE:

**Table 4.8:** N-CREATE Response Status.

<b>Service Status</b>	<b>Further Meaning</b>	<b>Error Code</b>	<b>Reason</b>
Success	Success	\$0000	Operation successfully completed
Failure	Invalid Attribute Value	\$0106	Invalid Attribute Value
Failure	Class-Instance Conflict	\$0119	Returned if the SOP Class Instance UID is not defined for the given SOP Class.
Failure	SOP Class Not Supported	\$0122	Returned if the SOP Class specified in the DIMSE command does not match the SOP Class UID.

The following status codes are returned in response to a N-SET:

**Table 4.9:** N-SET Response Status.

<b>Service Status</b>	<b>Further Meaning</b>	<b>Error Code</b>	<b>Reason</b>
Success	Success	\$0000	Operation successfully completed
Failure	Invalid Attribute Value	\$0106	Invalid Attribute Value
Failure	Class-Instance Conflict	\$0119	Returned if the SOP Class Instance UID is not defined for the given SOP Class.
Failure	SOP Class Not Supported	\$0122	Returned if the SOP Class specified in the DIMSE command does not match the SOP Class UID.

The following status codes are returned in response to a N-ACTION:

**Table 4.10:** N-ACTION Response Status.

<b>Service Status</b>	<b>Further Meaning</b>	<b>Error Code</b>	<b>Reason</b>
Success	Success	\$0000	Operation successfully completed
Failure	SOP Class Not Supported	\$0122	Returned if the SOP Class specified in the DIMSE command does not match the SOP Class UID.
Failure		\$C600	Film Session SOP Instance hierarchy does not contain Film Box SOP Instances
Failure		\$C601	Returned if the support for the Print Job Class was negotiated and the creation of the Print Job Instance failed.

The following status codes are returned in response to a N-DELETE:

**Table 4.11:** N-DELETE Response Status.

<b>Service Status</b>	<b>Further Meaning</b>	<b>Error Code</b>	<b>Reason</b>
Success	Success	\$0000	Operation successfully completed
Failure	SOP Class Not Supported	\$0122	Returned if the SOP Class specified in the DIMSE command does not match the SOP Class UID.



### 4.1.7.2.2 Specific Conformance for Basic Film Box SOP Class

The Hipax Print Broker accepts all DICOM Print related attributes. Any DICOM Print related attributes that are received by Hipax Print Broker are sent to the associated printer without change. Attributes not related to DICOM Print are not sent to the associated printers.

The Hipax Print Broker provides support for the following DIMSE Services:

- N-CREATE
- N-SET
- N-ACTION
- N-DELETE

#### 4.1.7.2.2.1 Basic Film Box SOP Class Operations for N-CREATE

The Hipax Print Broker provides the following support for the Basic Film Box attributes sent by the N-CREATE DIMSE service:

**Table 4.12:** Basic Film Box SOP Class Operations for N-CREATE.

<i>Attribute</i>	<i>Tag</i>	<i>Valid Range</i>	<i>Default Value</i>
Image Display Format	(2010,0010)	STANDARD ROW COL	
Annotation Display Format ID	(2010,0030)	ANNOTATION (Supported if the Annotation SOP Class is accepted)	
Film Orientation	(2010,0040)		
Film Size ID	(2010,0050)		
Magnification Type	(2010,0060)		
Smoothing Type	(2010,0080)		
Border Density	(2010,0100)		
Empty Image Density	(2010,0110)		
Minimum Density	(2010,0120)		
Maximum Density	(2010,0130)		
Trim	(2010,0140)		
Configuration Information	(2010,0150)		
Illumination	(2010,015E)		
Reflective Ambient Light	(2010,0160)		



The following status codes are returned in response to a N-CREATE:

**Table 4.13:** N-CREATE Response Status.

<b>Service Status</b>	<b>Further Meaning</b>	<b>Error Code</b>	<b>Reason</b>
Success	Success	\$0000	Operation successfully completed
Failure	Invalid Attribute Value	\$0106	Invalid Attribute Value
Failure	Class-Instance Conflict	\$0119	Returned if the SOP Class Instance UID is not defined for the given SOP Class.
Failure	SOP Class Not Supported	\$0122	Returned if the SOP Class specified in the DIMSE command does not match the SOP Class UID.

The following status codes are returned in response to a N-SET:

**Table 4.14:** N-SET Response Status.

<b>Service Status</b>	<b>Further Meaning</b>	<b>Error Code</b>	<b>Reason</b>
Success	Success	\$0000	Operation successfully completed
Failure	Invalid Attribute Value	\$0106	Invalid Attribute Value
Failure	Class-Instance Conflict	\$0119	Returned if the SOP Class Instance UID is not defined for the given SOP Class.
Failure	SOP Class Not Supported	\$0122	Returned if the SOP Class specified in the DIMSE command does not match the SOP Class UID.

The following status codes are returned in response to a N-ACTION:

**Table 4.15:** N-ACTION Response Status.

<b>Service Status</b>	<b>Further Meaning</b>	<b>Error Code</b>	<b>Reason</b>
Success	Success	\$0000	Operation successfully completed
Failure	SOP Class Not Supported	\$0122	Returned if the SOP Class specified in the DIMSE command does not match the SOP Class UID.
Failure		\$C600	Film Session SOP Instance hierarchy does not contain Film Box SOP Instances
Failure		\$C602	Returned if the support for the Print Job Class was negotiated and the creation of the Print Job Instance failed.

The following status codes are returned in response to a N-DELETE:

**Table 4.16:** N-DELETE Response Status.

<b>Service Status</b>	<b>Further Meaning</b>	<b>Error Code</b>	<b>Reason</b>
Success	Success	\$0000	Operation successfully completed
Failure	SOP Class Not Supported	\$0122	Returned if the SOP Class specified in the DIMSE command does not match the SOP Class UID.



### 4.1.7.2.3 Specific Conformance for Basic Grayscale Image Box SOP Class

The Hipax Print Broker accepts all DICOM Print related attributes. Any DICOM Print related attributes that are received by Hipax Print Broker are sent to the associated printer without change. Attributes not related to DICOM Print are not sent to the associated printers.

The Hipax Print Broker provides support for the following DIMSE Services:

- N-SET

#### 4.1.7.2.3.1 Basic Grayscale Image Box SOP Class Operations for N-SET

The Hipax Print Broker provides the following support for the Basic Grayscale Image Box attributes sent by the N-SET DIMSE service:

**Table 4.17:** Basic Grayscale Image Box SOP Class Operations for N-SET.

<i>Attribute</i>	<i>Tag</i>	<i>Valid Range</i>	<i>Default Value If not sent by SCU or invalid value received</i>
Magnification Type	(2010,0060)		
Smoothing Type	(2010,0080)		
Minimum Density	(2010,0120)		
Maximum Density	(2010,0130)		
Image Position	(2020,0010)		
Polarity	(2020,0020)		
Requested Image Size	(2020,0030)		
Basic Grayscale Image Sequence	(2020,0110)		
>Samples Per Pixel	(0028,0002)		
>Photometric Interpretation	(0028,0004)		
>Rows	(0028,0010)		
>Columns	(0028,0011)		
>Pixel Aspect Ratio	(0028,0034)		
>Bits Allocated	(0028,0100)		
>Bits Stored	(0028,0101)		
>High Bit	(0028,0102)		
>Pixel Representation	(0028,0103)		
>Pixel Data	(7FE0,0010)		
Ref. Presentation LUT Seq.	(2050,0500)		
>Ref. SOP Class UID	(0008,1150)		
>Ref. SOP Instance UID	(0008,1155)		



The following status codes are returned in response to a N-SET:

**Table 4.18:** N-SET Response Status.

<b>Service Status</b>	<b>Further Meaning</b>	<b>Error Code</b>	<b>Reason</b>
Success	Success	\$0000	Operation successfully completed
Failure	Invalid Attribute Value	\$0106	Invalid Attribute Value
Failure	Class-Instance Conflict	\$0119	Returned if the SOP Class Instance UID is not defined for the given SOP Class.
Failure	SOP Class Not Supported	\$0122	Returned if the SOP Class specified in the DIMSE command does not match the SOP Class UID.

#### 4.1.7.2.4 Specific Conformance for Printer SOP Class

The Hipax Print Broker provides support for the following DIMSE Services:

- N-EVENT-REPORT
- N-GET

##### 4.1.7.2.4.1 Printer SOP Class Operations for N-EVENT-REPORT

At any time during the Association, Hipax Print Broker may send an N-EVENT-REPORT to the SCU to report the printer status (which is a combined status of all associated printers). Sending a N-EVENT-REPORT is disabled by default. It is also possible to configure Hipax Print Broker so that the printer status attribute WARNING and FAILURE is sent along with more detailed Printer Status Information. Sending WARNINGS and FAILURES is disabled by default.

The ADPM provides the following support for the Printer attributes sent by the N-EVENT-REPORT DIMSE service

**Table 4.19:** N-EVENT-REPORT Request Attributes.

<b>Printer Status (2110,0010)</b>	<b>Printer Status Info (2110,0020)</b>	<b>Meaning</b>
NORMAL	NORMAL	Normal printer operation
FAILURE (*)	UNKNOWN	
WARNING (*)	BAD SUPPLY MGZ	Film supply tray open.
WARNING (*)	COVER OPEN	The printer top cover is open.
WARNING (*)	SUPPLY EMPTY	The printer film tray is empty. Spooling of print jobs is still possible.
WARNING (*)	FILM JAM	Film jam.

(\*) Warnings and failures are not returned by default. Enabling warnings and failures is explained in the user Manual.

Any other attributes provided by any associated printer will also be forwarded to the associated SCU if failures and warnings are enabled.



### 4.1.7.2.4.2 Printer SOP Class Operations for N-GET

At any time during the Association, Hipax Print Broker application may receive an N-GET request asking for the printer status. It is also possible to configure Hipax Print Broker so that the printer status attributes WARNING and FAILURE is sent along with more detailed Printer Status Information. The sending of WARNINGS and FAILURES is disabled by default.

Hipax Print Broker sends the following attributes in response to an N-GET request:

**Table 4.20:** Printer SOP Class N-GET Request Attributes.

<b>Attribute</b>	<b>Tag</b>	<b>Valid Range</b>
Printer Status	(2110,0010)	NORMAL FAILURE WARNING
Printer Status Info	(2110,0020)	NORMAL UNKNOWN
Printer Name	(2110,0030)	
Manufacturer	(0008,0070)	
Manufacturer Model Name	(0008,1090)	
Device Serial Number	(0018,1000)	
Software Versions	(0018,1020)	
Date Last Calibration	(0018,1200)	
Time Last Calibration	(0018,1201)	

**Table 4.21:** N-GET Response Status.

<b>Service Status</b>	<b>Further Meaning</b>	<b>Error Code</b>	<b>Reason</b>
Success	Success	\$0000	Operation successfully completed
Failure	Class-Instance Conflict	\$0119	Returned if the SOP Class Instance UID is not defined for the given SOP Class.
Failure	SOP Class Not Supported	\$0122	Returned if the SOP Class specified in the DIMSE command does not match the SOP Class UID.

### 4.1.7.3 Specific Conformance to Color Print Management Meta SOP Class:

The Hipax Print Broker supports the following mandatory SOP classes as defined by the Basic Color Print Management Meta SOP Class:

**Table 4.22:** SOP classes for Basic Color Print Management Meta SOP Class.

<b>Name</b>	<b>UID</b>	<b>SCU</b>	<b>SCP</b>
Basic Film Session	1.2.840.10008.5.1.1.1	Yes	Yes
Basic Film Box	1.2.840.10008.5.1.1.2	Yes	Yes
Basic Color Image Box	1.2.840.10008.5.1.1	Yes	Yes
Printer	1.2.840.10008.5.1.1.16	Yes	Yes



### 4.1.7.3.1 Specific Conformance for Basic Film Session SOP Class

Refer to the 'Basic Film Session SOP Class' for 'Basic Grayscale Print Management'.

### 4.1.7.3.2 Specific Conformance for Basic Film Box SOP Class

Refer to the 'Basic Film Box SOP Class' for 'Basic Grayscale Print Management'.

### 4.1.7.3.3 Specific Conformance for Basic Color Image Box SOP Class

The Hipax Print Broker accepts all DICOM Print related attributes. Any DICOM Print related attributes that are received by Hipax Print Broker are sent to the associated printer without change. Attributes not related to DICOM Print are not sent to the associated printers.

The Hipax Print Broker provides support for the following DIMSE Services:

- N-SET

#### 4.1.7.3.3.1 Basic Color Image Box SOP Class Operations for N-SET

The Hipax Print Broker provides the following support for the Basic Color Image Box attributes sent by the N-SET DIMSE service:

**Table 4.23:** Basic Color Image Box SOP Class Operations for N-SET.

<b>Attribute</b>	<b>Tag</b>	<b>Valid Range</b>
Magnification Type	(2010,0060)	
Smoothing Type	(2010,0080)	
Minimum Density	(2010,0120)	
Maximum Density	(2010,0130)	
Image Position	(2020,0010)	
Polarity	(2020,0020)	
Requested Image Size	(2020,0030)	
Basic Color Image Sequence	(2020,0111)	
>Samples Per Pixel	(0028,0002)	
>Photometric Interpretation	(0028,0004)	
>Planar Configuration	(0028,0006)	
>Rows	(0028,0010)	
>Columns	(0028,0011)	
>Pixel Aspect Ratio	(0028,0034)	
>Bits Allocated	(0028,0100)	
>Bits Stored	(0028,0101)	
>High Bit	(0028,0102)	
>Pixel Representation	(0028,0103)	
>Pixel Data	(7FE0,0010)	



The following status codes are returned in response to a N-SET:

**Table 4.24:** N-SET Response Status.

<b>Service Status</b>	<b>Further Meaning</b>	<b>Error Code</b>	<b>Reason</b>
Success	Success	\$0000	Operation successfully completed
Failure	Invalid Attribute Value	\$0106	Invalid Attribute Value
Failure	Class-Instance Conflict	\$0119	Returned if the SOP Class Instance UID is not defined for the given SOP Class.
Failure	SOP Class Not Supported	\$0122	Returned if the SOP Class specified in the DIMSE command does not match the SOP Class UID.

#### 4.1.7.3.4 Specific Conformance for Basic Color Printer SOP Class

Refer to the 'Basic Film Box SOP Class' for 'Basic Grayscale Print Management'.

#### 4.1.7.4 Specific Conformance to Basic Annotation Box SOP Class

The Hipax Print Broker provides standard conformance to the DICOM Basic Annotation Box SOP Class.

The Hipax Print Broker accepts all DICOM Print related attributes. Any DICOM Print related attributes that are received by Hipax Print Broker are sent to the associated printer without change. Attributes not related to DICOM Print are not sent to the associated printers.

The Basic Annotation Box SOP Instance is created by a N-CREATE of the Film Box SOP Class, if the "Annotation Display Format ID" attribute (2010,0030) has the value "ANNOTATION".

The Hipax Print Broker provides support for the following DIMSE Services:

- N-SET

Hipax Print Broker sends the following attributes in response to an N-SET request:

**Table 4.25:** Basic Annotation Box SOP Class N-SET Request Attributes.

<b>Attribute</b>	<b>Tag</b>	<b>Valid Range</b>
Annotation Position	(2030,0010)	1-6 (One for each Annotation Box)
Text String	(2030,0020)	Refer to explanation below

The following status codes are returned in response to a N-SET:

**Table 4.26:** N-SET Response Status.

<b>Service Status</b>	<b>Further Meaning</b>	<b>Error Code</b>	<b>Reason</b>
Success	Success	\$0000	Operation successfully completed
Failure	Class-Instance Conflict	\$0119	Returned if the SOP Class Instance UID is not defined for the given SOP Class.
Failure	SOP Class Not Supported	\$0122	Returned if the SOP Class specified in the DIMSE command does not match the SOP Class UID.



### 4.1.7.5 Specific Conformance to Print Job SOP Class

The Hipax Print Broker provides standard conformance to the DICOM Print Job SOP Class.

The Print Job SOP Instance is created by a N-ACTION of the Film Session SOP Class or the N-ACTION of the Film Box SOP Class. The Print Job instance is deleted after all films are printed or in event of an error.

The following DIMSE Service are supported:

- N-EVENT REPORT
- N-GET

#### 4.1.7.5.1 Print Job SOP Class N-EVENT-REPORT

N-EVENT-REPORT is used to report execution status changes to the SCU in an asynchronous way. N-EVENT-REPORT is disabled by default. Following Event Types and ID's are supported:

**Table 4.27:** Print Job SOP Class N-EVENT-REPORT.

<i>Event Type Name</i>	<i>Event Type ID</i>	<i>Description</i>
PENDING	1	Print job is pending
PRINTING	2	Print job is being printed
DONE	3	Print job is printed
FAILURE	4	Print job failed

The following values are supported for "Execution Status Info" attribute:

**Table 4.28:** Supported Attributes.

<i>Attribute</i>	<i>Tag</i>	<i>Valid Range</i>
Execution Status Info	(2100,0030)	Refer to Table 4.7: Execution Status Info
Print Job ID	(2100,0010)	
Film Session Label	(2000,0050)	
Printer Name	(2110,0030)	HipaxDicomPrint

The following values are supported for "Execution Status Info" attribute:

**Table 4.29:** Supported Attributes.

<i>Execution Status Info</i>	<i>Execution Status Info</i>	<i>Meaning</i>
PRINTING	"NORMAL"	
DONE	"NORMAL"	
FAILURE	"INVALID PAGE DES"	The specified page layout cannot be printed or other page description errors have been detected.
FAILURE	„INSUFFIC MEMORY“	There is not enough memory available to complete this.



### 4.1.7.5.2 Print Job SOP Class N-GET Request

N-GET is used to retrieve an instance of the Print Job SOP Class. Hipax Print Broker sends the following attributes in response to a Print Job N-GET request:

**Table 4.30:** Print Job SOP Class N-GET Request Attributes.

Attribute	Tag	Valid Range
Execution Status	(2100,0020)	PENDING PRINTING DONE FAILURE
Execution Status Info	(2100,0030)	Refer to Table 4.7: Execution Status Info
Print Priority	(2000,0020)	LOW MED HIGH
Creation Date	(2100,0040)	Date of print job creation
Creation Time	(2100,0050)	Time of print job creation
Originator	(2100,0070)	Calling AE Title
Printer Name	(2110,0030)	

### 4.1.7.6 Specific Conformance to Presentation LUT SOP Class

The Hipax Print Broker provides standard conformance to the DICOM Presentation LUT Service Class.

The Hipax Print Broker provides support for the following DIMSE Services:

- N-CREATE
- N-DELETE

#### 4.1.7.6.1 Presentation LUT SOP Class N-CREATE Request

N-CREATE is used to create a Presentation LUT SOP Instance. The Hipax Print Broker provides the following support for the Presentation LUT SOP Class attributes sent by the N-CREATE DIMSE service:

**Table 4.31:** Presentation LUT SOP Class N-Create Request Attributes.

Attribute	Tag	Valid Range
Presentation LUT sequence	(2050,0010)	
>LUT Descriptor	(0028,3002)	
>LUT Explanation	(0028,3003)	
>LUT Data	(0028,3006)	
Presentation LUT Shape	(2050,0020)	IDENTITY LIN OD



The following status codes are returned in response to a N-CREATE:

**Table 4.32:** N-CREATE Response Status.

<i>Service Status</i>	<i>Further Meaning</i>	<i>Error Code</i>	<i>Reason</i>
Success	Success	\$0000	Operation successfully completed
Failure	Invalid Attribute Value	\$0106	Invalid Attribute Value

#### 4.1.7.6.2 Presentation LUT SOP Class N-DELETE Request

Used to delete a Presentation LUT SOP Instance. The following status codes are returned in response to a N-DELETE:

**Table 4.33:** N-DELETE Response Status.

<i>Service Status</i>	<i>Further Meaning</i>	<i>Error Code</i>	<i>Reason</i>
Success	Success	\$0000	Operation successfully completed
Failure	SOP Class Not Supported	\$0122	Returned if the SOP Class specified in the DIMSE command does not match the SOP Class UID.

#### 4.1.7.7 Specific Conformance to Print Queue Management SOP Class

The Hipax Print Broker provides standard conformance to the DICOM Queue Management Service Class.

The Printer Queue Management SOP Class is used to monitor and manipulate the print queue of the printer. The Printer Queue Management SOP instance is created by ADPM during the start-up of the device and has a SOP instance UID of: **1.2.840.10008.5.1.1.25**. The print queue is restored after application startup.

The Hipax Print Broker provides support for the following DIMSE Services:

- N-EVENT-REPORT
- N-GET
- N-ACTION

##### 4.1.7.7.1 Print Queue SOP Class N-EVENT-REPORT

N-EVENT-REPORT is used to report execution status changes to the SCU in an asynchronous way. N-EVENT-REPORT is disabled by default. The following Event Names and Event Type ID's are supported:

**Table 4.34:** N-EVENT-REPORT Event Type ID's.

<i>Event Type Name</i>	<i>Event Type ID</i>	<i>Description</i>
HALTED	1	Queue operation is halted
FULL	2	Queue is full
NORMAL	3	Queue is operational



### 4.1.7.7.2 Print Queue SOP Class N-GET Request

N-GET is used by the SCU to retrieve an instance of the Print Queue Management SOP Class from Hipax Print Broker. Hipax Print Broker provides the following support for the Print Queue Management SOP Class attributes sent by the N-GET DIMSE service:

**Table 4.35:** Print Job SOP Class N-GET Request Attributes.

<b>Attribute</b>	<b>Tag</b>	<b>Valid Range</b>
Queue Status	(2120,0010)	FULL HALTED NORMAL
Print Job Description Sequence	(2120,0050)	
>Print Job ID	(2100,0010)	
>Execution Status	(2100,0020)	PENDING PRINTING DONE FAILURE
>Execution Status Info	(2100,0030)	Refer to Table 4.6: Execution Status Info
>Creation Date	(2100,0040)	Date of print job creation
>Creation Time	(2100,0005)	Time of print job creation
>Print Priority	(2000,0020)	
>Origin AE	(2100,0070)	Calling AE Title
>Destination AE	(2100,0140)	Called AE Title
>Printer Name	(2110,0030)	
>Film Destination	(2000,0040)	
>Film Session Label	(2000,0050)	
>Medium Type	(2000,0030)	
>Number Of Films	(2100,0170)	
>Referenced Print Job Sequence	(2120,0070)	
>>Referenced SOP Class UID	(0008,1150)	
>>Referenced SOP Instance UID	(0008,1155)	

The following status codes are returned in response to a N-GET:

**Table 4.36:** N-GET Response Status.

<b>Service Status</b>	<b>Further Meaning</b>	<b>Error Code</b>	<b>Reason</b>
Success	Success	\$0000	Operation successfully completed
Failure	SOP Class Not Supported	\$0122	Returned if the SOP Class specified in the DIMSE command does not match the SOP Class UID.



### 4.1.7.7.3 Print Queue SOP Class N-ACTION Request

N-ACTION is used by the SCU to manipulate the Brokers queue content. Following "Event Types" are supported:

**Table 4.37:** "Event Types" supported.

<i>Event Type Name</i>	<i>Event Type ID</i>	<i>Description</i>
PRIORITIZE	1	Change priority of queue entry
DELETE	2	Delete queue entry

The Hipax Print Broker provides the following support for the Print Queue Management SOP Class attributes sent by the N-ACTION DIMSE service:

**Table 4.38:** Supported Attributes.

<i>Attribute</i>	<i>Tag</i>	<i>Valid Range</i>
Print Job ID	(2100,0010)	
Print Priority	(2000,0020)	HIGH LOW
Owner ID	(2100,0160)	

The following status codes are returned in response to a N-ACTION:

**Table 4.39:** N-ACTION Response Status.

<i>Service Status</i>	<i>Further Meaning</i>	<i>Error Code</i>	<i>Reason</i>
Success	Success	\$0000	Operation successfully completed
Failure	SOP Class Not Supported	\$0122	Returned if the SOP Class specified in the DIMSE command does not match the SOP Class UID.

## 4.2 Network Interfaces

The Hipax Print Broker provides DICOM V3.0 TCP/IP Network Communication Support as defined in PS 3.8 of the DICOM Standard.

## 5. Media Interchange

Not supported.

### Steinhart Medizinsysteme GmbH

Grubstraße 6-8 – D-79279 Vörstetten

Tel.: +49(0)7666-9007-0 – Fax: +49(0)7666-9007-11 – E-mail: info@hipax.de – Internet: www.hipax.de