

# MAMORAY HDR-C FILM

HDR-C film is an integral part of the Agfa HD Mammography System.

- Exceptionally high contrast
- High dynamic range
- RP (standard) processing cycle
- Anti-halation layer for artifact resistance

## HDR-C film provides optimum contrast throughout all density ranges

HDR-C Film has excellent diagnostic capabilities in dense breast tissue. This unique film delivers high contrast in the shoulder of the H&D curve and latitude for dense structures within the toe portion of the H&D curve.

## Split Emulsion Layering Technology Provides Superior Detail


HDR-C uses Split Emulsion Layering, the latest in emulsion technology. These layers provide the viewer with the latitude necessary to visualize clinically significant information in dense breast tissue and the high contrast necessary for perception of microcalcifications. The cool, blue image tone reduces eye fatigue and maximizes the ability to visualize small details and fine structures within the breast.

## Consistent Processing Results

Due to Cubic Plus crystal characteristics, the HDR-C sensitometric results are minimally affected by normal fluctuations in processing conditions. Low inherent film grain reduces visible noise on the final image. The anti-halation layer and protective topcoating minimize film handling and processing artifacts.

## HDR-C Technical Data

HDR-C is a single-sided, high contrast, green light sensitive, orthochromatic film for mammography. This film is designed for use with single sided, green-light emitting intensifying screens. HDR-C can be combined with the Agfa Mammography Screen that provides the image quality and system speed that best meets the user's needs.



## Technical Specifications

### GENERAL

#### Safelight requirements

- Use a safelight filter appropriate for Orthochromatic Film, such as a GBX-2
- Use a 7.5w frosted bulb placed at least 4 feet from the work surface.

#### Processing

- RP (Standard) processing cycles.

#### Chemistry

- Agfa Mammography chemistry is strongly recommended.
- Equivalent mammography chemistries may be used.
- **Developer temperature:** should be maintained between 33°C - 35°C

#### Storage and operating conditions


- **Temperature:** 4 - 25° C (39 - 77° F)
- **Humidity:** 30 - 50% Rh
- Shield film from sources of heat and all penetrating radiation
- Observe expiration dates

#### Agfa's HD System Speeds


- **HD screen:** 100 Speed Class
- **HD-S screen:** 150 Speed Class

### CURVES


#### Speed Dynamic Curve



#### Contrast Dynamic Curve




#### Sensitometric and contrast curve



Machine processing MAMORAY HT-300 with RP processing cycle. Agfa Mammography Chemistry. X-Rite Sensitometer.

Speed and Contrast are relative data; measurements are influenced by the sensitometer and densitometer used.

#### Darkroom Sensitivity



For more information on Agfa, please visit our website on [www.agfa.com](http://www.agfa.com) ■

Agfa and the Agfa rhombus are trademarks of Agfa-Gevaert NV, Belgium, or its affiliates. All rights reserved. All information contained herein is intended for guidance purposes only, and characteristics of the products and services described in this publication can be changed at any time without notice. Products and services may not be available for your local area. Please contact your local sales representative for availability information. Agfa-Gevaert NV diligently strives to provide as accurate information as possible, but shall not be responsible for any typographical error.

© 2018 Agfa NV  
All rights reserved  
Published by Agfa NV  
Septestraat 27 - 2640 Mortsel  
Belgium

GB 00201803