

# AGFA HEALTHCARE DICOM Conformance Statement

→ **Workstation NX 3.0.8950**

Document No. 001541  
Revision 3  
Livelink NodeID: 55300739

**When printed, this is NOT a controlled copy**

## Document Information

---

|  |   |   |
|--|---|---|
| <b>Service-related contact information worldwide</b> | All service-related contact information is available on this URL→ | <a href="http://www.agfahealthcare.com/global/en/main/contact/index.jsp">http://www.agfahealthcare.com/global/en/main/contact/index.jsp</a> |
|--|---|---|

Issued by:  
Agfa HealthCare  
V&V Connectivity  
Septestraat 27  
B-2640 Mortsel  
Belgium

tel: 32 3 444 7588  
email: [connectivity@agfa.com](mailto:connectivity@agfa.com)

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

**Copyright © 2016**  
**Agfa HealthCare**  
**All rights reserved**

## Conformance Statement Overview

NX 3.0.8950 implements the necessary DICOM services to facilitate receiving of unprocessed images from an AGFA CR/DR Modality, performing image processing and sending the processed images to a remote storage or print device over a Medical Imaging Systems network. NX 3.0.8950 may also acquire patient information from a Radiology Information System (RIS) for use in identifying processed images.

The table below provides an overview of the network services supported by NX 3.0.8950.

**Table 1.1-1: Network Services Supported**

| SOP Classes   | User of Service (SCU) | Provider of Service (SCP) |
|---|-----------------------|---------------------------|
| <b>Transfer</b>   |                       |                           |
| Computed Radiography Image Storage                      | Yes                   | No                        |
| Digital X-Ray Image Storage – For Presentation          | Option                | No                        |
| Digital X-Ray Image Storage – For Processing            | Option                | No                        |
| Grayscale Softcopy Presentation State Storage SOP Class | Yes                   | No                        |
| Digital Mammography Image Storage – For Presentation    | Option                | No                        |
| Digital Mammography Image Storage – For Processing      | Option                | No                        |
| X-Ray Radiation Dose SR Storage                         | Yes                   | No                        |
| <b>Workflow Management</b>                              |                       |                           |
| Storage Commitment Push Model SOP Class                 | Yes                   | No                        |
| Modality Performed Procedure Step SOP Class             | Option                | No                        |
| Modality Worklist Information Model – FIND              | Option                | No                        |
| <b>Print Management</b>                                 |                       |                           |
| Basic Grayscale Print Management Meta SOP Class         | Yes                   | No                        |
| Basic Film Session SOP Class                            | Yes                   | No                        |
| Basic Film Box SOP Class                                | Yes                   | No                        |
| Basic Grayscale Image Box SOP Class                     | Yes                   | No                        |
| Printer SOP Class                                       | Yes                   | No                        |
| Print Job SOP Class                                     | Yes                   | No                        |
| Presentation LUT SOP Class                              | Yes                   | No                        |

NX 3.0.8950 supports Media Services.

**Table 1.1-2: Media Services Supported**

| Media Storage Application Profile | Write Files (FSC or FSU) | Read Files (FSR) |
|-----------------------------------|--------------------------|------------------|
| Compact Disk - Recordable         |                          |                  |
| General Purpose CD-R/DVD          | Yes                      | No               |

## Table of Contents

|             |   |    |
|-------------|---|----|
| 1           | Introduction .....  | 7  |
| 1.1         | Revision Record .....   | 7  |
| 1.2         | Purpose and Intended Audience of this Document .....                                    | 7  |
| 1.3         | General Remarks .....   | 7  |
| 1.3.1       | Integration and Validation Activities .....   | 7  |
| 1.3.2       | Future Evolution .....  | 7  |
| 1.4         | Acronyms and Abbreviations .....  | 8  |
| 1.5         | Related Documents .....   | 8  |
| 2           | Networking .....  | 9  |
| 2.1         | Implementation Model .....  | 9  |
| 2.1.1       | Application Data Flow Diagram .....   | 9  |
| 2.1.2       | Functional Definitions of AE's .....  | 10 |
| 2.1.2.1     | Functional Definition of DICOM Store Application Entity .....                           | 10 |
| 2.1.2.2     | Functional Definition of DICOM Print Application Entity .....                           | 10 |
| 2.1.2.3     | Functional Definition of RIS Application Entity .....                                   | 10 |
| 2.1.2.3.1   | DICOM Modality Worklist .....   | 10 |
| 2.1.2.3.2   | DICOM Modality Worklist query on Accession Number .....                                 | 10 |
| 2.1.2.4     | Functional Definition of MPPS Application Entity .....                                  | 10 |
| 2.1.3       | Sequencing of Real World Activities .....   | 11 |
| 2.2         | AE Specifications .....   | 12 |
| 2.2.1       | DICOM Store Application Entity Specification .....                                      | 12 |
| 2.2.1.1     | SOP Classes Supported .....   | 12 |
| 2.2.1.2     | Association Establishment Policies .....  | 12 |
| 2.2.1.2.1   | General .....   | 12 |
| 2.2.1.2.2   | Number of Associations .....  | 12 |
| 2.2.1.2.3   | Asynchronous Nature .....   | 13 |
| 2.2.1.2.4   | Implementation Identifying Information .....  | 13 |
| 2.2.1.3     | Association Initiation Policies .....   | 13 |
| 2.2.1.3.1   | Activity – Send Images .....  | 13 |
| 2.2.1.3.1.1 | Description and Sequencing of Activity .....  | 13 |
| 2.2.1.3.1.2 | Proposed Presentation Contexts .....  | 16 |
| 2.2.1.3.1.3 | SOP Specific Conformance .....  | 17 |
| 2.2.1.4     | Association Acceptance Policies .....   | 28 |
| 2.2.1.4.1   | Receive Storage Commitment Response .....   | 28 |
| 2.2.1.4.1.1 | Description and Sequencing of Activity .....  | 28 |
| 2.2.1.4.1.2 | Accepted Presentation Contexts .....  | 28 |
| 2.2.1.4.1.3 | SOP Specific Conformance – Storage Commitment SOP Class<br>(1.2.840.10008.1.20.1) ..... | 29 |
| 2.2.1.4.1.4 | SOP Specific Conformance – Verification SOP Class<br>(1.2.840.10008.1.1) .....          | 30 |
| 2.2.2       | DICOM Print Application Entity Specification .....                                      | 31 |
| 2.2.2.1     | SOP Classes Supported .....   | 31 |
| 2.2.2.2     | Association Establishment Policies .....  | 31 |
| 2.2.2.2.1   | General .....   | 31 |
| 2.2.2.2.2   | Number of Associations .....  | 31 |
| 2.2.2.2.3   | Asynchronous Nature .....   | 31 |
| 2.2.2.2.4   | Implementation Identifying Information .....  | 32 |
| 2.2.2.3     | Association Initiation Policies .....   | 32 |
| 2.2.2.3.1   | Activity – Print Images .....   | 32 |

|             |   |    |
|-------------|---|----|
| 2.2.2.3.1.1 | Description and Sequencing of Activity.....   | 32 |
| 2.2.2.3.1.2 | Proposed Presentation Contexts.....   | 34 |
| 2.2.2.3.1.3 | SOP Specific Conformance - Basic Grayscale Print Management<br>Meta SOP Class .....     | 34 |
| 2.2.2.3.1.4 | SOP Specific Conformance Print Job SOP Class<br>(1.2.840.10008.5.1.1.14).....           | 39 |
| 2.2.2.3.1.5 | SOP Specific Conformance Presentation LUT SOP Class<br>(1.2.840.10008.5.1.1.23).....    | 41 |
| 2.2.2.3.1.6 | SOP Specific Conformance Verification SOP Class<br>(1.2.840.10008.1.1).....             | 41 |
| 2.2.3       | RIS Application Entity Specification.....   | 42 |
| 2.2.3.1     | SOP Classes Supported .....   | 42 |
| 2.2.3.2     | Association Establishment Policies .....  | 42 |
| 2.2.3.2.1   | General .....   | 42 |
| 2.2.3.2.2   | Number of Associations.....   | 42 |
| 2.2.3.2.3   | Asynchronous Nature .....   | 42 |
| 2.2.3.2.4   | Implementation Identifying Information.....   | 42 |
| 2.2.3.3     | Association Initiation Policies .....   | 43 |
| 2.2.3.3.1   | Activity – Query RIS.....   | 43 |
| 2.2.3.3.1.1 | Description and Sequencing of Activity.....   | 43 |
| 2.2.3.3.1.2 | Proposed Presentation Contexts.....   | 43 |
| 2.2.3.3.1.3 | SOP Specific Conformance – Modality Worklist SOP Class<br>(1.2.840.10008.5.1.4.31)..... | 43 |
| 2.2.3.4     | User Defined fields .....   | 47 |
| 2.2.3.5     | Custom Patient Field .....  | 47 |
| 2.2.4       | MPPS Application Entity Specification .....   | 48 |
| 2.2.4.1     | SOP Classes Supported .....   | 48 |
| 2.2.4.2     | Association Establishment Policies .....  | 48 |
| 2.2.4.2.1   | General .....   | 48 |
| 2.2.4.2.2   | Number of Associations.....   | 48 |
| 2.2.4.2.3   | Asynchronous Nature .....   | 48 |
| 2.2.4.2.4   | Implementation Identifying Information.....   | 48 |
| 2.2.4.3     | Association Initiation Policies .....   | 49 |
| 2.2.4.3.1   | Activity – Acquire Images.....  | 49 |
| 2.2.4.3.1.1 | Description and Sequencing of Activity.....   | 49 |
| 2.2.4.3.1.2 | Proposed Presentation Contexts.....   | 50 |
| 2.2.4.3.1.3 | SOP Specific Conformance – MPPS (1.2.840.10008.3.1.2.3.3).....                          | 50 |
| 2.3         | Network Interfaces.....   | 54 |
| 2.3.1       | Physical Medium Support .....   | 54 |
| 2.3.2       | Additional Protocols.....   | 54 |
| 2.4         | Configuration.....  | 55 |
| 2.4.1       | AE Title/ Presentation Address Mapping .....  | 55 |
| 2.4.1.1     | Local AE Titles .....   | 55 |
| 2.4.1.2     | Remote AE Title/ Presentation Address Mapping.....                                      | 55 |
| 2.4.1.2.1   | Storage.....  | 55 |
| 2.4.1.2.1.1 | Archive options .....   | 55 |
| 2.4.1.2.2   | RIS.....  | 58 |
| 2.4.1.2.2.1 | Configuration of a RIS.....   | 58 |
| 2.4.1.2.2.2 | Configuration of RIS Mapping .....  | 59 |
| 2.4.1.2.3   | MPPS.....   | 59 |
| 2.4.1.2.4   | Printing.....   | 59 |
| 2.4.2       | Queue Management & Job description .....  | 59 |
| 2.4.2.1     | DICOM Store AE .....  | 59 |
| 2.4.2.1.1   | Queue description.....  | 59 |
| 2.4.2.1.2   | Job description .....   | 60 |
| 2.4.2.2     | DICOM Print AE .....  | 60 |

|             |  |    |
|-------------|--|----|
| 2.4.3       | Parameters.....  | 60 |
| 3           | Media Interchange.....   | 63 |
| 3.1         | Implementation Model .....                                       | 63 |
| 3.1.1       | Application Data Flow Diagram.....                               | 63 |
| 3.1.2       | Functional Definition of AEs.....                                | 63 |
| 3.1.3       | Sequencing of Real World Activities.....                         | 63 |
| 3.1.4       | File Meta Information for Implementation Class and Version ..... | 64 |
| 3.2         | AE Specifications.....   | 64 |
| 3.2.1       | Offline-Media Application Entity Specification.....              | 64 |
| 3.2.1.1     | Real World Activities.....                                       | 64 |
| 3.2.1.1.1   | Real World Activity - Export to CD-R.....                        | 64 |
| 3.2.1.1.1.1 | Media Storage Application Profile.....                           | 65 |
| 3.3         | Media Configuration .....  | 66 |
| 4           | Support for Extended Character Sets.....                         | 67 |
| 5           | Security .....   | 68 |
| 5.1         | Security Profiles .....  | 68 |
| 5.2         | Association Level Security .....                                 | 68 |
| 5.3         | Application Level Security.....                                  | 68 |
| 6           | Annexes .....  | 69 |
| 6.1         | IOD Contents.....  | 69 |
| 6.1.1       | Created SOP Instances .....                                      | 69 |
| 6.1.1.1     | Common Modules .....   | 69 |
| 6.1.1.2     | CR.....  | 74 |
| 6.1.1.2.1   | CR Image IOD.....  | 74 |
| 6.1.1.2.2   | CR Modules.....  | 74 |
| 6.1.1.3     | DX.....  | 75 |
| 6.1.1.3.1   | DX Image IOD .....   | 75 |
| 6.1.1.3.2   | DX Modules.....  | 76 |
| 6.1.1.4     | MG .....   | 80 |
| 6.1.1.4.1   | MG Image IOD .....   | 80 |
| 6.1.1.4.2   | MG Modules .....   | 81 |
| 6.1.1.5     | GSPS.....  | 82 |
| 6.1.1.5.1   | GSPS IOD.....  | 82 |
| 6.1.1.5.2   | GSPS Modules.....  | 82 |
| 6.1.1.6     | X-Ray Radiation Dose SR Storage .....                            | 85 |
| 6.1.1.6.1   | X-Ray Radiation Dose SR IOD .....                                | 85 |
| 6.1.1.6.2   | X-Ray Radiation Dose SR Modules .....                            | 85 |
| 6.2         | Data dictionary of Private Attributes .....                      | 94 |
| 6.3         | Attribute Mapping .....  | 97 |
| 6.4         | Grayscale Image Consistency .....                                | 98 |
| 6.5         | Veterinary Extensions.....                                       | 99 |

# 1 INTRODUCTION

## 1.1 Revision Record

For detailed version history and version numbers, see Livelink.

| DICOM Conformance Statement NX 3.0.8950 |            |                                  |
|---|------------|----------------------------------|
| Version                                 | Date       | Reason for Change                |
| 3                                       | 2016-06-24 | Reviewed version for publication |

## 1.2 Purpose and Intended Audience of this Document

This document is a DICOM Conformance Statement for the DICOM Services of the *NX 3.0.8950* product. It is written according to part PS 3.2 of Digital Imaging and Communications in Medicine (DICOM) 3.0, NEMA PS 3.1-3.16, 2008.

The user of this document is involved with system integration and/or software design. We assume that the reader is familiar with the terminology and concepts that are used in the DICOM 3.0 standard and the IHE Technical Framework.

Readers not familiar with DICOM 3.0 terminology should first read the appropriate parts of the DICOM standard itself, prior to reading this conformance statement.

Although the use of this conformance statement in conjunction with the DICOM 3.0 standard is intended to facilitate communication with Agfa ADC imaging equipment, it is not sufficient to guarantee, by itself, the inter-operation of the connection.

## 1.3 General Remarks

### 1.3.1 Integration and Validation Activities

The integration of any device into a system of interconnected devices goes beyond the scope of the DICOM 3.0 standard and this conformance statement when interoperability is desired. The responsibility for analyzing the applications requirements and developing a solution that integrates the Agfa equipment with other vendors' systems is the user's responsibility and should not be underestimated.

In some circumstances it might be necessary to perform a validation to make sure that functional interoperability between the Agfa equipment and non-Agfa devices works as expected. The user should ensure that any non-Agfa provider accepts responsibility for any validation required for their connection with the Agfa equipment.

### 1.3.2 Future Evolution

As the DICOM 3.0 standard evolves to meet the user's growing requirements and to incorporate new features and technologies, Agfa will follow the evolution of the standard. This evolution of the standard may require changes to devices that have implemented DICOM 3.0. The user should ensure that any non-Agfa provider, who connects with Agfa devices, also plans for future evolution of the DICOM standard. A refusal to do so may result in the loss of functionality and/or connectivity between the different products.

## 1.4 Acronyms and Abbreviations

Definitions, terms and abbreviations used in this document are defined within the different parts of the DICOM standard. Abbreviations and terms are as follows:

|       |  |
|-------|--|
| AE    | DICOM Application Entity                       |
| AET   | Application Entity Title                       |
| ACSE  | Association Control Service Element            |
| CD-R  | Compact Disk Recordable                        |
| DICOM | Digital Imaging and Communications in Medicine |
| FSC   | File-Set Creator                               |
| FSU   | File-Set Updater                               |
| FSR   | File-Set Reader                                |
| GSDF  | Grayscale Standard Display Function            |
| GSPS  | Grayscale Softcopy Presentation State          |
| IE    | Information Entity                             |
| IOD   | (DICOM) Information Object Definition          |
| ISO   | International Standard Organization            |
| LUT   | Lookup Table                                   |
| MPPS  | Modality Performed Procedure Step              |
| MSPS  | Modality Scheduled Procedure Step              |
| PDU   | DICOM Protocol Data Unit                       |
| SCU   | DICOM Service Class User (DICOM client)        |
| SCP   | DICOM Service Class Provider (DICOM server)    |
| SOP   | DICOM Service-Object Pair                      |
| SR    | Structured Report                              |
| UID   | Unique Identifier                              |
| VR    | Value Representation                           |

## 1.5 Related Documents

- ACR-NEMA Digital Imaging and Communications in Medicine (DICOM) V3.0.2016b.
- IHE Radiology Technical Framework Revision 14.0 – Final Text, July 2015.



## 2 NETWORKING

### 2.1 Implementation Model

#### 2.1.1 Application Data Flow Diagram




Figure 2.1-1: Functional Overview – Application Data Flow

## 2.1.2 Functional Definitions of AE's

### 2.1.2.1 Functional Definition of DICOM Store Application Entity

The DICOM Store Application Entity sends images and Presentation States to a remote AE (this can be a configured Archiving station, a diagnostic workstation ...). It is associated with the local real-world activity "Send Images, GSPS & Dose SR". If Storage Commitment is configured for the archive destination, the DICOM Store AE will request Storage Commitment and a job will only be marked as successful when the commitment is successfully obtained. An image that has been successfully sent to an archive cannot be sent to that archive again.

### 2.1.2.2 Functional Definition of DICOM Print Application Entity

The DICOM Print Application Entity prints images on a remote AE (a configured Printer). It is associated with the local real-world activity "Request to print film sheets". A job will only be marked as "successful" when the sheet is successfully printed.

### 2.1.2.3 Functional Definition of RIS Application Entity

The RIS Application Entity receives Worklist information from a remote AE. It is associated with the local real-world activity "Query RIS". This can be triggered manually by clicking the "Query RIS" button or triggered automatically at a specified interval. When properly configured, the RIS can also be queried by Accession Number.

The RIS query can be configured in different ways (for the possible configurations, see § [2.4.1.2.2.1](#)). The two that are relevant for DICOM are described in the following paragraphs:

#### 2.1.2.3.1 DICOM Modality Worklist

The use of DICOM Modality Worklist (DMWL) queries the RIS by means of query keys. The query keys can be changed by the user in the configuration tool (see § [2.4.1.2.2.1](#)). Based on these keys, the worklist will be populated during the next RIS query.


#### 2.1.2.3.2 DICOM Modality Worklist query on Accession Number

By querying by Accession Number, only the worklist items that contain the specified Accession Number are returned. These results are parsed, mapped and shown as a Worklist.

### 2.1.2.4 Functional Definition of MPPS Application Entity

The MPPS Application Entity sends MPPS information to a remote AE when MPPS reporting is configured. The local real-world activities "Acquire Images", "Archive Images" and "Print Images" can trigger the creation of an MPPS instance. This can also be done by adding a new SOP Instance UID to an MPPS, by adding a Dose to an MPPS or by adding a Print Sheet to a MPPS. The local real-world activity "Close Session" or the modification of the session's Accession Number will complete the MPPS. This happens completely automatically and no user intervention is possible.

### 2.1.3 Sequencing of Real World Activities



**Figure 2.1-2: sequencing constraints**

Under normal scheduled workflow conditions the sequencing constraints illustrated in Figure 2.1-2 apply:

1. Query RIS
2. Receive Worklist
3. Select Work item
4. Create MPPS
5. Acquire Images
6. Print Acquired Images
7. Store Acquired Images (+ GSPS and Dose SR)
8. Commit Acquired Images
9. Finalize MPPS

## 2.2 AE Specifications

### 2.2.1 DICOM Store Application Entity Specification

#### 2.2.1.1 SOP Classes Supported

This Application Entity provides Standard Conformance to the following SOP Class(es):

**Table 2.2-1: SOP Class(es) for the DICOM Store Application Entity**

| SOP Class Name  | SOP Class UID                 | SCU | SCP |
|---|-------------------------------|-----|-----|
| Computed Radiography Image Storage                      | 1.2.840.10008.5.1.4.1.1.1     | Yes | No  |
| Digital X-ray Image Storage – for presentation          | 1.2.840.10008.5.1.4.1.1.1.1   | Yes | No  |
| Digital X-ray Image Storage – for processing            | 1.2.840.10008.5.1.4.1.1.1.1.1 | Yes | No  |
| Digital Mammography Image Storage – For Presentation    | 1.2.840.10008.5.1.4.1.1.1.2   | Yes | No  |
| Digital Mammography Image Storage – For Processing      | 1.2.840.10008.5.1.4.1.1.1.2.1 | Yes | No  |
| Grayscale Softcopy Presentation State Storage SOP Class | 1.2.840.10008.5.1.4.1.1.11.1  | Yes | No  |
| X-Ray Radiation Dose SR Storage                         | 1.2.840.10008.5.1.4.1.1.88.67 | Yes | No  |
| Storage Commitment Push Model SOP Class                 | 1.2.840.10008.1.20.1          | Yes | No  |
| Verification SOP Class                                  | 1.2.840.10008.1.1             | Yes | Yes |

#### 2.2.1.2 Association Establishment Policies

##### 2.2.1.2.1 General

The DICOM standard Application context is always proposed:

**Table 2.2-2: DICOM Application Context**

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

##### 2.2.1.2.2 Number of Associations

NX 3.0.8950 initiates one association at a time for each destination to which a transfer request is being processed in the active job queue list. Only one job per destination will be active at a time, the other remains pending until the active job for that destination is completed or failed. There can however be several simultaneous associations to different destinations.

**Table 2.2-3: Number of Associations as an Association Initiator for DICOM Store AE**

|   |                                       |
|---|---------------------------------------|
| Maximum number of simultaneous associations initiated | <b>1 per destination<br/>(32 max)</b> |
|---|---------------------------------------|

**Table 2.2-4: Number of Associations as an Association Acceptor for DICOM Store AE**

|  |  |
|--|--|
| Maximum number of simultaneous associations accepted | <b>1 for storage<br/>commit response</b> |
|--|--|

### 2.2.1.2.3 Asynchronous Nature

**Table 2.2-5: Asynchronous Nature as an Association Initiator for DICOM Store AE**

|   |   |
|---|---|
| Maximum number of outstanding asynchronous transactions | 1 |
|---|---|

The DICOM Store AE allows a single outstanding operation on any association. Therefore, it does not support asynchronous operations window negotiation, other than the default as specified by the DICOM specification.

### 2.2.1.2.4 Implementation Identifying Information

**Table 2.2-6: DICOM implementation Class and Version for DICOM Store AE**

|                             |                                   |
|-----------------------------|-----------------------------------|
| Implementation Class UID    | <b>1.3.51.0.1.3</b>               |
| Implementation Version Name | <b>AGFA DTF1.0.XX<sup>1</sup></b> |

## 2.2.1.3 Association Initiation Policies

### 2.2.1.3.1 Activity – Send Images

#### 2.2.1.3.1.1 Description and Sequencing of Activity

A user can select an image and request it to be sent to a destination. The request is forwarded to the job queue and then processed. An image can also be sent automatically when closing an exam.

---

**Note:**

An image can only be sent to an archiving destination if it hasn't been archived there yet.

---

The DICOM Store AE is then invoked by the queue that is responsible for processing network archival tasks for a specific destination. The DICOM Store AE will attempt to initiate a new Association in order to issue a C-STORE request. If the job contains multiple images, then multiple C-STORE requests will be issued over the same Association. The association will be closed when the last image (or GSPS) is sent.

If the association cannot be opened, the job is set to a retry state. If after a configurable number of retries the job still fails, it is set to an error state ("Failed"). It can then be restarted by the user through the job control interface. If three successive jobs to a same destination fail, the queue for that destination will be stalled. It will retry to process the job three (3) more times. If this fails, the queue for that destination will be stalled for a longer time, meaning that it will only retry jobs to that destination every five (5) minutes.

When the association was rejected by the device due to a configuration issue, the queue for that device will be stalled when three (3) successive jobs experience a device failure.

If the Remote AE is configured to support Storage Commit, the DICOM Store AE will send a Storage Commit request (N-ACTION) over a new association and will then wait for an N-EVENT-REPORT. If the N-EVENT-REPORT does not arrive within the waiting period, the AE closes the association and assumes that a separate association will be set up for the N-


---

<sup>1</sup> XX is the build version number.

EVENT-REPORT. If this report does not arrive within a (configurable) amount of time, the job will be marked as FAILED.

NX 3.0.8950 does not foresee additional logic when the system is shut down, regarding storage commitment time outs. In other words: when a storage commit reply does not reach NX 3.0.8950 successfully, because NX 3.0.8950 was/is down at that point in time, the job will time out and go to FAILED.

A possible sequence of iterations between the DICOM Store AE and an Image Manager is illustrated in Figure 2.2-1.



**Figure 2.2-1: Example of a Storage AE Sequencing Diagram (with Storage Commit)**

The DICOM Store AE may reject the association attempts as shown in the table below:

**Table 2.2-7: Association Reject Reasons**

| Result                                    | Response | Reason/Diag  |
|---|----------|--|
| The DICOM library refuses the association |          | SCU Device not known   |
| The DICOM library refuses the association |          | Maximum number of association processing threads exceeded  |
| The DICOM library refuses the association | A799     | DICOM Library is not allowed to accept C-STORE commands from this device and therefore refuses the association.                                      |
| The DICOM library refuses the association | A0102    | SOP Class is not found back in the configuration.  |
| The DICOM library refuses the association | A702     | The DICOM library is unable to create the Dicom media file due to the fact that the disk is full.  |
| The DICOM library refuses the association | A703     | The DICOM library is unable to create the DICOM media file due to resource problems other than 'disk is full'.                                       |
| The DICOM library refuses the association | A703     | The DICOM library is unable to acquire the complete C-STORE request due to network problems.   |
| The DICOM library refuses the association | A701     | The DICOM library is unable to acquire the complete C-STORE request due to network problems (time-out while reading data from socket) <sup>2</sup> . |

<sup>2</sup> The time-out value is hard-coded in the DICOM library as being 60 seconds.

### 2.2.1.3.1.2 Proposed Presentation Contexts

The DICOM Store AE is capable of proposing the Presentation Contexts shown in the following table:

**Table 2.2-8: Presentation Contexts Proposed by DICOM Store AE**

| Presentation Context Table                           |                               |   |  |      |                      |
|--|-------------------------------|---|--|------|----------------------|
| Abstract Syntax                                      |                               | Transfer Syntax   |  | Role | Extended Negotiation |
| Name   | UID                           | Name List   | UID List   |      |                      |
| CR Image Storage                                     | 1.2.840.10008.5.1.4.1.1.1     | JPEG LLNH1 Encoded (i.e. lossless)<br>JPEG LLNHF Encoded (i.e. lossless)<br>JPEG Lossy 8 bit<br>JPEG Lossy 12 bit<br>Implicit VR Little Endian<br>Explicit VR Little Endian | 1.2.840.10008.1.2.4.57<br>1.2.840.10008.1.2.4.70<br>1.2.840.10008.1.2.4.50<br>1.2.840.10008.1.2.4.51<br>1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCU  | None                 |
| Digital X-ray Image Storage – For Presentation       | 1.2.840.10008.5.1.4.1.1.1.1   | JPEG LLNH1 Encoded (i.e. lossless)<br>JPEG LLNHF Encoded (i.e. lossless)<br>JPEG Lossy 8 bit<br>JPEG Lossy 12 bit<br>Implicit VR Little Endian<br>Explicit VR Little Endian | 1.2.840.10008.1.2.4.57<br>1.2.840.10008.1.2.4.70<br>1.2.840.10008.1.2.4.50<br>1.2.840.10008.1.2.4.51<br>1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCU  | None                 |
| Digital X-ray Image Storage – For Processing         | 1.2.840.10008.5.1.4.1.1.1.1.1 | JPEG LLNH1 Encoded (i.e. lossless)<br>JPEG LLNHF Encoded (i.e. lossless)<br>JPEG Lossy 8 bit<br>JPEG Lossy 12 bit<br>Implicit VR Little Endian<br>Explicit VR Little Endian | 1.2.840.10008.1.2.4.57<br>1.2.840.10008.1.2.4.70<br>1.2.840.10008.1.2.4.50<br>1.2.840.10008.1.2.4.51<br>1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCU  | None                 |
| Digital Mammography Image Storage – For Presentation | 1.2.840.10008.5.1.4.1.1.1.2   | JPEG LLNH1 Encoded (i.e. lossless)<br>JPEG LLNHF Encoded (i.e. lossless)<br>JPEG Lossy 8 bit<br>JPEG Lossy 12 bit<br>Implicit VR Little Endian<br>Explicit VR Little Endian | 1.2.840.10008.1.2.4.57<br>1.2.840.10008.1.2.4.70<br>1.2.840.10008.1.2.4.50<br>1.2.840.10008.1.2.4.51<br>1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCU  | None                 |
| Digital Mammography Image Storage – For Processing   | 1.2.840.10008.5.1.4.1.1.1.2.1 | JPEG LLNH1 Encoded (i.e. lossless)<br>JPEG LLNHF Encoded (i.e. lossless)<br>JPEG Lossy 8 bit<br>JPEG Lossy 12 bit<br>Implicit VR Little Endian<br>Explicit VR Little Endian | 1.2.840.10008.1.2.4.57<br>1.2.840.10008.1.2.4.70<br>1.2.840.10008.1.2.4.50<br>1.2.840.10008.1.2.4.51<br>1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCU  | None                 |
| X-Ray Radiation Dose SR Storage                      | 1.2.840.10008.5.1.4.1.1.88.67 | Implicit VR Little Endian<br>Explicit VR Little Endian  | 1.2.840.10008.1.2<br>1.2.840.10008.1.2.1   | SCU  | None                 |



| Presentation Context Table                              |                              |  |  |      |                      |
|---|------------------------------|--|--|------|----------------------|
| Abstract Syntax   |                              | Transfer Syntax  |  | Role | Extended Negotiation |
| Name  | UID                          | Name List  | UID List                                 |      |                      |
| Grayscale Softcopy Presentation State Storage SOP Class | 1.2.840.10008.5.1.4.1.1.11.1 | Implicit VR Little Endian<br>Explicit VR Little Endian | 1.2.840.10008.1.2<br>1.2.840.10008.1.2.1 | SCU  | None                 |
| Storage Commitment Push Model SOP Class                 | 1.2.840.10008.1.20.1         | Implicit VR Little Endian                              | 1.2.840.10008.1.2                        | SCU  | None                 |
| Verification SOP Class                                  | 1.2.840.10008.1.1            | Implicit VR Little Endian                              | 1.2.840.10008.1.2                        | SCU  | None                 |

Depending on the configuration of NX 3.0.8950, either the CR Image Storage, the Digital X-ray Image Storage – for processing or the Digital X-ray Image Storage – for presentation, the Digital Mammography Image Storage – For Presentation or the Digital Mammography Image Storage – For Processing will be proposed.

The Grayscale Softcopy Presentation State Storage SOP Class and the Storage Commitment Push Model SOP Class will only be proposed if configured.

### 2.2.1.3.1.3 SOP Specific Conformance

#### 2.2.1.3.1.3.1 Image & Presentation State Storage

##### 2.2.1.3.1.3.1.1 Computed Radiography Image Storage SOP Class (1.2.840.10008.5.1.4.1.1.1)

The Computed Radiography Image Storage SOP class is a Storage Standard SOP Class that uses the CR IOD (§ [6.1.1.2.1](#)).

##### 2.2.1.3.1.3.1.2 Digital X-ray Image Storage – for presentation SOP Class (1.2.840.10008.5.1.4.1.1.1.1)

The Digital X-Ray Image Storage - For Presentation SOP Class uses the DX IOD (§[6.1.1.3.1](#)) with an Enumerated Value of FOR PRESENTATION for Presentation Intent Type (0008, 0068).

##### 2.2.1.3.1.3.1.3 Digital X-ray Image Storage – for processing SOP Class (1.2.840.10008.5.1.4.1.1.1.1)

The Digital X-Ray Image Storage - For Processing SOP Class uses the DX IOD (§[6.1.1.3.1](#)) with an Enumerated Value of FOR PROCESSING for Presentation Intent Type (0008, 0068).

As a SCU of the Digital X-Ray Image Storage - For Processing SOP Class, it also supports the Digital X-Ray Image Storage - For Presentation SOP Class.

##### 2.2.1.3.1.3.1.4 Digital Mammography Image Storage – for presentation SOP Class (1.2.840.10008.5.1.4.1.1.1.2)

The Digital Mammography Image Storage - For Presentation SOP Class uses the MG IOD (§[6.1.1.4.1](#)) with an Enumerated Value of FOR PRESENTATION for Presentation Intent Type (0008, 0068).

### 2.2.1.3.1.3.1.5 Digital Mammography Image Storage – for processing SOP Class (1.2.840.10008.5.1.4.1.1.1.2.1)

The Digital Mammography Image Storage - For Processing SOP Class uses the MG IOD (§6.1.1.4.1) with an Enumerated Value of FOR PROCESSING for Presentation Intent Type (0008, 0068).

### 2.2.1.3.1.3.1.6 Grayscale Softcopy Presentation State Storage SOP Class (1.2.840.10008.5.1.4.1.1.11.1)

The Grayscale Softcopy Presentation State Storage SOP Class extends the functionality of the Storage Service class to add the ability to convey an intended presentation state or record an existing presentation state.

It includes capabilities for specifying:

- the output grayscale space in P-Values
- grayscale contrast transformations including modality and VOI LUT
- selection of the area of the image to display
- image and display relative annotations, including graphics, text and overlays

#### Note:

A GSPS always refers to exactly one (1) image. Since re-sending an image is prohibited in NX 3.0.8950, an archived image can also have maximally one (1) GSPS that refers to it.

Depending on the configured SOP class to be used (CR (§ 2.2.1.3.1.3.1.1), DX for presentation (§2.2.1.3.1.3.1.2) or MG for presentation (§ 2.2.1.3.1.3.1.4) annotations and shutters are stored in the GSPS or burned in the image as described in the following table:

**Table 2.2-9: Use of GSPS vs. Burning in the image**

|  | CR<br>DX for presentation<br>MG for presentation  | DX for processing |
|--|---|-------------------|
| <b>Image</b><br>(processing,W/L,collimation) | Processed pixels<br>+<br>LUTs   | RAW pixels        |
| <b>Annotations</b>                           | GSPS or Burned in   | Ignored           |
| <b>Shutters</b>                              | GSPS or<br>when GSPS is not used:<br>In case of CR the shutter is<br>burned into the Pixel data.<br>In case of DX and MG the<br>shutter is described in the<br>shutter module of the image<br>header. |                   |
| <b>Zoom</b>                                  | GSPS or discarded   |                   |
| <b>Markers</b>                               | Always burned in  |                   |

The following paragraphs describe into detail how each annotation is embedded in the GSPS:

**2.2.1.3.1.3.1.6.1 Text**

**Text text text text**

**Graphic Layer Module**

| (0070,0060) | Graphic Layer Sequence     |       |
|-------------|----------------------------|-------|
| (0070,0002) | >Graphic Layer             | ID123 |
| (0070,0068) | >Graphic Layer Description | Text  |

**Graphic Annotation Module**

| (0070,0001) | Graphic Annotation Sequence                   |         |
|-------------|---|---------|
| (0070,0002) | > Graphic Layer                               | ID123   |
| (0070,0008) | > Text Object Sequence                        |         |
| (0070,0003) | >> Bounding box annotation units              | PIXEL   |
| (0070,0004) | >> Anchor Point Annotation Units              | -       |
| (0070,0006) | >> Unformatted Text Value                     | [value] |
| (0070,0010) | >> Bounding Box Top Left Hand Corner          | [value] |
| (0070,0011) | >> Bounding Box Top Right Hand Corner         | [value] |
| (0070,0012) | >> Bounding Box Text Horizontal Justification | LEFT    |
| (0070,0014) | >> Anchor Point                               | -       |
| (0070,0015) | >> Anchor Point Visibility                    | -       |

**2.2.1.3.1.3.1.6.2 Arrow**

A text is foreseen for each arrow, but it may be an empty string. The arrow consists of two polylines.

**Graphic Layer Module**

| (0070,0060) | Graphic Layer Sequence     |       |
|-------------|----------------------------|-------|
| (0070,0002) | >Graphic Layer             | ID    |
| (0070,0068) | >Graphic Layer Description | Arrow |

**Graphic Annotation Module**

| (0070,0001)        | Graphic Annotation Sequence                   |                                 |
|--------------------|---|---------------------------------|
| (0070,0002)        | > Graphic Layer                               | [ID of the corresponding layer] |
| (0070,0008)        | > Text Object Sequence                        |                                 |
| (0070,0003)        | >> Bounding box annotation units              | PIXEL                           |
| (0070,0004)        | >> Anchor Point Annotation Units              | PIXEL                           |
| (0070,0006)        | >> Unformatted Text Value                     | [value]                         |
| (0070,0010)        | >> Bounding Box Top Left Hand Corner          | [value]                         |
| (0070,0011)        | >> Bounding Box Top Right Hand Corner         | [value]                         |
| (0070,0012)        | >> Bounding Box Text Horizontal Justification | LEFT                            |
| (0070,0014)        | >> Anchor Point                               | [value of arrow point]          |
| (0070,0015)        | >> Anchor Point Visibility                    | N                               |
| (0070,0009)        | > Graphic Object Sequence                     |                                 |
| <i>Arrow line</i>  |   |                                 |
| (0070,0005)        | >> Graphic Annotation Units                   | PIXEL                           |
| (0070,0020)        | >> Graphic Dimensions                         | 2                               |
| (0070,0021)        | >> Number of Graphic Points                   | 2                               |
| (0070,0022)        | >> Graphic Data                               | [values]                        |
| (0070,0023)        | >> Graphic Type                               | POLYLINE                        |
| (0070,0024)        | >> Graphic Filled                             | -                               |
| <i>Arrow point</i> |   |                                 |
| (0070,0005)        | >> Graphic Annotation Units                   | PIXEL                           |
| (0070,0020)        | >> Graphic Dimensions                         | 2                               |
| (0070,0021)        | >> Number of Graphic Points                   | 3                               |
| (0070,0022)        | >> Graphic Data                               | [values]                        |
| (0070,0023)        | >> Graphic Type                               | POLYLINE                        |
| (0070,0024)        | >> Graphic Filled                             | -                               |

**2.2.1.3.1.3.1.6.3 Rectangle**

A rectangle is always provided with a corresponding text (for measurements). If the text is moved by the operator, an extra text object is added to the GSPS.

**Graphic Layer Module**

| (0070,0060) | Graphic Layer Sequence     |                            |
|-------------|----------------------------|----------------------------|
| (0070,0002) | >Graphic Layer             | ID                         |
| (0070,0068) | >Graphic Layer Description | Rectangle or Rectangle_SAL |

**Graphic Annotation Module**

| (0070,0001) | Graphic Annotation Sequence                   |                                 |
|-------------|---|---------------------------------|
| (0070,0002) | > Graphic Layer                               | [ID of the corresponding layer] |
| (0070,0008) | > Text Object Sequence                        |                                 |
| (0070,0003) | >> Bounding box annotation units              | PIXEL                           |
| (0070,0004) | >> Anchor Point Annotation Units              | -                               |
| (0070,0006) | >> Unformatted Text Value                     | [value]                         |
| (0070,0010) | >> Bounding Box Top Left Hand Corner          | [value]                         |
| (0070,0011) | >> Bounding Box Top Right Hand Corner         | [value]                         |
| (0070,0012) | >> Bounding Box Text Horizontal Justification | LEFT                            |
| (0070,0014) | >> Anchor Point                               | -                               |
| (0070,0015) | >> Anchor Point Visibility                    | -                               |
| (0070,0009) | > Graphic Object Sequence                     |                                 |
| (0070,0005) | >> Graphic Annotation Units                   | PIXEL                           |
| (0070,0020) | >> Graphic Dimensions                         | 2                               |
| (0070,0021) | >> Number of Graphic Points                   | 5                               |
| (0070,0022) | >> Graphic Data                               | [values]                        |
| (0070,0023) | >> Graphic Type                               | POLYLINE                        |
| (0070,0024) | >> Graphic Filled                             | N                               |

**2.2.1.3.1.3.1.6.4 Circle**

A circle is always provided with a corresponding text (for measurements). If the text is moved by the operator, an extra text object is added to the GSPS.

**Graphic Layer Module**

| (0070,0060) | Graphic Layer Sequence     |        |
|-------------|----------------------------|--------|
| (0070,0002) | >Graphic Layer             | ID     |
| (0070,0068) | >Graphic Layer Description | Circle |

**Graphic Annotation Module**

| (0070,0001) | Graphic Annotation Sequence                   |                                 |
|-------------|---|---------------------------------|
| (0070,0002) | > Graphic Layer                               | [ID of the corresponding layer] |
| (0070,0008) | > Text Object Sequence                        |                                 |
| (0070,0003) | >> Bounding box annotation units              | PIXEL                           |
| (0070,0004) | >> Anchor Point Annotation Units              | -                               |
| (0070,0006) | >> Unformatted Text Value                     | [value]                         |
| (0070,0010) | >> Bounding Box Top Left Hand Corner          | [value]                         |
| (0070,0011) | >> Bounding Box Top Right Hand Corner         | [value]                         |
| (0070,0012) | >> Bounding Box Text Horizontal Justification | LEFT                            |
| (0070,0014) | >> Anchor Point                               | -                               |
| (0070,0015) | >> Anchor Point Visibility                    | -                               |
| (0070,0009) | > Graphic Object Sequence                     |                                 |
| (0070,0005) | >> Graphic Annotation Units                   | PIXEL                           |
| (0070,0020) | >> Graphic Dimensions                         | 2                               |
| (0070,0021) | >> Number of Graphic Points                   | 2                               |
| (0070,0022) | >> Graphic Data                               | [values]                        |
| (0070,0023) | >> Graphic Type                               | CIRCLE                          |
| (0070,0024) | >> Graphic Filled                             | N                               |

**2.2.1.3.1.3.1.6.5 Polygon**

A polygon is always provided with a corresponding text (for measurements). If the text is moved by the operator, an extra text object is added to the GSPS.

**Graphic Layer Module**

| (0070,0060) | Graphic Layer Sequence     |         |
|-------------|----------------------------|---------|
| (0070,0002) | >Graphic Layer             | ID      |
| (0070,0068) | >Graphic Layer Description | Polygon |

**Graphic Annotation Module**

| (0070,0001) | Graphic Annotation Sequence                   |                                 |
|-------------|---|---------------------------------|
| (0070,0002) | > Graphic Layer                               | [ID of the corresponding layer] |
| (0070,0008) | > Text Object Sequence                        |                                 |
| (0070,0003) | >> Bounding box annotation units              | PIXEL                           |
| (0070,0004) | >> Anchor Point Annotation Units              | -                               |
| (0070,0006) | >> Unformatted Text Value                     | [value]                         |
| (0070,0010) | >> Bounding Box Top Left Hand Corner          | [value]                         |
| (0070,0011) | >> Bounding Box Top Right Hand Corner         | [value]                         |
| (0070,0012) | >> Bounding Box Text Horizontal Justification | LEFT                            |
| (0070,0014) | >> Anchor Point                               | -                               |
| (0070,0015) | >> Anchor Point Visibility                    | -                               |
| (0070,0009) | > Graphic Object Sequence                     |                                 |
| (0070,0005) | >> Graphic Annotation Units                   | PIXEL                           |
| (0070,0020) | >> Graphic Dimensions                         | 2                               |
| (0070,0021) | >> Number of Graphic Points                   | [value]                         |
| (0070,0022) | >> Graphic Data                               | [values]                        |
| (0070,0023) | >> Graphic Type                               | POLYLINE                        |
| (0070,0024) | >> Graphic Filled                             | N                               |

**2.2.1.3.1.3.1.6.6 Freehand****Graphic Layer Module**


| (0070,0060) | Graphic Layer Sequence     |          |
|-------------|----------------------------|----------|
| (0070,0002) | >Graphic Layer             | ID       |
| (0070,0068) | >Graphic Layer Description | Freehand |

**Graphic Annotation Module**

| (0070,0001) | Graphic Annotation Sequence                   |                                 |
|-------------|---|---------------------------------|
| (0070,0002) | > Graphic Layer                               | [ID of the corresponding layer] |
| (0070,0008) | > Text Object Sequence                        |                                 |
| (0070,0003) | >> Bounding box annotation units              | PIXEL                           |
| (0070,0004) | >> Anchor Point Annotation Units              | -                               |
| (0070,0006) | >> Unformatted Text Value                     | [value]                         |
| (0070,0010) | >> Bounding Box Top Left Hand Corner          | [value]                         |
| (0070,0011) | >> Bounding Box Top Right Hand Corner         | [value]                         |
| (0070,0012) | >> Bounding Box Text Horizontal Justification | LEFT                            |
| (0070,0014) | >> Anchor Point                               | -                               |
| (0070,0015) | >> Anchor Point Visibility                    | -                               |
| (0070,0009) | > Graphic Object Sequence                     |                                 |
| (0070,0005) | >> Graphic Annotation Units                   | PIXEL                           |
| (0070,0020) | >> Graphic Dimensions                         | 2                               |
| (0070,0021) | >> Number of Graphic Points                   | [value]                         |
| (0070,0022) | >> Graphic Data                               | [values]                        |
| (0070,0023) | >> Graphic Type                               | INTERPOLATED                    |
| (0070,0024) | >> Graphic Filled                             | N                               |

**2.2.1.3.1.3.1.6.7 Line**

A line consists of two parts: the actual line, and a small line indicating the middle point of the line:

**Graphic Layer Module**


| (0070,0060) | Graphic Layer Sequence     |      |
|-------------|----------------------------|------|
| (0070,0002) | >Graphic Layer             | ID   |
| (0070,0068) | >Graphic Layer Description | Line |

**Graphic Annotation Module**

| (0070,0001)                                  | Graphic Annotation Sequence |                                 |
|--|-----------------------------|---------------------------------|
| (0070,0002)                                  | > Graphic Layer             | [ID of the corresponding layer] |
| (0070,0009)                                  | > Graphic Object Sequence   |                                 |
| <i>Item for the actual line</i>              |                             |                                 |
| (0070,0005)                                  | >> Graphic Annotation Units | PIXEL                           |
| (0070,0020)                                  | >> Graphic Dimensions       | 2                               |
| (0070,0021)                                  | >> Number of Graphic Points | 2                               |
| (0070,0022)                                  | >> Graphic Data             | [values]                        |
| (0070,0023)                                  | >> Graphic Type             | POLYLINE                        |
| (0070,0024)                                  | >> Graphic Filled           | -                               |
| <i>Item for the indication of the middle</i> |                             |                                 |
| (0070,0005)                                  | >> Graphic Annotation Units | PIXEL                           |
| (0070,0020)                                  | >> Graphic Dimensions       | 2                               |
| (0070,0021)                                  | >> Number of Graphic Points | 2                               |
| (0070,0022)                                  | >> Graphic Data             | [values]                        |
| (0070,0023)                                  | >> Graphic Type             | POLYLINE                        |
| (0070,0024)                                  | >> Graphic Filled           | -                               |

**2.2.1.3.1.3.1.6.8 Perpendicular**

A perpendicular also consists of two lines A and B, as illustrated in the example below:

**Graphic Layer Module**

| (0070,0060) | Graphic Layer Sequence     |               |
|-------------|----------------------------|---------------|
| (0070,0002) | >Graphic Layer             | ID            |
| (0070,0068) | >Graphic Layer Description | Perpendicular |

**Graphic Annotation Module**

| (0070,0001)                        | Graphic Annotation Sequence |                                 |
|------------------------------------|-----------------------------|---------------------------------|
| (0070,0002)                        | > Graphic Layer             | [ID of the corresponding layer] |
| (0070,0009)                        | > Graphic Object Sequence   |                                 |
| <i>Item for the first line (A)</i> |                             |                                 |
| (0070,0005)                        | >> Graphic Annotation Units | PIXEL                           |
| (0070,0020)                        | >> Graphic Dimensions       | 2                               |
| (0070,0021)                        | >> Number of Graphic Points | 2                               |
| (0070,0022)                        | >> Graphic Data             | [values]                        |
| (0070,0023)                        | >> Graphic Type             | POLYLINE                        |
| (0070,0024)                        | >> Graphic Filled           | -                               |

| (0070,0001)                         | Graphic Annotation Sequence |          |
|-------------------------------------|-----------------------------|----------|
| <i>Item for the second line (B)</i> |                             |          |
| (0070,0005)                         | >> Graphic Annotation Units | PIXEL    |
| (0070,0020)                         | >> Graphic Dimensions       | 2        |
| (0070,0021)                         | >> Number of Graphic Points | 2        |
| (0070,0022)                         | >> Graphic Data             | [values] |
| (0070,0023)                         | >> Graphic Type             | POLYLINE |
| (0070,0024)                         | >> Graphic Filled           | -        |

#### 2.2.1.3.1.3.1.6.9 Calibration ruler

The calibration ruler is displayed at the side of the image and is the result of a calibration action of the operator.

It consists out of a ruler (constructed using several lines) and some text.

#### Graphic Layer Module


| (0070,0060) | Graphic Layer Sequence     |                   |
|-------------|----------------------------|-------------------|
| (0070,0002) | >Graphic Layer             | ID                |
| (0070,0068) | >Graphic Layer Description | Calibration Ruler |

#### Graphic Annotation Module

| (0070,0001)             | Graphic Annotation Sequence                   |                                 |
|-------------------------|---|---------------------------------|
| (0070,0002)             | > Graphic Layer                               | [ID of the corresponding layer] |
| (0070,0008)             | > Text Object Sequence                        |                                 |
| (0070,0003)             | >> Bounding box annotation units              | PIXEL                           |
| (0070,0004)             | >> Anchor Point Annotation Units              | -                               |
| (0070,0006)             | >> Unformatted Text Value                     | [value]                         |
| (0070,0010)             | >> Bounding Box Top Left Hand Corner          | [value]                         |
| (0070,0011)             | >> Bounding Box Top Right Hand Corner         | [value]                         |
| (0070,0012)             | >> Bounding Box Text Horizontal Justification | LEFT                            |
| (0070,0014)             | >> Anchor Point                               | -                               |
| (0070,0015)             | >> Anchor Point Visibility                    | -                               |
| (0070,0009)             | > Graphic Object Sequence                     |                                 |
| <i>First ruler line</i> |   |                                 |
| (0070,0005)             | >> Graphic Annotation Units                   | PIXEL                           |
| (0070,0020)             | >> Graphic Dimensions                         | 2                               |
| (0070,0021)             | >> Number of Graphic Points                   | 2                               |
| (0070,0022)             | >> Graphic Data                               | [values]                        |
| (0070,0023)             | >> Graphic Type                               | CIRCLE                          |
| (0070,0024)             | >> Graphic Filled                             | N                               |
| ...                     |   |                                 |

#### 2.2.1.3.1.3.1.6.10 Angle

An angle measurement is a combination of 2 lines and 2 short lines, with 2 angle texts.

**Graphic Layer Module**

| (0070,0060) | Graphic Layer Sequence     |       |
|-------------|----------------------------|-------|
| (0070,0002) | >Graphic Layer             | ID    |
| (0070,0068) | >Graphic Layer Description | Angle |

**Graphic Annotation Module**

| (0070,0001) | Graphic Annotation Sequence |                                 |
|-------------|-----------------------------|---------------------------------|
| (0070,0002) | > Graphic Layer             | [ID of the corresponding layer] |
| (0070,0008) | > Text Object Sequence      | -> 2 texts for the degrees      |
| (0070,0009) | > Graphic Object Sequence   | Long line 1                     |
|             |                             | Long line 2                     |
|             |                             | Small middle line 1             |
|             |                             | Small middle line 2             |

**2.2.1.3.1.3.1.6.11 Distance**

A distance is composed of a line with a text value.

**Graphic Layer Module**


| (0070,0060) | Graphic Layer Sequence     |          |
|-------------|----------------------------|----------|
| (0070,0002) | >Graphic Layer             | ID       |
| (0070,0068) | >Graphic Layer Description | Distance |

**Graphic Annotation Module**

| (0070,0001)                                  | Graphic Annotation Sequence |                                 |
|--|-----------------------------|---------------------------------|
| (0070,0002)                                  | > Graphic Layer             | [ID of the corresponding layer] |
| (0070,0008)                                  | > Text Object Sequence      | -> used for the distance text   |
| (0070,0009)                                  | > Graphic Object Sequence   |                                 |
| Item for the actual line                     |                             |                                 |
| (0070,0005)                                  | >> Graphic Annotation Units | PIXEL                           |
| (0070,0020)                                  | >> Graphic Dimensions       | 2                               |
| (0070,0021)                                  | >> Number of Graphic Points | 2                               |
| (0070,0022)                                  | >> Graphic Data             | [values]                        |
| (0070,0023)                                  | >> Graphic Type             | POLYLINE                        |
| (0070,0024)                                  | >> Graphic Filled           | -                               |
| Items for the ruler (cfr. Calibration ruler) |                             |                                 |
| (0070,0005)                                  | >> Graphic Annotation Units | PIXEL                           |
| (0070,0020)                                  | >> Graphic Dimensions       | 2                               |
| (0070,0021)                                  | >> Number of Graphic Points | 2                               |
| (0070,0022)                                  | >> Graphic Data             | [values]                        |
| (0070,0023)                                  | >> Graphic Type             | POLYLINE                        |
| (0070,0024)                                  | >> Graphic Filled           | -                               |



**2.2.1.3.1.3.1.6.12 Leg Length Difference**




**Graphic Layer Module**

| (0070,0060) | Graphic Layer Sequence     |                       |
|-------------|----------------------------|-----------------------|
| (0070,0002) | >Graphic Layer             | ID                    |
| (0070,0068) | >Graphic Layer Description | Leg Length Difference |

**Graphic Annotation Module**

| (0070,0001) | Graphic Annotation Sequence |                                 |
|-------------|-----------------------------|---------------------------------|
| (0070,0002) | > Graphic Layer             | [ID of the corresponding layer] |
| (0070,0008) | > Text Object Sequence      | Measurement text                |
| (0070,0009) | > Graphic Object Sequence   | First horizontal line           |
|             |                             | Second horizontal line          |
|             |                             | Vertical line                   |
|             |                             | Arrow point 1                   |
|             |                             | Arrow point 2                   |
|             |                             | Arrow line                      |

**2.2.1.3.1.3.1.6.13 Scoliosis**



**Graphic Layer Module**

| (0070,0060) | Graphic Layer Sequence     |                       |
|-------------|----------------------------|-----------------------|
| (0070,0002) | >Graphic Layer             | ID                    |
| (0070,0068) | >Graphic Layer Description | Scoliosis Measurement |

**Graphic Annotation Module**

| (0070,0001) | Graphic Annotation Sequence |                                 |
|-------------|-----------------------------|---------------------------------|
| (0070,0002) | > Graphic Layer             | [ID of the corresponding layer] |
| (0070,0008) | > Text Object Sequence      | Measurement text                |
| (0070,0009) | > Graphic Object Sequence   | Upper line                      |
|             |                             | Upper line perpendicular        |
|             |                             | Lower line                      |
|             |                             | Lower line perpendicular        |

**2.2.1.3.1.3.1.6.14 Display Shutter****Display Shutter Module**

|             |                               |             |
|-------------|-------------------------------|-------------|
| (0018,1600) | Shutter Shape                 | RECTANGULAR |
| (0018,1602) | Shutter Left Vertical Edge    | [value]     |
| (0018,1604) | Shutter Right Vertical Edge   | [value]     |
| (0018,1606) | Shutter Upper Horizontal Edge | [value]     |
| (0018,1608) | Shutter Lower Horizontal Edge | [value]     |

The manner in which the display area occluded by the shutter is neutralized (black-out, gray, or other means) is defined by the Attribute Shutter Presentation Value (0018,1622). This attribute present in the Presentation LUT module is mandatory when a display shutter is present in the GSPS. The value can go from 0000xH (black) to FFFFxH (white).

**2.2.1.3.1.3.1.6.15 Flipping/ rotating/ zooming**

Flipping and rotating is always done on pixel level. The relevant DICOM tags can be found in the following table:

**Table 2.2-10: DICOM spatial transformation module attributes.**

| attribute             | tag         | Actions  |   |                         |
|-----------------------|-------------|--|---|-------------------------|
|                       |             | Rotate 90° clockwise   | Rotate 90° counter clockwise  | Flip horizontal         |
| Image Rotation        | (0070,0042) | If image is not flipped : increased with 90 (modulo 360)<br>Otherwise : decreased with 90 (modulo 360) | If image is not flipped: decreased with 90 (modulo 360)<br>Otherwise : increased with 90 (modulo 360) |                         |
| Image horizontal flip | (0070,0041) |  |   | inverted (true ↔ false) |

Zoom/Pan is integrated in the Displayed Area module of the GSPS.

**2.2.1.3.1.3.2 Storage Commitment Push Model SOP Class (1.2.840.10008.1.20.1)**

When the Storage Commitment Push Model has been configured, the DICOM Store AE will request storage commitment for instances of the Image Storage SOP Class and Grayscale Softcopy Presentation State Storage SOP Class with each successfully completed “sent” job. The DICOM Store AE transmits the SOP Instances to the Remote AE. The request for storage commitment is transmitted to the Remote AE together with a list of references to one or more SOP Instances. If the Provider accepts the Storage Commitment with Success Status, the generated Transaction UID, together with study identification data and a time-stamp, is kept. Success or failure of storage commitment is subsequently indicated by a notification from the Remote AE to NX 3.0.8950.

The DIMSE-N Services applicable to the Storage Commitment Push Model SOP Class are:

- N-EVENT\_REPORT
- N-ACTION

The Storage Commitment Request operation allows a DICOM Store AE to request an SCP to commit to the safekeeping of a set of SOP Instances as described above. This operation is invoked through the N-ACTION primitive.

The N-ACTION is invoked by NX 3.0.8950 and is sent by creating a new association.

The behavior of Storage AE when encountering status codes in an N-ACTION response is summarized in the Table below:

**Table 2.2-11: Storage Commitment N-ACTION Information**

| Action Type Name            | Action Type ID | Attribute                    | Tag         | Requirement Type SCU |
|-----------------------------|----------------|------------------------------|-------------|----------------------|
| Request Storage Commitment  | 1              | Transaction UID              | (0008,1195) | 1                    |
|                             |                | Storage Media File-Set ID    | (0088,0130) | 3                    |
|                             |                | Storage Media File-Set UID   | (0088,0140) | 3                    |
|                             |                | Referenced SOP Sequence      | (0008,1199) | 1                    |
|                             |                | >Referenced SOP Class UID    | (0008,1150) | 1                    |
|                             |                | >Referenced SOP Instance UID | (0008,1155) | 1                    |
|                             |                | >Storage Media File-Set ID   | (0088,0130) | 3                    |
| >Storage Media File-Set UID | (0088,0140)    | 3                            |             |                      |

**Table 2.2-12: Storage Commitment N-ACTION Response Status Handling Behavior**

| Service Status | Further Meaning | Error Code | Behavior                |
|----------------|-----------------|------------|-------------------------|
| Success        |                 | 0000       | successful Notification |
| Warning        |                 |            |                         |
| Error          |                 |            |                         |

The behavior of the AE during communication failure is summarized in a table as follows:

**Table 2.2-13: DICOM Command Communication Failure Behavior**

| Exception           | Behavior   |
|---------------------|--|
| Timeout             | e.g. The Association is aborted using A-ABORT and command marked as failed. The reason is logged and reported to the user. |
| Association aborted | e.g. The command is marked as failed. The reason is logged and reported to the user.                                       |

The DICOM Store AE does not wait for an N-EVENT-REPORT. It closes the association as soon as it receives the N-ACTION-RP from the Remote AE. So, NX does not support N-EVENT-REPORT within the same association as the N-ACTION.

### 2.2.1.3.1.3.3 Verification SOP Class (1.2.840.10008.1.1)

The Storage AE provides standard conformance to the Verification SOP Class as an SCU. This verification is accomplished on an established Association using the C-ECHO DIMSE-C service.

These tests can be executed in the "Workstation Service & Configuration Tool".

The Configuration Tool opens an association when testing of a remote application is requested during a configuration session. This can be done when entering new data for remote application configuration or to verify existing configuration data using the C-ECHO DIMSE-C service.

## 2.2.1.4 Association Acceptance Policies

### 2.2.1.4.1 Receive Storage Commitment Response

#### 2.2.1.4.1.1 Description and Sequencing of Activity

Each Storage Commitment Request that NX sends, is uniquely identified by the Transaction UID Attribute (0008,1195) value that is generated by NX. After sending a Storage Commitment Request, NX expects an N-EVENT-REPORT from the SCP. NX will then respond with an N-EVENT-REPORT response primitive with a status code.

The NX DICOM Store AE will accept associations in order to receive responses to a Storage Commitment Request.




Figure 2.2-1: Sequencing of Receive Storage Commitment response

#### 2.2.1.4.1.2 Accepted Presentation Contexts

Table 2.2-14: Acceptable Presentation Contexts for Receive Storage Commitment Response

| Presentation Context Table    |                      |                           |                     |      |                      |
|-------------------------------|----------------------|---------------------------|---------------------|------|----------------------|
| Abstract Syntax               |                      | Transfer Syntax           |                     | Role | Extended Negotiation |
| Name                          | UID                  | Name List                 | UID List            |      |                      |
| Storage Commitment Push Model | 1.2.840.10008.1.20.1 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None                 |
|                               |                      | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |                      |
| Verification                  | 1.2.840.10008.1.1    | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCP  | None                 |
|                               |                      | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |                      |

### 2.2.1.4.1.3 SOP Specific Conformance – Storage Commitment SOP Class (1.2.840.10008.1.20.1)

The AE will consider Storage Commitment FAILED if no N-EVENT-REPORT is received for a Transaction UID within a configurable time period.

**Table 2.2-15: Storage Commitment N-EVENT-REPORT expected Information**

| Action Type Name                                   | Action Type ID | Attribute                    | Tag         | Requirement SCU | Type |
|--|----------------|------------------------------|-------------|-----------------|------|
| Storage Commitment Request Successful              | 1              | Transaction UID              | (0008,1195) | 1               |      |
|  |                | Storage Media File-Set ID    | (0088,0130) | 3               |      |
|  |                | Storage Media File-Set UID   | (0088,0140) | 3               |      |
|  |                | Referenced SOP Sequence      | (0008,1199) | 1               |      |
|  |                | >Referenced SOP Class UID    | (0008,1150) | 1               |      |
|  |                | >Referenced SOP Instance UID | (0008,1155) | 1               |      |
|  |                | >Storage Media File-Set ID   | (0088,0130) | 3               |      |
|  |                | >Storage Media File-Set UID  | (0088,0140) | 3               |      |
| Storage Commitment Request Complete-Failures Exist | 2              | Transaction UID              | (0008,1195) | 1               |      |
|  |                | Referenced SOP Sequence      | (0008,1199) | 1               |      |
|  |                | >Referenced SOP Class UID    | (0008,1150) | 1               |      |
|  |                | >Referenced SOP Instance UID | (0008,1155) | 1               |      |
|  |                | Failed SOP sequence          | (0008,1198) | 1               |      |
|  |                | >Referenced SOP Class UID    | (0008,1150) | 1               |      |
|  |                | >Referenced SOP Instance UID | (0008,1155) | 1               |      |
|  |                | >Failure Reason              | (0008,1197) | 1               |      |

The behavior of Storage AE when receiving Event Types within the N-EVENT-REPORT is summarized in the Table below:

**Table 2.2-16: Storage Commitment N-EVENT-REPORT Behavior**

| Event Type Name                                      | Event Type ID | Behavior                             |
|--|---------------|--------------------------------------|
| Storage Commitment Request Successful                | 1             | The job will be marked as SUCCESSFUL |
| Storage Commitment Request Complete – Failure exists | 2             | The job will be marked as FAILED.    |

The reasons for returning specific status codes in an N-EVENT-REPORT response are summarized in the table below:

**Table 2.2-17: Storage Commitment N-EVENT-REPORT Response Status Reasons**

| Service Status | Further Meaning | status Code | Reasons   |
|----------------|-----------------|-------------|---|
| Success        |                 | 0000        | The SCP has successfully returned all matching information. |

The behavior of the DICOM Store AE when receiving Event types over this association is the same as when receiving them over the same association as documented in section 2.2.1.3.1.3.2 .

**2.2.1.4.1.4      SOP Specific Conformance – Verification SOP Class (1.2.840.10008.1.1)**

The Storage AE provides standard conformance to the Verification SOP Class as an SCP. This verification is accomplished on an established Association using the C-ECHO DIMSE-C service.

## 2.2.2 DICOM Print Application Entity Specification

### 2.2.2.1 SOP Classes Supported

This Application Entity provides Standard Conformance to the following SOP Class(es):

**Table 2.2-18: SOP Class(es) for the DICOM Print Application Entity**

| SOP Class                                       | SOP Class UID          | SCU | SCP |
|---|------------------------|-----|-----|
| Verification                                    | 1.2.840.10008.1.1      | Yes | No  |
| Basic Grayscale Print Management Meta SOP Class | 1.2.840.10008.5.1.1.9  | Yes | No  |
| Basic Film Session SOP Class                    | 1.2.840.10008.5.1.1.1  | Yes | No  |
| Basic Film Box SOP Class                        | 1.2.840.10008.5.1.1.2  | Yes | No  |
| Basic Grayscale Image Box SOP Class             | 1.2.840.10008.5.1.1.4  | Yes | No  |
| Printer SOP Class                               | 1.2.840.10008.5.1.1.16 | Yes | No  |
| Print Job SOP Class                             | 1.2.840.10008.5.1.1.14 | Yes | No  |
| Presentation LUT SOP Class                      | 1.2.840.10008.5.1.1.23 | Yes | No  |

### 2.2.2.2 Association Establishment Policies

#### 2.2.2.2.1 General

The DICOM standard Application context is always proposed:

**Table 2.2-19: DICOM Application Context**

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

#### 2.2.2.2.2 Number of Associations

NX 3.0.8950 initiates one association at a time for each destination to which a print request is being processed in the active job queue list. Only one job per destination will be active at a time, the others remain pending until the active job for that destination is completed or failed. There can however be several simultaneous associations to different destinations.

**Table 2.2-20: Number of Associations as an Association Initiator for DICOM Print AE**

|   |                               |
|---|-------------------------------|
| Maximum number of simultaneous associations initiated | 1 per destination<br>(32 max) |
|---|-------------------------------|

#### 2.2.2.2.3 Asynchronous Nature

**Table 2.2-21: Asynchronous Nature as an Association Initiator for DICOM Print AE**

|   |      |
|---|------|
| Maximum number of outstanding asynchronous transactions | None |
|---|------|

DICOM Print AE allows a single outstanding operation on any association. Therefore it does not support asynchronous operations window negotiation, other than the default as specified by the DICOM specification.

### 2.2.2.2.4 Implementation Identifying Information

**Table 2.2-22: DICOM implementation Class and Version for DICOM Print AE**

|                             |                            |
|-----------------------------|----------------------------|
| Implementation Class UID    | <b>1.3.51.0.1.3</b>        |
| Implementation Version Name | <b>DPM1.XX<sup>3</sup></b> |

### 2.2.2.3 Association Initiation Policies

#### 2.2.2.3.1 Activity – Print Images

##### 2.2.2.3.1.1 Description and Sequencing of Activity

The user composes images into film sheets and requests them to be sent to a specific hardcopy device. A priori, the desired film format can be selected. Each sheet is internally processed, converted to a STANDARD/1,1 page and then one print job is forwarded to the job queue of the destination and processed individually.

The DICOM Print AE or NX 3.0.8950 will initiate a separate Association for each Print Session.

If the Printer rejects the Association, then NX 3.0.8950 issues a warning message. In case of a time-out (e.g. no answer from the Printer) or a warning message, the request will be retried after at least 20 seconds. In the meantime requests to other destinations will be handled.

After an association is established, the NX 3.0.8950 will send one film session to the Printer. Each film session will contain one film box, which in turn contains one image box.

The N-ACTION DIMSE service on Film Session SOP class instructs the printer to print the film session.

The print job has finished printing when the job is transferred to the printer or when the printer has sent the N-EVENT-REPORT "Done" (in case print job sop class is supported by the printer).

The NX 3.0.8950 releases the association. In case of N-EVENT-REPORT it will not wait for print job status "Done".

The default PDU size negotiated by the NX 3.0.8950 is 65542 bytes.

It is possible to print up to 16 bit (see §2.4.1.2.4).

<sup>3</sup> XX is the build version number






Figure 2.2-2: Sequencing of Print Images

### 2.2.2.3.1.2 Proposed Presentation Contexts

The DICOM Print AE is capable of proposing the Presentation Contexts shown in the following table:

**Table 2.2-23: Presentation Contexts Proposed by DICOM Print AE**

| SOP Class                                       |                        | Transfer Syntax                 |                   | Role | Extended Negotiation |
|---|------------------------|---------------------------------|-------------------|------|----------------------|
| SOP Class                                       | SOP Class UID          | Name                            | UID               |      |                      |
| Verification                                    | 1.2.840.10008.1.1      | DICOM Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None                 |
| Basic Grayscale Print Management Meta SOP Class | 1.2.840.10008.5.1.1.9  | DICOM Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None                 |
| Print Job SOP Class                             | 1.2.840.10008.5.1.1.14 | DICOM Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None                 |
| Basic Annotation Box SOP Class                  | 1.2.840.10008.5.1.1.15 | DICOM Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None                 |
| Presentation LUT SOP Class                      | 1.2.840.10008.5.1.1.23 | DICOM Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None                 |
| Print Queue Management SOP Class                | 1.2.840.10008.5.1.1.26 | DICOM Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None                 |

#### **Note:**

The Presentation Context shall use Abstract Syntax IDs which correspond to the SOP Classes UID of the Meta SOP Class specified in the first column of the Transfer Syntax Table or included SOP Classes. None of the included SOP Classes supports extended negotiation.

### 2.2.2.3.1.3 SOP Specific Conformance - Basic Grayscale Print Management Meta SOP Class

NX 3.0.8950 provides Standard conformance to the DICOM Basic Grayscale Print Management Class (1.2.840.10008.1.1.9) as SCU.

Support for Basic Grayscale Print Management as SCU also implies support for the following SOP Classes as SCU:

1. Basic Film Session SOP Class
2. Basic Film Box SOP Class
3. Basic Grayscale Image Box SOP Class
4. Printer SOP Class

#### 2.2.2.3.1.3.1 Basic Film Session SOP Class (1.2.840.10008.5.1.1.1)

The Basic Film Session IOD describes the presentation parameters which are common for all the films of a film session (e.g. number of films, film destination, ...)

The Basic Film Session SOP Instance refers to one or more Basic Film Box SOP Instance (§2.2.2.3.1.3.2).

NX 3.0.8950 (SCU) can send the following DIMSE services:

- N-CREATE
- N-ACTION

**N-CREATE** is issued by NX 3.0.8950 (SCU) to create a Basic Film Session SOP instance, when an Association has been established. The N-CREATE causes the Basic Film Session (root element) to be created by the SCP AE and its attributes initialized. The following attributes are supported:

**Table 2.2-24: Supported N-CREATE Attributes for a Basic Film Session**

| Attribute Name   | Tag         | Value                                       |
|------------------|-------------|---|
| Number of Copies | (2000,0010) | Always set to 1                             |
| Print priority   | (2000,0020) | Configurable                                |
| Medium Type      | (2000,0030) | BLUE FILM ,CLEAR FILM, PAPER (configurable) |
| Film destination | (2000,0040) | MAGAZINE, BIN-i , PROCESSOR (configurable)  |

NX 3.0.8950 (SCU) will process the N-CREATE confirmation and response Status codes. The following status codes are recognized:

**Table 2.2-25: N-CREATE Status Codes**

| Code  | Status  | Meaning   |
|-------|---|---|
| 0000  | Success   | Film Session successfully created   |
| 0116H | Warning (or Failure, depends on the selected Printer Profile) | The SCP AE returns the attribute "Value Out of Range". This may result in Image Quality Degradation.<br>NX 3.0.8950 will continue when its destination is an AGFA printer, but stops (ABORT) when its destination is a non-AGFA printer (might lead to unacceptable image quality). |
| 0106H | Error   | The SCP AE returns the attribute "Invalid Attribute Value" . This indicates that the requested memory allocation cannot be provided.  |
| 0213H |   | Resource Limitation is returned by the SCP AE for the Basic Film Session SOP Class to indicate that the requested allocation can temporarily not be provided.   |
| *     | Other Status codes  | Other DICOM error codes result in the failure of the job. Other warnings are not communicated to the user.  |

**N-ACTION** is issued by NX 3.0.8950 (SCU) to print a Film Session. This means that all subordinate Basic Film Boxes will be assembled into a print job for printing (the job can therefore contain more than one film).

NX 3.0.8950 (SCU) will process the N-ACTION confirmation and response Status codes. The following status codes in Table 2.2-26 are recognized:

**Table 2.2-26: N-ACTION Status Codes**

| Code | Status             | Meaning  |
|------|--------------------|--|
| 0000 | Success            | Normally returned by the SCP AE. Film(s) belonging to the film session are accepted for printing.          |
| B601 | Warning            | Film session printing (collation) is not supported.  |
| B602 |                    | Film Session SOP Instance hierarchy does not contain Image Box SOP Instances (empty page).                 |
| B604 |                    | Image size is larger than Image Box size, the image has been magnified.                                    |
| C600 | Error              | Film Session SOP Instance hierarchy does not contain Film Box SOP instances.                               |
| C601 |                    | Unable to create Film SOP Instance; The print queue is full (device failure).                              |
| C603 |                    | Image size is larger than Image box size   |
| *    | Other Status codes | Other DICOM error codes result in the failure of the job. Other warnings are not communicated to the user. |

**N-DELETE** is issued by NX 3.0.8950 (SCU) to delete a Film Session. This means that the complete Film Session SOP Instance hierarchy will be deleted.

### 2.2.2.3.1.3.2 Basic Film Box SOP Class (1.2.840.10008.5.1.1.2)

The Basic Film Box IOD is an abstraction of the presentation of one film of the film session. The Basic Film Box IOD describes the presentation parameters which are common for all images on a given sheet of film.

The Basic Film Box SOP Instance refers to one or more Image Box SOP Instances, zero or more film related Annotation Box SOP Instances, and zero or one Presentation LUT SOP Instance.

NX 3.0.8950 (SCU) can send the following DIMSE services:

➤ N-CREATE

**N-CREATE** is issued by NX 3.0.8950 (SCU) to create a Basic Film Box under the created Film Session and initialize its attributes. (The creation of a Basic Film Box also causes the subordinate Basic Image Boxes to be created for each location in the film format.)

The supported Film Box N-CREATE attributes are listed in the table below:

**Table 2.2-27: Supported N-CREATE Attributes for a Basic Film Box**

| Attribute Name            | Tag         | Value  |
|---------------------------|-------------|--|
| Image Display Format      | (2010,0010) | STANDARD\1,1   |
| Film Orientation          | (2010,0040) | PORTRAIT,LANDSCAPE   |
| Film Size ID              | (2010,0050) | 8INX10IN, 10INX12IN, 10INX14IN, 11INX14IN, 14INX14IN, 14INX17IN, A4, A3, 14INx36IN, 14INx51IN (configurable) |
| Magnification type        | (2010,0060) | NONE, REPLICATE  |
| Border Density            | (2010,0100) | WHITE, BLACK or ODx100 (configurable)  |
| Empty Image Density       | (2010,0110) | WHITE, BLACK or ODx100 (configurable)  |
| Min Density               | (2010,0120) | Configurable   |
| Max Density               | (2010,0130) | Configurable   |
| Trim                      | (2010,0140) | NO   |
| Configuration Information | (2010,0150) | Configurable   |

| Attribute Name                       | Tag         | Value   |
|--------------------------------------|-------------|---|
| Illumination                         | (2010,015E) | Configurable, only applicable for printing with P-Values  |
| Reflective Ambient Light             | (2010,0160) | Configurable, only applicable for printing with P-Values  |
| Referenced Film Session Sequence     | (2010,0500) | A sequence which provides references to a Film Session SOP Class/Instance pairs. Only a single item is permitted in this Sequence.                  |
| >Referenced SOP Class UID            | (0008,1150) | Uniquely identifies the referenced SOP Class  |
| >Referenced SOP Instance UID         | (0008,1155) | Uniquely identifies the referenced SOP Instance.  |
| Referenced Presentation LUT Sequence | (2050,0500) | A sequence which provides references to a Presentation LUT related SOP Class/Instance pairs. Only a single item shall be included in this sequence. |
| >Referenced SOP Class UID            | (0008,1150) | Uniquely identifies the referenced SOP Class  |
| >Referenced SOP Instance UID         | (0008,1155) | Uniquely identifies the referenced SOP Instance.  |

NX 3.0.8950 (SCU) will process the N-CREATE confirmation and response Status codes. The status codes below are recognized:

**Table 2.2-28: N-CREATE Status Codes**

| Code | Status             | Meaning  |
|------|--------------------|--|
| 0000 | Success            | Normally returned by the SCP AE. Film Box successfully created.  |
| B605 | Warning            | Requested Min Density or Max Density outside of printer's operating range. The printer will use its respective minimum or maximum density value instead. |
| *    | Other Status codes | Other DICOM error codes result in the failure of the job. Other warnings are not communicated to the user.   |

### 2.2.2.3.1.3.3 Basic Grayscale Image Box SOP Class (1.2.840.10008.5.1.1.4)

The Basic Image Box IOD is an abstraction of the presentation of an image and image related data in the image area of a film. The Basic Image Box IOD describes the presentation parameters and image pixel data which apply to a single image of a sheet of film.

The Basic Grayscale Image Box SOP Instance is created by the SCP at the time the Basic Film Box SOP Instance (§ 2.2.2.3.1.3.2) is created, based on the value of the Basic Film Box Attribute Image Display Format (2010, 0010).

The Basic Grayscale Image Box SOP Instance refers to zero or one Image Overlay Box SOP Instance and zero or one Presentation LUT SOP Instance.

NX 3.0.8950 (SCU) can send the following DIMSE services:

- N-SET

**N-SET** is issued by NX 3.0.8950 (SCU) to update an instance of the Grayscale Image Box SOP Class.

When all needed Basic Grayscale Image Boxes have been set, NX 3.0.8950 (SCU) issues a print command. There can be empty image positions. By using N-SET, NX 3.0.8950 (SCU) can instruct the SCP to erase the image in the image position by setting a zero length and no value in the attribute Basic Grayscale Image Sequence. The N-Set attributes for the Basic Grayscale Image Box are listed below:

**Table 2.2-29: Supported N-SET Attributes for a Basic Grayscale Image Box**

| Attribute Name                 | Tag         | Value   |
|--------------------------------|-------------|---|
| Image Position                 | (2020,0010) | 1   |
| Polarity                       | (2020,0020) | NORMAL  |
| Basic Grayscale Image Sequence | (2020,0110) |   |
| >Samples Per Pixel             | (0028,0002) | 1   |
| >Photometric Interpretation    | (0028,0004) | MONOCHROME2   |
| >Rows                          | (0028,0010) | larger than 0   |
| >Columns                       | (0028,0011) | larger than 0   |
| >Bits Allocated                | (0028,0100) | 8 – 16  |
| >Bits Stored                   | (0028,0101) | 8 – 16  |
| >High Bit                      | (0028,0102) | 7 to 15   |
| >Pixel Representation          | (0028,0103) | 0   |
| >Pixel Data                    | (7FE0,0010) | A data stream of the pixel samples that comprise the Image. |

NX 3.0.8950 (SCU) will process the N-SET confirmation and response Status codes. The status codes listed below in Table 2.2-30 are recognized:

**Table 2.2-30: N-SET Status Codes**

| Code | Status             | Meaning  |
|------|--------------------|--|
| 0000 | Success            | Normally returned by the SCP AE. Image successfully stored in Image Box  |
| B604 | Warning            | Image size is larger than Image Box size, the image has been magnified.  |
| B605 |                    | Requested Min Density or Max Density outside of printer's operating range. The printer will use its respective minimum or maximum density value instead. |
| C603 | Error              | Image size is larger than Image box size   |
| C605 |                    | Insufficient memory in printer to store the image.   |
| *    | Other Status codes | Other DICOM error codes result in the failure of the job. Other warnings are not communicated to the user.   |

#### 2.2.2.3.1.3.4 Printer SOP Class (1.2.840.10008.5.1.1.16)

The Printer IOD is an abstraction of the hard copy printer and is the basic Information Entity to monitor the status of the printer. The Printer SOP Instance is created by the SCP during start-up of the hard copy printer and has a well-known SOP Instance UID.

The Printer SOP Class is used to monitor the status of the printer.

NX 3.0.8950 (SCU) will accept the following DIMSE services:

- N-EVENT-REPORT

NX 3.0.8950 (SCU) can send the following DIMSE services:

- N-GET

**N-EVENT-REPORT** is used to report the changes of the printer status in an asynchronous way. The SCP uses the N-EVENT-REPORT to inform NX 3.0.8950 about each execution change. NX 3.0.8950 will return the confirmation of the N-EVENT-REPORT operation.

**Table 2.2-31: Notification Event Information**

| Event Type Name | Event Type ID | Attribute           | Tag         |
|-----------------|---------------|---------------------|-------------|
| Normal          | 1             | Printer Status Info | (2110,0020) |
|                 |               | Film Destination    | (2000,0040) |
|                 |               | Printer Name        | (2110,0030) |
| Warning         | 2             | Printer Status Info | (2110,0020) |
|                 |               | Film Destination    | (2000,0040) |
|                 |               | Printer Name        | (2110,0030) |
| Failure         | 3             | Printer Status Info | (2110,0020) |
|                 |               | Film Destination    | (2000,0040) |
|                 |               | Printer Name        | (2110,0030) |

**Note:**

If the Event Type Name = Warning or Failure then the warning/failure condition can be stored by the SCP Printer Status Info argument.

**N-GET** is issued by NX 3.0.8950 (SCU) to retrieve an instance of the Printer SOP class. NX 3.0.8950 specifies the UID of the SOP Instance to be retrieved. The supported N-GET attributes are listed in the table below:

**Table 2.2-32: Supported Attributes for N-GET on a Printer**

| Attribute Name          | Tag         | Value   |
|-------------------------|-------------|---|
| Printer Status          | (2110,0010) | NORMAL, WARNING, FAILURE                                      |
| Printer Status Info     | (2110,0020) | Printer dependent   |
| Printer Name            | (2110,0030) | User defined name identifying the printer                     |
| Manufacturer            | (0008,0070) | Manufacturer of the printer                                   |
| Manufacturer Model Name | (0008,1090) | Manufacturer's model number of the printer                    |
| Device Serial Number    | (0018,1000) | Manufacturer's serial number of the printer                   |
| Software Versions       | (0018,1020) | Manufacturer's designation of software version of the printer |
| Date Last Calibration   | (0018,1200) | Date when the printer was last calibrated                     |
| Time Last Calibration   | (0018,1201) | Time when the printer was last calibrated                     |

**2.2.2.3.1.4 SOP Specific Conformance Print Job SOP Class (1.2.840.10008.5.1.1.14)**

The Print Job IOD is an abstraction of the Print Job transaction and is the basic information entity to monitor the execution of the Print Process. A Print Job contains one film or multiple films, all belonging to the same film session.

The Print Job SOP Class is created by N-ACTION operation of the Film Session SOP Class (§2.2.2.3.1.3.1), Film Box SOP Class (§2.2.2.3.1.3.2), or Pull Print Request SOP Class. The Print Job SOP Instance is deleted after the films are printed or after a failure condition.

NX 3.0.8950 (SCU) will accept the following DIMSE services:

- N-EVENT-REPORT

NX 3.0.8950 (SCU) can send the following DIMSE services:

- N-GET

**N-EVENT-REPORT** is used by the SCP to report execution status changes to NX 3.0.8950 (SCU) in an asynchronous way.

N-EVENT-REPORT has the following arguments

**Table 2.2-33: Notification Event Information**

| Event Type Name | Event Type ID | Attribute             | Tag         |
|-----------------|---------------|-----------------------|-------------|
| Pending         | 1             | Execution Status Info | (2100,0030) |
|                 |               | Print Job ID          | (2100,0010) |
|                 |               | Film Session Label    | (2000,0050) |
|                 |               | Printer Name          | (2110,0030) |
| Printing        | 2             | Execution Status Info | (2100,0030) |
|                 |               | Print Job ID          | (2100,0010) |
|                 |               | Film Session Label    | (2000,0050) |
|                 |               | Printer Name          | (2110,0030) |
| Done            | 3             | Execution Status Info | (2100,0030) |
|                 |               | Print Job ID          | (2100,0010) |
|                 |               | Film Session Label    | (2000,0050) |
|                 |               | Printer Name          | (2110,0030) |
| Failure         | 4             | Execution Status Info | (2100,0030) |
|                 |               | Print Job ID          | (2100,0010) |
|                 |               | Film Session Label    | (2000,0050) |
|                 |               | Printer Name          | (2110,0030) |

**Note:**

The SCU only releases the Association after the receipt of the event type Done or Failure, if the print job sop class is supported.

NX 3.0.8950 (SCU) returns the confirmation from the N-EVENT-REPORT operation.

**N-GET** is used to retrieve an instance of the Print Job SOP Class. NX 3.0.8950 (SCU) uses the N-GET to request the SCP to get a Print Job SOP Instance. NX 3.0.8950 specifies the UID of the SOP Instance to be retrieved.

**Table 2.2-34: Supported N-GET Attributes for a Print Job**

| Attribute Name       | Tag         | Value  |
|----------------------|-------------|--|
| Execution Status     | (2100,0020) | PENDING, PRINTING, DONE, FAILURE               |
| Execution State Info | (2100,0030) | Printer dependent                              |
| Print Priority       | (2000,0020) | HIGH, MED, LOW                                 |
| Creation Date        | (2100,0040) | Date of print job creation                     |
| Creation Time        | (2100,0050) | Time of print job creation                     |
| Printer Name         | (2110,0030) | User defined name identifying the printer      |
| Originator           | (2100,0070) | DICOM AE title that issued the print operation |



### 2.2.2.3.1.5 SOP Specific Conformance Presentation LUT SOP Class (1.2.840.10008.5.1.1.23)

The Presentation LUT Information Object is an abstraction of a Presentation LUT. The objective of the Presentation LUT is to realize image display tailored for specific modalities, applications, and user preferences. It is used to prepare image pixel data for display on devices that conform to the Grayscale Standard Display Function.

The output of the Presentation LUT is Presentation Values (P-Values).

An **N-CREATE** is issued by NX 3.0.8950 (SCU) to create a Presentation LUT SOP Instance. The supported Presentation LUT attributes are listed below:

**Table 2.2-35: Supported Attributes for Presentation LUT**

| Tag         | Name                   | Supported       | Default  |
|-------------|------------------------|-----------------|----------|
| (2050,0020) | Presentation LUT Shape | IDENTITY (Note) | IDENTITY |

**Note:**

'Presentation LUT Sequence' is not supported.

NX 3.0.8950 (SCU) will process the N-CREATE confirmation and response Status codes. The status codes listed below are recognized:

**Table 2.2-36: N-CREATE confirmation and response Status codes**

| Code | Status  | Meaning  |
|------|---------|--|
| 0000 | Success | Presentation LUT successfully created  |
| B605 | Warning | Requested Min or Max Density outside of printer's operating range. The printer will use its respective minimum or maximum value instead. |

**Note:**

NX 3.0.8950 (SCU) uses the N-CREATE Service Element to request the SCP to create a Presentation LUT SOP Instance. NX 3.0.8950 shall initialize the Attributes of the SOP Class.

The Presentation LUT persists in the SCP as long as the Association in which it was created is open or an explicit N-DELETE is issued by the SCU.

### 2.2.2.3.1.6 SOP Specific Conformance Verification SOP Class (1.2.840.10008.1.1)

See § [2.2.1.4.1.4](#).

## 2.2.3 RIS Application Entity Specification

### 2.2.3.1 SOP Classes Supported

This Application Entity provides Standard Conformance to the following SOP Class(es):

**Table 2.2-37: SOP Class(es) for the DICOM Store Application Entity**

| SOP Class Name                             | SOP Class UID          | SCU | SCP |
|--|------------------------|-----|-----|
| Modality Worklist Information Model – FIND | 1.2.840.10008.5.1.4.31 | Yes | No  |

### 2.2.3.2 Association Establishment Policies

#### 2.2.3.2.1 General

The DICOM standard Application context is always proposed:

**Table 2.2-38: DICOM Application Context**

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

#### 2.2.3.2.2 Number of Associations

NX 3.0.8950 initiates one association at a time to query the worklist.

**Table 2.2-39: Number of Associations as an Association Initiator for RIS AE**

|   |   |
|---|---|
| Maximum number of simultaneous associations initiated | 1 |
|---|---|

#### 2.2.3.2.3 Asynchronous Nature

**Table 2.2-40: Asynchronous Nature as an Association Initiator for RIS AE**

|   |   |
|---|---|
| Maximum number of outstanding asynchronous transactions | 1 |
|---|---|

NX 3.0.8950 does not support asynchronous communication (multiple outstanding transactions over a single connection).

#### 2.2.3.2.4 Implementation Identifying Information

**Table 2.2-41: DICOM implementation Class and Version for DICOM RIS AE**

|                             |                             |
|-----------------------------|-----------------------------|
| Implementation Class UID    | 1.3.51.0.1.3                |
| Implementation Version Name | AGFA DTF1.0.XX <sup>4</sup> |

<sup>4</sup> XX is the build version number

### 2.2.3.3 Association Initiation Policies

#### 2.2.3.3.1 Activity – Query RIS

##### 2.2.3.3.1.1 Description and Sequencing of Activity

The request for Query RIS is initiated by user interaction (pressing the “Query RIS” button) or automatically at specific time intervals (configurable by the user). Depending on the configuration this can either be a query based on a user provided accession number or be a complete worklist update based on date, modality and Scheduled Station AE title.




Figure 2.2-3: Sample Sequencing Diagram for Refresh Worklist

##### 2.2.3.3.1.2 Proposed Presentation Contexts

The RIS Application Entity is capable of proposing the Presentation Contexts shown in the following table:

Table 2.2-42: Presentation Contexts Proposed by DICOM Store AE

| Presentation Context Table                 |                        |                           |                   |      |                      |
|--|------------------------|---------------------------|-------------------|------|----------------------|
| Abstract Syntax                            |                        | Transfer Syntax           |                   | Role | Extended Negotiation |
| Name                                       | UID                    | Name List                 | UID List          |      |                      |
| Modality Worklist Information Model – FIND | 1.2.840.10008.5.1.4.31 | Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None                 |

##### 2.2.3.3.1.3 SOP Specific Conformance – Modality Worklist SOP Class (1.2.840.10008.5.1.4.31)

The Modality Worklist SOP class, defined within the Basic Worklist Management Service Class, defines an application-level class of service which facilitates the communication of information to the imaging modality about Scheduled Procedure Steps, and entities related to the Scheduled Procedure Steps.

This worklist is structured according to Scheduled Procedure Steps. A procedure step is a unit of service in the context of a requested imaging procedure.

The behavior of a RIS AE when encountering status codes in a C-FIND response is summarized in the Table below:

**Table 2.2-43: C-Find Response Status Handling Behavior**

| Service Status | Status Code        | Status Code  |
|----------------|--------------------|--|
| Success        | 0000               | Matching is complete   |
| Refused        | A700               | Out of resources   |
| Failed         | A900               | Identifier does not match SOP Class  |
|                | C000-CFFF          | Unable to Process  |
| Cancel         | FE00               | Matching terminated due to Cancel request  |
| Pending        | FF00               | Matches are continuing – Current Match is supplied and any Optional Keys were supported in the same manner as Required Keys. |
|                | FF01               | Matches are continuing – Warning that one or more Optional Keys were not supported   |
| *              | Other Status codes | Other DICOM error codes result in the failure of the job. Other warnings are not communicated to the user.                   |

The behavior of the AE during communication failure is summarized in a table as follows:

**Table 2.2-44: DICOM Command Communication Failure Behavior**

| Exception           | Behavior |
|---------------------|----------|
| Timeout             |          |
| Association aborted |          |

The table below provides a description of the Worklist Request Identifier. Unexpected attributes returned in a C-FIND response are ignored.

**Table 2.2-45: Worklist request identifiers**

| Attribute name                               | Tag         | Broad Query           | Accession number query |
|--|-------------|-----------------------|------------------------|
| <b>Scheduled Procedure Step</b>              |             |                       |                        |
| Scheduled Station AE Title                   | (0040,0001) | Single Value Matching |                        |
| Scheduled Procedure Step Start Date          | (0040,0002) | Range Matching        |                        |
| Scheduled Procedure Step Start Time          | (0040,0003) |                       |                        |
| Scheduled Procedure Step Location            | (0040,0011) |                       |                        |
| Modality                                     | (0008,0060) | Single Value Matching |                        |
| Scheduled Procedure Step Description         | (0040,0007) |                       |                        |
| Scheduled Performing Physician's Name        | (0040,0006) |                       |                        |
| Scheduled Performing Physician's ID Sequence | (0040,000B) |                       |                        |
| Scheduled Station Name                       | (0040,0010) |                       |                        |
| Scheduled Protocol Code Sequence             | (0040,0008) |                       |                        |
| Pre-Medication                               | (0040,0012) |                       |                        |
| Scheduled Procedure Step ID                  | (0040,0009) |                       |                        |
| Request Contrast Agent                       | (0032,1070) |                       |                        |
| Scheduled Procedure Step Status              | (0040,0020) |                       |                        |
| Scheduled Procedure Step Comment             | (0040,0400) |                       |                        |

| Attribute name   | Tag         | Broad Query | Accession number query |
|--|-------------|-------------|------------------------|
| <b>Requested Procedure</b>                             |             |             |                        |
| Requested Procedure Description                        | (0032,1060) |             |                        |
| Requested Procedure Code Sequence                      | (0032,1064) |             |                        |
| Requested Procedure ID                                 | (0040,1001) |             |                        |
| Study Instance UID                                     | (0020,000D) |             |                        |
| Referenced Study Sequence                              | (0008,1110) |             |                        |
| Reason for the Requested Procedure                     | (0040,1002) |             |                        |
| Requested Procedure Priority                           | (0040,1003) |             |                        |
| Patient Transport Arrangement                          | (0040,1004) |             |                        |
| Requested Procedure Location                           | (0040,1005) |             |                        |
| Requested Procedure Comments                           | (0040,1400) |             |                        |
| Confidentiality Code                                   | (0040,1008) |             |                        |
| Reporting Priority                                     | (0040,1009) |             |                        |
| Names of Intended Recipients of Results                | (0040,1010) |             |                        |
| Intended Recipients of Results Identification Sequence | (0040,1011) |             |                        |
| <b>Imaging Service Request</b>                         |             |             |                        |
| Accession Number                                       | (0008,0050) |             | Wildcard matching      |
| Requesting Physician                                   | (0032,1032) |             |                        |
| Requesting Physician Identification Sequence           | (0032,1031) |             |                        |
| Referring Physician's Name                             | (0008,0090) |             |                        |
| Referring Physician Identification Sequence            | (0008,0096) |             |                        |
| Requesting Service                                     | (0032,1033) |             |                        |
| Placer Order Number                                    | (0040,2016) |             |                        |
| Filler Order Number                                    | (0040,2017) |             |                        |
| Imaging Service Request Comments                       | (0040,2400) |             |                        |
| Issue Date of Imaging Service Request                  | (0040,2004) |             |                        |
| Issue Time of Imaging Service Request                  | (0040,2005) |             |                        |
| Order Entered by...                                    | (0040,2008) |             |                        |
| Order Enterer's Location                               | (0040,2009) |             |                        |
| Order Callback Phone Number                            | (0040,2010) |             |                        |
| Reason for imaging service                             | (0040,2001) |             |                        |
| <b>Visit Identification</b>                            |             |             |                        |
| Admission ID   | (0038,0010) |             |                        |
| Institution Name                                       | (0008,0080) |             |                        |
| Institution Address                                    | (0008,0081) |             |                        |
| Institution Code Sequence                              | (0008,0082) |             |                        |
| Issuer of Admission ID                                 | (0038,0011) |             |                        |
| <b>Visit Status</b>                                    |             |             |                        |
| Current Patient Location                               | (0038,0300) |             |                        |
| Visit Status ID  | (0038,0008) |             |                        |
| Patient's Institution Residence                        | (0038,0400) |             |                        |
| Visit Comments   | (0038,4000) |             |                        |
| <b>Visit Relationship</b>                              |             |             |                        |
| Referenced Patient Sequence                            | (0008,1120) |             |                        |
| <b>Visit Admission</b>                                 |             |             |                        |
| Route of Admissions                                    | (0038,0016) |             |                        |
| Admitting Date   | (0038,0020) |             |                        |
| Admitting Time   | (0038,0021) |             |                        |

| Attribute name                                     | Tag         | Broad Query | Accession number query |
|--|-------------|-------------|------------------------|
| Scheduled Admission Date                           | (0038,001A) |             |                        |
| Scheduled Admission Time                           | (0038,001B) |             |                        |
| Admitting Diagnosis Description                    | (0008,1080) |             |                        |
| Admitting Diagnosis Code Sequence                  | (0008,1084) |             |                        |
| <b>Patient Identification</b>                      |             |             |                        |
| Patient's Name                                     | (0010,0010) |             |                        |
| Patient ID   | (0010,0020) |             |                        |
| Other Patient IDs                                  | (0010,1000) |             |                        |
| Issuer of Patient ID                               | (0010,0021) |             |                        |
| Other Patient Names                                | (0010,1001) |             |                        |
| Patient's Birth Name                               | (0010,1005) |             |                        |
| Patient's Mother's Birth Name                      | (0010,1060) |             |                        |
| Medical Record Locator                             | (0010,1090) |             |                        |
| Referenced Patient Alias Sequence                  | (0038,0004) |             |                        |
| <b>Patient Demographic</b>                         |             |             |                        |
| Patients Birth Date                                | (0010,0030) |             |                        |
| Patient's Sex                                      | (0010,0040) |             |                        |
| Patient's Birth Time                               | (0010,0032) |             |                        |
| Patient's Primary Language Code Sequence           | (0010,0101) |             |                        |
| >Patient's Primary Language Code Modifier Sequence | (0010,0102) |             |                        |
| Patient's Weight                                   | (0010,1030) |             |                        |
| Patient's Size                                     | (0010,1020) |             |                        |
| Patient's Age                                      | (0010,1010) |             |                        |
| Military Rank                                      | (0010,1080) |             |                        |
| Branch of Service                                  | (0010,1081) |             |                        |
| Ethnic Group                                       | (0010,2160) |             |                        |
| Occupation   | (0010,2180) |             |                        |
| Patient Comment                                    | (0010,4000) |             |                        |
| Confidentiality constraint on patient data         | (0040,3001) |             |                        |
| Patient's Insurance Plan Code Sequence             | (0010,0050) |             |                        |
| Patient's Address                                  | (0010,1040) |             |                        |
| Country of Residence                               | (0010,2150) |             |                        |
| Region of Residence                                | (0010,2152) |             |                        |
| Patient's Telephone Numbers                        | (0010,2154) |             |                        |
| Patient's Religious Preference                     | (0010,21F0) |             |                        |
| <b>Patient Medical</b>                             |             |             |                        |
| Patient State                                      | (0038,0500) |             |                        |
| Pregnancy Status                                   | (0010,21C0) |             |                        |
| Medical Alerts                                     | (0010,2000) |             |                        |
| Contrast Allergies                                 | (0010,2110) |             |                        |
| Special Needs                                      | (0038,0050) |             |                        |
| Additional Patient History                         | (0010,21B0) |             |                        |
| Last Menstrual Date                                | (0010,21D0) |             |                        |
| Smoking Status                                     | (0010,21A0) |             |                        |

### 2.2.3.4 User Defined fields

There are five (5) optional (type 3) User Defined fields that can be mapped from the RIS and which can be sent out during archiving and exporting.

These fields can be found in the following Table:

**Table 2.2-46: List of User defined DICOM fields**

| Tag         | Description          | VR |
|-------------|----------------------|----|
| (0019,10F0) | User Defined field 1 | LO |
| (0019,10F1) | User Defined field 2 | LO |
| (0019,10F2) | User Defined field 3 | LO |
| (0019,10F3) | User Defined field 4 | LO |
| (0019,10F4) | User Defined field 5 | LO |

### 2.2.3.5 Custom Patient Field

There are five (5) optional (type 3) Custom Patient fields that can be mapped on the RIS and which can be sent out during archiving and exporting.

These fields can be found in the following Table:

**Table 2.2-47: List of Custom Patient DICOM fields**

| Tag         | Description            | VR |
|-------------|------------------------|----|
| (0019,10D0) | Custom Patient Field 1 | LO |
| (0019,10D1) | Custom Patient Field 2 | LO |
| (0019,10D2) | Custom Patient Field 3 | LO |
| (0019,10D3) | Custom Patient Field 4 | LO |
| (0019,10D4) | Custom Patient Field 5 | LO |

## 2.2.4 MPPS Application Entity Specification

### 2.2.4.1 SOP Classes Supported

This Application Entity provides Standard Conformance to the following SOP Class(es):

**Table 2.2-48: SOP Class(es) for the DICOM Store Application Entity**

| SOP Class Name                    | SOP Class UID           | SCU | SCP |
|-----------------------------------|-------------------------|-----|-----|
| Modality Performed Procedure Step | 1.2.840.10008.3.1.2.3.3 | Yes | No  |

### 2.2.4.2 Association Establishment Policies

#### 2.2.4.2.1 General

The DICOM standard Application context is always proposed:

**Table 2.2-49: DICOM Application Context**

|                          |                       |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |
|--------------------------|-----------------------|

#### 2.2.4.2.2 Number of Associations

NX 3.0.8950 initiates one association at a time to query the worklist.

**Table 2.2-50: Number of Associations as an Association Initiator for MPPS AE**

|   |   |
|---|---|
| Maximum number of simultaneous associations initiated | 1 |
|---|---|

#### 2.2.4.2.3 Asynchronous Nature

**Table 2.2-51: Asynchronous Nature as an Association Initiator for MPPS AE**

|   |   |
|---|---|
| Maximum number of outstanding asynchronous transactions | 1 |
|---|---|

NX 3.0.8950 does not support asynchronous communication (multiple outstanding transactions over a single connection)

#### 2.2.4.2.4 Implementation Identifying Information

**Table 2.2-52: DICOM implementation Class and Version for DICOM MPPS AE**

|                             |                             |
|-----------------------------|-----------------------------|
| Implementation Class UID    | 1.3.51.0.1.3                |
| Implementation Version Name | AGFA DTF1.0.XX <sup>5</sup> |

<sup>5</sup> XX is the build version number



## 2.2.4.3 Association Initiation Policies

### 2.2.4.3.1 Activity – Acquire Images

#### 2.2.4.3.1.1 Description and Sequencing of Activity

When the first image of a session arrives, an association will be opened to create an MPPS instance (via N-CREATE) on the MPPS manager. NX 3.0.8950 will wait for an N-CREATE response from the Department Scheduler (status success). When the N-CREATE response is received, the association with the Department Scheduler will be closed. When the session is closed on NX 3.0.8950, an association will be opened to complete the MPPS instance (via N-SET). Again NX 3.0.8950 will wait for an N-SET response from the Department Scheduler (status success). Once the response is received, the association will be closed.

When all images of a study are transferred to another study and no new images are acquired before the session is closed, an N-SET MPPS discontinued will be sent instead of the normal N-SET response.

When no images are available for 1 of 2 selected SPS, N-SET MPPS discontinued will be sent with discontinue reason code: "Incorrect Worklist Entry Selected".




Figure 2.2-4: Sample Sequencing Diagram for MPPS

### 2.2.4.3.1.2 Proposed Presentation Contexts

The MPPS Application Entity is capable of proposing the Presentation Contexts shown in the following table:

**Table 2.2-53: Presentation Contexts Proposed by DICOM Store AE**

| Presentation Context Table        |                         |                           |                     |      |                      |
|-----------------------------------|-------------------------|---------------------------|---------------------|------|----------------------|
| Abstract Syntax                   |                         | Transfer Syntax           |                     | Role | Extended Negotiation |
| Name                              | UID                     | Name List                 | UID List            |      |                      |
| Modality Performed Procedure Step | 1.2.840.10008.3.1.2.3.3 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None                 |
|                                   |                         | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |                      |

### 2.2.4.3.1.3 SOP Specific Conformance – MPPS (1.2.840.10008.3.1.2.3.3)

NX 3.0.8950 (SCU) can send the following DIMSE services:

- N-CREATE
- N-SET

An **N-CREATE** allows NX 3.0.8950 to create an instance of the Modality Performed Procedure Step SOP Class and provide information about a specific real-world Performed Procedure Step that is under control of NX 3.0.8950.

An **N-SET** allows NX 3.0.8950 to set Attribute Values of an instance of the Modality Performed Procedure Step SOP Class and provide information about a specific real-world Modality Performed Procedure Step that is under control of NX 3.0.8950.

#### Note:

NX 3.0.8950 informs the Information System as soon as possible that the performance of the Procedure Step has been started by sending the N-CREATE Service Request. This allows an SCP of the Modality Worklist SOP Class (if supported) to update the Modality Worklist. Some of the attribute values are already known at the beginning of the Procedure Step, they are sent in the N-CREATE command. Other mandatory attributes are known only at the end of the Performed Procedure Step, they are assigned a value in the N-SET command.

The behavior of MPPS AE when encountering status codes in an MPPS N-CREATE or N-SET response is summarized in the Table below:

**Table 2.2-54: N-CREATE/N-SET Response Status Handling Behavior**

| Service Status | Further Meaning   | Status Code | Behavior   |
|----------------|---|-------------|--|
| Success        | Matching is complete  | 0000        | The SCP has completed the operation successfully   |
| Failure        | Processing Failure – Performed Procedure Step Object may no longer be updated | 0110H       |  |
| *              | Other Status codes  | *           | Other DICOM error codes result in the failure of the job. Other warnings are not communicated to the user. |

The behavior of the AE during communication failure is summarized in a table as follows:

**Table 2.2-55: DICOM Command Communication Failure Behavior**

| Exception           | Behavior   |
|---------------------|--|
| Timeout             | If the sending fails, the system will retry automatically until the session is cleaned up. |
| Association aborted | If the sending fails, the system will retry automatically until the session is cleaned up. |

The following table provides a description of the MPPS N-CREATE and N-SET request identifiers sent by the MPPS AE. Empty cells in the N-CREATE and N-SET columns indicate that the attribute is not sent. An "X" indicates that an appropriate value will be sent.

**Table 2.2-56: MPPS N-CREATE/N-SET Request Identifier**

| Attribute  | Tag         | N-CREATE               | N-SET |
|--|-------------|------------------------|-------|
| Specific Character Set                                       | (0008,0005) | X                      | X     |
| <b>Performed Procedure Step Relationship</b>                 |             |                        |       |
| Scheduled Step Attribute Sequence                            | (0040,0270) |                        |       |
| >Study Instance UID  | (0020,000D) | X                      |       |
| >Referenced Study Sequence                                   | (0008,1110) | From MWL or user input |       |
| >>Referenced SOP Class UID                                   | (0008,1150) | From MWL or user input |       |
| >>Referenced SOP Instance UID                                | (0008,1155) | From MWL or user input |       |
| >Accession Number  | (0008,0050) | From MWL or user input | N/A   |
| >Placer Order Number/Imaging Service Request                 | (0040,2016) | From MWL or user input |       |
| >Filler Order Number/Imaging Service Request                 | (0040,2017) | From MWL or user input |       |
| >Requested Procedure ID                                      | (0040,1001) | From MWL or user input |       |
| >Requested Procedure Description                             | (0032,1060) | From MWL or user input |       |
| >Scheduled Procedure Step ID                                 | (0040,0009) | From MWL or user input |       |
| >Scheduled Procedure Step Description                        | (0040,0007) | From MWL or user input |       |
| >Scheduled Protocol Code Sequence                            | (0040,0008) | From MWL or user input |       |
| >>Code Value   | (0008,0100) | From MWL or user input |       |
| >>Coding Scheme designator                                   | (0008,0102) | From MWL or user input |       |
| >>Coding Scheme Version                                      | (0008,0103) | From MWL or user input |       |
| >>Code Meaning   | (0008,0104) | From MWL or user input |       |
| >>All other Attributes from Scheduled Protocol Code Sequence |             |                        |       |
| Patient's Name   | (0010,0010) | From MWL or user input |       |
| Patient ID   | (0010,0020) | From MWL or user input |       |
| Patient's Birth Date   | (0010,0030) | From MWL or user input |       |
| Patient's Sex  | (0010,0040) | From MWL or user input |       |
| Referenced Patient Sequence                                  | (0008,1120) | From MWL or user input |       |
| >Referenced SOP Class UID                                    | (0008,1150) | From MWL or user input |       |
| >Referenced Instance UID                                     | (0008,1155) | From MWL or user input |       |
| <b>Performed Procedure Step Information</b>                  |             |                        |       |
| Performed Procedure Step ID                                  | (0040,0253) | X                      |       |
| Performed Station AE Title                                   | (0040,0241) | X                      |       |
| Performed Station Name                                       | (0040,0242) | X                      |       |
| Performed Location   | (0040,0243) | X                      |       |
| Performed Procedure Step Start Date                          | (0040,0244) | X                      |       |
| Performed Procedure Step Start Time                          | (0040,0245) | X                      |       |

| Attribute   | Tag         | N-CREATE  | N-SET   |
|---|-------------|---|---|
| Performed Procedure Step Status                               | (0040,0252) | "IN PROGRESS"   | "COMPLETED"   |
| Performed Procedure Step Description                          | (0040,0254) | Exam group names of the exposures in the study                                | X   |
| Performed Procedure Type Description                          | (0040,0255) | set to the Requested Procedure Description from the SPS                       | X   |
| Procedure Code Sequence                                       | (0008,1032) | Zero length   | Zero length   |
| Performed Procedure Step End Date                             | (0040,0250) | Empty   | close session date  |
| Performed Procedure Step End Time                             | (0040,0251) | Empty   | close session time  |
| Comments on the Performed Procedure Step                      | (0040,0280) | Empty   | only if user provided this information  |
| Performed Procedure Step Discontinuation Reason Code Sequence | (0040,0281) | Empty   | Available if discontinued   |
| >Code Value   | (0008,0100) | Empty   | ""  |
| >Coding Scheme Designator                                     | (0008,0102) | Empty   | ""  |
| >Coding Scheme Version  | (0008,0103) | Empty   | ""  |
| >Code Meaning   | (0008,0104) | Empty   | ""  |
| <b>Image Acquisition Results</b>                              |             |   |   |
| Modality  | (0008,0060) | CR or DX  |   |
| Study ID  | (0020,0010) | set to Requested Procedure ID or automatically generated for unscheduled case |   |
| Performed Protocol Code Sequence                              | (0040,0260) | Empty   | Empty   |
| >Code Value   | (0008,0100) | Empty   | Empty   |
| >Coding Scheme Designator                                     | (0008,0102) | Empty   | Empty   |
| >Coding Scheme Version  | (0008,0103) | Empty   | Empty   |
| >Code Meaning   | (0008,0104) | Empty   | Empty   |
| >All other Attributes from Performed Protocol Code Sequence   |             |   |   |
| Performed Series Sequence                                     | (0040,0340) | Empty   |   |
| >Performing Physician's Name                                  | (0008,1050) | Empty   | If user provided this information. Otherwise empty.                                       |
| >Protocol Name  | (0018,1030) | Empty   | concatenation of "CR " or "DX " plus all exam groups of the exposures, taken in this MPPS |
| >Operator's Name  | (0008,1070) | Empty   | Operator who was logged in at the Create of this MPPS                                     |
| >Series Instance UID  | (0020,000E) | Empty   | X   |
| >Series Description   | (0008,103E) | Empty   | BODY PART [space] VIEW POSITION   |
| >Retrieve AE Title  | (0008,0054) | Empty   | Blank   |
| >Referenced Image Sequence                                    | (0008,1140) | Empty   | X   |
| >>Referenced SOP Class UID                                    | (0008,1150) | Empty   | X   |
| >>Referenced SOP Instance UID                                 | (0008,1155) | Empty   | X   |
| >Referenced Non-Image Composite SOP Instance Sequence         | (0040,0220) | Empty   | X   |
| >>Referenced SOP Class UID                                    | (0008,1150) | Empty   | X   |
| >>Referenced SOP Instance UID                                 | (0008,1155) | Empty   | X   |

| Attribute  | Tag         | N-CREATE | N-SET |
|--|-------------|----------|-------|
| <b>Billing And Material Management Code Module</b> |             |          |       |
| Film Consumption Sequence                          | (0040,0321) | Empty    | X     |
| >Number of Films                                   | (2100,0170) | Empty    | X     |
| >Medium Type                                       | (2000,0030) | Empty    | X     |
| >Film Size ID                                      | (2010,0050) | Empty    | X     |
| <b>Radiation Dose Module</b>                       |             |          |       |
| Total Number of Exposures                          | (0040,0301) | Empty    | X     |
| Distance Source to Detector                        | (0018,1110) | Empty    | X     |
| Distance Source to Entrance                        | (0040,0306) | Empty    | X     |
| Entrance Dose                                      | (0040,0302) | Empty    | X     |
| Entrance Dose in mGy                               | (0040,8302) | Empty    | X     |
| Exposed Area                                       | (0040,0303) | Empty    | X     |
| Image and Fluoroscopy Area Dose Product            | (0018,115E) | Empty    | X     |
| Exposure Dose Sequence                             | (0040,030E) | Empty    | X     |
| >KVp   | (0018,0060) | Empty    | X     |
| >X-ray Tube Current in $\mu$ A                     | (0018,8151) | Empty    | X     |
| >Exposure Time                                     | (0018,1150) | Empty    | X     |
| >Filter Type <sup>6</sup>                          | (0018,1160) | Empty    | X     |
| >Filter Material                                   | (0018,7050) | Empty    | X     |
| >Comments on Radiation Dose <sup>7</sup>           | (0040,0310) | Empty    | X     |

**Note:**

1. The requirement for the final state is that which applies at the time that the Performed Procedure Step Status (0040,0252) is N-SET to a value of COMPLETED or DISCONTINUED. It is only described if it is different from the SCP requirement for the N-CREATE.
2. The Performed Series Sequence (0040,0340) may not be empty (zero length) at the time that the Performed Procedure Step Status (0040,0252) is N-SET to a value of COMPLETED or DISCONTINUED. In other words a Series must exist for every Performed Procedure Step, though it may contain no Images or Non-Image Composite objects, if none were created.
3. Only attributes that are specified in a SOP Instance at N-CREATE may later be updated through the N-SET. If NX 3.0.8950 wishes to use the PPS Discontinuation Reason Code Sequence (0040,0281), it creates that attribute (zero-length) during MPPS N-CREATE.

<sup>6</sup> The attribute 'Filter Type' (0018,1160) is implemented but no information is passed yet.

<sup>7</sup> The attribute 'Comments on Radiation Dose'(0040,0310) is implemented but no information is passed yet.

## 2.3 Network Interfaces

NX 3.0.8950 provides DICOM V3.0 TCP/IP Network Communication Support as defined in PS 3.8 of the DICOM Standard (2011). NX 3.0.8950 inherits its TCP/IP stack from the computer system upon which it executes.

### 2.3.1 Physical Medium Support

NX 3.0.8950 is indifferent to the physical medium over which TCP/IP executes; it inherits the medium from the computer system upon which it is being executed.

### 2.3.2 Additional Protocols

NX 3.0.8950 can use DNS to resolve hostnames. It will use the TCP/IP stack from the Windows System it runs on.

For audit trailing the NX 3.0.8950 will query an NTP server as time reference. This NTP Server can be configured in the "Workstation Service & Configuration Tool". For all other cases the local system clock will be used as a time reference.

## 2.4 Configuration

### 2.4.1 AE Title/ Presentation Address Mapping

#### 2.4.1.1 Local AE Titles

NX 3.0.8950 uses the AE Titles and TCP/IP Ports configured by means of the "Workstation Service & Configuration Tool". The Field Service Engineer can configure the TCP Port via the "Workstation Service & Configuration Tool". No Default AE Titles are provided. The AE Titles must be configured during installation. The local AE Title used by each individual application can be configured independently of the AE Title used by other local applications.

**Table 2.4-1: AE Title Configuration Table**

| Application Entity | Default AE Title | Default TCP/IP Port |
|--------------------|------------------|---------------------|
| Storage            | No Default       | 104                 |
| Hardcopy           | No Default       | Not Applicable      |
| RIS                | No Default       | Not Applicable      |
| MPPS               | No Default       | Not Applicable      |

#### 2.4.1.2 Remote AE Title/ Presentation Address Mapping

The AE title, host names and port numbers of remote applications are configured using the "Workstation Service & Configuration Tool".

##### 2.4.1.2.1 Storage

The "Workstation Service & Configuration Tool" must be used to set the AE Titles, port-numbers, host-names and capabilities for the remote Storage SCPs. Associations will only be accepted from known AE Titles while associations from unknown AE Titles will be rejected (an AE Title is known if it can be selected within the "Workstation Service & Configuration Tool"). Multiple remote Storage SCPs can be defined. Any Storage SCP can be configured to be an "Archive" device causing storage commitment to be requested for images or presentation states transmitted to the device.

Storage Commit can be enabled or disabled.

The Archive Connection can be tested by means of a DICOM Ping (C-echo) test to the archive.

##### 2.4.1.2.1.1 Archive options

Sending images to an archive can be done by means of the following SOP classes: CR, DX and MG for presentation or DX and MG for processing. The latter is a RAW image, without any annotations or measurements burned in.

Image pixels can be sent in the requested output type (configurable) depending on the used SOP class and whether GSPS is supported or not<sup>8</sup>.

<sup>8</sup> anonymous archiving (i.e. without patient demographics) is not supported

**CR image, WITHOUT GSPS:**

- 8 bit OD
- 8 bit OD Gamma correction
- 8 bit P-Value
- 12-bit OD
- 12-bit OD Gamma correction
- 12 bit P-Value
- 12-bit OD REL
- 12-bit Image + VOI Lut OD REL
- 12-bit Image + VOI Lut Gamma
- 15-bit Image + VOI Lut P-Value

Note 1: '15-bit Image VOI Lut P-Value' is the AGFA recommended post-processing when it is supported by the SCP.

Note 2: it is possible to send P-values even in CR images, using the (2050,0020) Presentation Lut Shape = "Identity" in the General Image Module.

**CR image, GSPS ENABLED:**

→ Only P-value output formats are allowed, since the output of GSPS has to be P-values anyway:

- 8 bit P-Value
- 12 bit P-Value
- 15-bit Image + VOI Lut P-Value

Note 1: '15-bit Image VOI Lut P-Value' is the AGFA recommended post-processing when it is supported by the SCP.

Note 2: when GSPS is enabled, the LUT tables in the GSPS are a copy of the LUTs in the image object.

**DX image "for presentation" and/or MG image "for presentation", both with or without GSPS:**

- 8 bit P-Value
- 12 bit P-Value
- 15-bit Image + VOI Lut P-Value

Note: '15-bit Image VOI Lut P-Value' is the AGFA recommended post-processing when it is supported by the SCP.

**DX image "for processing" or MG image "for processing"**

→ this is always without GSPS, as it is the raw image.

→ the output type cannot be chosen either.



**Table 2.4-2: Overview of the PACS DICOM CR-Image Types supported by NX 3.0.8950.**

|                                       | PI      | M-Lut Rescale |          |           | VOI-Lut       |           |          | P-LUT     | MU1 | MU2 |
|---------------------------------------|---------|---------------|----------|-----------|---------------|-----------|----------|-----------|-----|-----|
|                                       | (1)     | Intercept     | Slope    | Type      | Explanation   | WinCenter | WinWidth | LUT-Shape | (2) | (3) |
| Description                           | 28,0004 | 28,1052       | 28,1053  | 28,1054   | 28,3003       | 28,1050   | 28,1051  | 2050,0020 |     |     |
| 8-bit OD                              | M1      | 200           | 10.9804  | OD        | -             | 1600      | 2800     | -         | X   | X   |
|                                       | M2      | 3000          | -10.9804 | OD        | -             | 1600      | 2800     | -         | X   | X   |
| 8-bit ODGAMMA                         | M1      | 200           | 10.9804  | OD        | -             | 1600      | 2800     | -         | X   | X   |
|                                       | M2      | 3000          | -10.9804 | OD        | -             | 1600      | 2800     | -         | X   | X   |
| 8-bit P-Value                         | M1      | 0.0           | 1.0      | P-VALUES  | -             | 128       | 256      | INVERSE   | X   | X   |
|                                       | M2      | 0.0           | 1.0      | P-VALUES  | -             | 128       | 256      | IDENTITY  | X   | X   |
| 12-bit OD                             | M1      | 200           | 0.684    | OD        | -             | 1600      | 2800     | -         | X   | X   |
|                                       | M2      | 3000          | -0.684   | OD        | -             | 1600      | 2800     | -         | X   | X   |
| 12-bit ODGAMMA                        | M1      | 200           | 0.684    | OD        | -             | 1600      | 2800     | -         | X   | X   |
|                                       | M2      | 3000          | -0.684   | OD        | -             | 1600      | 2800     | -         | X   | X   |
| 12-bit P-Value                        | M1      | 0.0           | 1.0      | P-VALUES  | -             | 2048      | 4096     | INVERSE   | X   | X   |
|                                       | M2      | 0.0           | 1.0      | P-VALUES  | -             | 2048      | 4096     | IDENTITY  | X   | X   |
| 12-bit OD REL                         | M1      | 0.0           | 1.0      | OD REL    | -             | 2048      | 4096     | -         | X   | X   |
|                                       | M2      | 0.0           | 1.0      | OD REL    | -             | 2048      | 4096     | -         | X   | X   |
| 12-bit Image --> VOI Lut OD REL (12)  | M1      | 0.0           | 1.0      | LOG_E REL | E25, NK5. ... | - (7)     | - (7)    | -         | X   |     |
|                                       | M2      | 0.0           | 1.0      | LOG_E REL | E25, NK5. ... | - (7)     | - (7)    | -         | X   |     |
| 12-bit image --> VOI Lut ODGamma (12) | M1      | 0.0           | 1.0      | LOG_E REL | E25, NK5. ... | - (7)     | - (7)    | -         | X   |     |
|                                       | M2      | 0.0           | 1.0      | LOG_E REL | E25, NK5. ... | - (7)     | - (7)    | -         | X   |     |
| 15-bit Image --> VOI Lut P-Value (15) | M1      | 0.0           | 1.0      | P-VALUES  | E25, NK5. ... | - (7)     | - (7)    | INVERSE   | X   |     |
|                                       | M2      | 0.0           | 1.0      | P-VALUES  | E25, NK5. ... | - (7)     | - (7)    | IDENTITY  | X   |     |
| 12-bit Image --> VOI Lut OD REL (12)  | M1      | 0.0           | 1.0      | LOG_E REL | -             | - (7)     | - (7)    | -         |     | X   |
|                                       | M2      | 0.0           | 1.0      | LOG_E REL | -             | - (7)     | - (7)    | -         |     | X   |
| 12-bit Image --> VOI Lut ODGamma (12) | M1      | 0.0           | 1.0      | LOG_E REL | -             | - (7)     | - (7)    | -         |     | X   |
|                                       | M2      | 0.0           | 1.0      | LOG_E REL | -             | - (7)     | - (7)    | -         |     | X   |
| 15-bit Image --> VOI Lut P-Value Lin  | M1      | 0.0           | 1.0      | P-VALUES  | -             | WC (9)    | WW (9)   | INVERSE   |     | X   |
|                                       | M2      | 0.0           | 1.0      | P-VALUES  | -             | WC (9)    | WW (9)   | IDENTITY  |     | X   |
| 15-bit Image --> VOI Lut P-Value (15) | M1      | 0.0           | 1.0      | P-VALUES  | -             | WC (7)    | WW (7)   | INVERSE   |     | X   |
|                                       | M2      | 0.0           | 1.0      | P-VALUES  | -             | WC (7)    | WW (7)   | IDENTITY  |     | X   |

(1) Photometric Interpretation: MONOCHROME1 (M1) or MONOCHROME2 (M2)

(2) Formats supported by Musica1 processing

(3) Formats supported by Musica2 processing

(7) WinCenter/WinWidth are not present. Instead LUT Data (0028,3006) in VOI LUT sequence is present.

(9) Values are in pixel values (e.g. 16384, 16384). LUT Data (0028,3006) is not present (no VOI LUT sequence).

(12) After application of the VOI Lut, the bit depth will be equal to the bit depth before application of the VOI Lut, i.e. 12 bit

(15) After application of the VOI Lut, the bit depth will be equal to the bit depth before application of the VOI Lut, i.e. 15 bit

(15) '15-bit Image – VOI Lut P-Value' is the AGFA recommended post-processing (when supported by the SCP)

The above table should be read as follows:

- This table contains the properties for each DICOM CR-Image Type. The relevant DICOM tags are shown in the column headers and the values that define a specific DICOM CR-Image Type can be found below the corresponding tag.
- The PACS DICOM CR-Image Types supported by NX 3.0.8950 are marked with an "X".

### 2.4.1.2.2 RIS

The "Workstation Service & Configuration Tool" must be used to set the AE Title, port-number, host-name and capabilities of the remote Modality Worklist SCP.

Only one single remote Modality Worklist SCP can be defined at the same time.

#### 2.4.1.2.2.1 Configuration of a RIS

##### 2.4.1.2.2.1.1 Query keys

Query keys are used when the RIS connection is configured for DICOM Modality Worklist. In the configuration Tool, the user is allowed to enter and/or modify certain query keys. Based on these keys, the worklist will be populated at the next RIS query.

See Table 2.2-45.

##### 2.4.1.2.2.1.2 Protocol Codes

###### Note on Japanese Protocol Codes

In Japan, the use of protocol codes is slightly different. See the guidelines of the JIRA and JAHIS [JAPAN].

Each Scheduled *Protocol Code Sequence* consists of the following information:

**Table 2.4-3: Information in a Japanese Scheduled Protocol Code Sequence.**

| Item              | Coding Scheme Designator | No. of items | Code Value/Code Meaning               |
|-------------------|--------------------------|--------------|---------------------------------------|
| Procedure content | JJ1017T                  | 1            | Specified in Section 5.3 <sup>9</sup> |
| Target region     | JJ1017P                  | 1            | Specified in Section 5.4 <sup>9</sup> |
| Imaging direction | JJ1017D                  | 0-N          | Specified in Section 5.5 <sup>9</sup> |

(Table copied from [JAPAN], p.9)

Each item mentioned in the table is wrapped in a scheduled protocol code sequence item. As shown in the table, for each item a different coding scheme designator is used.

In practice, a triplet (with the imaging direction possibly omitted) refers to a single Exposure Type in the NX 3.0.8950 exposure tree.

###### Example of a Japanese SPS

Example of a Scheduled Protocol Code Sequence in a Japanese SPS:

| Coding Scheme Designator | Coding Scheme Version | Code Value | Code Meaning                            |
|--------------------------|-----------------------|------------|---|
| JJ1017T                  | 1.0                   | GX.01.00   | Radiography – General Radiography – NOS |
| JJ1017P                  | 1.0                   | 25.6.201   | Chest – Respiratory System – Lung       |
| JJ1017D                  | 1.0                   | G-5200     | Anterior to posterior                   |
| JJ1017D                  | 1.0                   | G-A101     | Left lateral                            |

This SPS schedules 2 exposures: Chest AP and Chest LL

NX 3.0.8950 can be configured for the use of the Japanese system of Protocol Codes (JJ1017) by means of the "Workstation Service & Configuration Tool".

<sup>9</sup> See guidelines of the JIRA and JAHIS [JAPAN].

### 2.4.1.2.2 Configuration of RIS Mapping

RIS mapping defines how incoming worklist attributes from the RIS are mapped to the SPS elements.

With RIS mapping, the incoming SPS attributes of the RIS are mapped to the internal data structure of NX 3.0.8950. Standard the system will provide a default (1:1) mapping. However in cases where the RIS sends out the RIS data in a non DICOM conformant way, the customization of the mapping needs to be done by the operator. This can be done in the "Workstation Service & Configuration Tool".

### 2.4.1.2.3 MPPS

MPPS reporting can be enabled/ disabled by the operator.

The "Workstation Service & Configuration Tool" must be used to set the AE Title, port-number, host-name and capabilities of the remote MPPS SCP. Only one single remote MPPS SCP can be defined at a time.

### 2.4.1.2.4 Printing

The "Workstation Service & Configuration Tool" must be used to set the AE Titles, port-numbers, host-names and capabilities for the printers. Only the supported printers can be configured. Multiple printers can be defined.

The DICOM attributes per printer that are configurable are defined in Table 2.2-24, Table 2.2-27 and Table 2.2-29.

Additionally, the following parameters can be configured:

**Table 2.4-4: Configurable Printer parameters.**

| Parameter    | Configurable (yes/no) |
|--------------|-----------------------|
| Name         | Yes                   |
| Description  | Yes                   |
| Use N-Events | Yes                   |
| IP Address   | Yes                   |
| Port Number  | Yes                   |
| AE Title     | Yes                   |
| SSL Enabled  | Yes                   |

## 2.4.2 Queue Management & Job description

### 2.4.2.1 DICOM Store AE

#### 2.4.2.1.1 Queue description

Each Remote Store AE destination has its own configuration (as described in § [2.4.1.2.1](#)) and queue. For each destination, a different retry policy can be configured.

It is possible to configure the rerouting of a DICOM Store queue. A typical example for when this is needed is a PACS that is temporarily down. The queue can be rerouted to a web server. This web server will then later on forward the images to the PACS. A DICOM Store queue can only be rerouted to another DICOM Store queue with:

- The same output type configured (this includes the SOP class to be used)
- The same GSPS support

If the new destination is configured for Storage Commitment, the storage commit request will be sent to the new destination when processing the job.

By means of a SENT-flag, an image that has already been successfully sent to a specific archive cannot be sent to that archive twice. The SENT-flag can be applied to all DICOM store destinations, including both the archive and the other softcopy destinations (e.g. viewing stations).

#### 2.4.2.1.2 Job description

One archive job can contain multiple images. All these images and their GSPS's will be sent to the archive through a single association.

The user can perform the following actions on existing jobs: "Delete job" and "Expedite job".

When an image in the archive job does not reach its destination, or when the storage commit replies "time out", the job is FAILED.

#### 2.4.2.2 DICOM Print AE

There is only one entry per print job (with possible multiple sheets) in the queue.

All queue actions are done on one print job: retry, reprint, expedite ...

#### 2.4.3 Parameters

The specification of important operational parameters, and if configurable, their default value and range, are specified in the table below.

The parameters that apply to all Application Entities are specified in the "General Parameters" section. Those specific to particular Applications are specified in separate sections specific to each AE.

**Table 2.4-5: Configurable Parameters for NX 3.0.8950.**

| Parameter   | Configurable (yes/no) | Default value |
|---|-----------------------|---------------|
| <b>General Parameters</b>   |                       |               |
| Max PDU Receive Size  | No                    | 65542         |
| Max PDU Send Size (larger PDUs will never be sent, even if the receiver supports a larger Max PDU Receive Size. If the receiver supports a smaller Max PDU Receive Size then the Max PDU Send Size will be reduced accordingly for the duration of the Association. Max PDU Receive Size information is exchanged during DICOM Association Negotiation in the Maximum Length Sub-Item of the A-ASSOCIATION-RQ and A-ASSOCIATE-AC) | No                    | 65542         |
| Time-out waiting for an acceptance or rejection response to an Association Request (Application Level Timeout)  | No                    | 15 minutes    |
| Time-out waiting for a response to an Association release request (Application Level Timeout)   | No                    | 15 minutes    |
| Time-out waiting for completion of a TCP/IP connect request (Low-level timeout)   | No                    | 100 seconds   |
| Time-out awaiting a Response to a DIMSE Request (Low-Level Timeout)   | No                    | 100 seconds   |

| Parameter   | Configurable (yes/no) | Default value  |
|---|-----------------------|--|
| Time-out for waiting for data between TCP/IP-packets (Low Level Timeout)  | No                    | 100 seconds  |
| <b>Storage parameters</b>   |                       |  |
| Storage SCU time-out waiting for a response to a C-STORE-RQ   | No                    | 15 minutes   |
| Number of times a failed send job may be retried  | No                    | Manually: there's no restriction of the number of retries of failed jobs as long as they aren't deleted. |
| Delay between retrying failed send jobs   | No                    | N/A  |
| Supported Transfer Syntaxes (separately configurable for each remote AE)  | Yes                   |  |
| Secure DICOM (SSL)  | Yes                   |  |
| Is Archive? <sup>10</sup>   | Yes                   |  |
| <b>Storage Commit Parameters</b>  |                       |  |
| Enable Storage Commit   | Yes                   |  |
| Enable secure DICOM Connection  | Yes                   |  |
| Timeout waiting for a Storage Commitment Notification (maximum duration of applicability for a Storage Commitment Transaction UID). | Yes                   |  |
| <b>Modality Worklist Parameters</b>   |                       |  |
| Modality Worklist SCU time-out waiting for the final response to a C-FIND-RQ  | No                    | 15 minutes   |
| Maximum number of Worklist Items  | No                    | 400  |
| Supported Transfer Syntaxes for Modality Worklist   | No                    | Implicit VR Little Endian<br>Explicit VR Little Endian   |
| Delay between automatic Worklist Updates  | Yes                   |  |
| <b>MPPS Parameters</b>  |                       |  |
| Enable MPPS Reporting   | Yes                   |  |
| MPPS SCU time-out waiting for a response to a N-CREATE-RQ   | No                    | 15 minutes   |
| MPPS SCU time-out waiting for a response to a N-SET-RQ  | No                    | 15 minutes   |
| Supported Transfer Syntaxes for MPPS  | No                    | Implicit VR Little Endian<br>Explicit VR Little Endian   |
| <b>Print Parameters</b>   |                       |  |
| Print SCU time-out waiting for a response to a N-CREATE-RQ  | No                    | 10 minutes   |
| Print SCU time-out waiting for a response to a N-SET-RQ   | No                    | 10 minutes   |
| Print SCU time-out waiting for a response to a N-ACTION-RQ  | No                    | 10 minutes   |
| Supported Transfer Syntaxes (separately configurable for each remote printer)   | No                    | Taken from printer device model provided by Agfa.  |
| Number of times a failed print-job may be retried   | No                    | Automatically : 8 times  |

<sup>10</sup> If yes, images sent to this destination will be referenced in the MPPS.


| Parameter  | Configurable (yes/no) | Default value   |
|--|-----------------------|---|
| Delay between retrying failed print-jobs                   | No                    | 2 times immediately, then the job gets parked for 3 minutes and retried 3 times again. In case the job still fails, it is parked for 5 minutes and retried 3 times. |
| Print SCU time-out waiting for a response to a N-CREATE-RQ | No                    | 10 minutes  |

## 3 MEDIA INTERCHANGE

NX 3.0.8950 is able to create or read DICOM Interchange media. The related capabilities are described in the following sections.

### 3.1 Implementation Model

#### 3.1.1 Application Data Flow Diagram



**Figure 3.1-1: Application Data Flow Diagram for Media Storage**

The Offline-Media Application Entity exports DICOM images and Presentation States to a CD-R Storage medium. It is associated with the local real-world activity “Export Images”. “Export Images” is performed upon user request.

#### 3.1.2 Functional Definition of AEs

#### 3.1.3 Sequencing of Real World Activities

At least one image or presentation state must exist and be selected before the Offline-Media Application Entity can be invoked. The Offline-Media Application Entity is invoked through the local real-world activity “Export Images”.

The operator must insert a new (blank) CD-R media before invocation of the Offline-Media Application Entity. If no CD-R is inserted, the Offline-Media Application Entity will prompt for a media to be inserted before starting to write to the CD-R device. The export job can be canceled from the job queue by clicking “Cancel” on this prompt.

Processed Images are exported in standard DICOM format. A Dicom Viewer is burned with the images on the CD, to view these images later on. The images can be read by any Dicom Compatible application.

Export to a hard disk is NOT supported by NX 3.0.8950.

### 3.1.4 File Meta Information for Implementation Class and Version

The implementation information written to the File Meta Header in each file is:

**Table 3.1-1: File Meta Implementation Class and Version Name**

|                             |                              |
|-----------------------------|------------------------------|
| Implementation Class UID    | 1.3.51.0.1.3                 |
| Implementation Version Name | AGFA DTF1.0.XX <sup>11</sup> |

## 3.2 AE Specifications

### 3.2.1 Offline-Media Application Entity Specification

The Offline-Media Application Entity provides standard conformance to the DICOM Interchange Option of the Media Storage Service Class. The Application Profiles and roles are listed below:

**Table 3.2-1: AE Related Application Profiles, Real World Activities and Roles**

| Application Profile Supported | Real-World Activity | Roles                           | SC Option   |
|-------------------------------|---------------------|---------------------------------|-------------|
| Portable data for imaging     | Export to CD-R      | FSC<br>(creation of a File-Set) | Interchange |

#### 3.2.1.1 Real World Activities

##### 3.2.1.1.1 Real World Activity - Export to CD-R

The Offline-Media Application Entity acts as an FSC using the interchange option when requested to export SOP Instances from the local database to a CD-R medium. If the current contents selection does not fit on a single media, an error message is displayed and the export will fail. The user will be prompted to insert an empty CD-R for each export job. The contents of the export job will be written together with a corresponding DICOMDIR, which is placed in the root directory, to a single-session CD-R. Writing in multi-session mode is supported. The user can cancel an export job in the job queue.

The file names and directory names may reflect the patient name or ID if necessary, but in case the operator selects “anonymous export” these names and IDs should not refer to the actual patient. The aim of anonymous export is that any data, which could be used to derive the identity of the patient, is removed from the image. This includes all dates, UID and other ID's, patient related information and all related physician's names. All attributes from the following modules are therefore blanked:

- Patient
- Patient Study
- Patient Medical
- General Study
- General Series
- DX Series
- SOP Common
- General Image

<sup>11</sup> XX is the build version number.



Exceptions exist for attributes which are necessary to correctly import or interpret the image. Therefore the following values are kept or replaced by generated values:

**Table 3.2-2: Necessary Values when performing an Anonymous Export.**

| Attribute                      | Tag         | Set to ...  |
|--------------------------------|-------------|-------------|
| Patient's Name                 | (0010,0010) | "Anonymous" |
| Study Instance UID             | (0020,000D) | Keep value  |
| Modality                       | (0008,0060) | Keep value  |
| Series Instance UID            | (0020,000E) | Keep value  |
| Laterality                     | (0020,0060) | Keep value  |
| Series description             | (0008,103E) | Keep value  |
| Protocol name                  | (0018,1030) | Keep value  |
| Body Part Examined             | (0018,0015) | Keep value  |
| Smallest Pixel Value in Series | (0028,0108) | Keep value  |
| Largest Pixel Value in Series  | (0028,0109) | Keep value  |
| Presentation Intent Type       | (0008,0068) | Keep value  |
| SOP Instance UID               | (0008,0018) | Keep value  |
| SOP Class UID                  | (0008,0016) | Keep value  |
| Specific character set         | (0008,0005) | Keep value  |
| Patient Orientation            | (0020,0020) | Keep value  |
| Image type                     | (0008,0008) | Keep value  |
| Image comments                 | (0020,4000) | Keep value  |
| Presentation LUT shape         | (2050,0020) | Keep value  |

### 3.2.1.1.1.1 Media Storage Application Profile

The Offline-Media Application Entity supports the STD-GEN-CD Application Profile.

#### 3.2.1.1.1.1.1 Options

The Offline-Media Application Entity supports the SOP Classes and Transfer Syntaxes listed in the Table below:

**Table 3.2-3: IOD'S, SOP Classes and Transfer syntaxes for offline media**

| IOD  | SOP Class UID                 | Transfer Syntax           | Transfer Syntax UID |
|--|-------------------------------|---------------------------|---------------------|
| CR Image Storage                                     | 1.2.840.10008.5.1.4.1.1.1     | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |
| Digital X-Ray Image Storage – For Presentation       | 1.2.840.10008.5.1.4.1.1.1.1   | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |
| Digital X-Ray Image Storage – For Processing         | 1.2.840.10008.5.1.4.1.1.1.1.1 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |
| Grayscale Softcopy Presentation State Storage        | 1.2.840.10008.5.1.4.1.1.11.1  | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |
| Digital Mammography Image Storage – For Presentation | 1.2.840.10008.5.1.4.1.1.1.2   | Implicit VR Little Endian | 1.2.840.10008.1.2   |
| Digital Mammography Image Storage – For Processing   | 1.2.840.10008.5.1.4.1.1.1.2.1 | Implicit VR Little Endian | 1.2.840.10008.1.2   |

### 3.3 Media Configuration

The Source AE Title in the File Meta Information is the AE Title that can be configured in the General Settings in the "Workstation Service & Configuration Tool".

## 4 SUPPORT FOR EXTENDED CHARACTER SETS

NX 3.0.8950 supports the following character sets:

|                                    |                           |
|------------------------------------|---------------------------|
| • ISO-IR 6 (default)               | Basic G0 Set              |
| • ISO-IR 100                       | Latin Alphabet No. 1      |
| • ISO-IR 101                       | Latin Alphabet No. 2      |
| • ISO-IR 109                       | Latin Alphabet No. 3      |
| • ISO-IR 110                       | Latin Alphabet No. 4      |
| • ISO-IR 148                       | Latin Alphabet No. 5      |
| • ISO-IR 126                       | Greek                     |
| • ISO-IR 144                       | Cyrillic                  |
| • ISO-IR 127                       | Arabic                    |
| • ISO-IR 13                        | Japanese                  |
| • ISO-IR 138                       | Hebrew                    |
| <b>Chinese Character Sets</b>      |                           |
| • GB18030                          | GB18030-2000              |
| • ISO 2022 B5                      | Big 5                     |
| • ISO 2022 GBK                     | GB 2313-80                |
| <b>Japanese Character Sets</b>     |                           |
| • ISO 2022 IR 87                   | ASCII and Kanji           |
| • ISO 2022 IR 13<br>ISO 2022 IR 87 | Katakana, Roman and Kanji |
| • ISO 2022 IR 159                  | Japanese suppl. Kanji     |
| <b>Other Character Sets</b>        |                           |
| • ISO 2022 IR 126                  | Greek                     |
| • ISO 2022 IR 144                  | Cyrillic                  |
| • ISO 2022 IR 138                  | Hebrew                    |
| • ISO 2022 IR 149                  | Korean                    |
| • ISO 2022 IR 127                  | Arabic                    |
| • ISO 2022 IR 100                  | Latin alphabet No 1       |
| • ISO 2022 IR 101                  | Latin alphabet No 2       |
| • ISO 2022 IR 109                  | Latin alphabet No 3       |
| • ISO 2022 IR 110                  | Latin alphabet No 4       |
| • ISO 2022 IR 148                  | Latin alphabet No 5       |
| • UTF-8                            | Unicode                   |

## 5 SECURITY

### 5.1 Security Profiles

NX 3.0.8950 supports the following profiles:

- Basic TLS Secure Transport Connection Profile
- Basic Network Address Management Profile
- Basic Time synchronization Profile

### 5.2 Association Level Security

NX 3.0.8950 supports secure associations using TLS. Associations are only allowed to be opened when they are received from an AE title that is known to the Workstation (i.e. the specific AE title has been configured in the "Workstation Service & Configuration Tool").

### 5.3 Application Level Security

NX 3.0.8950 supports application level security by means of role based access control. These user roles are associated to the Windows OS users. In the "Workstation Service & Configuration Tool", user roles can be assigned to certain accounts. Some of the user roles can be modified or new ones can be created by a user with the appropriate rights.

By default the workstation will be configured to automatically log the current user out after a certain period of non-activity.

Audit logging to an external repository is supported.

NX 3.0.8950 can be configured to get its time from an NTP-server.

## 6 ANNEXES

### 6.1 IOD Contents

#### 6.1.1 Created SOP Instances

The following tables use a number of abbreviations. The abbreviations used in the “Presence of Value”- column, are:

- VNAP Value Not Always Present (attribute sent zero length if no value is present)
- ANAP Attribute Not Always Present
- ALWAYS Always Present with a value
- EMPTY Attribute is sent without a value

The abbreviations used in the “Source”-column:

- USER the attribute value source is from User input
- AUTO the attribute value is generated automatically
- MWL the attribute value is the value received from the Modality Worklist
- CONFIG the attribute value source is a configurable parameter

#### 6.1.1.1 Common Modules

**Table 6.1-1: Common Modules**

| Attribute Name   | Tag         | VR | Value | Presence of Value | Source     |
|--|-------------|----|-------|-------------------|------------|
| <b>Patient</b>   |             |    |       |                   |            |
| Patient Identification                                 |             |    |       |                   |            |
| Patient's Name   | (0010,0010) | PN |       | VNAP              | User / MWL |
| Patient ID   | (0010,0020) | LO |       | VNAP              | User / MWL |
| Issuer of Patient ID                                   | (0010,0021) | LO |       | ANAP              | MWL        |
| Other Patient IDs                                      | (0010,1000) | LO |       | ANAP              | MWL / MWL  |
| Other Patient Names                                    | (0010,1001) | PN |       | ANAP              | MWL        |
| Patient's Birth Name                                   | (0010,1005) | PN |       | ANAP              | MWL        |
| Patient's Mother's Birth Name                          | (0010,1060) | PN |       | ANAP              | MWL        |
| Medical Record Locator                                 | (0010,1090) | LO |       | ANAP              | MWL        |
| Patient Demographic                                    |             |    |       |                   |            |
| Patient's Age  | (0010,1010) | AS |       | ANAP              | User / MWL |
| Occupation   | (0010,2180) | SH |       | ANAP              | MWL        |
| Confidentiality Constraint on Patient Data Description | (0040,3001) | LO |       | ANAP              | MWL        |
| Patient's Birth Date                                   | (0010,0030) | DA |       | VNAP              | User / MWL |
| Patient's Birth Time                                   | (0010,0032) | TM |       | VNAP              | User / MWL |
| Patient's Sex  | (0010,0040) | CS |       | VNAP              | User / MWL |
| Patient's Insurance Plan Code Sequence                 | (0010,0050) | SQ |       | ANAP              | MWL        |
| > Code Sequence  |             |    |       |                   |            |
| > Patient's Primary Language Code                      | (0010,0102) | SQ |       | ANAP              | MWL        |

| Attribute Name                                     | Tag         | VR | Value  | Presence of Value | Source     |
|--|-------------|----|--|-------------------|------------|
| Modifier Sequence                                  |             |    |  |                   |            |
| >> Code Sequence                                   |             |    |  |                   |            |
| Patient's Size                                     | (0010,1020) | DS |  | ANAP              | User / MWL |
| Patient's Weight                                   | (0010,1030) | DS |  | ANAP              | User / MWL |
| Patient's Address                                  | (0010,1040) | LO |  | ANAP              | MWL        |
| Military Rank                                      | (0010,1080) | LO |  | ANAP              | User / MWL |
| Branch of Service                                  | (0010,1081) | LO |  | ANAP              | User / MWL |
| Country of Residence                               | (0010,2150) | LO |  | ANAP              | MWL        |
| Region of Residence                                | (0010,2152) | LO |  | ANAP              | MWL        |
| Patient's Telephone Numbers                        | (0010,2154) | SH |  | ANAP              | MWL        |
| Ethnic Group                                       | (0010,2160) | SH |  | ANAP              | User / MWL |
| Patient's Religious Preference                     | (0010,21F0) | LO |  | ANAP              | MWL        |
| Patient Comments                                   | (0010,4000) | LT |  | ANAP              | User / MWL |
| <b>Visit Status</b>                                |             |    |  |                   |            |
| Referenced Patient Sequence                        | (0008,1120) | SQ |  | ANAP              | MWL        |
| >Referenced SOP Class UID                          | (0008,1150) | UI |  | ANAP              | MWL        |
| >Referenced SOP Instance UID                       | (0008,1155) | UI |  | ANAP              | MWL        |
| Visit Status ID                                    | (0038,0008) | CS |  | ANAP              | MWL        |
| Current Patient Location                           | (0038,0300) | LO |  | ANAP              | MWL        |
| Patient's Institution Residence                    | (0038,0400) | LO |  | ANAP              | MWL        |
| Visit Comments                                     | (0038,4000) | LT |  | ANAP              | Auto       |
| <b>Patient Medical</b>                             |             |    |  |                   |            |
| Medical Alerts                                     | (0010,2000) | LO |  | ANAP              | MWL        |
| Contrast Allergies                                 | (0010,2110) | LO |  | ANAP              | MWL        |
| Smoking Status                                     | (0010,21A0) | CS |  | ANAP              | User / MWL |
| Pregnancy Status                                   | (0010,21C0) | US |  | ANAP              | User / MWL |
| Last Menstrual Date                                | (0010,21D0) | DA |  | ANAP              | User / MWL |
| Special Needs                                      | (0038,0050) | LO |  | ANAP              | MWL        |
| Patient State                                      | (0038,0500) | LO |  | ANAP              | MWL        |
| Additional Patient's History                       | (0010,21B0) | LT |  | ANAP              | MWL        |
| <b>Study</b>                                       |             |    |  |                   |            |
| <b>General Study</b>                               |             |    |  |                   |            |
| Study Instance UID                                 | (0020,000D) | UI |  | ALWAYS            | MWL / Auto |
| Study Date   | (0008,0020) | DA |  | ALWAYS            | Auto       |
| Study Time   | (0008,0030) | TM |  | ALWAYS            | Auto       |
| Referring Physician's Name                         | (0008,0090) | PN |  | VNAP              | User / MWL |
| Referring Physician Identification Sequence        | (0008,0096) | SQ |  | ANAP              | MWL        |
| Study ID   | (0020,0010) | SH | Requested Procedure ID or generated for unscheduled (app A IHE note 5) | ALWAYS            | Auto / MWL |
| Accession Number                                   | (0008,0050) | SH |  | VNAP              | User / MWL |
| Study Description                                  | (0008,1030) | LO | Performed Procedure step Description                                   | ANAP              | Auto       |
| Physician(s) of Record                             | (0008,1048) | PN |  | ANAP              | MWL        |
| Physician(s) of Record Identification Sequence     | (0008,1049) | SQ |  | ANAP              | MWL        |
| Name of Physician(s) Reading Study                 | (0008,1060) | PN |  | ANAP              | MWL        |
| Physician(s) Reading Study Identification Sequence | (0008,1062) |    |  | ANAP              | User       |

| Attribute Name                                 | Tag         | VR | Value  | Presence of Value | Source                                    |
|--|-------------|----|--|-------------------|---|
| Referenced Patient Sequence                    | (0008,1120) | SQ |  | ANAP              | MWL                                       |
| >Referenced SOP Class UID                      | (0008,1150) | UI |  | ANAP              | MWL                                       |
| >Referenced SOP Instance UID                   | (0008,1155) | UI |  | ANAP              | MWL                                       |
| Procedure Code Sequence                        | (0008,1032) | SQ |  | ANAP              | MWL                                       |
| <b>Patient Study</b>                           |             |    |  |                   |   |
| Admitting Diagnoses Description                | (0008,1080) | LO |  | ANAP              | MWL                                       |
| Admitting Diagnoses Code Sequence              | (0008,1084) | SQ |  | ANAP              | MWL                                       |
| Patient's Age                                  | (0010,1010) | AS |  | ANAP              | User / MWL                                |
| Patient's Size                                 | (0010,1020) | DS |  | ANAP              | User / MWL                                |
| Patient's Weight                               | (0010,1030) | DS |  | ANAP              | User / MWL                                |
| Occupation                                     | (0010,2480) | SH |  | ANAP              | MWL                                       |
| Additional Patient's History                   | (0010,21B0) | LT |  | ANAP              | MWL                                       |
| <b>Series</b>                                  |             |    |  |                   |   |
| <b>General Series</b>                          |             |    |  |                   |   |
| Modality                                       | (0008,0060) | CS |  | ALWAYS            | Config,<br>>Archive settings<br>SOP Class |
| Series Instance UID                            | (0020,000E) | UI | Different for each image                                       | ALWAYS            | Auto                                      |
| Series Number                                  | (0020,0011) | IS |  | ALWAYS            | Auto                                      |
| Laterality                                     | (0020,0060) | CS | Based on protocol code   | ALWAYS            | User / Auto                               |
| Performing Physicians' Name                    | (0008,1050) | PN |  | ANAP              | User                                      |
| Performing Physician's Identification Sequence | (0008,1050) | SQ |  | ANAP              | MWL                                       |
| Protocol Name                                  | (0018,1030) | LO |  | ANAP              | Auto                                      |
| Series Description                             | (0008,103E) | LO | Exposure Type Name   | ALWAYS            | Auto                                      |
| Operators' Name                                | (0008,1070) | PN | Depending on security setting either auto filled in with login | ALWAYS            | User / Auto                               |
| Referenced Performed Procedure Step Sequence   | (0008,1111) | SQ |  | ANAP              | Auto                                      |
| >Referenced SOP Class UID                      | (0008,1150) | UI |  | ANAP              | Auto                                      |
| >Referenced SOP Instance UID                   | (0008,1155) | UI |  | ANAP              | Auto                                      |
| Patient Position                               | (0018,5100) | CS | Empty  | ANAP              |   |
| Smallest Pixel Value in Series                 | (0028,0108) | US |  | NEVER             |   |
| Largest Pixel Value in Series                  | (0028,0109) | US |  | NEVER             |   |
| Request Attributes Sequence                    | (0040,0275) | SQ |  | ANAP              | MWL                                       |
| >Requested Procedure ID                        | (0040,1001) | SH |  | ANAP              | MWL                                       |
| >Reason for the Requested Procedure            | (0040,1002) | LO |  |                   |   |
| >Reason for Requested Procedure Code Sequence  | (0040,100A) | SQ | Screening/Diagnostic   |                   |   |
| >Scheduled Procedure Step ID                   | (0040,0009) | SH |  | ANAP              | MWL                                       |
| >Scheduled Procedure Step Description          | (0040,0007) | LO |  | ANAP              | MWL                                       |
| >Scheduled Protocol Code Sequence              | (0040,0008) | SQ | <sup>12</sup>  | ANAP              | MWL                                       |

<sup>12</sup> The following rules are applicable for DICOM:

- When 0008,102 was not supplied by the RIS, its value will be "UNKNOWN".
- When 0008,103 was not supplied by the RIS, its value will be "UNKNOWN".
- When 0008,104 was not supplied by the RIS, its value will be "UNKNOWN".

| Attribute Name                           | Tag         | VR | Value | Presence of Value                          | Source                |
|--|-------------|----|-------|--|-----------------------|
| >> Protocol Context Sequence             | (0040,0440) | SQ |       |  |                       |
| >>>Content Item Modifier Sequence        | (0040,0441) | SQ |       |  |                       |
| Requested Procedure Description          | (0032,1060) | LO |       | ANAP                                       | Manual / MWL          |
| Performed Procedure Step ID              | (0040,0253) | SH |       | ANAP                                       | Auto                  |
| Performed Procedure Step Description     | (0040,0254) | LO |       | Exam group names of exposures in the study | Auto                  |
| Performed Procedure Step Start Date      | (0040,0244) | DA |       |  | Auto                  |
| Performed Procedure Step Start Time      | (0040,0245) | TM |       |  | Auto                  |
| Comments on the Performed Procedure Step | (0040,0280) | LO |       | ANAP                                       | Auto                  |
| <b>Equipment</b>                         |             |    |       |  |                       |
| General Equipment                        |             |    |       |  |                       |
| Station Name                             | (0008,1010) | SH |       |  | Auto                  |
| Institution Name                         | (0008,0080) | LO |       | ANAP                                       | MWL/ Config           |
| Institution Address                      | (0008,0081) | ST |       | ANAP                                       | MWL /Config           |
| Institutional Department Name            | (0008,1040) | LO |       | ANAP                                       | MWL/ Config           |
| Manufacturer                             | (0008,0070) | LO |       | ALWAYS                                     | Auto (from digitizer) |
| Manufacturer's Model Name                | (0008,1090) | LO |       | ALWAYS                                     | Auto (from digitizer) |
| Device Serial Number                     | (0018,1000) | LO |       | ALWAYS                                     | Auto                  |
| Gantry ID                                | (0018,1008) | LO |       | ALWAYS                                     | Auto                  |
| Software Versions                        | (0018,1020) | LO |       | ANAP                                       | Auto                  |
| <b>Image</b>                             |             |    |       |  |                       |
| General Image                            |             |    |       |  |                       |
| Instance Number                          | (0020,0013) | CS |       | ALWAYS                                     | Auto                  |
| Patient Orientation                      | (0020,0020) | IS |       | ALWAYS                                     | Auto                  |
| Content Date                             | (0008,0023) | DT |       | VNAP                                       | Auto                  |
| Content Time                             | (0008,0033) | TM |       | VNAP                                       | Auto                  |
| Image Type                               | (0008,0008) | CS |       | ALWAYS                                     | Auto                  |
| Acquisition Number                       | (0020,0012) | IS |       | ANAP                                       | Auto                  |
| Acquisition Date                         | (0008,0022) | DA |       | ALWAYS                                     | Auto / User           |
| Acquisition Time                         | (0008,0032) | TM |       | ALWAYS                                     | Auto / User           |
| Acquisition Datetime                     | (0008,002A) | DT |       | ALWAYS                                     | Auto                  |
| Referenced Image Sequence                | (0008,1140) | SQ |       | ANAP                                       |                       |
| > Referenced SOP Class UID               | (0008,1150) | UI |       | ANAP                                       |                       |
| > Referenced SOP Instance UID            | (0008,1155) | UI |       | ANAP                                       |                       |
| > Referenced Frame Number                | (0008,1160) | IS |       | ANAP                                       |                       |
| > Purpose of Reference Code Sequence     | (0040,A170) | SQ |       | ANAP                                       |                       |
| Derivation description                   | (0008,2111) | ST |       | ANAP                                       | Auto                  |
| Source Image Sequence                    | (0008,2112) | SQ |       | ANAP                                       | Auto                  |
| > Referenced SOP Class UID               | (0008,1150) | UI |       |  |                       |
| > Referenced SOP Instance UID            | (0008,1155) | UI |       |  |                       |
| > Referenced Frame Number                | (0008,1160) | IS |       |  |                       |
| > Purpose of Reference Code Sequence     | (0040,A170) | SQ |       |  |                       |
| > Purpose of Reference Code Sequence     | (0040,A170) | SQ |       | ANAP                                       | Auto                  |



| Attribute Name                           | Tag         | VR | Value  | Presence of Value    | Source |
|--|-------------|----|--|----------------------|--------|
| Image Comments                           | (0020,4000) | LT |  | ANAP                 | User   |
| Quality Control Image                    | (0028,0300) | CS | When exposure is defined as QC in study tree   | ANAP                 | Auto   |
| Burned in Annotation                     | (0028,0301) | CS | No   | ALWAYS               | Fixed  |
| Lossy Image Compression                  | (0028,2110) | CS |  | EMPTY                | Empty  |
| Lossy Image Compression Ratio            | (0028,2112) | DS |  | ANAP                 | Auto   |
| Lossy Image Compression Method           | (0028,2114) | CS |  | ANAP                 | Auto   |
| Presentation LUT Shape                   | (2050,0020) | CS |  | ALWAYS               | Config |
| Irradiation Event UID                    | (0008,3010) | CS |  | ALWAYS               | Auto   |
| <b>Image Pixel</b>                       |             |    |  |                      |        |
| Samples per Pixel                        | (0028,0002) | US | 1  | ALWAYS               | Fixed  |
| Photometric Interpretation <sup>13</sup> | (0028,0004) | CS | M1 or M2   | ALWAYS               | Auto   |
| Rows                                     | (0028,0010) | US |  | ALWAYS               | Auto   |
| Columns                                  | (0028,0011) | US |  | ALWAYS               | Auto   |
| Pixel Aspect Ratio                       | (0028,0034) | IS |  | ANAP                 | Auto   |
| Bits Allocated                           | (0028,0100) | US |  | ALWAYS               | Auto   |
| Bits Stored                              | (0028,0101) | US |  | ALWAYS               | Auto   |
| High Bit                                 | (0028,0102) | US |  | ALWAYS               | Auto   |
| Pixel Representation <sup>14</sup>       | (0028,0103) | US |  | ALWAYS               | Auto   |
| Pixel Spacing Calibration Type           | (0028,0A02) | CS | FIDUCIAL GEOMETRY  | ANAP                 | Auto   |
| Pixel Spacing Calibration Description    | (0028,0A04) | LO | <ul style="list-style-type: none"> <li>▪ In case of FLFS stitched image: "Stitching grid, 5x5cm"</li> <li>▪ In case of ERMF tool: "Manually entered distance in cm"</li> <li>▪ In case of calibration tool: "Manually entered size in cm of an object in the image"</li> </ul> | ALWAYS <sup>15</sup> | Auto   |
| Pixel Data                               | (7FE0,0010) | OB |  | ALWAYS               | Auto   |
| <b>SOP Common</b>                        |             |    |  |                      |        |
| SOP Class UID                            | (0008,0016) | UI |  | ALWAYS               | Fixed  |
| SOP Instance UID                         | (0008,0018) | UI |  | ALWAYS               | Auto   |
| Specific Character Set                   | (0008,0005) | CS |  | ANAP                 | Config |
| Instance Creation Date                   | (0008,0012) | DA |  | ALWAYS               | Auto   |
| Instance Creation Time                   | (0008,0013) | TM |  | ALWAYS               | Auto   |
| Instance Number                          | (0020,0013) | IS |  | ALWAYS               | Auto   |

<sup>13</sup> Fixed M1 when using DX for Processing<sup>14</sup> Depends on digitizer type. In case of DX for Processing, the pixel representation is fixed 0<sup>15</sup> Attribute always present when (0028,0A02) is present

**6.1.1.2 CR****6.1.1.2.1 CR Image IOD****Table 6.1-2: IOD of Created CR Image SOP Instances**

| IE        | Module                 | Reference   | Presence of Module   |
|-----------|------------------------|-------------|----------------------|
| Patient   | Patient Identification | Table 6.1-1 | ALWAYS               |
|           | Extended Patient       | Table 6.1-1 | If received from MWL |
| Study     | General Study          | Table 6.1-1 | ALWAYS               |
|           | Patient Study          | Table 6.1-1 | ALWAYS               |
| Series    | General Series         | Table 6.1-1 | ALWAYS               |
|           | CR Series              | Table 6.1-3 | ALWAYS               |
| Equipment | General Equipment      | Table 6.1-1 | ALWAYS               |
| Image     | General Image          | Table 6.1-1 | ALWAYS               |
|           | Image Pixel            | Table 6.1-1 | ALWAYS               |
|           | CR Image               | Table 6.1-3 | ALWAYS               |
|           | Modality LUT           | Table 6.1-3 | ALWAYS               |
|           | VOI LUT                | Table 6.1-3 | ALWAYS               |
|           | SOP Common             | Table 6.1-1 | ALWAYS               |

**6.1.1.2.2 CR Modules****Table 6.1-3: CR Modules of Created SOP Instances**

| Attribute Name              | Tag         | VR | Value   | Presence of Value   | Source                   |
|-----------------------------|-------------|----|---|---------------------|--------------------------|
| <b>Series</b>               |             |    |   |                     |                          |
| CR Series                   |             |    |   |                     |                          |
| View Position               | (0018,5101) | CS | Based on protocol code                        | ALWAYS              | Auto / User              |
| Focal Spot                  | (0018,1190) | DS | 0 to 1 (mm), 0 = small focus, 1 = large focus |                     | Auto/XRDI <sup>(1)</sup> |
| Filter Type                 | (0018,1160) | SH |   |                     | Auto                     |
| Body Part Examined          | (0018,0015) | CS | Based on protocol code                        | ALWAYS              | Auto / User              |
| Collimator/ Grid Name       | (0018,1180) | SH |   |                     |                          |
| Plate Type                  | (0018,1260) | SH |   | ALWAYS              | Auto                     |
| Phosphor Type               | (0018,1261) | LO |   | ANAP                | Auto                     |
| <b>Image</b>                |             |    |   |                     |                          |
| CR Image                    |             |    |   |                     |                          |
| Photometric Interpretation  | (0028,0004) | CS |   | ALWAYS              | Fixed per archive model  |
| KVP                         | (0018,0060) | DS |   | ANAP <sup>(1)</sup> | Auto                     |
| Plate ID                    | (0018,1004) | LO |   | ALWAYS              | Auto                     |
| Distance Source to Detector | (0018,1110) | DS |   | ANAP <sup>(1)</sup> | Auto                     |
| Distance Source to Patient  | (0018,1111) | DS |   | ANAP <sup>(1)</sup> | Auto                     |
| Exposure Time               | (0018,1150) | IS |   | ANAP <sup>(1)</sup> | Auto                     |
| X-Ray Tube Current          | (0018,1151) | IS |   | ANAP <sup>(1)</sup> | Auto                     |
| Exposure                    | (0018,1152) | IS |   | ANAP <sup>(1)</sup> | Auto                     |
| Exposure in $\mu$ As        | (0018,1153) | IS |   | ANAP <sup>(1)</sup> | Auto                     |
| Imager Pixel Spacing        | (0018,1164) | DS |   | ALWAYS              | Auto                     |
| Pixel Spacing               | (0028,0030) | DS |   | ALWAYS              | Auto                     |
| Generator Power             | (0018,1170) | IS |   | ANAP <sup>(1)</sup> | Auto                     |

| Attribute Name                            | Tag         | VR | Value                       | Presence of Value  | Source        |
|---|-------------|----|-----------------------------|--------------------|---------------|
| Acquisition Device Processing Description | (0018,1400) | LO |                             | ANAP               | Auto          |
| Acquisition Device Processing Code        | (0018,1401) | LO |                             | ANAP               | Auto          |
| Cassette Orientation                      | (0018,1402) | CS |                             | ALWAYS             | Auto / User   |
| Cassette Size                             | (0018,1403) | CS |                             | ALWAYS             | Auto          |
| Relative X-Ray Exposure                   | (0018,1405) | IS | Lgm value multiplied by 100 | ANAP               | Auto          |
| Exposure Index                            | (0018,1411) | DS |                             | ANAP <sup>16</sup> | Auto          |
| Target Exposure Index                     | (0018,1412) | DS |                             | ANAP <sup>17</sup> | Auto          |
| Deviation Index                           | (0018,1413) | DS |                             | ANAP <sup>18</sup> | Auto          |
| Sensitivity                               | (0018,6000) | DS |                             | ALWAYS             | Config / User |
| <b>Modality LUT</b>                       |             |    |                             |                    |               |
| Modality LUT Sequence                     | (0028,3000) | SQ |                             | ANAP               | Auto          |
| > LUT Descriptor                          | (0028,3002) | SS |                             | ANAP               | Auto          |
| > LUT Explanation                         | (0028,3003) | LO |                             | ANAP               | Auto          |
| > Modality LUT Type                       | (0028,3004) | LO |                             |                    |               |
| > LUT Data                                | (0028,3006) | SS |                             | ANAP               | Auto          |
| Rescale Intercept                         | (0028,1052) | DS |                             | ANAP               | Auto          |
| Rescale Slope                             | (0028,1053) | DS |                             | ANAP               | Auto          |
| Rescale Type                              | (0028,1054) | LO |                             | ANAP               | Auto          |
| <b>VOI LUT</b>                            |             |    |                             |                    |               |
| VOI LUT Sequence                          | (0028,3010) | SQ |                             | ANAP               | Config        |
| > LUT Descriptor                          | (0028,3002) | SS |                             | ANAP               | Auto          |
| > LUT Explanation                         | (0028,3003) | LO |                             | ANAP               | Auto          |
| > LUT Data                                | (0028,3006) | SS |                             | ANAP               | Auto          |
| Window Center                             | (0028,1050) | DS |                             | ANAP               | Auto          |
| Window Width                              | (0028,1051) | DS |                             | ANAP               | Auto          |
| Window Center & Width Explanation         | (0028,1055) | LO |                             | ANAP               | Auto          |

(1): Can only be available in case of XRDI.

### 6.1.1.3 DX

#### 6.1.1.3.1 DX Image IOD

**Table 6.1-4: IOD of Created DX Image SOP Instances**

| IE      | Module                 | Reference   | Presence of Module   |
|---------|------------------------|-------------|----------------------|
| Patient | Patient Identification | Table 6.1-1 | ALWAYS               |
|         | Extended Patient       | Table 6.1-1 | If received from MWL |
| Study   | General Study          | Table 6.1-1 | ALWAYS               |
|         | Patient Study          | Table 6.1-1 | ALWAYS               |
| Series  | General Series         | Table 6.1-1 | ALWAYS               |
|         | DX Series              | Table 6.1-5 | ALWAYS               |

<sup>16</sup> Present when NX is configured for EI and not for lgm

<sup>17</sup> Present when NX is configured for EI and not for lgm and when the dose statistic is complete or fixed

<sup>18</sup> Present when NX is configured for EI and not for lgm and when the dose statistic is complete or fixed

| IE        | Module                 | Reference   | Presence of Module                |
|-----------|------------------------|-------------|-----------------------------------|
| Equipment | General Equipment      | Table 6.1-1 | ALWAYS                            |
| Image     | General Image          | Table 6.1-1 | ALWAYS                            |
|           | Image Pixel            | Table 6.1-1 | ALWAYS                            |
|           | Display Shutter        | Table 6.1-5 | When shutter is applied           |
|           | DX Anatomy             | Table 6.1-5 | ALWAYS                            |
|           | DX Image               | Table 6.1-5 | ALWAYS                            |
|           | DX Detector            | Table 6.1-5 | ALWAYS                            |
|           | X-Ray acquisition dose | Table 6.1-5 | When the XRDI Module is installed |
|           | X-Ray Collimator       | Table 6.1-5 | When the XRDI Module is installed |
|           | DX Positioning         | Table 6.1-5 |                                   |
|           | VOI LUT                | Table 6.1-5 | ALWAYS                            |
|           | Acquisition Context    | Table 6.1-5 | ALWAYS                            |
|           | SOP Common             | Table 6.1-1 | ALWAYS                            |

### 6.1.1.3.2 DX Modules

Table 6.1-5: DX Module of Created SOP Instances

| Attribute Name                               | Tag         | VR | Value   | Presence of Value | Source        |
|--|-------------|----|---|-------------------|---------------|
| <b>Series</b>                                |             |    |   |                   |               |
| DX Series                                    |             |    |   |                   |               |
| Modality                                     | (0008,0060) | CS | DX  | ALWAYS            | Config        |
| Referenced Performed Procedure Step Sequence | (0008,1111) | SQ |   | ANAP              | Auto          |
| > Referenced SOP Class UID                   | (0008,1150) | UI |   | ANAP              | Auto          |
| > Referenced SOP Instance UID                | (0008,1155) | UI |   | ANAP              | Auto          |
| Presentation Intent Type                     | (0008,0068) | CS | For Processing/For Presentation   | ALWAYS            | Config        |
| <b>Image</b>                                 |             |    |   |                   |               |
| Display Shutter                              |             |    |   |                   |               |
| Shutter Shape                                | (0018,1600) | CS |   | ANAP              | User          |
| Shutter Left Vertical Edge                   | (0018,1602) | IS |   | ANAP              | User          |
| Shutter Right Vertical Edge                  | (0018,1604) | IS |   | ANAP              | User          |
| Shutter Upper Horizontal Edge                | (0018,1606) | IS |   | ANAP              | User          |
| Shutter Lower Horizontal Edge                | (0018,1608) | IS |   | ANAP              | User          |
| Shutter Presentation Value                   | (0018,1622) | US |   | ANAP              | User          |
| DX Anatomy                                   |             |    |   |                   |               |
| Anatomic Region Sequence                     | (0008,2218) | SQ |   | ANAP              | Config / User |
| > Anatomic Region Modifier Sequence          | (0008,2220) | SQ |   | ANAP              | Config        |
| Image Laterality                             | (0020,0062) | CS |   | ALWAYS            | Config / User |
| DX Image                                     |             |    |   |                   |               |
| Image Type                                   | (0008,0008) | CS | ORIGINAL when coming from a digitizer, DERIVED when the image comes from a Save as new.<br><br>The sub-attribute IOD Specific Characteristics is filled in with the exposure type name for QA images (used by Auto-QC2) | ALWAYS            | Auto          |

| Attribute Name                            | Tag         | VR | Value   | Presence of Value   | Source |
|---|-------------|----|---|---------------------|--------|
| Samples per Pixel                         | (0028,0002) | US |   | ALWAYS              | Auto   |
| Photometric Interpretation                | (0028,0004) | CS |   | ALWAYS              | Config |
| Bits Allocated                            | (0028,0100) | US |   | ALWAYS              | Auto   |
| Bits Stored                               | (0028,0101) | US |   | ALWAYS              | Auto   |
| High Bit                                  | (0028,0102) | US |   | ALWAYS              | Auto   |
| Pixel Representation                      | (0028,0103) | US |   | ALWAYS              | Auto   |
| Pixel Intensity Relationship              | (0028,1040) | CS |   | ALWAYS              | Auto   |
| Pixel Intensity Relationship Sign         | (0028,1041) | SS |   | ALWAYS              | Auto   |
| Rescale Intercept <sup>19</sup>           | (0028,1052) | DS | 0 and 1   | ALWAYS              | Config |
| Rescale Slope <sup>19</sup>               | (0028,1053) | DS | 0 and 1   | ALWAYS              | Config |
| Rescale Type                              | (0028,1054) | LO |   | ALWAYS              | Config |
| Presentation LUT Shape                    | (2050,0020) | CS | IDENTITY  | ALWAYS              | Fixed  |
| Lossy Image Compression                   | (0028,2110) | CS | 00  | ALWAYS              | Fixed  |
| Lossy Image Compression Ratio             | (0028,2112) | DS |   | ANAP                | Auto   |
| Derivation Description                    | (0008,2111) | ST |   | ANAP                | Auto   |
| Acquisition Device Processing Description | (0018,1400) | LO |   | ANAP                | Auto   |
| Acquisition Device Processing Code        | (0018,1401) | LO |   | ANAP                | Auto   |
| Patient Orientation                       | (0020,0020) | CS |   | ALWAYS              | Auto   |
| Calibration Image                         | (0050,0004) | CS |   | ANAP                | Config |
| Burned in Annotation                      | (0028,0301) | CS | YES   | ALWAYS              | Fixed  |
| VOI LUT Sequence                          | (0028,3010) | SQ |   | ANAP                | Auto   |
| >LUT Descriptor                           | (0028,3002) | SS |   | ANAP                | Auto   |
| >LUT Explanation                          | (0028,3003) | LO |   | ANAP                | Auto   |
| >LUT Data                                 | (0028,3006) | SS |   | ANAP                | Auto   |
| Window Center                             | (0028,1050) | DS |   | ANAP                | Auto   |
| Window Width                              | (0028,1051) | DS |   | ANAP                | Auto   |
| Window Center & Width Explanation         | (0028,1055) | LO |   | ANAP                | Auto   |
| <b>DX Detector</b>                        |             |    |   |                     |        |
| Detector Type                             | (0018,7004) | CS |   | ANAP                | Auto   |
| Detector Configuration                    | (0018,7005) | CS |   | ANAP <sup>(1)</sup> | Auto   |
| Detector Description                      | (0018,7006) | LT | For the DX/S family of digitizers:<br>0 = Phosphor IP<br>1 = Needle IP<br><br>For the ADC Compact family of digitizers:<br>- Label of the type of storage phosphor plates (e.g. MD10) | ALWAYS              | Auto   |
| Detector Mode                             | (0018,7008) | LT |   | ANAP <sup>(1)</sup> | Auto   |
| Detector ID                               | (0018,700A) | SH | For the DX/S family of digitizers, this is the Image plate ID.<br><br>For the ADC Compact family of digitizers, this is the cassette ID.  | ALWAYS              | Auto   |
| Date of Last Detector Calibration         | (0018,700C) | DA |   | ANAP                | Auto   |

<sup>19</sup> Fixed when using DX for Processing

| Attribute Name                               | Tag         | VR | Value                            | Presence of Value   | Source        |
|--|-------------|----|----------------------------------|---------------------|---------------|
| Time of Last Detector Calibration            | (0018,700E) | TM |                                  | ANAP                | Auto          |
| Exposures on Detector Since Last Calibration | (0018,7010) | IS |                                  | ANAP                | Auto          |
| Exposures on Detector Since Manufactured     | (0018,1404) | IS |                                  | ANAP                | Auto          |
| Detector Time Since Last Exposure            | (0018,7012) | DS |                                  | ANAP <sup>(1)</sup> | Auto          |
| Detector Active Time                         | (0018,7014) | DS |                                  | ANAP <sup>(1)</sup> | Auto          |
| Detector Activation Offset From Exposure     | (0018,7016) | DS |                                  | ANAP <sup>(1)</sup> | Auto          |
| Detector Binning                             | (0018,701A) | DS |                                  | ANAP <sup>(1)</sup> | Auto          |
| Detector Manufacturer Name                   | (0018,702A) | LO |                                  | ANAP <sup>(1)</sup> | Auto          |
| Detector Manufacturer's Model Name           | (0018,702B) | LO |                                  | ANAP <sup>(1)</sup> | Auto          |
| Detector Conditions Nominal Flag             | (0018,7000) | CS |                                  | ANAP <sup>(1)</sup> | Auto          |
| Detector temperature                         | (0018,7001) | DS |                                  | ANAP <sup>(1)</sup> | Auto          |
| Sensitivity                                  | (0018,6000) | DS |                                  | ALWAYS              | Config / User |
| Field of View Shape                          | (0018,1147) | CS |                                  | ANAP <sup>(1)</sup> | Auto          |
| Field of view Dimension(s)                   | (0018,1149) | IS |                                  | ANAP <sup>(1)</sup> | Auto          |
| Field of View Origin                         | (0018,7030) | DS |                                  | ANAP <sup>(1)</sup> | Auto          |
| Field of View Rotation                       | (0018,7032) | DS |                                  | ANAP                | User          |
| Field of View Horizontal Flip                | (0018,7034) | CS |                                  | ANAP                | User          |
| Imager Pixel Spacing                         | (0018,1164) | DS |                                  | ALWAYS              | Auto          |
| Pixel Spacing                                | (0028,0030) | DS |                                  | ALWAYS              | Auto          |
| Detector Element Physical Size               | (0018,7020) | DS |                                  | ANAP                | Auto          |
| Detector Element Spacing                     | (0018,7022) | DS |                                  | ANAP                | Auto          |
| Detector Active Shape                        | (0018,7024) | CS | RECTANGLE                        | ALWAYS              | Fixed         |
| Detector Active Dimension(s)                 | (0018,7026) | DS | Cassette size translated into MM | ALWAYS              | Auto          |
| Detector Active Origin                       | (0018,7028) | DS |                                  | ANAP                | Auto          |
| <b>X-Ray Collimator</b>                      |             |    |                                  |                     |               |
| Collimator Shape                             | (0018,1700) | CS |                                  | ANAP                | Auto / User   |
| Collimator Left Vertical Edge                | (0018,1702) | IS |                                  | ANAP                | Auto / User   |
| Collimator Right Vertical Edge               | (0018,1704) | IS |                                  | ANAP                | Auto / User   |
| Collimator Upper Horizontal Edge             | (0018,1706) | IS |                                  | ANAP                | Auto / User   |
| Collimator Lower Horizontal Edge             | (0018,1708) | IS |                                  | ANAP                | Auto / User   |
| Center of Circular Collimator                | (0018,1710) | IS |                                  | ANAP                | Auto / User   |
| Radius of Circular Collimator                | (0018,1712) | IS |                                  | ANAP                | Auto / User   |
| Vertices of the Polygonal Collimator         | (0018,1720) | IS |                                  | ANAP                | Auto / User   |
| <b>DX Positioning</b>                        |             |    |                                  |                     |               |
| Projection Eponymous Name Code Sequence      | (0018,5104) | SQ |                                  | ANAP                | Auto          |
| Patient Position                             | (0018,5100) | CS |                                  | ANAP                | User          |
| View Position                                | (0018,5101) | CS | Based on protocol code           | ALWAYS              | Auto / User   |
| View Code Sequence                           | (0054,0220) | SQ |                                  | ANAP                | Auto          |
| Patient Orientation Code Sequence            | (0054,0410) | SQ |                                  | ANAP <sup>(1)</sup> | Auto          |
| Patient Orientation Modifier Code Sequence   | (0054,0412) | SQ |                                  | ANAP <sup>(1)</sup> | Auto          |
| Patient Gantry Relationship Code Sequence    | (0054,0414) | SQ |                                  | ANAP <sup>(1)</sup> | Auto          |
| Distance Source to Patient                   | (0018,1111) | DS |                                  | ANAP <sup>(1)</sup> | Auto          |

| Attribute Name                              | Tag         | VR | Value   | Presence of Value   | Source                   |
|---|-------------|----|---|---------------------|--------------------------|
| Distance Source to Detector                 | (0018,1110) | DS |   | ANAP <sup>(1)</sup> | Auto                     |
| Estimated Radiographic Magnification Factor | (0018,1114) | DS | 0.00 to 9.99  | ANAP <sup>(1)</sup> | User/XRDI <sup>(1)</sup> |
| Positioner Type                             | (0018,1508) | CS |   | ANAP <sup>(1)</sup> | Auto                     |
| Positioner Primary Angle                    | (0018,1510) | DS | -180 to +180 (Degree)   | ANAP <sup>(1)</sup> | Auto/XRDI <sup>(1)</sup> |
| Positioner Secondary Angle                  | (0018,1511) | DS |   | ANAP <sup>(1)</sup> | Auto                     |
| Detector Primary Angle                      | (0018,1530) | DS |   | ANAP <sup>(1)</sup> | Auto                     |
| Detector Secondary Angle                    | (0018,1531) | DS |   | ANAP <sup>(1)</sup> | Auto                     |
| Body Part Thickness                         | (0018,11A0) | DS | 0 to 999 (mm)   | ANAP                | User/XRDI <sup>(1)</sup> |
| Compression Force                           | (0018,11A2) | DS | 0 to 99 (kg)  | ANAP                | User/XRDI <sup>(1)</sup> |
| <b>X-Ray Acquisition dose</b>               |             |    |   |                     |                          |
| KVP   | (0018,0060) | DS | 0 to 99   | ANAP <sup>(1)</sup> | User/XRDI <sup>(1)</sup> |
| X-Ray Tube Current                          | (0018,1151) | IS |   | ANAP <sup>(1)</sup> | User/XRDI <sup>(1)</sup> |
| X-Ray Tube Current in µA                    | (0018,8151) | IS |   | ANAP <sup>(1)</sup> | User/XRDI <sup>(1)</sup> |
| Exposure Time                               | (0018,1150) | IS | 0 to 99999 (ms)   | ANAP <sup>(1)</sup> | User/XRDI <sup>(1)</sup> |
| Exposure Time in µs                         | (0018,8150) | IS |   | ANAP <sup>(1)</sup> | User/XRDI <sup>(1)</sup> |
| Exposure                                    | (0018,1152) | IS | 0 to 999  | ANAP <sup>(1)</sup> | User/XRDI <sup>(1)</sup> |
| Exposure in µAs                             | (0018,1153) | IS | 0 to 999  | ANAP <sup>(1)</sup> | Auto/XRDI <sup>(1)</sup> |
| Distance Source to Detector                 | (0018,1110) | DS |   | ANAP <sup>(1)</sup> | User/XRDI <sup>(1)</sup> |
| Distance Source to Patient                  | (0018,1111) | DS |   | ANAP <sup>(1)</sup> | Auto                     |
| Image Area Dose Product                     | (0018,115E) | DS |   | ANAP <sup>(1)</sup> | Auto                     |
| Body Part Thickness                         | (0018,11A0) | DS | 0 to 999  | ANAP <sup>(1)</sup> | User/XRDI <sup>(1)</sup> |
| Relative X-Ray Exposure                     | (0018,1405) | IS | Lgm value multiplied by 100   | ANAP                | Auto                     |
| Entrance Dose                               | (0040,0302) | US |   | ANAP <sup>(1)</sup> | User/XRDI <sup>(1)</sup> |
| Entrance Dose in mGy                        | (0040,8302) | US |   | ANAP <sup>(1)</sup> | User/XRDI <sup>(1)</sup> |
| Exposed Area                                | (0040,0303) | US |   | ANAP <sup>(1)</sup> | Auto                     |
| Distance Source to Entrance                 | (0040,0306) | DS |   | ANAP <sup>(1)</sup> | Auto                     |
| Comments on Radiation Dose                  | (0040,0310) | ST |   | ANAP <sup>(1)</sup> | Auto                     |
| X-Ray Output                                | (0040,0312) | DS |   | ANAP <sup>(1)</sup> | Auto                     |
| Half Value Layer                            | (0040,0314) | DS |   | ANAP <sup>(1)</sup> | Auto                     |
| Organ Dose                                  | (0040,0316) | DS | 0.0 to 99.9, 100 to 999, Calculation error = "----" (4 blanks), >999mGy = "----" (4 blanks) | ANAP <sup>(1)</sup> | User/XRDI <sup>(1)</sup> |
| Organ Exposed                               | (0040,0318) | CS |   | ANAP <sup>(1)</sup> | Auto                     |
| Exposure Index                              | (0018,1411) | DS |   | ANAP <sup>20</sup>  | Auto                     |
| Target Exposure Index                       | (0018,1412) | DS |   | ANAP <sup>21</sup>  | Auto                     |
| Deviation Index                             | (0018,1413) | DS |   | ANAP <sup>22</sup>  | Auto                     |
| Anode target Material                       | (0018,1191) | CS | 0 to 1, 0 = W, 1 = Mo   | ANAP <sup>(1)</sup> | User/XRDI <sup>(1)</sup> |
| Filter Material                             | (0018,7050) | CS | 0 to 1, 0 = Rh, 1 = Mo  | ANAP <sup>(1)</sup> | User/XRDI <sup>(1)</sup> |
| Filter Thickness Minimum                    | (0018,7052) | DS |   | ANAP <sup>(1)</sup> | Auto                     |
| Filter Thickness Maximum                    | (0018,7054) | DS |   | ANAP <sup>(1)</sup> | Auto                     |
| Grid  | (0018,1166) | CS | 0 to 1, 0 = Grid, 1 = No Grid   | ANAP <sup>(1)</sup> | User/XRDI <sup>(1)</sup> |
| Grid Absorbing Material                     | (0018,7040) | LT |   | ANAP <sup>(1)</sup> | User/XRDI <sup>(1)</sup> |
| Grid Pitch                                  | (0018,7044) | DS |   | ANAP <sup>(1)</sup> | User/XRDI <sup>(1)</sup> |
| Grid Thickness                              | (0018,7042) | DS |   | ANAP <sup>(1)</sup> | User/XRDI <sup>(1)</sup> |
| Exposure Control Mode                       | (0018,7060) | CS | - Manual  | ANAP <sup>(1)</sup> | User/XRDI <sup>(1)</sup> |

<sup>20</sup> Present when NX is configured for EI and not for lgm<sup>21</sup> Present when NX is configured for EI and not for lgm and when the dose statistic is complete or fixed<sup>22</sup> Present when NX is configured for EI and not for lgm and when the dose statistic is complete or fixed



| Attribute Name                          | Tag         | VR | Value       | Presence of Value   | Source                   |
|---|-------------|----|-------------|---------------------|--------------------------|
|   |             |    | - Automatic |                     |                          |
| Exposure Control Mode Description       | (0018,7062) | LT |             | ANAP <sup>(1)</sup> | User/XRDI <sup>(1)</sup> |
| Phototimer Setting                      | (0018,7065) | DS |             | ANAP <sup>(1)</sup> | User/XRDI <sup>(1)</sup> |
| Image and Fluoroscopy Area Dose Product | (0018,115E) | DS |             | ANAP <sup>(1)</sup> | User/XRDI <sup>(1)</sup> |
| Grid Focal Distance                     | (0018,704C) | DS |             | ANAP <sup>(1)</sup> | User/XRDI <sup>(1)</sup> |
| Rectification Type                      | (0018,1156) | CS |             | ANAP <sup>(1)</sup> | Auto                     |
| <b>VOI LUT</b>                          |             |    |             |                     |                          |
| VOI LUT Sequence                        | (0028,3010) | SQ |             | ANAP                | Auto                     |
| >LUT Descriptor                         | (0028,3002) | SS |             | ANAP                | Auto                     |
| >LUT Explanation                        | (0028,3003) | LO |             | ANAP                | Auto                     |
| >LUT Data                               | (0028,3006) | SS |             | ANAP                | Auto                     |
| Window Center                           | (0028,1050) | DS |             | ANAP                | Auto                     |
| Window Width                            | (0028,1051) | DS |             | ANAP                | Auto                     |
| Window Center & Width Explanation       | (0028,1055) | LO |             | ANAP                | Auto                     |
| <b>Acquisition Context</b>              |             |    |             |                     |                          |
| Acquisition Context Sequence            | (0040,0555) | SQ |             | EMPTY               | Fixed                    |
| > Value Type                            | (0040,A040) | CS |             | EMPTY               | Fixed                    |
| > Concept Name Code Sequence            | (0040,A043) | SQ |             | EMPTY               | Fixed                    |
| > Referenced Frame numbers              | (0040,A136) | US |             | EMPTY               | Fixed                    |
| > Numeric Value                         | (0040,A30A) | DS |             | EMPTY               | Fixed                    |
| > Measurement Units Code Sequence       | (0040,08EA) | SQ |             | EMPTY               | Fixed                    |
| > Date                                  | (0040,A121) | DT |             | EMPTY               | Fixed                    |
| > Time                                  | (0040,A122) | TM |             | EMPTY               | Fixed                    |
| > Person Name                           | (0040,A123) | PN |             | EMPTY               | Fixed                    |
| > Text Value                            | (0040,A160) | UT |             | EMPTY               | Fixed                    |
| > Concept Code Sequence                 | (0040,A168) | SQ |             | EMPTY               | Fixed                    |
| Acquisition Context Description         | (0040,0556) | ST |             | ANAP <sup>(2)</sup> | Auto                     |

<sup>(1)</sup> : Can only be available in case of XRDI.

<sup>(2)</sup> : Never present

#### 6.1.1.4 MG

##### 6.1.1.4.1 MG Image IOD

**Table 6.1-6: IOD of Created MG Image SOP Instances**

| IE      | Module                 | Reference   | Presence of Module   |
|---------|------------------------|-------------|----------------------|
| Patient | Patient Identification | Table 6.1-1 | ALWAYS               |
|         | Extended Patient       | Table 6.1-1 | If received from MWL |
| Study   | General Study          | Table 6.1-1 | ALWAYS               |
|         | Patient Study          | Table 6.1-1 | ALWAYS               |
| Series  | General Series         | Table 6.1-1 | ALWAYS               |
|         | DX Series              | Table 6.1-5 | ALWAYS               |
|         | Mammography Series     | Table 6.1-7 | ALWAYS               |



| IE        | Module                 | Reference   | Presence of Module                |
|-----------|------------------------|-------------|-----------------------------------|
| Equipment | General Equipment      | Table 6.1-1 | ALWAYS                            |
| Image     | General Image          | Table 6.1-1 | ALWAYS                            |
|           | Image Pixel            | Table 6.1-1 | ALWAYS                            |
|           | Display Shutter        | Table 6.1-5 | When shutter is applied           |
|           | DX Anatomy             | Table 6.1-5 | ALWAYS                            |
|           | DX Image               | Table 6.1-5 | ALWAYS                            |
|           | DX Detector            | Table 6.1-5 | ALWAYS                            |
|           | X-Ray acquisition dose | Table 6.1-5 | When the XRDI Module is installed |
|           | X-Ray Collimator       | Table 6.1-5 | When the XRDI Module is installed |
|           | DX Positioning         | Table 6.1-5 |                                   |
|           | Mammography Image      | Table 6.1-7 | ALWAYS                            |
|           | VOI LUT                | Table 6.1-5 | ALWAYS                            |
|           | Acquisition Context    | Table 6.1-5 | ALWAYS                            |
|           | SOP Common             | Table 6.1-1 | ALWAYS                            |

#### 6.1.1.4.2 MG Modules

Table 6.1-7: MG Module of Created SOP Instances

| Attribute Name              | Tag         | VR       | Value  | Presence of Value | Source                  |
|-----------------------------|-------------|----------|--|-------------------|-------------------------|
| <b>Series</b>               |             |          |  |                   |                         |
| Mammography Series          |             |          |  |                   |                         |
| Modality                    | (0008,0060) | CS       | MG   | ALWAYS            | Config                  |
| Request Attributes Sequence | (0040,0275) | DS       | Sequence that contains attributes from the Imaging Service Request | ANAP              | MWL                     |
| <b>Equipment</b>            |             |          |  |                   |                         |
| General Equipment           |             |          |  |                   |                         |
| Pixel Padding Value         | (0028,0120) | US or SS |  | ANAP              | Auto                    |
| <b>Image</b>                |             |          |  |                   |                         |
| Mammography Image           |             |          |  |                   |                         |
| Positioner Type             | (0018,1508) | CS       | Mammographic<br>None   | ALWAYS            | Auto                    |
| Distance Source to detector | (0018,1110) | DS       |  | ANAP              | Auto                    |
| Distance Source to patient  | (0018,1111) | DS       |  | ANAP              | Auto                    |
| Positioner Primary Angle    | (0018,1510) | DS       |  | ANAP              | User/XRDI <sup>23</sup> |
| Positioner Secondary Angle  | (0018,1511) | DS       |  | ANAP              | User/XRDI               |
| Image Laterality            | (0020,0062) | CS       | R = right<br>L = left<br>B = both (e.g. cleavage)                  | ALWAYS            | Auto                    |
| Organ Exposed               | (0040,0318) | CS       | BREAST   | ALWAYS            | Auto                    |
| Implant Present             | (0028,1300) | DS       | YES<br>NO  | ANAP              | Config                  |
| View Code Sequence          | (0054,0220) | CS       |  | Always            | Config / User           |
| View Modifier Code Sequence | (0054,0222) | SQ       |  | ANAP              | Config / User           |
| <b>Image Pixel</b>          |             |          |  |                   |                         |
| Pixel Padding Range Limit   | (0028,0121) | US or SS |  | ANAP              | Auto                    |

<sup>23</sup> Value is supplied by user of XRDI component in case of XRDI configuration

**6.1.1.5 GSPS****6.1.1.5.1 GSPS IOD****Table 6.1-8: IOD of Created GSPS SOP Instances**

| IE        | Module                    | Reference   | Presence of Module                   |
|-----------|---------------------------|-------------|--------------------------------------|
| Patient   | Patient Identification    | Table 6.1-1 | ALWAYS                               |
|           | Extended Patient          | Table 6.1-1 | If received from MWL                 |
| Study     | General Study             | Table 6.1-1 | ALWAYS                               |
|           | Patient Study             | Table 6.1-1 | ALWAYS                               |
| Series    | General Series            | Table 6.1-1 | ALWAYS                               |
|           | Presentation Series       | Table 6.1-9 | ALWAYS                               |
| Equipment | General Equipment         | Table 6.1-1 | ALWAYS                               |
|           | Image                     |             |                                      |
| Image     | Presentation State        | Table 6.1-9 | ALWAYS                               |
|           | Display Shutter           | Table 6.1-9 | When shutter is applied              |
|           | Bitmap Display Shutter    | Table 6.1-9 | When shutter is applied              |
|           | Overlay Plane             | Table 6.1-9 | When shutter is applied              |
|           | Displayed Area            | Table 6.1-9 | ALWAYS                               |
|           | Graphic Annotation        | Table 6.1-9 | When graphic annotations are present |
|           | Graphic Layer             | Table 6.1-9 | When graphic annotations are present |
|           | Modality LUT              | Table 6.1-9 | ALWAYS                               |
|           | Softcopy VOI LUT          | Table 6.1-9 | ALWAYS                               |
|           | Softcopy Presentation LUT | Table 6.1-9 | ALWAYS                               |
|           | Spatial Transformation    | Table 6.1-9 | ALWAYS                               |
|           | SOP Common                | Table 6.1-1 | ALWAYS                               |

**6.1.1.5.2 GSPS Modules****Table 6.1-9: GSPS Modules of Created SOP Instances**

| Attribute Name                           | Tag         | VR | Value | Presence of Value | Source |
|--|-------------|----|-------|-------------------|--------|
| <b>Series</b>                            |             |    |       |                   |        |
| Presentation Series                      |             |    |       |                   |        |
| Modality                                 | (0008,0060) | CS | PR    | ALWAYS            | Fixed  |
| <b>Presentation State</b>                |             |    |       |                   |        |
| Presentation State Identification Module |             |    |       |                   |        |
| Instance Number                          | (0020,0013) | IS |       | ALWAYS            | Auto   |
| Content Label                            | (0070,0080) | CS | GSPS  | ALWAYS            | Auto   |
| Content Description                      | (0070,0081) | LO |       | Empty             | Empty  |
| Presentation Creation Date               | (0070,0082) | DA |       | ALWAYS            | Auto   |
| Presentation Creation Time               | (0070,0083) | TM |       | ALWAYS            | Auto   |
| Content Creator's Name                   | (0070,0084) | PN |       | VNAP              | Auto   |
| Presentation State Relationship Module   |             |    |       |                   |        |
| Referenced Series Sequence               | (0008,1115) | SQ |       | ALWAYS            | Auto   |
| >Series Instance UID                     | (0020,000E) | UI |       | ALWAYS            | Auto   |
| >Referenced Image Sequence               | (0008,1140) | SQ |       | ALWAYS            | Auto   |
| >>Referenced SOP Class UID               | (0008,1150) | UI |       | ALWAYS            | Auto   |
| >>Referenced SOP Instance UID            | (0008,1155) | UI |       | ALWAYS            | Auto   |

| Attribute Name                           | Tag         | VR       | Value                | Presence of Value | Source |
|--|-------------|----------|----------------------|-------------------|--------|
| >> Referenced Frame Number               | (0008,1160) | IS       |                      | ALWAYS            | Auto   |
| <b>Presentation State Mask</b>           |             |          |                      |                   |        |
| Mask Subtraction Sequence                | (0028,6100) | SQ       |                      | ANAP              | Auto   |
| >Mask Operation                          | (0028,6101) | CS       |                      | ANAP              | Auto   |
| > Contrast Frame Averaging               | (0028,6112) | US       |                      | ANAP              | Auto   |
| Recommended Viewing Mode                 | (0028,1090) | CS       |                      | ANAP              | Auto   |
| <b>Display Shutter</b>                   |             |          |                      |                   |        |
| Shutter Shape                            | (0018,1600) | CS       |                      | ANAP              | User   |
| Shutter Left Vertical Edge               | (0018,1602) | IS       |                      | ANAP              | User   |
| Shutter Right Vertical Edge              | (0018,1604) | IS       |                      | ANAP              | User   |
| Shutter Upper Horizontal Edge            | (0018,1606) | IS       |                      | ANAP              | User   |
| Shutter Lower Horizontal Edge            | (0018,1608) | IS       |                      | ANAP              | User   |
| Vertices of the Polygonal Shutter        | (0018,1620) | IS       |                      | ANAP              | User   |
| <b>Bitmap Display Shutter</b>            |             |          |                      |                   |        |
| Shutter Shape                            | (0018,1600) | CS       | BITMAP               | ANAP              | User   |
| Shutter Presentation Value               | (0018,1622) | US       | 0                    | ANAP              | User   |
| Shutter Overlay Group                    | (0018,1623) | US       | 6000                 | ANAP              | Auto   |
| <b>Overlay Plane</b>                     |             |          |                      |                   |        |
| Overlay Rows                             | (60xx,0010) | US       |                      | ANAP              | Auto   |
| Overlay Columns                          | (60xx,0011) | US       |                      | ANAP              | Auto   |
| Overlay Description                      | (60xx,0022) | LO       | Background darkening | ANAP              | Auto   |
| Overlay Type                             | (60xx,0040) | CS       | G                    | ANAP              | Auto   |
| Overlay Subtype                          | (60xx,0045) | LO       | AUTOMATED            | ANAP              | Auto   |
| Overlay Origin                           | (60xx,0050) | SS       | 1                    | ANAP              | Auto   |
| Overlay Bits Allocated                   | (60xx,0100) | US       | 1                    | ANAP              | Auto   |
| Overlay Bit Position                     | (60xx,0102) | US       | 0                    | ANAP              | Auto   |
| Overlay Activation Layer                 | (60xx,1001) | CS       |                      | ANAP              | Auto   |
| Overlay Label                            | (60xx,1500) | LO       | Mask image           | ANAP              | Auto   |
| Overlay Data                             | (60xx,3000) | OB or OW |                      | ANAP              | Auto   |
| <b>Displayed Area</b>                    |             |          |                      |                   |        |
| Displayed Area Selection Sequence        | (0070,005A) | SQ       |                      | ALWAYS            | Auto   |
| >Referenced Image Sequence               | (0008,1140) | SQ       |                      | ANAP              | Auto   |
| >>Referenced SOP Class UID               | (0008,1150) | UI       |                      | ANAP              | Auto   |
| >>Referenced SOP Instance UID            | (0008,1155) | UI       |                      | ANAP              | Auto   |
| >Displayed Area Top Left Hand Corner     | (0070,0052) | SL       |                      | ALWAYS            | Auto   |
| >Displayed Area Bottom Right Hand Corner | (0070,0053) | SL       |                      | ALWAYS            | Auto   |
| >Presentation Size Mode                  | (0070,0100) | CS       | SCALE TO FIT         | ALWAYS            | Auto   |
| >Presentation Pixel Spacing              | (0070,0101) | DS       |                      | ALWAYS            | Auto   |
| > Presentation Pixel Aspect Ratio        | (0070,0102) | IS       |                      | ANAP              | Auto   |
| > Presentation Pixel Magnification Ratio | (0070,0103) | FL       |                      | ANAP              | Auto   |
| <b>Spatial Transformation</b>            |             |          |                      |                   |        |
| Image Rotation                           | (0070,0042) | US       |                      | ANAP              | User   |
| Image Horizontal Flip                    | (0070,0041) | CS       |                      | ANAP              | User   |
| <b>Graphic Annotation</b>                |             |          |                      |                   |        |
| Graphic Annotation Sequence              | (0070,0001) | SQ       |                      | ALWAYS            | User   |

| Attribute Name                                     | Tag         | VR | Value    | Presence of Value | Source |
|--|-------------|----|----------|-------------------|--------|
| >Referenced Image Sequence                         | (0008,1140) | SQ |          | ALWAYS            | Auto   |
| >>Referenced SOP Class UID                         | (0008,1150) | UI |          | ALWAYS            | Auto   |
| >>Referenced SOP Instance UID                      | (0008,1155) | UI |          | ALWAYS            | Auto   |
| >>Referenced Frame number                          | (0008,1155) | IS |          |                   |        |
| >Graphic Layer                                     | (0070,0002) | CS |          | ALWAYS            | Auto   |
| >Text Object Sequence                              | (0070,0008) | SQ |          | ANAP              | Auto   |
| >>Bounding Box Annotation Units                    | (0070,0003) | CS | PIXEL    | ANAP              | Auto   |
| >>Anchor Point Annotation Units                    | (0070,0004) | CS | PIXEL    | ANAP              | Auto   |
| >>Unformatted Text Value                           | (0070,0006) | ST |          | ANAP              | Auto   |
| >>Bounding Box Top Left Hand Corner                | (0070,0010) | FL |          | ANAP              | Auto   |
| >>Bounding Box Bottom Right Hand Corner            | (0070,0011) | FL |          | ANAP              | Auto   |
| >>Bounding Box Text Horizontal Justification       | (0070,0012) | CS |          | ANAP              | Auto   |
| >>Anchor Point                                     | (0070,0014) | FL |          | ANAP              | Auto   |
| >>Anchor Point Visibility                          | (0070,0015) | CS |          | ANAP              | Auto   |
| >Graphic Object Sequence                           | (0070,0009) | SQ |          | ANAP              | Auto   |
| >>Graphic Annotation Units                         | (0070,0005) | CS |          | ANAP              | Auto   |
| >>Graphic Dimensions                               | (0070,0020) | US |          | ANAP              | Auto   |
| >>Number of Graphic Points                         | (0070,0021) | US |          | ANAP              | Auto   |
| >>Graphic Data                                     | (0070,0022) | FL |          | ANAP              | Auto   |
| >>Graphic Type                                     | (0070,0023) | CS |          | ANAP              | Auto   |
| >>Graphic Filled                                   | (0070,0024) | CS |          | ANAP              | Auto   |
| <b>Softcopy Presentation LUT</b>                   |             |    |          |                   |        |
| Presentation LUT Shape                             | (2050,0020) | CS | IDENTITY | ALWAYS            | Auto   |
| <b>Graphic Layer</b>                               |             |    |          |                   |        |
| Graphic Layer Sequence                             | (0070,0060) | SQ |          | ALWAYS            | Auto   |
| >Graphic Layer                                     | (0070,0002) | CS |          | ALWAYS            | Auto   |
| >Graphic Layer Order                               | (0070,0062) | IS |          | ALWAYS            | Auto   |
| >Graphic Layer Recommended Display Grayscale Value | (0070,0066) | US |          | ANAP              | Auto   |
| >Graphic Layer Recommended Display RGB Value       | (0070,0067) | US |          | ANAP              | Auto   |
| >Graphic Layer Description                         | (0070,0068) | LO |          | ANAP              | Auto   |
| <b>Modality LUT</b>                                |             |    |          |                   |        |
| Rescale Intercept                                  | (0028,1052) | DS |          | ANAP              | Auto   |
| Rescale Slope                                      | (0028,1053) | DS |          | ANAP              | Auto   |
| Rescale Type                                       | (0028,1054) | LO |          | ANAP              | Auto   |
| <b>VOI LUT</b>                                     |             |    |          |                   |        |
| VOI LUT Sequence                                   | (0028,3010) | SQ |          | ANAP              | Auto   |
| >LUT Descriptor                                    | (0028,3002) | SS |          | ANAP              | Auto   |
| >LUT Explanation                                   | (0028,3003) | LO |          | ANAP              | Auto   |
| >LUT Data  | (0028,3006) | SS |          | ANAP              | Auto   |
| Window Center                                      | (0028,1050) | DS |          | ANAP              | Auto   |
| Window Width                                       | (0028,1051) | DS |          | ANAP              | Auto   |
| Window Center & Width Explanation                  | (0028,1055) | LO |          | ANAP              | Auto   |

**6.1.1.6 X-Ray Radiation Dose SR Storage****6.1.1.6.1 X-Ray Radiation Dose SR IOD****Table 6.1-10: IOD of Created X-Ray Radiation Dose SR SOP Instances**

| IE        | Module                     | Reference    | Presence of Module |
|-----------|----------------------------|--------------|--------------------|
| Patient   | Patient                    | Table 6.1-11 |                    |
|           | Patient Identification     | Table 6.1-11 |                    |
|           | Patient Demographic        | Table 6.1-11 |                    |
|           | Patient Medical            | Table 6.1-11 |                    |
| Study     | General Study              | Table 6.1-11 |                    |
|           | Patient Study              | Table 6.1-11 |                    |
| Series    | SR Document Series         | Table 6.1-11 |                    |
| Equipment | General Equipment          | Table 6.1-11 |                    |
|           | Enhanced General Equipment | Table 6.1-11 |                    |
| Document  | SR Document General        | Table 6.1-11 |                    |
|           | SR Document Content        | Table 6.1-12 |                    |
|           | SOP Common                 | Table 6.1-11 |                    |

**6.1.1.6.2 X-Ray Radiation Dose SR Modules****Table 6.1-11: X-Ray Radiation Dose SR Modules of Created SOP Instances**

| Attribute Name   | Tag         | VR | Value | Presence of Value | Source     |
|--|-------------|----|-------|-------------------|------------|
| <b>Patient</b>   |             |    |       |                   |            |
| Patient Identification                                 |             |    |       |                   |            |
| Patient's Name   | (0010,0010) | PN |       | VNAP              | User / MWL |
| Patient ID   | (0010,0020) | LO |       | VNAP              | User / MWL |
| Issuer of Patient ID                                   | (0010,0021) | LO |       | ANAP              | MWL        |
| Other Patient IDs                                      | (0010,1000) | LO |       | ANAP              | MWL / MWL  |
| Other Patient Names                                    | (0010,1001) | PN |       | ANAP              | MWL        |
| Patient's Birth Name                                   | (0010,1005) | PN |       | ANAP              | MWL        |
| Patient's Mother's Birth Name                          | (0010,1060) | PN |       | ANAP              | MWL        |
| Medical Record Locator                                 | (0010,1090) | LO |       | ANAP              | MWL        |
| Patient Demographic                                    |             |    |       |                   |            |
| Patient's Age  | (0010,1010) | AS |       | ANAP              | User / MWL |
| Occupation   | (0010,2180) | SH |       | ANAP              | MWL        |
| Confidentiality Constraint on Patient Data Description | (0040,3001) | LO |       | ANAP              | MWL        |
| Patient's Birth Date                                   | (0010,0030) | DA |       | VNAP              | User / MWL |
| Patient's Birth Time                                   | (0010,0032) | TM |       | VNAP              | User / MWL |
| Patient's Sex  | (0010,0040) | CS |       | VNAP              | User / MWL |
| Patient's Insurance Plan Code Sequence                 | (0010,0050) | SQ |       | ANAP              | MWL        |
| > Code Sequence  |             |    |       |                   |            |
| > Patient's Primary Language Code Modifier Sequence    | (0010,0102) | SQ |       | ANAP              | MWL        |
| >> Code Sequence                                       |             |    |       |                   |            |
| Patient's Size   | (0010,1020) | DS |       | ANAP              | User / MWL |

| Attribute Name                                     | Tag         | VR | Value | Presence of Value | Source     |
|--|-------------|----|-------|-------------------|------------|
| Patient's Weight                                   | (0010,1030) | DS |       | ANAP              | User / MWL |
| Patient's Address                                  | (0010,1040) | LO |       | ANAP              | MWL        |
| Military Rank                                      | (0010,1080) | LO |       | ANAP              | User / MWL |
| Branch of Service                                  | (0010,1081) | LO |       | ANAP              | User / MWL |
| Country of Residence                               | (0010,2150) | LO |       | ANAP              | MWL        |
| Region of Residence                                | (0010,2152) | LO |       | ANAP              | MWL        |
| Patient's Telephone Numbers                        | (0010,2154) | SH |       | ANAP              | MWL        |
| Ethnic Group                                       | (0010,2160) | SH |       | ANAP              | User / MWL |
| Patient's Religious Preference                     | (0010,21F0) | LO |       | ANAP              | MWL        |
| Patient Comments                                   | (0010,4000) | LT |       | ANAP              | User / MWL |
| Current Patient Location                           | (0038,0300) | LO |       | ANAP              | MWL        |
| Patient's Institution Residence                    | (0038,0400) | LO |       | ANAP              | MWL        |
| Referenced Patient Sequence                        | (0008,1120) | SQ |       | ANAP              | MWL        |
| >Referenced SOP Class UID                          | (0008,1150) | UI |       | ANAP              | MWL        |
| >Referenced SOP Instance UID                       | (0008,1155) | UI |       | ANAP              | MWL        |
| <b>Patient Medical</b>                             |             |    |       |                   |            |
| Medical Alerts                                     | (0010,2000) | LO |       | ANAP              | MWL        |
| Contrast Allergies                                 | (0010,2110) | LO |       | ANAP              | MWL        |
| Smoking Status                                     | (0010,21A0) | CS |       | ANAP              | User / MWL |
| Pregnancy Status                                   | (0010,21C0) | US |       | ANAP              | User / MWL |
| Last Menstrual Date                                | (0010,21D0) | DA |       | ANAP              | User / MWL |
| Special Needs                                      | (0038,0050) | LO |       | ANAP              | MWL        |
| Patient State                                      | (0038,0500) | LO |       | ANAP              | MWL        |
| Additional Patient's History                       | (0010,21B0) | LT |       | ANAP              | MWL        |
| Sex neutered                                       | (0010,2203) | CS |       | VNAP              |            |
| <b>Study</b>                                       |             |    |       |                   |            |
| <b>General Study</b>                               |             |    |       |                   |            |
| Study Instance UID                                 | (0020,000D) | UI |       | ALWAYS            | MWL / Auto |
| Study Date   | (0008,0020) | DA |       | VNAP              | Auto       |
| Study Time   | (0008,0030) | TM |       | VNAP              | Auto       |
| Referring Physician's Name                         | (0008,0090) | PN |       | VNAP              | User / MWL |
| Referring Physician Identification Sequence        | (0008,0096) | SQ |       | ANAP              | MWL        |
| Study ID   | (0020,0010) | SH |       | VNAP              | Auto / MWL |
| Accession Number                                   | (0008,0050) | SH |       | VNAP              | User / MWL |
| Study Description                                  | (0008,1030) | LO |       | ANAP              | Auto       |
| Physician(s) of Record                             | (0008,1048) | PN |       | ANAP              | MWL        |
| Physician(s) of Record Identification Sequence     | (0008,1049) | SQ |       | ANAP              | MWL        |
| Name of Physician(s) Reading Study                 | (0008,1060) | PN |       | ANAP              | MWL        |
| Physician(s) Reading Study Identification Sequence | (0008,1062) | SQ |       | ANAP              | User       |
| Referenced Study Sequence                          | (0008,1110) | SQ |       | ANAP              | MWL        |
| >Referenced SOP Class UID                          | (0008,1150) | UI |       | ANAP              | MWL        |
| >Referenced SOP Instance UID                       | (0008,1155) | UI |       | ANAP              | MWL        |
| Procedure Code Sequence                            | (0008,1032) | SQ |       | ANAP              | MWL        |
| <b>Patient Study</b>                               |             |    |       |                   |            |
| Admitting Diagnoses Description                    | (0008,1080) | LO |       | ALWAYS            | MWL        |

| Attribute Name                                   | Tag         | VR | Value                      | Presence of Value | Source      |
|--|-------------|----|----------------------------|-------------------|-------------|
| Admitting Diagnoses Code Sequence                | (0008,1084) | SQ |                            | ALWAYS            | MWL         |
| <b>Series</b>                                    |             |    |                            |                   |             |
| SR Document Series                               |             |    |                            |                   |             |
| Modality   | (0008,0060) | CS |                            | ALWAYS            |             |
| Series Instance UID                              | (0020,000E) | UI | Different for each image   | ALWAYS            | Auto        |
| Series Number                                    | (0020,0011) | IS |                            | ALWAYS            | Auto        |
| Series Date                                      | (0008,0021) | DA |                            | ANAP              |             |
| Series Time                                      | (0008,0031) | TM |                            | ANAP              |             |
| Series Description                               | (0008,103E) | LO | Radiation Dose Information | ALWAYS            |             |
| Series Description Code Sequence                 | (0008,103F) | SQ |                            | ANAP              |             |
| Referenced Performed Procedure Step Sequence     | (0008,1111) | SQ |                            | ALWAYS            |             |
| >Referenced SOP Class UID                        | (0008,1150) | UI |                            | ALWAYS            |             |
| >Referenced SOP Instance UID                     | (0008,1155) | UI |                            | ALWAYS            |             |
| <b>Equipment</b>                                 |             |    |                            |                   |             |
| General Equipment                                |             |    |                            |                   |             |
| Institution Name                                 | (0008,0080) | LO |                            | ANAP              | MWL/ Config |
| Institution Address                              | (0008,0081) | ST |                            | ANAP              | MWL /Config |
| Station Name                                     | (0008,1010) | SH |                            | ANAP              | Auto        |
| Institutional Department Name                    | (0008,1040) | LO |                            | ANAP              | MWL/ Config |
| Gantry ID  | (0018,1008) | LO |                            | ANAP              | Auto        |
| Spatial Resolution                               | (0018,1050) | DS |                            | ANAP              |             |
| Date of Last Calibration                         | (0018,1200) | DA |                            | ANAP              |             |
| Time of Last Calibration                         | (0018,1201) | TM |                            | ANAP              |             |
| Date of Last Detector Calibration                | (0018,700C) | DA |                            | ANAP              |             |
| Time of Last Detector Calibration                | (0018,700E) | TM |                            | ANAP              |             |
| Pixel Padding Value                              | (0028,0120) |    |                            | ALWAYS            |             |
| Enhanced General Equipment                       |             |    |                            |                   |             |
| Manufacturer                                     | (0008,0070) | LO |                            | ALWAYS            |             |
| Manufacturer's Model Name                        | (0008,1090) | LO |                            | ALWAYS            |             |
| Device Serial Number                             | (0018,1000) | LO |                            | ALWAYS            |             |
| Software Versions                                | (0018,1020) | LO |                            | ALWAYS            |             |
| <b>Document</b>                                  |             |    |                            |                   |             |
| SR Document General                              |             |    |                            |                   |             |
| Instance Number                                  | (0020,0013) | CS |                            | ALWAYS            | Auto        |
| Preliminary Flag                                 | (0040,A496) | CS |                            | ANAP              |             |
| Completion Flag                                  | (0040,A491) | CS |                            | ALWAYS            |             |
| Completion Flag Description                      | (0040,A492) | LO |                            | ANAP              |             |
| Verification Flag                                | (0040,A493) | CS |                            | ALWAYS            |             |
| Content Date                                     | (0008,0023) | DA |                            | ALWAYS            |             |
| Content Time                                     | (0008,0033) | TM |                            | ALWAYS            |             |
| Verifying Observer Sequence                      | (0040,A073) | SQ |                            | ALWAYS            |             |
| >Verifying Observer Name                         | (0040,A075) | PN |                            | ALWAYS            |             |
| >Verifying Observer Identification Code Sequence | (0040,A088) | SQ |                            | VNAP              |             |
| >>Include 'Code Sequence Macro'                  |             |    |                            |                   |             |
| >Verifying Organization                          | (0040,A027) | LO |                            | ALWAYS            |             |
| >Verification DateTime                           | (0040,A030) | DT |                            | ALWAYS            |             |



| Attribute Name  | Tag         | VR | Value | Presence of Value | Source |
|---|-------------|----|-------|-------------------|--------|
| Author Observer Sequence                              | (0040,A078) | SQ |       | ANAP              |        |
| > Code sequence                                       |             |    |       |                   |        |
| Participant Sequence                                  | (0040,A07A) | SQ |       | ANAP              |        |
| >Participation Type                                   | (0040,A080) | CS |       | ALWAYS            |        |
| >Participation DateTime                               | (0040,A082) | DT |       | VNAP              |        |
| >Include 'Identified Person or Device Macro'          |             |    |       |                   |        |
| Custodial Organization Sequence                       | (0040,A07C) | SQ |       | ANAP              |        |
| >Institution Name                                     | (0008,0080) | LO |       | VNAP              |        |
| >Institution Code Sequence                            | (0008,0082) | SQ |       |                   |        |
| >>Include 'Code Sequence Macro'                       |             |    |       |                   |        |
| Predecessor Documents Sequence [1..n]                 | (0040,A360) | SQ |       | ALWAYS            |        |
| > Study Instance UID                                  | (0020,000D) | UI |       | ALWAYS            |        |
| > Referenced Series Sequence [1..n]                   | (0008,1115) | SQ |       | ALWAYS            |        |
| >> Series Instance UID                                | (0020,000E) | UI |       | ALWAYS            |        |
| >> Retrieve AE Title                                  | (0008,0054) | AE |       | ANAP              |        |
| >> Retrieve Location UID                              | (0040,E011) | UI |       | ANAP              |        |
| >> Storage Media File-Set ID                          | (0088,0130) | SH |       | ANAP              |        |
| >> Storage Media File-Set UID                         | (0088,0140) | UI |       | ANAP              |        |
| >> Referenced SOP Sequence                            | (0008,1199) | SQ |       | ALWAYS            |        |
| >>> Referenced SOP Class UID                          | (0008,1150) | UI |       | ALWAYS            |        |
| >>> Referenced SOP Instance UID                       | (0008,1155) | UI |       | ALWAYS            |        |
| >>> Purpose of Reference Code Sequence                | (0040,A170) | SQ |       | ANAP              |        |
| >>>> Include 'Code Sequence Macro'                    |             |    |       |                   |        |
| >>>> Referenced Digital Signature Sequence            | (0400,0402) | SQ |       | ANAP              |        |
| >>>> Digital Signature UID                            | (0400,0100) | UI |       | ALWAYS            |        |
| >>>> Signature  | (0400,0120) | OB |       | ALWAYS            |        |
| >>>> Referenced SOP Instance MAC Sequence             | (0400,0403) | SQ |       | ANAP              |        |
| >>>> MAC Calculation Transfer Syntax UID              | (0400,0010) | UI |       | ALWAYS            |        |
| >>>> MAC Algorithm                                    | (0400,0015) | CS |       | ALWAYS            |        |
| >>>> Data Elements Signed                             | (0400,0020) | AT |       | ALWAYS            |        |
| >>>> MAC  | (0400,0404) | OB |       | ALWAYS            |        |
| Identical Documents Sequence                          | (0040,A525) | SQ |       | ALWAYS            |        |
| >Include ' Hierarchical SOP Instance Reference Macro' |             |    |       |                   |        |
| Referenced Request Sequence [1..n]                    | (0040,A370) | SQ |       | ALWAYS            |        |
| >Study Instance UID                                   | (0020,000D) | UI |       | ALWAYS            |        |
| >Referenced Study Sequence                            | (0008,1110) | SQ |       | VNAP              |        |
| >> Include 'SOP Instance Reference Macro'             |             |    |       |                   |        |
| >Accession Number                                     | (0008,0050) | SH |       | VNAP              |        |
| >Issuer of Accession Number Sequence                  | (0008,0051) | SQ |       | ANAP              |        |
| >>Include HL7v2 Hierarchic Designator Macro           |             |    |       |                   |        |
| >Placer Order Number/Imaging                          | (0040,2016) | LO |       | VNAP              |        |



| Attribute Name  | Tag         | VR | Value | Presence of Value | Source |
|---|-------------|----|-------|-------------------|--------|
| Service Request                                       |             |    |       |                   |        |
| >Order Placer Identifier Sequence                     | (0040,0026) | SQ |       | ANAP              |        |
| >>Include HL7v2 Hierarchic Designator Macro           |             |    |       |                   |        |
| >Filler Order Number/Imaging Service Request          | (0040,2017) | LO |       | VNAP              |        |
| >Order Filler Identifier Sequence                     | (0040,0027) | SQ |       | ANAP              |        |
| >>Include HL7v2 Hierarchic Designator Macro           |             |    |       |                   |        |
| >Requested Procedure ID                               | (0040,1001) | SH |       | VNAP              |        |
| >Requested Procedure Description                      | (0032,1060) | LO |       | ALWAYS            |        |
| >Requested Procedure Code Sequence                    | (0032,1064) | SQ |       | VNAP              |        |
| >>Include 'Code Sequence Macro'                       |             |    |       |                   |        |
| >Reason for the Requested Procedure                   | (0040,1002) | LO |       | ALWAYS            |        |
| >Reason for Requested Procedure Code Sequence         | (0040,100A) | SQ |       | ALWAYS            |        |
| >>Include 'Code Sequence Macro'                       |             |    |       |                   |        |
| Performed Procedure Code Sequence [0..n]              | (0040,A372) | SQ |       | ALWAYS            |        |
| > Code Value  | (0008,0100) | SH |       | ALWAYS            |        |
| > Coding Scheme Designator                            | (0008,0102) | SH |       | ALWAYS            |        |
| > Coding Scheme Version                               | (0008,0103) | SH |       | ALWAYS            |        |
| > Code Meaning  | (0008,0104) | LO |       | ALWAYS            |        |
| Current Requested Procedure Evidence Sequence         | (0040,A375) | SQ |       | ALWAYS            |        |
| >Include ' Hierarchical SOP Instance Reference Macro' |             |    |       |                   |        |
| Pertinent Other Evidence Sequence                     | (0040,A385) | SQ |       | ALWAYS            |        |
| >Include ' Hierarchical SOP Instance Reference Macro' |             |    |       |                   |        |
| Referenced Instance Sequence                          | (0008,114A) | SQ |       | ALWAYS            |        |
| >Include 'SOP Instance Reference Macro'               |             |    |       |                   |        |
| >Purpose of Reference Code Sequence                   | (0040,A170) | SQ |       | ALWAYS            |        |
| <b>SOP Common</b>                                     |             |    |       |                   |        |
| SOP Class UID   | (0008,0016) | UI |       | ALWAYS            |        |
| SOP Instance UID                                      | (0008,0018) | UI |       | ALWAYS            |        |
| Specific Character Set                                | (0008,0005) | CS |       | ALWAYS            |        |
| Instance Creation Date                                | (0008,0012) | DA |       | ANAP              |        |
| Instance Creation Time                                | (0008,0013) | TM |       | ANAP              |        |
| Instance Creator UID                                  | (0008,0014) | UI |       | ANAP              |        |
| Related General SOP Class UID                         | (0008,001A) | UI |       | ANAP              |        |
| Original Specialized SOP Class UID                    | (0008,001B) | UI |       | ANAP              |        |
| Coding Scheme Identification Sequence                 | (0008,0110) | SQ |       | ANAP              |        |
| >Coding Scheme Designator                             | (0008,0102) | SH |       | ALWAYS            |        |
| >Coding Scheme Registry                               | (0008,0112) | LO |       | ALWAYS            |        |
| >Coding Scheme UID                                    | (0008,010C) | UI |       | ALWAYS            |        |
| >Coding Scheme External ID                            | (0008,0114) | ST |       | VNAP              |        |

| Attribute Name                               | Tag         | VR | Value | Presence of Value | Source |
|--|-------------|----|-------|-------------------|--------|
| >Coding Scheme Name                          | (0008,0115) | ST |       | ANAP              |        |
| >Coding Scheme Version                       | (0008,0103) | SH |       | ANAP              |        |
| >Coding Scheme Responsible Organization      | (0008,0116) | ST |       | ANAP              |        |
| Timezone Offset From UTC                     | (0008,0201) | SH |       | ANAP              |        |
| Contributing Equipment Sequence              | (0018,A001) | SQ |       | ALWAYS            |        |
| > Purpose of Reference Code Sequence         | (0040,A170) | SQ |       | ALWAYS            |        |
| >> Code Value, Code Meaning, ...             |             |    |       |                   |        |
| > Manufacturer                               | (0008,0070) | LO |       | ALWAYS            |        |
| >Institution Name                            | (0008,0080) | LO |       | ANAP              |        |
| >Institution Address                         | (0008,0081) | ST |       | ANAP              |        |
| >Station Name                                | (0008,1010) | SH |       | ANAP              |        |
| >Institutional Department Name               | (0008,1040) | LO |       | ANAP              |        |
| >Operators' Name                             | (0008,1070) | PN |       | ANAP              |        |
| >Operator Identification Sequence            | (0008,1072) | SQ |       | ANAP              |        |
| >>Include 'Person Identification Macro'      |             |    |       | ANAP              |        |
| >Manufacturer's Model Name                   | (0008,1090) | LO |       | ANAP              |        |
| >Device Serial Number                        | (0018,1000) | LO |       | ANAP              |        |
| >Software Versions                           | (0018,1020) | LO |       | ANAP              |        |
| >Spatial Resolution                          | (0018,1050) | DS |       | ANAP              |        |
| >Date of Last Calibration                    | (0018,1200) | DA |       | ANAP              |        |
| >Time of Last Calibration                    | (0018,1201) | TM |       | ANAP              |        |
| >Contribution DateTime                       | (0018,A002) | DT |       | ANAP              |        |
| >Contribution Description                    | (0018,A003) | ST |       | ANAP              |        |
| Instance Number                              | (0020,0013) | IS |       | ANAP              |        |
| SOP Instance Status                          | (0100,0410) | CS |       | ANAP              |        |
| SOP Authorization DateTime                   | (0100,0420) | DT |       | ANAP              |        |
| SOP Authorization Comment                    | (0100,0424) | LT |       | ANAP              |        |
| Authorization Equipment Certification Number | (0100,0426) | LO |       | ANAP              |        |
| > Include 'Digital Signatures Macro'         |             |    |       |                   |        |
| Encrypted Attributes Sequence                | (0400,0500) | SQ |       | ALWAYS            |        |
| > Code sequence                              |             |    |       |                   |        |
| Original Attributes Sequence                 | (0400,0561) | SQ |       | ANAP              |        |
| > Code sequence                              |             |    |       |                   |        |
| HL7 Structured Document Reference Sequence   | (0040,A390) | SQ |       | ALWAYS            |        |

Table 6.1-12: SR Document Content

| Report ID | Line Number | Reference ID | Attribute   | VM  | Req Type | Condition  | Value  | Units    |
|-----------|-------------|--------------|---|-----|----------|--|--|----------|
| 1         |             | TID 10001    | <b>X-Ray Radiation Dose</b>                             | 1   | M        |  |  |          |
| 1.1       | 1           | 113701       | X-Ray Radiation Dose Report Container                   |     |          |  |  |          |
| 1.2       | 2           | 121058       | Procedure reported                                      | 1   | M        |  | "Projection X-Ray" / "Mammography"                         |          |
| 1.2.1     | 3           | G-C0E8       | Has Intent  | 1   | M        |  | "Diagnostic Intent"  |          |
| 1.3       | 4           | TID 1002     | <b>Observer Context</b>                                 | 1-n | M        |  |  |          |
| 1.3.1     | 1           | 121005       | Observer Type   | 1   | MC       | IFF Observer type is device  | "Device"   |          |
| 1.3.2     | 3           | TID 1004     | <b>Device Observer Identifying Attributes Container</b> | 1   | MC       | IFF row 1 value is "Device"  |  |          |
| 1.3.2.1   | 1           | 121012       | Device Observer UID                                     | 1   | M        |  | "1.3.51.0.7.1104921769.4906.39497.36025.41322.41911.30231" |          |
| 1.3.2.2   | 2           | 121013       | Device Observer Name                                    | 1   | U        |  | NX workstation Name  |          |
| 1.3.2.3   | 3           | 121014       | Device Observer Manufacturer                            | 1   | U        |  | "Agfa"   |          |
| 1.3.2.4   | 4           | 121015       | Device Observer Model Name                              | 1   | U        |  | "NX"   |          |
| 1.3.2.5   | 5           | 121016       | Device Observer Serial Number                           | 1   | U        |  | xx-xx-xx-xx-xx-xx  |          |
| 1.4       | 5           | 113705       | Scope of Accumulation                                   | 1   | M        |  | "Performed Procedure Step"                                 |          |
| 1.4.1     | 6           | 121126       | Performed Procedure Step SOP instance UID               | 1   | M        |  | MPPSInstanceUID  |          |
| 1.5       | 7           | TID 10002    | <b>Accumulated X-Ray Dose Data</b>                      | 1   | MC       | IFF single plane system  |  |          |
| 1.5.1     | 1           | 113702       | Accumulated X-Ray Dose Data Container                   | 1   | M        |  |  |          |
| 1.5.2     | 2           | 113764       | Acquisition Plane                                       | 1   | M        |  | "Single Plane"   |          |
| 1.5.3     | 3           | 122505       | <b>Calibration Container</b>                            | 1-n | MC       | IFF Calibration Data is available  |  |          |
| 1.5.3.1   | 4           | 113794       | Dose Measurement Device                                 | 1   | M        |  | "Dosimeter"  |          |
| 1.5.3.2   | 5           | 113723       | Calibration Date  | 1   | M        |  |  |          |
| 1.5.3.3   | 6           | 122322       | Calibration Factor                                      | 1   | M        |  |  | No units |
| 1.5.3.4   | 7           | 113763       | Calibration Uncertainty                                 | 1   | M        |  |  | Percent  |
| 1.5.3.5   | 8           | 113724       | Calibration Responsible Party                           | 1   | M        |  |  |          |
| 1.5.4     | 10          | TID 10004    | <b>Accumulated Projection X-Ray Dose</b>                | 1   | MC       | XOR row 11, 12 ; IFF TID(10001) Row 2 = "Projection X-Ray"                                   |  |          |
| 1.5.4.1   | 1           | 113722       | Dose Area Product Total                                 | 1   | M        |  |  | Gym2     |
| 1.5.4.2   | 2           | 113725       | Dose (RP) Total   | 1   | MC       | IF any of the values of TID (10001) Row 14 are not "MPPS Content". May be present otherwise. |  | Gy       |

| Report ID | Line Number | Reference ID | Attribute                                 | VM  | Req Type | Condition  | Value  | Units |
|-----------|-------------|--------------|---|-----|----------|--|--|-------|
| 1.5.4.3   | 6           | 113727       | Acquisition Dose Area Product Total       | 1   | M        |  |  | Gym2  |
| 1.5.4.4   | 7           | 113729       | Acquisition Dose (RP) Total               | 1   | MC       | IF any of the values of TID (10001) Row 14 are not "MPPS Content". May be present otherwise.                 |  | Gy    |
| 1.5.4.5   | 8           | 113855       | Total Acquisition Time                    | 1   | M        |  |  | s     |
| 1.5.4.6   | 11          | 113780       | Reference Point Definition                | 1   | MC       | IF Row 2 or Row 7 is present AND Row 10 is not present   |  |       |
| 1.5.5     | 11          | TID 10005    | <b>Accumulated Mammography X-Ray Dose</b> | 1   | MC       | XOR row 10, 12 ; IFF TID(10001) Row 2 = "Mammography"  |  |       |
| 1.5.5.1   | 1           | 111637       | Accumulated Average Glandular Dose        | 1-2 | M        |  |  | dGy   |
| 1.5.5.2   | 2           | G-C171       | Laterality                                | 1   | M        |  | "Left breast" / "Right breast" / "Both breasts"                    |       |
| 1.6       | 10          | TID 10003    | <b>Irradiation Event X-Ray Data</b>       | 1-n | M        |  |  |       |
| 1.6.1     | 1           | 113706       | Irradiation Event X-Ray Data Container    | 1   | M        |  |  |       |
| 1.6.2     | 2           | 113764       | Acquisition Plane                         | 1   | M        |  | "Single Plane"   |       |
| 1.6.3     | 3           | 111526       | DateTime Started                          | 1   | M        |  |  |       |
| 1.6.4     | 4           | 113721       | Irradiation Event Type                    | 1   | M        |  | "Stationary Acquisition"   |       |
| 1.6.4.1   | 7           | G-C171       | Laterality                                | 1   | UC       | IF anatomy is bi-lateral   | "Right" / "Left" / "Right and left" / "Unilateral"                 |       |
| 1.6.5     | 8           | 113780       | Reference Point Definition                | 1   | MC       | IF Row 13 or Row 14 is present AND Row 9 is not present  |  |       |
| 1.6.6     | 10          | 113769       | Irradiation Event UID                     | 1   | M        |  |  |       |
| 1.6.7     | 11          | 122130       | Dose Area Product                         | 1   | MC       | IFF TID(10001) Row 2 = "Projection X-Ray"  |  | Gym2  |
| 1.6.8     | 12          | 111631       | Average Glandular Dose                    | 1   | MC       | IFF TID(10001) Row 2 = "Mammography"   |  | dGy   |
| 1.6.9     | 13          | 113738       | Dose (RP)                                 | 1   | MC       | IFF TID(10001) Row 2 = "Projection X-Ray" AND any of the values of TID (10001) Row 14 are not "MPPS Content" |  | Gy    |
| 1.6.10    | 14          | 111636       | Entrance Exposure at RP                   | 1   | MC       | IFF TID(10001) Row 2 = "Mammography"   |  | mGy   |
| 1.6.11    | 15          | 112011       | Positioner Primary Angle                  | 1   | UC       | XOR row 19   |  | °     |
| 1.6.12    | 16          | 112012       | Positioner Secondary Angle                | 1   | UC       | XOR row 19   |  | °     |
| 1.6.13    | 20          | 113790       | Collimated Field Area                     | 1   | U        |  |  | m2    |
| 1.6.14    | 21          | 113771       | <b>X-Ray Filters Container</b>            | 1-n | U        |  |  |       |
| 1.6.14.1  | 22          | 113772       | X-Ray Filter Type                         | 1   | U        |  | "Strip filter" / "Wedge filter" / "Butterfly filter" / "No filter" |       |

| Report ID | Line Number | Reference ID | Attribute                      | VM  | Req Type | Condition   | Value                       | Units    |
|-----------|-------------|--------------|--------------------------------|-----|----------|---|-----------------------------|----------|
| 1.6.14.2  | 23          | 113757       | X-Ray Filter Material          | 1   | U        |   |                             |          |
| 1.6.14.3  | 24          | 113758       | X-Ray Filter Thickness Minimum | 1   | U        |   |                             | mm       |
| 1.6.14.4  | 25          | 113773       | X-Ray Filter Thickness Maximum | 1   | U        |   |                             | mm       |
| 1.6.15    | 30          | 113733       | KVP                            | 1-n | U        |   |                             | kV       |
| 1.6.16    | 31          | 113734       | X-Ray Tube Current             | 1-n | U        |   |                             | mA       |
| 1.6.17    | 32          | 113735       | Exposure Time                  | 1   | U        |   |                             | ms       |
| 1.6.18    | 34          | 113736       | Exposure                       | 1   | U        |   |                             | uAs      |
| 1.6.19    | 35          | 113766       | Focal Spot Size                | 1   | U        |   |                             | mm       |
| 1.6.20    | 36          | 113742       | Irradiation Duration           | 1   | U        |   |                             | mA       |
| 1.6.21    | 37          | 113767       | Average X-Ray Tube Current     | 1   | U        |   |                             |          |
| 1.6.22    | 45          | 123014       | Target Region                  | 1   | M        |   |                             |          |
| 1.6.23    | 46          | 111632       | Anode Target Material          | 1   | U        |   |                             |          |
| 1.6.24    | 47          | 111633       | Compression Thickness          | 1   | U        |   |                             | mm       |
| 1.6.25    | 53          | TID 1021     | <b>Device Participant</b>      | 1   | MC       | Required if the irradiating device is not the recording device<br>XOR row 54 ; IF TID (10001) Row 2 is not "Cassette-based X-Ray" |                             |          |
| 1.6.25.1  | 1           | 113876       | Device Role in Procedure       | 1   | M        |   | "Irradiating Device"        |          |
| 1.6.25.2  | 2           | 113877       | Device Name                    | 1   | U        |   |                             |          |
| 1.6.25.3  | 3           | 113878       | Device Manufacturer            | 1   | M        |   |                             |          |
| 1.6.25.4  | 4           | 113879       | Device Model Name              | 1   | M        |   |                             |          |
| 1.6.25.5  | 5           | 113880       | Device Serial Number           | 1   | M        |   |                             |          |
| 1.6.25.6  | 6           | 121012       | Device Observer UID            | 1   | M        |   |                             |          |
| 1.6.26    | 55          | 113795       | Acquired Image                 | 1-n | MC       | IFF image object is created for this irradiation event  |                             |          |
| 1.6.27    | 56          | 113845       | Exposure Index                 | 1   | U        |   |                             | No units |
| 1.6.28    | 57          | 113846       | Target Exposure Index          | 1   | U        |   |                             | No units |
| 1.6.29    | 58          | 113847       | Deviation Index                | 1   | U        |   |                             | No units |
| 1.7       | 14          | 113854       | Source of Dose Information     | 1-n | M        |   | "Automated Data Collection" |          |

## 6.2 Data dictionary of Private Attributes

**Table 6.2-1: Private Attributes**

| Tag           | VR | VM    | Attribute Name                       | Meaning / Values   | Presence of Value                                      | SOP Class                  |
|---------------|----|-------|--------------------------------------|--|--|----------------------------|
| (0019,0010)   | LO | 1     | Private Creator                      | Agfa ADC NX  | ALWAYS   | -                          |
| (0019,1001)   | FL | n     | Contrast                             |  | Musica on Pacs : enabled                               | GSPS                       |
| (0019,1002)   | FL | n     | Brightness                           |  | Musica on Pacs : enabled                               | GSPS                       |
| (0019,1003)   | FL | 1     | Sharpness                            |  | Musica on Pacs : enabled                               | GSPS                       |
| (0019,1004)   | LO | 1     | Package                              |  | Musica on Pacs : enabled                               | GSPS                       |
| (0019,1005)   | IS | 1     | Processing Version                   | Indicates the version of the package used for processing the image | Musica on Pacs : enabled<br>Only present if version >1 | GSPS                       |
| (0019,1006)   | FL | 1     | Border Contrast Density              |  | Musica on Pacs : enabled                               | GSPS                       |
| (0019,1007)   | CS | 1     | Is Speedclass Free                   | Required for Auto-QC2<br>Enum: YES / NO                            | ALWAYS   | CR,DX,MG                   |
| (0019,1008)   | CS | 1     | Applicationtype                      |  | Musica on Pacs : enabled                               | CR,DX,MG                   |
| (0019,1009)   | SQ |       | Collimator shape sequence            | Collimator Presentation Selection List                             | Musica on Pacs : enabled                               | GSPS                       |
| > (0019,0010) | LO | 1     | Identification Code                  | "Agfa ADC NX"  | Musica on Pacs : enabled                               | GSPS                       |
| > (0019,100D) | CS | 1     | Is Diagnostic Area Inside            | Enum : YES/NO  | Musica on Pacs : enabled                               | GSPS                       |
| > (0018,1700) | CS | 1(-3) | Collimator Shape                     | "RECTANGULAR", "POLYGONAL"   | Musica on Pacs : enabled                               | GSPS                       |
| > (0018,1702) | IS | 1     | Collimator Left Vertical Edge        | If "RECTANGULAR"   | Musica on Pacs : enabled                               | GSPS                       |
| > (0018,1704) | IS | 1     | Collimator Right Vertical Edge       | If "RECTANGULAR"   | Musica on Pacs : enabled                               | GSPS                       |
| > (0018,1706) | IS | 1     | Collimator Upper Horizontal Edge     | If "RECTANGULAR"   | Musica on Pacs : enabled                               | GSPS                       |
| > (0018,1708) | IS | 1     | Collimator Lower Horizontal Edge     | If "RECTANGULAR"   | Musica on Pacs : enabled                               | GSPS                       |
| > (0018,1720) | IS | 2-2n  | Vertices of the Polygonal Collimator | List of vertices, if "POLYGONAL"                                   | Musica on Pacs : enabled                               | GSPS                       |
| (0019,100A)   | CS | 1     | Background darkening calculate       | Enum : YES/NO  | Musica on Pacs : enabled                               | GSPS                       |
| (0019,100B)   | CS | 1     | Background darkening padding smooth  | Enum : YES/NO  | Musica on Pacs : enabled                               | GSPS                       |
| (0019,100C)   | CS | 1     | Is Invert Enabled                    | Enum : YES/NO  | Musica on Pacs : enabled                               | GSPS                       |
| (0019,1010)   | ST | 1     | Image Processing Parameters          | Required for Auto-QC2 & Musica 1 (& 2)                             | When importing AutoQC2 Exam Tree                       | DX for Proc<br>MG for Proc |
| (0019,1020)   | CS | 1     | Is Border Contrast Enabled           | Enum: YES / NO   | Musica on Pacs : enabled                               | GSPS                       |
| (0019,1021)   | FL | 1     | Calibration Factor                   | Required for QC-Mammo  | ALWAYS   | CR,DX,MG                   |

| Tag         | VR | VM | Attribute Name         | Meaning / Values   | Presence of Value                 | SOP Class |
|-------------|----|----|------------------------|--|-----------------------------------|-----------|
| (0019,1026) | LO | n  | Stitch Zone Parameters | Used in NX for Sticking up to 4 images <sup>24</sup><br>Quality label for every stitch zone :<br>Vxx.x Hxx.x M Y<br>G N<br>A <sup>25</sup> | Case of image Sticking            | CR,DX,MG  |
| (0019,1028) | CS | 1  | Is IPSR enabled        | Enum: YES / NO   | If IPSR license is enabled        | CR,DX,MG  |
| (0019,10B0) | LO | 1  | Target AETitle         | Used for MusicaMicro: Target AETitle   | If Musica Micro License is active | MG        |
| (0019,10B1) | LO | 1  | Target IP address      | Used for MusicaMicro: Target IP address  | If Musica Micro License is active | MG        |
| (0019,10B2) | LO | 1  | Target Port number     | Used for MusicaMicro: Target Port number   | If Musica Micro License is active | MG        |
| (0019,10B3) | LO | 1  | Session UID            | Used for MusicaMicro: Session UID  | If Musica Micro License is active | MG        |
| (0019,10D0) | LO | 1  | Custom Patient Field 1 | Used to map unknown DICOM attributes from the RIS (Worklist). Stored on ID Session level   | Configurable                      | CR,DX,MG  |
| (0019,10D1) | LO | 1  | Custom Patient Field 2 |  |                                   | CR,DX,MG  |
| (0019,10D2) | LO | 1  | Custom Patient Field 3 |  |                                   | CR,DX,MG  |
| (0019,10D3) | LO | 1  | Custom Patient Field 4 |  |                                   | CR,DX,MG  |
| (0019,10D4) | LO | 1  | Custom Patient Field 5 |  |                                   | CR,DX,MG  |
| (0019,10E0) | LO | 1  | Custom Image Field 1   | Used to map unknown DICOM attributes from the RIS (Worklist). Stored on Workflow X-RAY image level   | Configurable                      | CR,DX,MG  |
| (0019,10E1) | LO | 1  | Custom Image Field 2   |  |                                   | CR,DX,MG  |
| (0019,10E2) | LO | 1  | Custom Image Field 3   |  |                                   | CR,DX,MG  |
| (0019,10E3) | LO | 1  | Custom Image Field 4   |  |                                   | CR,DX,MG  |
| (0019,10E4) | LO | 1  | Custom Image Field 5   |  |                                   | CR,DX,MG  |
| (0019,10F0) | LO | 1  | User defined field 1   | Used to map unknown DICOM attributes from the RIS (Worklist). Stored on SPS level.   | Configurable                      | CR,DX,MG  |
| (0019,10F1) | LO | 1  | User defined field 2   |  |                                   | CR,DX,MG  |
| (0019,10F2) | LO | 1  | User defined field 3   |  |                                   | CR,DX,MG  |
| (0019,10F3) | LO | 1  | User defined field 4   |  |                                   | CR,DX,MG  |
| (0019,10F4) | LO | 1  | User defined field 5   |  |                                   | CR,DX,MG  |

<sup>24</sup> Stitch Zone Parameters (0019,1026) - VR=LO - VM=1-n - Type=VR=LO

Used in CR, DX or MG images together with the 'Basic Pixel Spacing Calibration Macro' :

- Pixel spacing (0028,0030) - VR=DS - VM=2 - Type=1C (C: if image is calibrated)
- Pixel spacing Calibration Type (0028,0A02) - VR=CS - VM=1 - Type=3
- Pixel spacing Calibration Description (0028,0A04) - VR=LO - VM=1 - Type=1C (C: if Type is present)

<sup>25</sup> Where:

- Vxx.x Hxx.x: the deviation in cm from grid only stitched position in Vertical / Horizontal direction.
- M G A: indicates the way the images are stitched: Manually corrected or manipulated / GRID auto stitched / Anatomical auto stitched.
- Y N: result of patient movement estimation: Yes or No.

| Tag         | VR | VM | Attribute Name            | Meaning / Values                                  | Presence of Value                   | SOP Class |
|-------------|----|----|---------------------------|---|-------------------------------------|-----------|
| (0019,10F5) | CS | 1  | Cassette orientation      | Required for Auto-QC2<br>LANDSCAPE or PORTRAIT    | ALWAYS                              | DX        |
| (0019,10F6) | DS | 1  | Image plate sensitivity   | Required for Auto-QC2                             | ALWAYS                              | DX        |
| (0019,10F7) | DS | 1  | Image plate erasability   | Required for Auto-QC2                             | ALWAYS                              | DX        |
| (0019,10F8) | IS | 1  | Breast density percentage | Is filled-in by IPD in case the NX                | If Breast Density license is active | MG        |
| (0019,10FA) | IS | 1  | Exposure index            | Related with Xray Exposure Dose <sup>26</sup>     | ALWAYS                              | CR,DX,MG  |
| (0019,10FB) | FL | 1  | Deviation index           | Related with Xray Exposure Dose <sup>26</sup>     | ALWAYS                              | CR,DX,MG  |
| (0019,10FC) | IS | 1  | Target exposure index     | Related with Xray Exposure Dose <sup>26</sup>     | ALWAYS                              | CR,DX,MG  |
| (0019,10FD) | CS | 1  | Is overexposed            | Related with Xray Exposure Dose Enum:<br>YES / NO | ALWAYS                              | CR,DX,MG  |
| (0019,10FE) | CS | 1  | Study Priority ID         | Enum: HIGH, MED (=default), LOW,<br>NULL          | ALWAYS                              | CR,DX,MG  |

<sup>26</sup> Official tags were created to replace the private tags.  
Now both are used for compatibility reason but contain the same information.

| Name                  | Tag                             | Type | VR           | Units   | VM |
|-----------------------|---------------------------------|------|--------------|---------|----|
| Exposure Index        | (0018, 1411) <i>(0019,10FA)</i> | 3    | DS <i>IS</i> | 100.uGy | 1  |
| Target Exposure Index | (0018, 1412) <i>(0019,10FC)</i> | 3    | DS <i>IS</i> | 100.uGy | 1  |
| Deviation Index       | (0018, 1413) <i>(0019,10FB)</i> | 3    | DS <i>FL</i> | -       | 1  |



## 6.3 Attribute Mapping

The relationships between attributes received via Modality Worklist, stored in acquired images and communicated via MPPS are summarized in Table 6.3-1.

**Table 6.3-1: Attribute Mapping between Modality Worklist, Image and MPPS**

| Modality Worklist                     | Image IOD                                | MPPS IOD                                 |
|---------------------------------------|--|--|
| Patient Name                          | Patient Name                             | Patient Name                             |
| Patient ID                            | Patient ID                               | Patient ID                               |
| Patient's Birth Date                  | Patient's Birth Date                     | Patient's Birth Date                     |
| Patient's Sex                         | Patient's Sex                            | Patient's Sex                            |
| Patient's Weight                      | Patient's Weight                         |  |
| Referring Physician's Name            | Referring Physician's Name               |  |
| ----                                  | ----                                     | Scheduled Step Attributes Sequence       |
| Study Instance UID                    | Study Instance UID                       | > Study Instance UID                     |
| Referenced Study Sequence             | Referenced Study Sequence                | > Referenced Study Sequence              |
| Accession Number                      | Accession Number                         | > Accession Number                       |
| ----                                  | Request Attributes Sequence              | ----                                     |
| Requested Procedure ID                | > Requested Procedure ID                 | > Requested Procedure ID                 |
| Requested Procedure Description       |  | > Requested Procedure Description        |
| Scheduled Procedure Step ID           | > Scheduled Procedure Step ID            | > Scheduled Procedure Step ID            |
| Scheduled Procedure Step Description  | > Scheduled Procedure Step Description   | > Scheduled Procedure Step Description   |
| Scheduled Protocol Code Sequence      | > Scheduled Protocol Code Sequence       | ----                                     |
| ----                                  | Performed Protocol Code Sequence         | Performed Protocol Code Sequence         |
| ----                                  | Study ID                                 | Study ID                                 |
| ----                                  | Performed Procedure Step ID              | Performed Procedure Step ID              |
| ----                                  | Performed Procedure Step Start Date      | Performed Procedure Step Start Date      |
| ----                                  | Performed Procedure Step Start Time      | Performed Procedure Step Start Time      |
| ----                                  | Performed Procedure Step Description     | Performed Procedure Step Description     |
| ----                                  | Comments on the Performed Procedure Step | Comments on the Performed Procedure Step |
| ----                                  | ----                                     | Performed Series Sequence                |
| Scheduled Performing Physician's Name | Performing Physician's Name              | > Performing Physician's Name            |
| Requested Procedure Code Sequence     | ----                                     | Procedure Code Sequence                  |
| ----                                  | Referenced Study Component Sequence      | ----                                     |
| ----                                  | > Referenced SOP Class UID               | SOP Class UID                            |
| ----                                  | > Referenced SOP Instance UID            | SOP Instance UID                         |
| ----                                  | Protocol Name                            | Protocol Name                            |

## 6.4 Grayscale Image Consistency

The display monitor attached to NX 3.0.8950 can be calibrated according to the Grayscale Standard Display Function (GSDf).

## 6.5 Veterinary Extensions

This section lists the DICOM attributes that are supported by NX 3.0.8950 for the veterinary market. The attributes have been incorporated based on DICOM correction item CP-643.

| Attribute Name                 | Tag         | VR | Value | Presence of Value | Source   |
|--------------------------------|-------------|----|-------|-------------------|----------|
| <b>Patient</b>                 |             |    |       |                   |          |
| Patient species description    | (0010,2201) | LO |       |                   | User/MWL |
| Patient breed description      | (0010,2292) | LO |       | ANAP              | User/MWL |
| Breed registration sequence    | (0010,2294) | SQ |       | ANAP              | User/MWL |
| > Breed registration number    | (0010,2295) | LO |       |                   | User/MWL |
| > Breed registry code sequence | (0010,2296) | SQ |       |                   | User/MWL |
| Responsible person             | (0010,2297) | PN |       | ANAP              | MWL      |
| Responsible person Role        | (0010,2298) | PN |       | ANAP              | MWL      |
| Responsible organization       | (0010,2299) | LO |       | ANAP              | MWL      |
| Patient's Sex Neutered         | (0010,2203) | CS |       | ANAP              | User/MWL |