

# AGFA HEALTHCARE HL7 Conformance Statement

## → Enterprise Imaging 8.3.x

Document No. 001663 - Revision 5

Livelihood NodeID: 81707306

**When printed, this is NOT a controlled copy**

## Document Information

---

<b>Service-related contact information worldwide</b>	All service-related contact information is available on this URL→	<a href="https://www.agfahealthcare.com/contact/">https://www.agfahealthcare.com/contact/</a>
--	---	---

Issued by:  
Agfa HealthCare  
V&V Connectivity  
Septestraat 27  
B-2640 Mortsel  
Belgium

email : [connectivity@agfa.com](mailto:connectivity@agfa.com)

Agfa HealthCare shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment by Agfa personnel.

**Copyright © 2023**  
**Agfa HealthCare**  
**All rights reserved**

## Table of Contents

---

1	Introduction .....	4
1.1	Revision Record .....	4
1.2	Purpose and Intended Audience of this Document .....	4
1.3	Acronyms and Abbreviations .....	5
1.4	Related Documents .....	6
2	Inbound Messages .....	7
2.1	Core Server supported Trigger Events .....	7
2.1.1	Supported ACK Events .....	7
2.1.1.1	Supported ACK segments .....	7
2.1.2	Supported ADT Events .....	7
2.1.2.1	Supported ADT segments .....	8
2.1.3	Supported MFN Events .....	8
2.1.3.1	Supported MFN segments .....	8
2.1.4	Supported ORM Events .....	9
2.1.4.1	Supported ORM segments .....	9
2.1.5	Supported ORU Events .....	9
2.1.5.1	Supported ORU segments .....	9
2.2	Core Server supported fields in Inbound segments .....	10
2.2.1	MSH Segment .....	10
2.2.2	MSA Segment .....	10
2.2.3	ERR Segment .....	11
2.2.4	EVN Segment .....	11
2.2.5	PID Segment .....	11
2.2.6	PD1 Segment .....	14
2.2.7	AL1 Segment .....	14
2.2.8	MRG Segment .....	15
2.2.9	PV1 Segment .....	15
2.2.10	ORC Segment .....	17
2.2.11	OBR Segment .....	20
2.2.12	OBX Segment .....	24
2.2.13	ZDS Segment .....	28
2.2.14	NTE Segment .....	29
2.2.15	MFI Segment .....	29
2.2.16	MFE Segment .....	29
2.2.17	STF Segment .....	29
2.2.18	PRA Segment .....	31
3	Outbound Messages .....	32
3.1	Core Server supported Trigger Events .....	32
3.1.1	Supported ACK Events .....	32
3.1.2	Supported DFT Events .....	32
3.1.3	Supported ORM Events .....	32
3.1.4	Supported ORU Events .....	32
4	Queries .....	34
4.1	Web Server Query .....	34
4.1.1	Supported Query events .....	34

# 1 INTRODUCTION

## 1.1 Revision Record

HL7 Conformance Statement Template, Livelihood Node ID: 19565451		
Revision Number	Date	Reason for Change
8	March 2019	New structure for better quality communication with external parties

HL7 Conformance Statement Enterprise Imaging 8.3.x		
Revision Number	Date	Reason for Change
2	June 2023	Initial Revision for EI version 8.3.0.000
3-4	September 2023	Update for EI version 8.3.1.000 -update comment for PID-7.1 (only the date is stored) -update comment for AL1-6 -update comment for OBX-14 (only date is supported for Pregnancy) -DFT outbound FT1-20 & 21: removal of superfluous sub-components in the profile
5	December 2023	Update for EI version 8.3.2.000 Added MSH-6 and MSH-6.1 as optional to match the Conformance Profiles

## 1.2 Purpose and Intended Audience of this Document

This document is a HL7 Conformance Statement for the HL7 Services of Enterprise Imaging 8.3.x.

The user of this document is involved with system integration and/or software design. We assume that the reader is familiar with the terminology and concepts that are used in the HL7 standard and the IHE Technical Framework.

Readers not familiar with HL7 terminology should first read the appropriate parts of the HL7 standard itself, prior to reading this conformance statement.

Although the use of this conformance statement in conjunction with the HL7 standard is intended to facilitate communication with Enterprise Imaging 8.3.x, it is not sufficient to guarantee, by itself, the inter-operation of the connection between Enterprise Imaging 8.3.x and the 3<sup>rd</sup> party HL7-based system.

To help this integration validation we provide HL7 Conformance Profiles bundled together in a zip file "[001664 Enterprise Imaging 8.3.x HL7 Conformance Profiles](#)", also available on <http://www.agfahealthcare.com/hl7>. These conformance profiles are related to the Core Server component and the Web Server component.

**Note:**

*Core Server might also support trigger events that are not included in this conformance statement, provided that the message content can be converted to one of the supported trigger events. For example, Core Server has limited support for SIU messages since some SIU trigger events can be converted to ORM trigger events (with potential loss of some SIU data which isn't supported in an ORM trigger event).*

*Another example is the limited support for patient linking since the A24 trigger event can be converted to an A40 trigger event.*

*This requires a thorough analysis by Agfa Healthcare Professional Services that ultimately determine whether or not Core Server can support a non-listed trigger event.*

The integration of any device into a system of interconnected devices goes beyond the scope of the HL7 standard and this conformance statement when *interoperability* is desired. The responsibility for analyzing the applications requirements and developing a solution that integrates the Agfa equipment with other vendors' systems is the user's responsibility and should not be underestimated.

### 1.3 Acronyms and Abbreviations

Definitions, terms and abbreviations used in this document, many of which are defined within the HL7 standard. Abbreviations and terms are as follows:

ACK	Acknowledgement message
ADT	Admission, Discharge, and Transfer message
AL1	Patient Allergy Information segment
DFT	Detail Financial Transaction message
ERR	Error segment
EVN	Event Type segment
HL7	Health Level 7
IHE	Integrating the Healthcare Enterprise
MFE	Master File Entry segment
MFI	Master Files Identification segment
MFN	Master Files Notification message
MRG	Merge Patient Information segment
MSA	Message Acknowledgment segment
MSH	Message Header segment
NTE	Notes and comments segment
OBR	Observation Request segment
OBX	Observation/Result segment
OF	Order Filler – Device that sends filled orders (ORM)
OP	Order Placer
ORC	Common Order segment
ORM	Order Request Message
ORU	Observation Results - Unsolicited message
PD1	Patient Additional Demographic segment
PID	Patient ID segment
PRA	Practitioner Detail segment
PV1	Patient Visit segment
QRD	Query Definition segment
RIS	Radiology Information System
SCH	Scheduling Activity Information segment
SIU	Scheduling Information Unsolicited message
STF	Staff Identification segment

ZDS            EI Study Instance Identification segment

## 1.4 Related Documents

- [HL7 Standard v2.3.1, 2.4, 2.5](#)
- [IHE Radiology Technical Framework](#)
- [IHE IT Infrastructure Technical Framework](#)

## 2 INBOUND MESSAGES

### 2.1 Core Server supported Trigger Events

#### 2.1.1 Supported ACK Events

For more information about the processed segments and supported fields, refer to the document specified in the 'Conformance profile document' column or to chapter 2.1.1.1 and 2.2.

**Table 2-1 Supported ACK Events**

Func Area	Event Code	ACK Trigger Event Description	Conformance profile document
ACK	ALL	General ACK Message	ENTERPRISE IMAGING INBOUND ACK.pdf

#### 2.1.1.1 Supported ACK segments

The following segments are supported when Enterprise imaging receives an ACK message:

Segments without brackets are mandatory. Segment with square brackets are optional and segments with curly brackets (braces) are repeating.

**Table 2-2 Supported ACK segments**

HL7 ACK Event Code	Supported ACK segments
All	MSH MSA [ERR]

#### 2.1.2 Supported ADT Events

For more information about the processed segments and supported fields, refer to the document specified in the 'Conformance profile document' column or to chapter 2.1.2.1 and 2.2.

**Table 2-3 Supported ADT Events**

Func Area	Event Code	ADT Trigger Event Description	Conformance profile document
ADT	A01	Admit / Visit notification	ENTERPRISE IMAGING INBOUND ADT_A01.pdf
ADT	A02	Transfer a patient	ENTERPRISE IMAGING INBOUND ADT_A02.pdf
ADT	A03	Discharge / End visit	ENTERPRISE IMAGING INBOUND ADT_A03.pdf
ADT	A04	Register a patient	ENTERPRISE IMAGING INBOUND ADT_A04.pdf
ADT	A05	Pre-admit a patient	ENTERPRISE IMAGING INBOUND ADT_A05.pdf
ADT	A06	Change an outpatient to an inpatient	ENTERPRISE IMAGING INBOUND ADT_A06.pdf
ADT	A07	Change an inpatient to an outpatient	ENTERPRISE IMAGING INBOUND ADT_A07.pdf
ADT	A08	Update patient information	ENTERPRISE IMAGING INBOUND ADT_A08.pdf
ADT	A11	Cancel admit/visit notification	ENTERPRISE IMAGING INBOUND ADT_A11.pdf
ADT	A12	Cancel transfer	ENTERPRISE IMAGING INBOUND ADT_A12.pdf
ADT	A13	Cancel discharge/end visit	ENTERPRISE IMAGING INBOUND ADT_A13.pdf
ADT	A28	Add person information	ENTERPRISE IMAGING INBOUND ADT_A28.pdf
ADT	A31	Update person information	ENTERPRISE IMAGING INBOUND ADT_A31.pdf

Func Area	Event Code	ADT Trigger Event Description	Conformance profile document
ADT	A38	Cancel pre-admit	ENTERPRISE IMAGING INBOUND ADT_A38.pdf
ADT	A40	Merge patient - patient identifier list	ENTERPRISE IMAGING INBOUND ADT_A40.pdf
ADT	A47	Change patient identifier list	ENTERPRISE IMAGING INBOUND ADT_A47.pdf

### 2.1.2.1 Supported ADT segments

The following segments are supported when Enterprise imaging receives an ADT message:

Segments without brackets are mandatory. Segments with square brackets are optional and segments with curly brackets (braces) are repeating.

**Table 2-4 Supported ADT segments**

HL7 ADT Event Code	Supported ADT segments
A01; A03; A05; A06; A07; A13	MSH EVN PID [PD1] PV1 {[OBX]} {[AL1]}
A02, A11, A12, A38	MSH EVN PID [PD1] PV1 {[OBX]}
A04; A08; A28, A31	MSH EVN PID [PD1] [PV1] {[OBX]} {[AL1]}
A40, A47	MSH EVN PID [PD1] MRG

### 2.1.3 Supported MFN Events

For more information about the processed segments and supported fields, refer to the document specified in the 'Conformance profile document' column or to chapter 2.1.3.1 and 2.2.

**Table 2-5 Supported MFN Events**

Func Area	Event Code	MFN Trigger Event Description	Conformance profile document
MFN	M02	Staff/Practitioner master file	ENTERPRISE IMAGING INBOUND MFN_M02.pdf

#### 2.1.3.1 Supported MFN segments

The following segments are supported when Enterprise imaging receives an MFN message:

Segments without brackets are mandatory. Segments with square brackets are optional and segments with curly brackets (braces) are repeating.

**Table 2-6 Supported MFN segments**

HL7 MFN Event Code	Supported MFN segments
M02	MSH MFI { MFE STF [PRA] }



## 2.1.4 Supported ORM Events

For more information about the processed segments and supported fields, refer to the document specified in the 'Conformance profile document' column or to chapter 2.1.4.1 and 2.2.

**Table 2-7 Supported ORM Events**

Func Area	Event Code	ORM Trigger Event Description	Conformance profile document
ORM	O01	Order message	ENTERPRISE IMAGING INBOUND ORM_O01.pdf

### 2.1.4.1 Supported ORM segments

The following segments are supported when Enterprise imaging receives an ORM message:

Segments without brackets are mandatory. Segments with square brackets are optional and segments with curly brackets (braces) are repeating.

**Table 2-8 Supported ORM segments**

HL7 ORM Event Code	Supported ORM segments
O01	MSH PID [PD1] [PV1] {[AL1]} ORC OBR {[NTE]} {[OBX]} [ZDS]

## 2.1.5 Supported ORU Events

For more information about the processed segments and supported fields, refer to the document specified in the 'Conformance profile document' column or to chapter 2.1.5.1 and 2.2.

**Table 2-9 Supported ORU Events**

Func Area	Event Code	ORU Trigger Event Description	Conformance profile document
ORU	R01	Unsolicited transmission of an observation Message	ENTERPRISE IMAGING INBOUND ORU_R01.pdf
ORU	R01	Unsolicited transmission of an observation Message	ENTERPRISE IMAGING INBOUND ORU_R01.pdf

### 2.1.5.1 Supported ORU segments

The following segments are supported when Enterprise imaging receives an ORU message:

Segments without brackets are mandatory. Segments with square brackets are optional and segments with curly brackets (braces) are repeating.

**Table 2-10 Supported ORU segments**

HL7 ORU Event Code	Supported ORU segments
R01 - default	MSH PID [PD1] [PV1] {[ORC] OBR [ZDS] [NTE]} { Only the segments of the last iteration ORDER_OBSERVATION group will be implemented. OBX } }

HL7 ORU Event Code	Supported ORU segments
R01 (Heartstation)	MSH PID [PD1] [PV1] {[ORC] OBR [ZDS] [NTE] {OBX}}

## 2.2 Core Server supported fields in Inbound segments

The following sub-chapters will indicate what are the Required and optional fields supported by Enterprise Imaging for each segment.

The supported segments per message types are listed in the chapter 2.1

### 2.2.1 MSH Segment

*Table 2-11 MSH Fields Support in Enterprise Imaging*

Seq	HL7 Field Name	Option	Comments
1	Field Separator	R	Usually " "
2	Encoding Characters	R	Usually "^~\&"
3	Sending Application	O	
3.1	>namespace ID	O	
4	Sending Facility	O	
4.1	>namespace ID	O	
6	Receiving Facility	O	
6.1	>namespace ID	O	
7	Date/time of Message	O	
7.1	>Date/Time	O	YYYYMMDD[HHMM[SS]][+-ZZZZ] If timezone information is present, it is taken into account for all date/time fields in the message. If timezone information is present also in other datetime fields, it overrules the timezone information in this field.
9	Message Type	R	
9.1	>Message type	R	ADT or ORM or ORU or MFN of ACK
9.2	>Trigger event	R	See supported "event codes" in chapter 2.1.1, 2.1.2, 2.3.2, 2.4.2, 2.5.2
9.3	>Message structure	O	
10	Message Control ID	R	
11	Processing ID	O	
11.1	>Processing ID	O	
12	Version ID	R	Supported versions: 2.3 - 2.3.1 - 2.4 - 2.5.
12.1	>Version ID	R	
18	Character Set	O	

### 2.2.2 MSA Segment

*Table 2-12 MSA Fields Support in Enterprise Imaging*

Seq	HL7 Field Name	Option	Comments
1	Acknowledgment Code	R	Supported values: AA, AR, AE
2	Message Control ID	R	The message control id of the original message

## 2.2.3 ERR Segment

**Table 2-13 ERR Fields Support in Enterprise Imaging**

Seq	HL7 Field Name	Option	Comments
1	Error Code and Location	O	Indicates the position in the message which caused the processing problem
1.1	>Segment ID	O	
1.2	>Segment sequence	O	
1.3	>Field position	O	
1.4	>Code identifying error	O	
1.4.1	>>Identifier	O	
1.4.5	>>Alternate text	O	Description of the error in the processing application
2	Error location	O	Indicates the position in the message which caused the processing problem. Applicable for HL7 version 2.5.
2.1	>Segment ID	O	
2.2	>Segment sequence	O	
2.3	>Field position	O	
2.4	>Field repetition	O	
2.5	>Component number	O	
2.6	>Sub-component number	O	
3	HL7 Error Code	O	Applicable for HL7 version 2.5
3.1	>Identifier	O	
3.2	>text	O	
4	<b>Severity</b>	<b>R</b>	<b>Applicable for HL7 version 2.5</b>
5	<b>Application Error Code</b>	<b>R</b>	<b>The code and description of the error in the processing application. Applicable for HL7 version 2.5</b>
5.1	>Identifier	R	
5.2	>Text	R	

## 2.2.4 EVN Segment

**Table 2-14 EVN Fields Support in Enterprise Imaging**

Seq	HL7 Field Name	Option	Comments
2	Recorded Date/Time	R	
2.1	>Date/Time	R	YYYYMMDD[HHMM[SS]][+ZZZZ] If timezone information is present it overrules the timezone information if available in MSH-7

## 2.2.5 PID Segment

**Table 2-15 PID Fields Support in Enterprise Imaging**

Seq	HL7 Field Name	Option	Comments
3	Patient Identifier List	R	
3.1	>ID	R	Identification number or code of the patient

Seq	HL7 Field Name	Option	Comments
3.4	>Assigning Authority	O	The "Assigning authority" aka "Issuer of Patient ID" can consist of the following three parts: - Namespace ID // Issuer of Patient ID - Universal ID // Universal Entity ID - Universal ID Type // Universal Entity ID Type Namespace ID, or combination of Universal ID and Universal ID Type, unique define the issuer. Either only Namespace is specified, or only Universal ID and Universal ID Type are specified, or all three components are specified. The combination of "Namespace ID" - "Universal ID" - "Universal ID Type" must be unique in EI! Though the component is not required, it is strongly advised to provide a value.
3.4.1	>>namespace ID	O	Identifier of the Assigning Authority that issued the patient id. If not available the default patient assigning authority for the system will be used
3.4.2	>> universal ID	O	Universal Entity ID. This ID must be unique in EI and must be in the correct format, i.e. containing only digits and dots.
3.4.3	>>universal ID type	O	Universal Entity ID Type. Only supported value : - ISO
3.5	> Identifier Type code (ID)	O	Supported values: PI = Patient Primary Identifier SS = Social Security number If empty PI will be used by default
<b>5</b>	<b>Patient Name</b>	<b>R</b>	
<b>5.1</b>	<b>&gt;family name</b>	<b>R</b>	
<b>5.1.1</b>	<b>&gt;&gt;surname</b>	<b>R</b>	<b>Patient's last name</b>
5.1.2	>>own surname prefix	O	For alphabetic name / maiden name, if PID.5.1.1 is empty, PID.5.1.2 and PID.5.1.3 are concatenated as last name
5.1.3	>>own surname	O	For alphabetic name / maiden name, if PID.5.1.1 is empty, PID.5.1.2 and PID.5.1.3 are concatenated as last name
5.2	>given name	O	Patient's first name
5.3	>second and further given name or initials thereof	O	Patient Middle Name
5.4	>suffix	O	Patient's suffix (eg: JR or III)
5.5	>prefix	O	Patient's prefix (eg: DR)
5.7	>name type code	O	M = maiden name, only support storing last name for maiden name.
5.8	>name representation code	O	Supported are : A or empty = Alphabetic I = Ideographic P = Phonetic
7	Date/Time of Birth	O	
7.1	Date/Time	O	YYYYMMDD[HHMM[SS]] Only the Date is stored
8	Administrative Sex	O	M (Male), F (Female), O (Other)

Seq	HL7 Field Name	Option	Comments
			The DICOM Patient's Sex (0010,0040) attribute can have only the values M, F or O (for other). So when EI receives other values than the one indicated above, it will convert it to O
11	Patient Address	O	
11.1	>street address (SAD)	O	This component specifies the street or mailing address of a person or institution.
11.1.1	>>street or mailing address	O	This is used when PID-11.1.2 (Street Name) is empty
11.1.2	>>street name	O	Name of the street
11.1.3	>>dwelling number	O	House number
11.2	>other designation	O	
11.3	>city	O	Municipality code of the city. If field is empty, the city will be listed as 'Unknown' (code UKW).
11.4	>state or province	O	communication_channel.Address_Line3
11.5	>zip or postal code	O	Zip or postal code of the city
11.6	>country	O	3-digit code defining the country. This is necessary in EI for the address segment to be processed. It defaults to UKW (Unknown country) when left empty
11.7	>address type	O	O (office), H (home address) are supported. L (Legal Address) defaults to office and the rest defaults to H (home address)
11.8	>other geographic designation	O	
11.9	>country/parish code	O	communication_channel.Address_Line4
13	Phone Number – Home	O	This field is used for all home communication. First one is the primary home communication channel
13.1	>telephone number	O	In case of telephone number/fax number/cellular number, the phone number is put in this field [(999)] 999-9999 [X99999][C any text]
13.2	>telecommunication use code	R	<b>If this is empty the communication cannot be registered.</b> <b>Supported are:</b> - Primary Residence Number (PRN) - Work Number (WPN) - Network (email) Address (NET)
13.3	>telecommunication equipment type (ID)	R	<b>If this is empty the communication cannot be registered.</b> <b>Supported are:</b> - Telephone (PH) - Fax (FX) - Cellular Phone(CP) - Internet Address (Internet)
13.4	>Email address	O	In case of Network Email Address, this needs to be filled in
14	Phone Number – Business	O	This field is used for all business communication. First one is the primary business channel
14.1	>telephone number	O	In case of telephone number/fax number/cellular number, the phone number is put in this field [(999)] 999-9999 [X99999][C any text]
14.2	>telecommunication use code	R	<b>Required. If this is empty the communication cannot be registered.</b> <b>Supported are:</b> - Primary Residence Number (PRN) - Work Number (WPN)

Seq	HL7 Field Name	Option	Comments
			- Other Residence Number (ORN) - Network (email) Address (NET)
14.3	>telecommunication equipment type (ID)	R	Required. If this is empty the communication cannot be registered. Supported are: - Telephone (PH) - Fax (FX) - Cellular Phone(CP) - Internet Address (Internet)
14.4	>Email address	O	In case of Network Email Address, this needs to be filled in
19	SSN Number – Patient	O	Social Security number of the patient
23	Birth Place	O	
29	Patient Death Date/Time	C	Condition Predicate:If PID.30 has a value Y, this field also needs a value.
29.1	Date/Time	R	YYYYMMDD[HHMM[SS]][+ZZZ Z] or "" If timezone information is present it overrules the timezone information available in MSH-7
30	Patient Death Indicator	O	Y or N

## 2.2.6 PD1 Segment

*Table 2-16 PD1 Fields Support in Enterprise Imaging*

Seq	HL7 Field Name	Option	Comments
4	Patient Primary Care Provider Name & ID No	O	
4.1	>ID number (ST)	O	Primary Care Physician code
4.2	>family name	O	Physician Family Name
4.2.1	>>surname	O	
4.3	>given name	O	Physician Given Name
4.9	>assigning authority	O	Assigning authority that issued the physicians code
4.9.1	>>namespace ID	O	

## 2.2.7 AL1 Segment

*Table 2-17 AL1 Fields Support in Enterprise Imaging*

Seq	HL7 Field Name	Option	Comments
3	Allergen Code/Mnemonic/Description	R	Allergy code, name
3.1	>identifier	R	
4	Allergy Severity Code	O	Severity code, supported values : SV - Severe, MO - Moderate, MI - Mild, U - Unknown.
4.1	>identifier	O	
5	Allergy Reaction Code	O	Reaction code
6	Identification Date	O	The start date when the allergy was identified. Format : YYYYMMDD

## 2.2.8 MRG Segment

Determines the source patient of a patient merge when the message is an ADT\_A40 event. When it's an ADT\_A47 event, then this segment determines the source patient identifier(s), which is going to be replaced by the target patient identifier(s) present in the PID-segment.

**Table 2-18 MRG Fields Support in Enterprise Imaging**

Seq	HL7 Field Name	Option	Comments
1	Prior Patient Identifier List	R	<p>For ADT_A40: Identifiers specified in the MRG segment identify the source patient that needs to be merged to the target patient identified in the PID-segment</p> <p>For ADT_A47: One to one field corresponding list to PID identifiers list. Identifiers specified in MRG segment are to be replaced by field corresponding identifier in PID</p>
1.1	ID number	R	Identification number or code for the patient
1.4	Assigning authority	O	<p>The 'Assigning authority' aka 'Issuer of Patient ID' consist of the following three parts:</p> <ul style="list-style-type: none"> <li>- Namespace ID // Issuer of Patient ID</li> <li>- Universal ID // Universal Entity ID</li> <li>- Universal ID Type // Universal Entity ID Type</li> </ul> <p>The combination of 'Namespace ID' - 'Universal ID' - 'Universal ID Type' must be unique in EI!</p>
1.4.1	Namespace ID	O	Identifier of the Assigning Authority that issued the patient id. If not available an internal patient assigning authority for the system will be used
1.4.2	>>universal ID	O	<p>Universal Entity ID.</p> <p>This ID must be unique in EI and must be in the correct format, i.e. containing only digits and dots.</p> <p>Its total length &lt;= 64 characters!</p>
1.4.3	>>universal ID type		<p>Universal Entity ID Type.</p> <p>Only supported value :</p> <ul style="list-style-type: none"> <li>- ISO</li> </ul>
1.5	Identifier type code	O	<p>Supported values:</p> <ul style="list-style-type: none"> <li>- PI = Patient Primary Identifier</li> </ul> <p>If empty PI will be used by default</p> <p>Value entered should match PID 3.5 for successful merge</p>

## 2.2.9 PV1 Segment

**Table 2-19 PV1 Fields Support in Enterprise Imaging**

Seq	HL7 Field Name	Option	Comments
1	Set ID – PV1	O	
2	Patient Class	R	<p>Required for: A01, A04, A05, A06, A07</p> <p>I (Inpatient), O (Outpatient), E (Emergency), P (Preadmit), R (Recurring), B (Obstetrics), N (Not applicable)</p>
		O	Optional for A02, A03, A08, A11, A12, A13, A28, A31, A38
3	Assigned Patient Location	R	<p>Required for: A02, A12. At least one of the component must be filled in (3.1 or 3.2 or 3.3 or 3.4 or 3.11)</p>

Seq	HL7 Field Name	Option	Comments
		O	Optional for: A01, A03, A04, A05, A06, A07, A08, A11, A13, A28, A31, A38, ORM, ORU
3.1	>Point of care	O	Patient's location department code. When an unexisting department code is used, the department will be autocreated
3.2	>Room	O	Room in the department
3.3	>Bed	O	Bed in the room
3.4	>Facility (HD)	O	
3.4.1	>>Namespace ID	O	Patient's location facility code. When an unexisting facility code is used, the facility will be autocreated
3.11	>Assigning Authority for location	O	
3.11.1	>>Namespace ID	O	Identifier of the Assigning Authority that issued the department/facility code. If empty, the default assigning authority will be taken.
4	Admission Type	O	Type of admission, the circumstances under which the patient was or will be admitted
7	Attending Doctor	O	Attending Physician information
7.1	>ID number (ST)	O	physician identifier
7.2	>Family name	O	
7.2.1	>>Surname	O	Family name of the physician
7.3	>Given name	O	Given name of the physician
7.4	>Second and further given names or initials thereof	O	Middle name of the physician.
7.5	>Suffix (e.g., JR or III)	O	Name suffix, like SR
7.6	>Prefix (e.g., DR)	O	Name prefix, like Prof
7.9	>Assigning authority	O	Assigning Authority of attending physician code.
7.9.1	>>Namespace ID	O	Identifier of the Assigning Authority that issued the attending physician code. If empty, the default assigning authority will be taken.
8	Referring Doctor	O	Referring physician information (see detailed component description in PV1-7)
8.1	>ID number (ST)	O	
8.2	>Family name	O	
8.2.1	>>Surname	O	
8.3	>Given name	O	
8.4	> Second and further given names or initials thereof	O	
8.5	>Suffix	O	e.g: JR or III
8.6	>Prefix	O	e.g: DR
8.9	>Assigning authority	O	Assigning Authority of referring physician code.
8.9.1	>>Namespace ID	O	Identifier of the Assigning Authority that issued the referring physician code. If empty, the default assigning authority will be taken.
15	Ambulatory Status	O	Code indicating any permanent or transient handicapped conditions
16	VIP Indicator	O	This field identifies the type of VIP. - 1, set VIP or personnel as yes - 0, set as no
17	Admitting Doctor	O	Admitting physician information (see detailed component description in PV1-7)



Seq	HL7 Field Name	Option	Comments
17.1	>ID number (ST)	O	
17.2	>Family name	O	
17.2.1	>>Surname	O	
17.3	>Given name	O	
17.4	> Second and further given names or initials thereof	O	
17.5	>Suffix	O	e.g: JR or III
17.6	>Prefix	O	e.g: DR
17.9	>Assigning authority	O	Assigning Authority of admitting physician code.
17.9.1	>>Namespace ID	O	Identifier of the Assigning Authority that issued the admitting physician code. If empty, the default assigning authority will be taken.
18	Patient Type	O	Site-specific values that identify the patient type
19	Visit Number	R	<b>Required for A01, A02, A03, A04, A05, A06, A07, A11, A12, A13, A38</b>
		O	Optional for: A08, A28, A31, ORM, ORU
19.1	>ID	R	<b>Unique number identifying the admission.</b>
19.4	>Assigning authority	O	
19.4.1	>>Namespace ID	O	Identifier of the Assigning Authority that issued the admission number. If empty, the default assigning authority will be taken.
44	Admit Date/Time	O	When no admit date/time is present in the message, a fallback is done to EVN-2.
44.1	>Date/time	O	YYYYMMDD[HHHMM[SS]][+-ZZZ Z] If timezone information is present it overrules the timezone information if available in MSH-7
45	Discharge Date/Time	O	When no discharge date/time is present in the message, a fallback is done to EVN-2
45.1	>Date/time	O	YYYYMMDD[HHHMM[SS]][+-ZZZ Z] If timezone information is present it overrules the timezone information if available in MSH-7

## 2.2.10 ORC Segment

**Table 2-20 ORC Fields Support in Enterprise Imaging**

Seq	HL7 Field Name	Option	Comments
1	Order Control	R	<b>Supported in ORM:</b> - NW (New order) - SC/XO/XX(Change Order) - CA (Cancel Order) - DC (Discontinue order)  <b>Supported in ORU</b> <b>SC (status changed)</b>
2	Placer Order number	O	
2.1	>Entity identifier	R	<b>Unique identifier for the order, defined by the order placer</b>
2.2	>Namespace ID	O	Identifier of the Assigning Authority that issued the placer order number.

Seq	HL7 Field Name	Option	Comments
			If empty, the default assigning authority will be taken
3	Filler Order number	R	<b>Required for ORM</b>
		O	Optional for ORU
3.1	>Entity identifier	R	<b>Unique identifier for the order, defined by the order filler</b>
3.2	>Namespace ID	O	Identifier of the Assigning Authority that issued the filler order number. If empty, the default assigning authority will be taken
5	Order Status	R	<b>Required for ORU.</b> <b>-IP (In Progress)</b> <b>-CM (Completed)</b> <b>-DC (Discontinued)</b> <b>-CA (Cancelled)</b> <b>-SC (Scheduled)</b>
		O	Optional for ORM - empty, when updating an existing order no status change when creating a new order set SPS status as Scheduled - IP: set SPS status as STARTED - CM: set SPS status as COMPLETED - DC: set SPS status as DISCONTINUED - CA: set SPS status as CANCELLED - SC: set SPS status as Ordered/Scheduled - PA: <ul style="list-style-type: none"> <li>o When an existing procedure, then keep the existing SPS status + Set the Patient Arrived flag to True.</li> <li>o When a new procedure, then set SPS status as Ordered/Scheduled + Set the Patient Arrived flag to True.</li> </ul>
7	Quantity/Timing	O	Identical to the Quantity/Timing in OBR-Segment. If values provided in these fields differ, EI will consider value provided in OBR-27.
7.4	>Start date/time	O	
7.4.1	>>Date/time	O	YYYYMMDD[HHMM[SS]][+ZZZ Z] If timezone information is present it overrules the timezone information if available in MSH-7
7.5	>End date/time	O	
7.5.1	>>Date/time	O	YYYYMMDD[HHMM[SS]][+ZZZ Z] If timezone information is present it overrules the timezone information if available in MSH-7
7.6	>Priority	O	1) For existing sites (migration) - High Priority codes: A, T - Normal Priority codes: S, Unknown* - Low Priority Codes: C, P, R  2) For new installations - STAT priority code: S - Urgent priority codes: A, T - High Priority code: C - Normal Priority codes: P, Unknown* - Routine Priority Code: R

Seq	HL7 Field Name	Option	Comments
			These values are configurable. * A value for which no mapping is applicable.
9	Date/Time of Transaction	O	Fall back for ORC-15 being empty. Used as order creation date/time. When empty, the Date/Time of processing the message is taken
9.1	>Date/time	O	YYYYMMDD[HHHMM[SS]][+-ZZZ Z] If timezone information is present it overrules the timezone information if available in MSH-7
12	Ordering Provider	O	Requesting Physician. If this is not filled in OBR-16 is processed instead
12.1	>ID number	O	Unique code identifying the requesting physician
12.2	>Family name	O	
12.2.1	>>Surname	O	Requesting physician last name
12.3	>Given name	O	Requesting physician first name
12.4	>Second and further given names or initials thereof	O	Middle name
12.5	>Suffix	O	e.g: JR or III
12.6	>Prefix	O	e.g: DR
12.9	>Assigning authority	O	
12.9.1	>>Namespace ID	O	Identifier of the Assigning Authority that issued the physician id
15	Order Effective Date/Time	O	Fall back for OBR-6 being empty. Used as order creation date/time. When empty, ORC.9 is taken.
15.1	>Date/Time	O	YYYYMMDD[HHHMM[SS[.SSSS]]][+-ZZZ Z]
17	Entering Organization	O	Requesting department
<b>17.1</b>	<b>&gt;Identifier</b>	<b>R</b>	<b>Code of the department</b>
17.2	>Text	O	Name of the department. If empty, the identifier will be taken as name.
17.3	>Name of Coding system	O	Identifier of the Assigning Authority that issued the department id. If empty, the default system assigning authority will be taken.
21	Ordering facility name	O	Requesting Hospital
21.1	>Organization name	O	Name of the hospital. If empty, the identifier of the requesting facility will be taken as name
<b>21.3</b>	<b>&gt;ID number (NM)</b>	<b>R</b>	<b>Identifier of the hospital</b>
21.6	>assigning authority	O	
<b>21.6.1</b>	<b>&gt;&gt;Namespace ID</b>	<b>R</b>	<b>Identifier of the Assigning Authority that issued the requesting hospital id.</b> <b>If empty, the default assigning authority will be taken.</b>
25	Order Status Modifier	O	Order Status Modifier only supported when property on ORM inbound pipeline. Delete requested procedure upon ORM cancel with ORC-25 = DELETE, is true.
25.1	>identifier (ST)	O	Only supported value : DELETE. For a procedure where no reading workflow or QC workflow is started and no images are available the following applies depending on the ORM inbound pipeline property setting : - property true : When ORM message is sent to cancel a procedure ( Order control, ORC.1=CA and Order status, ORC.5=CA)and value of this field is "DELETE", the procedure will be deleted.

Seq	HL7 Field Name	Option	Comments
			-property false : A status update to procedure canceled is performed. In all other cases only status update to procedure canceled will be performed.

## 2.2.11 OBR Segment

**Table 2-21 OBR Fields Support in Enterprise Imaging**

Seq	HL7 Field Name	Option	Comments
1	Set ID - Observation Request	O	Set to 1,2,3...
2	Placer Order number	O	Identical to ORC-2
2.1	>Entity identifier	R	
2.2	>Namespace ID	O	
3	Filler Order number	R	Identical to ORC-3
3.1	>Entity identifier	R	
3.2	>Namespace ID	O	
4	Universal Service Identifier	R	
4.1	>Identifier	R	The code of the requested procedure definition
4.2	>Text	O	The description of the requested procedure definition
4.3	>Name of coding system	O	Code of the performing department. Any value in OBR-43 will overrule this value. Department will be created if not existing in EI. The combination of OBR-4.3 and the default assigning authority is used in the lookup
4.4	>Alternate Identifier	O	The code of the scheduled procedure step definition
4.5	>Alternate text	O	The description of the scheduled procedure step definition
4.6	>Name of alternate coding system	O	Coding system providing the code. If empty use issuer of order filler number. If there is no issuer for order filler number, use the default assigning authority.
6	Requested Date/Time	O	Used as order creation date/time. When empty a fall back on ORC-15 and ORC-9 is done. If both empty the Date/Time of processing the message is taken
6.1	>Date/Time	O	YYYYMMDD[HHMM[SS]][+-ZZZ Z] If timezone information is present it overrules the timezone information if available in MSH-7
13	Relevant Clinical Info	O	
16	Ordering Provider	O	The ordering provider should be the same as ORC-12. If ORC-12 is filled in, this will not be processed.
16.1	>ID number (ID)	O	Physician Identifier
16.2	>Family name	O	
16.2.1	>>Surname	O	Physician lastname
16.3	>Given name	O	Physician firstname
16.9	>Assigning authority	O	

Seq	HL7 Field Name	Option	Comments
16.9.1	>>Namespace ID	O	Identifier of the Assigning Authority that issued the physician code. If empty, the default assigning authority will be taken.
18	Placer Field 1	R	<b>Required for ORU. Accession number of the requested procedure.</b>
		C	Conditional for ORM. Accession number of the requested procedure Condition Predicate: Field is required when the message has no StudyInstanceUid specified. (no ZDS segment)
19	Placer Field 2	R	<b>Required for ORM. Requested procedure ID</b>
		O	Optional for ORU. Requested procedure ID
20	Filler Field 1	R	<b>Required for ORM. Scheduled Procedure Step ID</b>
		O	Optional for ORU. Scheduled Procedure Step ID
22	Result Rpt/Status Change - Date/Time	O	ORU only. Report Modification Date or Addendum Modification Date
22.1	>Date/Time	R	<b>The modification date of the report</b>
24	Diagnostic Serv Sect Id	O	Modality Type (Used when Procedure Definition is created on the fly)
25	Result Status	O	ORU only. Report status  Default values (property 'ORU controlled reporting workflow' is set to false) - P preliminary report - F signed report - C correction to results - D report deleted - X no results available, order cancelled  ORU controlled reporting workflow (property 'ORU controlled reporting workflow' is set to true) - UP unapproved preliminary report - P approved preliminary report - F signed report - W wet read report - PA preliminary addendum - A final addendum - D report deleted
27	Quantity/Timing	O	Identical to the Quantity/Timing in ORC-7 Segment. If values provided in these fields differ, EI will consider value provided in OBR- 27
27.4	>Start date/time	O	Used as the requested procedure scheduled study date/ time and the scheduled procedure step start date/ time. When the value is empty, the value of OBR-36 is taken as requested procedure study date/time and scheduled procedure step start date/time
27.4.1	>>Date/Time	O	YYYYMMDD[HHMM[SS]][+ -ZZZ Z] If timezone information is present it overrules the timezone information if available in MSH-7
27.5	>End date/time	O	Used as scheduled procedure step end date/time

Seq	HL7 Field Name	Option	Comments
27.5.1	>>Date/Time	O	YYYYMMDD[HHMM[SS]][+-ZZZ Z] If timezone information is present it overrules the timezone information if available in MSH-7
27.6	>Priority	O	1) For existing sites (migration) - High Priority codes: A, T - Normal Priority codes: S, Unknown* - Low Priority Codes: C, P, R  2) For new installations - STAT priority code: S - Urgent priority codes: A, T - High Priority code: C - Normal Priority codes: P, Unknown* - Routine Priority Code: R  These values are configurable. * A value for which no mapping is applicable
28	Result Copies To	O	
28.1	>ID number (ST)	R	code of the physician
28.2	>family name	O	
28.2.1	>>surname	O	family name of the physician
28.3	>given name	O	given name of the physician (optional)
28.9	>assigning authority	O	
28.9.1	>>namespace ID	O	Identifier of the Assigning Authority that issued the physician code. If empty, the default assigning authority will be taken.
31	Reason For Study	O	
31.2	>Text	O	Text containing clinical information for the study.
32	Principal Result Interpreter	O	For ORU 1) Default behaviour (property 'ORU controlled reporting workflow' is set to false). Final report : Physician that validated the report - if first iteration of OBR-33 is empty, this physician will also handled as author of the report. Preliminary report: - physician will be handled as author of the report.  2)ORU controlled reporting workflow (property 'ORU controlled reporting workflow' is set to true) - The physician specified in this field, will always be the attending physician which will get assigned to the signoff task in case the report is preliminary. - If no entry is specified, the first physician in OBR-33 will get the task assigned.  For ORM

Seq	HL7 Field Name	Option	Comments
			This field can be used to fill in the "reading physician". Some RIS want to indicate already to which radiologist the reading task should be assigned to. So when they send the ORM, EI stores the reading physician and when a reading task is created, it can be assigned to the already filled in "reading physician".
<b>32.1</b>	<b>&gt;Name</b>	<b>R</b>	
<b>32.1.1</b>	<b>&gt;&gt;ID number (ST)</b>	<b>R</b>	<b>Code identifying the physician</b>
<b>32.1.2</b>	<b>&gt;&gt;Family name</b>	<b>R</b>	<b>Last name of the Physician</b>
<b>32.1.3</b>	<b>&gt;&gt;Given name</b>	<b>R</b>	<b>Firstname of the physician</b>
32.1.9	>>Assigning authority	O	Identifier of the Assigning Authority that issued the principle result interpreter code. If empty, the default assigning authority will be taken.
33	Assistant Result Interpreter	O	<p>ORU only.</p> <p>1) Default behaviour (property 'ORU controlled reporting workflow' is set to false)</p> <p>Final report:</p> <ul style="list-style-type: none"> <li>- the first iteration will be handled as the author of the report.</li> <li>- All the other iterations are reviewers.</li> </ul> <p>Preliminary report:</p> <ul style="list-style-type: none"> <li>- all the iterations are handled as reviewers.</li> </ul> <p>2)ORU controlled reporting workflow (property 'ORU controlled reporting workflow' is set to true)</p> <ul style="list-style-type: none"> <li>-The physician specified in the first field, will always be the author of the report. If no entry is specified in OBR-32, the first physician in OBR-33 is assigned the created signoff task in case the report is preliminary.</li> <li>- All the other iterations are reviewers.</li> </ul>
<b>33.1</b>	<b>&gt;Name</b>	<b>R</b>	
<b>33.1.1</b>	<b>&gt;&gt;ID number (ST)</b>	<b>R</b>	<b>Code identifying the physician</b>
<b>33.1.2</b>	<b>&gt;&gt;Family name</b>	<b>R</b>	<b>Last name of the Physician</b>
<b>33.1.3</b>	<b>&gt;&gt;Given name</b>	<b>R</b>	<b>Firstname of the physician</b>
33.1.9	>>Assigning authority	O	Identifier of the Assigning Authority that issued the assistant result interpreter code. If empty, the default assigning authority will be taken.
34	Technician	O	
34.1	>Name	O	
34.1.1	>>ID number (ST)	R	<b>Code of the technician</b>
34.1.2	>> family name	O	Family name of the technician.
34.1.3	>>given name	O	First name of the technician
34.1.9	>>Assigning authority	O	Identifier of the Assigning Authority that issued the technician code. If empty, the default assigning authority will be taken.
34.5	>Room	O	Acquisition room
35	Transcriptionist	O	ORU only. Transcriptionist of the report
<b>35.1</b>	<b>&gt;name</b>	<b>R</b>	
<b>35.1.1</b>	<b>&gt;&gt;ID number (ST)</b>	<b>R</b>	<b>Code identifying the transcriptionist.</b>
<b>35.1.2</b>	<b>&gt;&gt;Family name</b>	<b>R</b>	<b>Family name of the transcriptionist.</b>

Seq	HL7 Field Name	Option	Comments
35.1.3	>>Given name	O	Given name of the transcriptionist.
35.1.9	>>Assigning authority	O	Identifier of the Assigning Authority that issued the transcriptionist code. If empty, the default assigning authority will be taken.
36	Scheduled Date/Time	O	Used for requested procedure scheduled study date/time and scheduled procedure step date/time when OBR-27.4 is empty.
36.1	>Date/Time	O	YYYYMMDD[HHMM[SS]][+ZZZ Z] If timezone information is present it overrules the timezone information if available in MSH-7
43	Planned Patient Transport Comment	O	This field can be used to provide performing department. When the performing department code provided in OBR-43.1 (with or without an issuer in OBR-43.3) is not found in the system, it will be created. When no name is present, the department code is used as the department name. Another way to provide the performing department is by means of OBR-4.3 (universal service identifier) . If present the value will be used in combination with the default assigning authority in the lookup.
43.1	>Identifier	O	performing department code
43.2	>Text	O	performing department name
43.3	>Name of coding system	O	issuer of performing department code

## 2.2.12 OBX Segment

**Table 2-22 OBX Fields Support in Enterprise Imaging**

Seq	HL7 Field Name	Option	Comments
1	Set ID – OBX	O	Set to 1,2,3...
2	Value Type	R	<p><b>Supported values For ADT:</b></p> <ul style="list-style-type: none"> <li>- TX (For Text Attachment type or Pregnancy information)</li> <li>- FT (For Text Attachment type)</li> </ul> <p><b>Supported values for ORM:</b></p> <ul style="list-style-type: none"> <li>- RP (For PDF/Image Attachment type)</li> <li>- TX (For Text Attachment type or Pregnancy information)</li> <li>- ED (For Attachment type PDF/Image - Base64 encoded encapsulated data)</li> </ul> <p><b>Supported values for ORU:</b></p> <p>In case of reports:</p> <ul style="list-style-type: none"> <li>- FT</li> <li>- RP (For reports send as PDF with reference pointer)</li> <li>- ED (For reports send as Base64 encoded data)</li> </ul> <p>In case of attachments:</p> <ul style="list-style-type: none"> <li>- TX (For Text attachment type)</li> <li>- RP (For PDF/Image attachment type)</li> </ul> <p>In case of pregnancy information :</p> <ul style="list-style-type: none"> <li>- TX</li> </ul>



Seq	HL7 Field Name	Option	Comments
3	Observation identifier	R	
3.1	>Identifier	R	<p><b>For ADT:</b> Must correspond with an attachment code that is setup as a "Patient level" attachment as defined in the Administrator desktop. For Pregnancy : PREGNANT to indicate Pregnancy information</p> <p><b>For ORM:</b> 1) OBX segment with attachment information: Value of this field must correspond with an attachment code as defined in the Administrator desktop. The following attachment codes are already delivered in dedicated HL7 fields and therefore should not be processed as attachments in OBX-segments: - REASONFORSTUDY (OBR(31)) - PROCEDURECOMMENT (NTE(3)) - CLINICALINFO (OBR(13) or NTE(3))</p> <p><b>Remark :</b> If the attachment code is not setup as a 'Study level' or 'Order level' or 'Report level' attachment in the Administrator Desktop, the attachment will not be processed. The reception of multiple attachments with the same attachment code OBX(3.1) is supported: - PDF/Image type attachments: when multiple attachments with the same attachment code are present, a new document/image is added per attachment. - Single value Text attachment: when multiple attachments with the same attachment code are present, the text from all these attachments is concatenated. A newline character is used to separate the texts from the different attachments. - Message board Text attachment: when multiple attachments with the same attachment code are present, a separate text item is added to the message board per attachment.</p> <p>2)OBX segment with Patient Arrived information: Value of OBX-3 is PATIENT_ARRIVED. See following details: The combination (OBX-2 = TX, OBX-3 = PATIENT_ARRIVED, OBX-5 = YES) indicates the patient's arrived status is true; whereas the combination (OBX-2 = TX, OBX-3 = PATIENT_ARRIVED, OBX-5 = NO) indicates the patient's arrived status is false.</p> <p>3) OBX segment holding Pregnancy information - PREGNANT indicate pregnancy information</p>

Seq	HL7 Field Name	Option	Comments
			<p><b>For ORU:</b></p> <p><b>In case of reports:</b></p> <ul style="list-style-type: none"> <li>- This field should contain &amp;GDT for specifying general report.</li> <li>- It should contain &amp;ADT for specifying addendum report.</li> <li>- When report type is FT (OBX-2) only 1 OBX segment with general report is allowed and multiple addendum OBX are allowed.</li> <li>- When report type is RP or ED (OBX-2) only 1 OBX segment is allowed. The OBX should specify a general report. When addendum is present, it is supposed to be included in the general report.</li> </ul> <p><b>In case of attachments:</b></p> <ul style="list-style-type: none"> <li>-Must correspond with an attachment code as defined in the administrator desktop.</li> <li>-The following attachment codes are already delivered in dedicated HL7 fields and therefore should not be processed as attachments in OBX-segments:                             <ul style="list-style-type: none"> <li>- REASONFORSTUDY (OBR(31))</li> <li>- PROCEDURECOMMENT (NTE(3))</li> <li>- CLINICALINFO (OBR(13) or NTE(3))</li> </ul> </li> </ul> <p><b>Remark:</b></p> <p>If the attachment code is not setup as a 'Study level', 'Order level' or 'Report level' attachment in the Administrator Desktop, the attachment will not be processed. The reception of multiple attachments with the same attachment code (OBX-3.1) is supported:</p> <ul style="list-style-type: none"> <li>- PDF/Image type attachments: when multiple attachments with the same attachment code are present, a new document/image is added per attachment.</li> <li>- Single value Text attachment: when multiple attachments with the same attachment code are present, the text from all these attachments is concatenated. A newline character is used to separate the texts from the different attachments.</li> <li>- Message board Text attachment: when multiple attachments with the same attachment code are present, a separate text item is added to the message board per attachment.</li> </ul> <p><b>For Pregnancy : PREGNANT to indicate Pregnancy information</b></p>
3.2	>Text	O	Description of the observation

Seq	HL7 Field Name	Option	Comments
3.4	>Alternate identifier	O	<p>For ORU only:</p> <p>Used in case of attachments : Should contain the accession number where the attachment belongs to (in case of study level attachment) or empty (In case of order/report level attachment).</p> <p>When no accession number is passed in this component, the requested procedure/study level attachments are to be added to all requested procedures this report is linked to.</p> <p>Example:</p> <p>An ORU message is sent containing two requested procedures belonging to the same service request and it contains two OBX-segments with attachments:</p> <ul style="list-style-type: none"> <li>- one attachment with code STDAT20 is a 'Study-level' attachment for requested procedure with accession number ACC123.</li> <li>- the other attachment with code SCANNEDREQUEST is an 'Order-level' attachment which applies to both requested procedures.</li> </ul> <p>OBX 1 RP STDAT20^^ACC123  \attachmentserver\foldername \attachmentfilename1.pdf    OBX 2 RP SCANNEDREQUEST  \attachmentserver\foldername \attachmentfilename2.pdf  </p> <p>Attachment with code STDAT20 is linked to requested procedure with accession number ACC123 and attachment with code SCANNEDREQUEST is linked to the service request.</p>
4	Observation Sub-ID	O	sub identifier for this observation
5	Observation Value	R	<p><b>For ADT:</b></p> <ul style="list-style-type: none"> <li>- Plain text (in case of Text attachment type)</li> <li>- YES/NO/UNKNOWN (in case of Pregnancy information)</li> </ul> <p><b>For ORM:</b></p> <p><b>If OBX-2 = RP</b></p> <ul style="list-style-type: none"> <li>- OBX-5 = path + file name for PDF/Image type attachment with Reference pointer</li> </ul> <p><b>If OBX-2 = TX + OBX-3 = PATIENT_ARRIVED</b></p> <ul style="list-style-type: none"> <li>- OBX-5 = YES/NO</li> </ul> <p><b>If OBX-2 = TX + OBX-3 = PREGNANT</b></p> <ul style="list-style-type: none"> <li>- OBX-5 = YES/NO/UNKNOWN</li> </ul> <p><b>If OBX-2 = ED</b></p> <ul style="list-style-type: none"> <li>- OBX-5.1 = AGILITY</li> <li>- OBX-5.2 = TEXT (in case of PDF) / IM (in case of Image)</li> <li>- OBX-5.3 = PDF (in case of PDF) / JPEG, TIFF, GIF, PNG, BMP, JPG, TIF (in case of image: the type of the image)</li> <li>- OBX-5.4 = Base64</li> <li>- OBX-5.5 = The encoded attachment as one PDF in Base64 encoding</li> </ul>

Seq	HL7 Field Name	Option	Comments
			<p>For ORU:</p> <p>In case of reports:</p> <ul style="list-style-type: none"> <li>- plain text(OBX-2=FT)</li> <li>- path + file name (in case of report send as PDF with reference pointer OBX-2=RP)</li> <li>- ^TEXT^PDF^BASE64^base64-encoded string of binary PDF data (OBX-2=ED)</li> </ul> <p>In case of attachments:</p> <ul style="list-style-type: none"> <li>- plain text (in case of Text attachment type)</li> <li>- path + file name (in case of PDF/Image attachment type)</li> <li>- "" (delete all attachments with code present in OBX-3.1)</li> </ul> <p>If empty, this OBX segment will not be processed.</p> <p>In case of pregnancy information :</p> <ul style="list-style-type: none"> <li>- YES/NO/UNKNOWN</li> </ul>
11	Observation Result Status	R	<p>for ORU only:</p> <p>Only required when ORU controlled reporting workflow is set to False (disabled). If value is not present a fall back is done to OBR-25. supported values:</p> <ul style="list-style-type: none"> <li>- F (final results)</li> <li>- P (preliminary)</li> <li>- C (correction)</li> <li>- X (report body no longer found - overwrite) Report gets final status.</li> <li>- D (deleted)</li> </ul>
14	Date/Time of the Observation	O	<ul style="list-style-type: none"> <li>- The creation date of the report</li> <li>- YYYYMMDD (date of last menstrual period - in case of pregnancy information)</li> </ul>
14.1	>Date/Time	R	<p>YYYYMMDD[HHMM[SS]][+ZZZ Z]</p> <p>If timezone information is present it overrules the timezone information if available in MSH-7. If no date time is filled in, the report will be saved with the datetime of saving/updating.</p>

## 2.2.13 ZDS Segment

Table 2-23 ZDS Fields Support in Enterprise Imaging

Seq	HL7 Field Name	Option	Comments
1	Study Instance UID	R	ORM and ORU only
1.1	Reference pointer	R	DICOM Study Instance UID (0020,000D)

## 2.2.14 NTE Segment

*Table 2-24 NTE Fields Support in Enterprise Imaging*

Seq	HL7 Field Name	Option	Comments
1	Set ID - NTE		Identifier
		R	Required for ORU
		O	Optional for ORM
3	Comment	O	Plain text of Clinical info (Comment type CI) or Procedure comment(Comment type PC)
4	Comment type	R	Required for ORU
		O	Optional for ORM
4.1	>Identifier	R	Supported values: - CI (clinical info) - PC (procedure comment)

## 2.2.15 MFI Segment

*Table 2-25 MFI Fields Support in Enterprise Imaging*

Seq	HL7 Field Name	Option	Comments
1	Master File Identifier	R	
1.1	>Identifier	R	Only values PRA/STF supported
1.2	>Text	O	
3	File-Level Event Code	R	Only value UPD supported
6	Response Level Code	R	Only value NE supported

## 2.2.16 MFE Segment

*Table 2-26 MFE Fields Support in Enterprise Imaging*

Seq	HL7 Field Name	Option	Comments
1	Record-Level Event Code	R	
2	MFN Control ID	C	
3	Effective Date/Time	O	
3.1	>Date/Time	O	The date of the validate of the MFN Event
4	Primary Key Value - MFE	R	Same value as STF-1-1
5	Primary Key Value Type	R	Only value CE supported

## 2.2.17 STF Segment

*Table 2-27 STF Fields Support in Enterprise Imaging*

Seq	HL7 Field Name	Option	Comments
1	Primary Key Value - STF	R	
1.1	>Identifier	R	The same value as MFE-4 / Primary doctor id
1.3	>Name of coding system	O	The Primary doctor id issuer (AA)

Seq	HL7 Field Name	Option	Comments
2	Staff ID Code	O	
2.1	>ID	O	The other doctor id
2.4	>Assigning Authority	O	
2.4.1	>>Namespace ID	O	The other doctor id issuer (AA)
3	Staff name	O	
3.1	>Family name	O	
3.1.1	>>Surname	O	
3.2	>Given name		
3.3	>Second and further given names or initials thereof	O	Middle name
3.4	>Suffix (e.g., JR or III)	O	
3.5	>Prefix (e.g., DR)	O	
4	Staff Type	O	Used as professional-department relationship type. Default = EMPLOYEE
5	Administrative Sex	O	
8	Department	O	Will be added as department relationship to the professional
8.1	>Identifier	O	Department id
8.3	>Name of coding system	O	The department id issuer (AA)
9	Hospital service	O	Will be used as the facility in the department-relationship
9.1	>Identifier	O	Facility id
9.3	>Name of coding system	O	
10	Phone	O	
10.1	>Telephone number	O	In case of telephone number/fax number/cellular number, the phone number is put in this field. [(999)] 999-9999 [X99999][C any text]
10.2	>Telecommunication use code	R	<b>If this is empty the communication cannot be registered.</b> <b>Supported are:</b> - <b>Primary Residence Number (PRN)</b> - <b>Work Number (WPN)</b> - <b>Network (email) Address (NET)</b>
10.3	>Telecommunication equipment type (ID)	R	<b>If this is empty the communication cannot be registered.</b> <b>Supported are:</b> - <b>Telephone (PH)</b> - <b>Fax (FX)</b> - <b>Cellular Phone(CP)</b> - <b>Internet Address (Internet)</b>
10.4	>Email address	O	In case of Network Email Address, this needs to be filled in
11	Office/Home Address	O	
11.1	>Street Address (SAD)	O	This component specifies the street or mailing address of a person or institution.
11.1.1	>>Street or mailing address	O	communication_channel.Street (if STF-11.1.2 is empty)
11.1.2	>>Street name	O	communication_channel.Street
11.1.3	>>Dwelling number	O	communication_channel.Streetnr
11.2	>other designation	O	communication_channel.Address_Line1
11.3	>City	O	municipality.name (if empty -> Unknown is used)
11.4	>State or province	O	communication_channel.Address_Line3
11.5	>Zip or postal code	O	municipality.zip (if empty -> code UKW is used - mapping to name = Unknown)

Seq	HL7 Field Name	Option	Comments
11.6	>Country	O	communication_channel.country (3 digit code - if empty: code UKW - mapping to name = Unknown)
11.7	>Address type	O	communciation_channel.communication_type (H - for home address / O - for office address / L maps to O / other codes map to H)
11.8	>other geographic designation	O	communication_channel.Address_Line2
11.9	>county/parish code	O	communication_channel.Address_Line4
12	Institution Activation Date	O	Begin date of the professional/department relationship
12.1	>Date	O	
<b>12.1.1</b>	<b>&gt;&gt;Date/Time</b>	<b>R</b>	
13	Institution Inactivation date	O	End date of the professional/department relationship
13.1	>Date	O	
<b>13.1.1</b>	<b>&gt;&gt;Date/Time</b>	<b>R</b>	

## 2.2.18 PRA Segment

*Table 2-28 PRA Fields Support in Enterprise Imaging*

Seq	HL7 Field Name	Option	Comments
3	Practitioner Category	O	Corresponds with professional's position in EI. Possible values: - RAD / PHYSICIAN - PHYS / PHYSICIAN - TECH / TECHNOLOGIST - TRANS / TRANSCRIPTIONIST - SYSADM / SYSTEM ADMINISTRATOR - MAN / MANAGER - AGFA / AGFA SERVICE - REC / RECEPTIONIST If another value is present in this field this will not be added as a position of the professional
7	Privileges	O	
7.1	>Privilege	O	
<b>7.1.1</b>	<b>&gt;&gt;Identifier</b>	<b>R</b>	<b>Supported value(s): INTERNAL, INT</b>
9	Institution	O	Will be used if STF-9 is empty
9.1	>Identifier	O	
9.2	>Text	O	
9.3	>Name of coding system	O	

## 3 OUTBOUND MESSAGES

### 3.1 Core Server supported Trigger Events

#### 3.1.1 Supported ACK Events

For more information about the processed segments and supported fields, refer to the document specified in the 'Conformance profile document' column.

**Table 3-1 Supported ACK Events**

Func Area	Event Code	ACK Trigger Event Description	Conformance profile document
ACK	ALL	General Acknowledgement message	ENTERPRISE IMAGING OUTBOUND ACK.pdf

#### 3.1.2 Supported DFT Events

For more information about the processed segments and supported fields, refer to the document specified in the 'Conformance profile document' column.

**Table 3-2 Supported DFT Events**

Func Area	Event Code	DFT Trigger Event Description	Conformance profile document
DFT	P03	Order message	ENTERPRISE IMAGING OUTBOUND DFT_P03.pdf

#### 3.1.3 Supported ORM Events

For more information about the processed segments and supported fields, refer to the document specified in the 'Conformance profile document' column.

**Table 3-3 Supported ORM Events**

Func Area	Event Code	ORM Trigger Event Description	Conformance profile document
ORM	O01	Order message	ENTERPRISE IMAGING OUTBOUND ORM_O01.pdf

#### 3.1.4 Supported ORU Events

For more information about the processed segments and supported fields, refer to the document specified in the 'Conformance profile document' column.

**Table 3-4 Supported ORU Events**



Func Area	Event Code	ORU Trigger Event Description	Conformance profile document
ORU	R01	Unsolicited transmission of an observation Message	ENTERPRISE IMAGING OUTBOUND ORU_R01 - ASCII ENCAPSULATED DATA-ORU_R01.pdf
ORU	R01	Unsolicited transmission of an observation Message	ENTERPRISE IMAGING OUTBOUND ORU_R01 - ASCII LINE BY LINE.pdf
ORU	R01	Unsolicited transmission of an observation Message	ENTERPRISE IMAGING OUTBOUND ORU_R01 - ASCII NARRATIVE REPORTING.pdf
ORU	R01	Unsolicited transmission of an observation Message	ENTERPRISE IMAGING OUTBOUND ORU_R01 - Base64 ENCAPSULATED DATA PDF.pdf
ORU	R01	Unsolicited transmission of an observation Message	ENTERPRISE IMAGING OUTBOUND ORU_R01 - FORMATTED TEXT.pdf
ORU	R01	Unsolicited transmission of an observation Message	ENTERPRISE IMAGING OUTBOUND ORU_R01 - REFERENCE POINTER.pdf

## 4 QUERIES

### 4.1 Web Server Query

Web Server translates user based patient-level queries into HL7 PDQ messages.

Web Server uses the HL7 PIX services to combine patients together as well as to allow querying for longitudinal records against external non-HL7 systems which require specific patient identities.

Web Server can use v2.5 and will accept the related responses, as well as v3 SOAP based queries.

#### 4.1.1 Supported Query events

For more information about the processed segments and segment mappings, refer to the document specified in the Conformance profile document column.

**Table 4-1 Supported Query Events**

Func Area	Event Code	QRY Trigger Event	Conformance profile document <sup>1</sup>
QBP	Q23	Query for associated patient identities	EnterprisImaging_WS_PIX_Outbound_QPB_Q23.pdf EnterprisImaging_WS_PIX_Inbound_RSP_K23.pdf
QBP	Q22	Query for patient demographics	EnterprisImaging_WS_PDQ_Outbound_QBP_Q22.pdf EnterprisImaging_WS_PDQ_Inbound_RSP_K22.pdf EnterprisImaging_WS_PDQ_Inbound_RSP_K21.pdf

<sup>1</sup> These documents are included in the zip file "[001664 Enterprise Imaging 8.3.x HL7 Conformance Profiles](#)", also available on <http://global.agfahealthcare.com/hl7>.



### Details as of PDF Creation Date

#### Document Metadata

<b>Title:</b>	001663_Enterprise Imaging 8.3.x HL7 Conformance Statement
<b>Livelihood ID:</b>	81707306
<b>Version#:</b>	5
<b>Version Date:</b>	2023-12-14 12:07 PM CET
<b>Status:</b>	Approved on 2023-12-14 03:48 PM CET
<b>Owner:</b>	Bruno Laffin (nawyv)
<b>Created By:</b>	Bruno Laffin (nawyv)
<b>Created Date:</b>	2023-05-30 09:30 AM CET
<b>PDF Creation Date:</b>	2023-12-14 03:49 PM CET

#### This document was approved by:

#### Signatures:

1. Bruno Laffin (nawyv) on 2023-12-14 03:47 PM CET
2. Bart De Smedt (awxma) on 2023-12-14 03:48 PM CET

#### Detailed Approver History:

- **Approval Workflow started on 2023-12-14 02:54 PM CET**
  - Approval task originally assigned to and completed by Bruno Laffin (nawyv) on 2023-12-14 03:47 PM CET
  - Approval task originally assigned to and completed by Bart De Smedt (awxma) on 2023-12-14 03:48 PM CET

#### Version & Status History

Version#	Date Created	Status
5	2023-12-14 12:07 PM CET	Approved - 2023-12-14 Reviewed - 2023-12-14
4	2023-09-21 05:34 PM CET	Published - 2023-10-04 Approved - 2023-09-21
3	2023-09-18 12:36 PM CET	Reviewed - 2023-09-21

2	2023-05-30 05:57 PM CET	Unpublished - 2023-10-04 Published - 2023-06-12 Approved - 2023-06-12 Reviewed - 2023-06-07
1	2023-05-30 09:30 AM CET	

**Applied Categories and Attributes:**

<b>ITCo Library</b>	
Document Type:	Conformance Statement
Category > SubCategory > Item:	Connectivity > Connectivity & Application > To Agfa
Content Manager:	Bruno Laffin (nawyv)
Summary:	The document 001663 is the HL7 Conformance Statement for Enterprise Imaging 8.3.1.000. This is the second revision for EI version 8.3.x. This document will be updated (content and title) for all subsequent releases of EI 8.3.x until EI 8.3 is release for Sale.
Language:	English
Availability:	Internet
Language Master Document:	
Software Download URL:	