


CASRN: 152876-71-6

1. Name:

Cyan dye

2. Chemical Structure



3. Features:

Formula	C ₂₈ H ₃₅ N ₇ O ₅
Absorption max	638 nm (methanol)
E max	87853 mol/l.cm
Melting point	127 °C
Appearance	green crystals

4. Specifications:

	Specification
Assay	≥ 97 %
Solubility g/100ml at 25 °C:	
water	none
MEK	> 20
CH ₂ Cl ₂	> 20
EtOH	< 1.0

5. Registration:

No registration (only temporarily for Germany).

R43, R99, S24, S37

