

AGFA HEALTHCARE IHE Integration Statement



WEB 1000

Status: Released

Document no : 000948

Revision 1.1

NotelD Livelink : 9765766

When printed, this is NOT a controlled copy

Document Information

Service-related contact information worldwide	All service-related contact information is available on this URL→	http://www.agfa.com/en/he/support/support_service/index.jsp

Issued by:
Agfa HealthCare
<SIV Connectivity >
<Septestraat 27>
<B-2640 Mortsel>

tel: <+32.3.444.7338>
email: <connectivity@agfa.com>

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

Copyright © December, 05
Agfa HealthCare
All rights reserved

Agfa HealthCare

REVISION RECORD

Revision Number	Date	Reason for Change
1.0	22 Apr 2004	Initial version
1.1	16 Dec 2005	Update according new template

Agfa HealthCare

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	Web 1000	4.1	22 Apr 2004
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Image Display	none	
Access to Radiology Information (ARI)	Image Display	none	
Internet address for vendor's IHE information: http://www.agfa.com/healthcare/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7			
DICOM	http://www.agfa.com/healthcare/dicom		
Links to general information on IHE			
In North America: www.rsna.org/IHE	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	