

# AGFA HEALTHCARE IHE Integration Statement



## IMPAX 5.2

Status: Released

Document no : 001011

Revision 1.1

NotelD Livelink : 9765624

**When printed, this is NOT a controlled copy**

## Document Information

---

<b>Service-related contact information worldwide</b>	All service-related contact information is available on this URL→	<a href="http://www.agfa.com/en/he/support/support_service/index.jsp">http://www.agfa.com/en/he/support/support_service/index.jsp</a>

Issued by:  
Agfa HealthCare  
<SIV Connectivity >  
<Septestraat 27>  
<B-2640 Mortsel>

tel: <+32.3.444.7338>  
email: <[connectivity@agfa.com](mailto:connectivity@agfa.com)>

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

**Copyright © December, 05**  
**Agfa HealthCare**  
**All rights reserved**

Agfa HealthCare

---

## REVISION RECORD

Revision Number	Date	Reason for Change
1.0	06 Jan 2005	Initial version
1.1	16 Dec 2005	Update according new template

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
<b>AGFA HealthCare</b>	<b>IMPAX</b>	<b>5.2</b>	<b>06 Jan 2005</b>
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
<b>Portable Data for Imaging</b>	<b>Portable Media Importer</b>	<b>none</b>	
	<b>Image Display</b>	<b>none</b>	
<b>Internet address for vendor's IHE information:</b> <a href="http://www.agfa.com/healthcare/ihe">http://www.agfa.com/healthcare/ihe</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>			
<b>DICOM</b>	<a href="http://www.agfa.com/healthcare/dicom">http://www.agfa.com/healthcare/dicom</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.rsna.org/IHE">www.rsna.org/IHE</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>	