

AGFA HEALTHCARE IHE Integration Statement



IMPAX 6.0

Status: Released

Document no : 001035

Revision 1.2

NoteID Livelink : 9765943

When printed, this is NOT a controlled copy

Document Information

Service-related contact information worldwide	All service-related contact information is available on this URL→	http://www.agfa.com/en/he/support/support_service/index.jsp

Issued by:
Agfa HealthCare
<SIV Connectivity >
<Septestraat 27>
<B-2640 Mortsel>

tel: <+32.3.444.7338>
email: <connectivity@agfa.com>

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

Copyright © May, 06
Agfa HealthCare
All rights reserved

Agfa HealthCare

REVISION RECORD

Revision Number	Date	Reason for Change
1.0	22 Aug 2005	Initial version
1.1	16 Dec 2005	Update according new template
1.2	24 April 2006	Update on implemented integration profiles

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	IMPAX	6.0	24 Apr 2006
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Image Manager / Archive	None	
	Image Display	None	
	Evidence Creator	None	
	PPS Manager	None	
Patient information Reconciliation	Image Manager / Archive	None	
	PPS Manager	None	
Portable data for Imaging	Image Display	None	
	Portable Media importer	None	
Access to Radiology Information	Image Manager / Archive	None	
	Image Display	None	
Consistent Time	Time Client	None	
Internet address for vendor's IHE information: http://www.agfa.com/healthcare/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	http://www.agfa.com/healthcare/hl7		
DICOM	http://www.agfa.com/healthcare/dicom		
Links to general information on IHE			
In North America: www.rsna.org/IHE	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	