


# AGFA HEALTHCARE IHE Integration Statement



## Qdoc 5.4.1

Status: Released

Document number 001061

Revision 1.0

NotelD Livelink : 9801874

**When printed, this is NOT a controlled copy**

Agfa HealthCare

## Document Information

---

<b>Service-related contact information worldwide</b>	All service-related contact information is available on this URL→	<a href="http://www.agfa.com/en/he/support/support_service/index.jsp">http://www.agfa.com/en/he/support/support_service/index.jsp</a>

Issued by:  
Agfa HealthCare  
SIV Connectivity  
Septestraat 27  
B-2640 Mortsel

tel: <+32.3.444.7338>  
email: <[connectivity@agfa.com](mailto:connectivity@agfa.com)>

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

**Copyright © December, 05**  
**Agfa HealthCare**  
**All rights reserved**

## REVISION RECORD

Revision Number	Date	Reason for Change
1.0	15 dec 2005	Initial revision

Agfa HealthCare

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
<b>AGFA HealthCare</b>	<b>Qdoc</b>	<b>5.4.1</b>	<b>Dec 2005</b>
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
Scheduled Workflow	Department System Scheduler/ Order Filler	Departmental appointment Notification	
	Performed Procedure Step Manager		
Patient Information Reconciliation	Department System Scheduler/ Order Filler		
	Performed Procedure Step Manager		
Presentation of grouped Procedures	Department System Scheduler/ Order Filler		
Basic Security	Secure Node		
	Department System Scheduler/ Order Filler		
Charge Posting	Department System Scheduler/ Order Filler		
Patient Demographics Query	Patient Demographics consumer		
Consistent Time	Time client		
<b>Internet address for vendor's IHE information:</b> <a href="http://www.agfa.com/healthcare/ihe">http://www.agfa.com/healthcare/ihe</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>			
<b>DICOM</b>	<a href="http://www.agfa.com/healthcare/dicom">http://www.agfa.com/healthcare/dicom</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.rsna.org/IHE">www.rsna.org/IHE</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>	