

AGFA HEALTHCARE IHE Integration Statement

→ **Drystar 5503 & 5503 Mammo**

Status: Released
Document No. 001142
Revision: 1.1
NodeID Livelink: 15532150

When printed, this is NOT a controlled copy

Document Information

Service-related contact information worldwide	All service-related contact information is available on this URL→	http://www.agfa.com/en/he/support/support_service/index.jsp

Issued by:
Agfa HealthCare
SIV Connectivity
Septestraat 27
B-2640 Mortsel

tel: +32.3.444.7588
email: connectivity@agfa.com

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

Copyright © December, 07
Agfa HealthCare
All rights reserved

REVISION RECORD

IHE Integration Statement Template NodeID : 8818216		
Revision Number	Date	Reason for Change
1.0	28 th Sep 2005	Initial revision
1.1	01 st Feb 2007	Update on Agfa logo and Revision Record layout.

IHE Integration Statement Drystar 5503 & 5503 Mammo		
Revision Number	Date	Reason for Change
1.0	20 th March 2007	Initial revision
1.1	23 rd October 2007	Updated with new Profiles and Actors: - ATNA & SA (Secure Application) - MAMMO & PS (Print Server) - CT & Time Client Added table 2 for Drystar 5503 Mammo with version 5.0.0 C2 Added table 3 for Drystar 5503 with version 5.2.0 C1 Added table 4 for Drystar 5503 Mammo with version 5.2.0 C1

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	Drystar 5503	5.0.0 C2	15.05.2006
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Consistent Presentation of Images (CPI)	Print Server (PS)	None	
Audit Trail and Node Authentication (ATNA)	Secure Application (SA) Secure Node (SN)	None	
Consistent Time (CT)	Time Client	None	
Internet address for vendor's IHE information: http://www.agfa.com/healthcare/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7			
DICOM	http://www.agfa.com/healthcare/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	Drystar 5503 Mammo	5.0.0 C2	15.05.2006
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Mammography Images (MAMMO)	Print Server (PS)	None	
Consistent Presentation of Images (CPI)	Print Server (PS)	None	
Audit Trail and Node Authentication (ATNA)	Secure Application (SA) Secure Node (SN)	None	
Consistent Time (CT)	Time Client	None	
Internet address for vendor's IHE information: http://www.agfa.com/healthcare/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7			
DICOM	http://www.agfa.com/healthcare/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	Drystar 5503	5.2.0 C1	01.06.2007
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Consistent Presentation of Images (CPI)	Print Server (PS)	None	
Audit Trail and Node Authentication (ATNA)	Secure Application (SA) Secure Node (SN)	None	
Consistent Time (CT)	Time Client	None	
Internet address for vendor's IHE information: http://www.agfa.com/healthcare/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7			
DICOM	http://www.agfa.com/healthcare/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	Drystar 5503 Mammo	5.2.0 C1	01.06.2007
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Mammography Images (MAMMO)	Print Server (PS)	None	
Consistent Presentation of Images (CPI)	Print Server (PS)	None	
Audit Trail and Node Authentication (ATNA)	Secure Application (SA) Secure Node (SN)	None	
Consistent Time (CT)	Time Client	None	
Internet address for vendor's IHE information: http://www.agfa.com/healthcare/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7			
DICOM	http://www.agfa.com/healthcare/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	



Document Metadata

Title:	Drystar 5503
Livelink ID:	15532150
Version#:	4
Version Date:	2007/12/04 02:57 PM CET
Status:	Approved on 2007/12/05 01:01 PM CET
Owner:	Peter Buytaert (awabr)
Created By:	Linda De Wolf (agldw)
Created Date:	2007/03/21 01:51 PM CET

This document was approved by:

Signatures:

1. Bruno Laffin (NAWYV) on 2007/12/05 09:44 AM CET

Version & Status History

Version#	Date Created	Status
4	2007/12/04 02:57 PM CET	Approved - 2007/12/05
3	2007/03/28 11:51 AM CET	Approved - 2007/03/28
2	2007/03/21 01:53 PM CET	Approved - 2007/03/28
1	2007/03/21 01:51 PM CET	