

AGFA HEALTHCARE IHE Integration Statement

IMPAX Data Center 1.0

Status: Released
Document No. 001181
Revision: 1.2
NodeID LiveLink : 20857760

When printed, this is NOT a controlled copy

Document Information

Service-related contact information worldwide	All service-related contact information is available on this URL→	http://www.agfa.com/en/he/support/support_service/index.jsp

Issued by:
Agfa HealthCare
SIV Connectivity
Septestraat 27
B-2640 Mortsel

tel: +32 3 444 7588
email: connectivity@agfa.com

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

Copyright © January, 10
Agfa HealthCare
All rights reserved

REVISION RECORD

IHE Integration Statement Template NodeID : 8818216		
Revision Number	Date	Reason for Change
1.0	28 Sep 2005	Initial revision
1.1	01 Feb 2007	Update on Agfa logo and Revision Record layout.
1.2	19 Mar 2008	Insert LiveLink NodeID in header and frontpage, adjust links.

IHE Integration Statement IMPAX Data Center 1.0		
Revision Number	Date	Reason for Change
1.0	August 2007	Initial revision (Waterloo)
1.1	November 2007	Adapted to Agfa Healthcare Template 8818216
1.2	January 2009	Update layout and Remove non supported workflow profiles

Agfa HealthCare

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	IMPAX Data Center	1.0	August 3, 2007
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Mammography Image	Image Manager/Archive	none	
NM Image	Image Manager/Archive	none	
Consistent Presentation of Images	Image Manager/Archive	none	
Image Fusion	Image Manager/Archive	none	
Evidence Documents	Image Manager/Archive	none	
Eye Care Evidence Documents	Image Manager/Archive	none	
Key Image Note	Image Manager/Archive	none	
Simple Image and Numeric Report	Report Repository	none	
Access to Radiology Information	Report Repository	none	
	Image Manager/Archive	none	
Patient Identifier Cross-referencing	Patient Identifier Cross-reference Consumer	PIX Update Notification	
Consistent Time	Time Client	none	
Audit Trail and Node Authentication	Secure Node	Encryption Option	
	{any supported IHE actor}	none	
	Audit Record Repository	none	
Internet address for vendor's IHE information: www.agfa.com/healthcare/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	www.agfa.com/healthcare/hl7		
DICOM	www.agfa.com/healthcare/dicom		
Links to general information on IHE			
In North America: www.ihe.net		In Europe: http://www.ihe-europe.net/	
In Japan: www.jira-net.or.jp/ihe-i			



Details as of PDF Creation Date

Document Metadata

Title:	001181_IMPAX Data Center 1.0 IHE Integration Statement
Livelink ID:	20857760
Version#:	3
Version Date:	2010/01/06 11:08 AM CET
Status:	Approved on 2010/01/07 12:54 PM CET
Owner:	Peter Buytaert (awabr)
Created By:	Peter Buytaert (awabr)
Created Date:	2007/11/28 10:56 AM CET
PDF Creation Date:	2010/01/07 12:54 PM CET

This document was approved by:

Signatures:

1. Bruno Laffin (nawyv) on 2010/01/07 12:16 PM CET

Detailed Approver History:

- Approval task originally assigned to and completed by Bruno Laffin (nawyv)

Version & Status History

Version#	Date Created	Status
3	2010/01/06 11:08 AM CET	Approved - 2010/01/07
2	2007/11/28 11:01 AM CET	Obsolescence Check Event Completed - 2008/12/01 Approved - 2007/11/30
1	2007/11/28 10:56 AM CET	