

AGFA HEALTHCARE IHE Integration Statement



ORBIS Service Pack 8

Status: Released
Document No. 001216
Revision: 1.0
LiveLink NodeID : 22846278

When printed, this is NOT a controlled copy

Document Information

Service-related contact information worldwide	All service-related contact information is available on this URL→	http://www.agfa.com/en/he/support/support_service/index.jsp

Issued by:
Agfa HealthCare
SIV Connectivity
Septestraat 27
B-2640 Mortsel

tel: +32.3.444.7588
email: connectivity@agfa.com

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

Copyright © February, 08
Agfa HealthCare
All rights reserved

REVISION RECORD

IHE Integration Statement Template NodeID : 8818216		
Revision Number	Date	Reason for Change
1.0	28 Sep 2005	Initial revision
1.1	01 Feb 2007	Update on Agfa logo and Revision Record layout.

IHE Integration Statement ORBIS		
Revision Number	Date	Reason for Change
1.0	15 Feb 2008	Initial Version

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	ORBIS	SP08	15 Feb 2008
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Audit Trail and Node Authentication (ATNA)	(SA) Secure Application	Radiology Audit Trail Option	
	(SN) Secure Node		
Consistent Time (CT)	(TIME_CLIENT) Time Client		
Patient Administration Management (PAM)	(PDS) Patient Demographics Supplier	Merge	
	(PES) Patient Encounter Supplier		
Patient Information Reconciliation (PIR)	(ADT) ADT Patient Registration		
	(OP) Order Placer		
Patient Identifier Cross-referencing for MPI (PIX)	(PAT_IDENTITY_SRC) Patient Identity Source		
Scheduled Workflow (SWF)	(ADT) ADT Patient Registration		
	(OP) Order Placer		
Cross-Enterprise Clinical Documents Share (XDS.a)	(DOC_SOURCE) Document Source		
Internet address for vendor's IHE information: http://www.agfa.com/healthcare/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	http://www.agfa.com/healthcare/hl7		
DICOM	http://www.agfa.com/healthcare/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	

**Document Metadata**

Title:	ORBIS SP8
Livelink ID:	22846278
Version#:	4
Version Date:	2008/02/22 12:37 PM CET
Status:	Approved on 2008/02/22 02:36 PM CET
Owner:	Peter Buytaert (awabr)
Created By:	Peter Buytaert (awabr)
Created Date:	2008/02/20 03:00 PM CET

This document was approved by:

Signatures:

1. Bruno Laffin (NAWYV) on 2008/02/22 01:59 PM CET

Version & Status History

Version#	Date Created	Status
4	2008/02/22 12:37 PM CET	Approved - 2008/02/22
3	2008/02/22 11:49 AM CET	Rejected - 2008/02/22
2	2008/02/20 03:02 PM CET	Approved - 2008/02/22
1	2008/02/20 03:00 PM CET	