

AGFA HEALTHCARE IHE Integration Statement

→ **NX 2.0.8100 / 3.0.8100**

Status: Released
Document No. 001279
Revision: 1.1
NodeID LiveLink: 27647569

When printed, this is NOT a controlled copy

Agfa HealthCare

Document Information

Service-related contact information worldwide	All service-related contact information is available on this URL→	http://www.agfa.com/en/he/support/support_service/index.jsp

Issued by:
Agfa HealthCare
SIV Connectivity
Septestraat 27
B-2640 Mortsel

tel: +32 3 444 7588
email: connectivity@agfa.com

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

Copyright © December, 09
Agfa HealthCare
All rights reserved

REVISION RECORD

IHE Integration Statement Template NodeID : 8818216		
Revision Number	Date	Reason for Change
1.0	28 Sep 2005	Initial revision
1.1	01 Feb 2007	Update on Agfa logo and Revision Record layout.
1.2	19 Mar 2008	Insert LiveLink NodeID in header and frontpage, adjust links.

IHE Integration Statement NX 2.0.8100/3.0.8100		
Revision Number	Date	Reason for Change
1.0	08 June 2009	Initial version
1.1	10 december 2009	Add Profile Consistent Time – Actor Time Client

Agfa HealthCare

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	NX	2.0.8100 / 3.0.8100	June 08, 2009
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Acquisition Modality	Broad Worklist Query	
		PPS Exception Management	
		Billing and Material Management	
Patient Info Reconciliation	Acquisition Modality	none	
Consistent Presentation of Images	Acquisition Modality	none	
	Print Composer	none	
Audit Trail and Node Authentication	Secure Application	none	
Consistent Time	Time Client	none	
Mammo Integration Profile	Acquisition Modality	none	
Portable Data for Imaging	Portable Media Creator	none	
Internet address for vendor's IHE information: www.agfa.com/healthcare/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	www.agfa.com/healthcare/hl7		
DICOM	www.agfa.com/healthcare/dicom		
Links to general information on IHE			
In North America: www.ihe.net		In Europe: http://www.ihe-europe.net/	
		In Japan: www.jira-net.or.jp/ihe-j	



Details as of PDF Creation Date

Document Metadata

Title:	001279_NX 2.0.8100/3.0.8100 IHE Integration Statement
Livelihood ID:	27647569
Version#:	3
Version Date:	2009/12/10 11:20 AM CET
Status:	Approved on 2009/12/21 09:21 AM CET
Owner:	Peter Buytaert (awabr)
Created By:	Peter Buytaert (awabr)
Created Date:	2009/05/14 01:30 PM CET
PDF Creation Date:	2009/12/21 09:21 AM CET

This document was approved by:

Signatures:

1. Bruno Laffin (nawyv) on 2009/12/18 05:24 PM CET

Detailed Approver History:

- Approval task originally assigned to and completed by Bruno Laffin (nawyv)

Version & Status History

Version#	Date Created	Status
3	2009/12/10 11:20 AM CET	Approved - 2009/12/21
2	2009/06/08 03:13 PM CET	Approved - 2009/06/17 Reviewed - 2009/06/16
1	2009/05/14 01:30 PM CET	