

# AGFA HEALTHCARE IHE Integration Statement

→ **NX 2.0.8200 / 3.0.8200**

Status: Released  
Document No. 001309  
Revision: 1.1  
NodeID LiveLink: 28582948

**When printed, this is NOT a controlled copy**

Agfa HealthCare

## Document Information

---

|  |   |   |
|--|---|---|
| <b>Service-related contact information worldwide</b> | All service-related contact information is available on this URL→ | <a href="http://www.agfa.com/en/he/support/support_service/index.jsp">http://www.agfa.com/en/he/support/support_service/index.jsp</a> |
|  |   |   |

Issued by:  
Agfa HealthCare  
SIV Connectivity  
Septestraat 27  
B-2640 Mortsel

tel: +32 3 444 7588  
email: [connectivity@agfa.com](mailto:connectivity@agfa.com)

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

**Copyright © December, 09**  
**Agfa HealthCare**  
**All rights reserved**

## REVISION RECORD

---

| IHE Integration Statement Template NodeID : 8818216 |             |   |
|---|-------------|---|
| Revision Number                                     | Date        | Reason for Change   |
| 1.0   | 28 Sep 2005 | Initial revision  |
| 1.1   | 01 Feb 2007 | Update on Agfa logo and Revision Record layout.               |
| 1.2   | 19 Mar 2008 | Insert LiveLink NodeID in header and frontpage, adjust links. |

| IHE Integration Statement NX 2.0.8200/3.0.8200 |                  |   |
|--|------------------|---|
| Revision Number                                | Date             | Reason for Change                               |
| 1.0  | 01 december 2009 | Initial version                                 |
| 1.1  | 09 december 2009 | Add Profile Consistent Time – Actor Time Client |

Agfa HealthCare

| <b>IHE Integration Statement</b>   |  |  |                       |
|--|--|--|-----------------------|
| <b>Vendor</b>  | <b>Product Name</b>  | <b>Version</b>   | <b>Date</b>           |
| <b>AGFA HealthCare</b>   | <b>NX</b>  | <b>2.0.8200 / 3.0.8200</b>   | <b>December, 2009</b> |
| This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below: |  |  |                       |
| <b>Integration Profiles Implemented</b>  | <b>Actors Implemented</b>  | <b>Options Implemented</b>   |                       |
| <b>Scheduled Workflow</b>  | <b>Acquisition Modality</b>  | <b>Broad Worklist Query</b>  |                       |
|  |  | <b>PPS Exception Management</b>  |                       |
|  |  | <b>Billing and Material Management</b>   |                       |
| <b>Patient Info Reconciliation</b>   | <b>Acquisition Modality</b>  | none   |                       |
| <b>Consistent Presentation of Images</b>   | <b>Acquisition Modality</b>  | none   |                       |
|  | <b>Print Composer</b>  | none   |                       |
| <b>Audit Trail and Node Authentication</b>   | <b>Secure Application</b>  | none   |                       |
| <b>Consistent Time</b>   | <b>Time Client</b>   | none   |                       |
| <b>Mammo Integration Profile</b>   | <b>Acquisition Modality</b>  | none   |                       |
| <b>Portable Data for Imaging</b>   | <b>Portable Media Creator</b>  | none   |                       |
| <b>Internet address for vendor's IHE information:</b> <a href="http://www.agfa.com/healthcare/ihe">www.agfa.com/healthcare/ihe</a>                         |  |  |                       |
| <b>Links to Standards Conformance Statements for the Implementation</b>  |  |  |                       |
| <b>HL7</b>   | <a href="http://www.agfa.com/healthcare/hl7">www.agfa.com/healthcare/hl7</a>     |  |                       |
| <b>DICOM</b>   | <a href="http://www.agfa.com/healthcare/dicom">www.agfa.com/healthcare/dicom</a> |  |                       |
| <b>Links to general information on IHE</b>   |  |  |                       |
| In North America: <a href="http://www.ihe.net">www.ihe.net</a>   |  | In Europe: <a href="http://www.ihe-europe.net/">http://www.ihe-europe.net/</a>   |                       |
|  |  | In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a> |                       |



### Details as of PDF Creation Date

#### Document Metadata

|                           |   |
|---------------------------|---|
| <b>Title:</b>             | 001309_NX 2.0.8200/3.0.8200 IHE Integration Statement |
| <b>Livelihood ID:</b>     | 28582948  |
| <b>Version#:</b>          | 3   |
| <b>Version Date:</b>      | 2009/12/09 04:24 PM CET                               |
| <b>Status:</b>            | Approved on 2009/12/21 09:25 AM CET                   |
| <b>Owner:</b>             | Peter Buytaert (awabr)                                |
| <b>Created By:</b>        | Peter Buytaert (awabr)                                |
| <b>Created Date:</b>      | 2009/10/19 02:11 PM CET                               |
| <b>PDF Creation Date:</b> | 2009/12/21 09:27 AM CET                               |

**This document was approved by:**

#### Signatures:

1. Bruno Laffin (nawyv) on 2009/12/18 05:22 PM CET

#### Detailed Approver History:

- Approval task originally assigned to and completed by Bruno Laffin (nawyv)

#### Version & Status History

| Version# | Date Created            | Status   |
|----------|-------------------------|--|
| 3        | 2009/12/09 04:24 PM CET | Approved - 2009/12/21                                    |
| 2        | 2009/12/01 02:41 PM CET | Approval Cancelled - 2009/12/09<br>Reviewed - 2009/12/09 |
| 1        | 2009/10/19 02:11 PM CET |  |