

AGFA HEALTHCARE IHE Integration Statement



Document No. 001312
Revision: 1.0
NodeID Livelink: 28864796

When printed, this is NOT a controlled copy

Document Information

Service-related contact information worldwide	All service-related contact information is available on this URL→	http://www.agfa.com/en/he/support/support_service/index.jsp

Issued by:
Agfa HealthCare
SIV Connectivity
Septestraat 27
B-2640 Mortsel

tel: +32 3 444 7588
email: connectivity@agfa.com

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

Copyright © December, 09
Agfa HealthCare
All rights reserved

REVISION RECORD

IHE Integration Statement Template NodeID : 8818216		
Revision Number	Date	Reason for Change
1.0	28 Sep 2005	Initial revision
1.1	01 Feb 2007	Update on Agfa logo and Revision Record layout.

IHE Integration Statement IMPAX RIS 5.7		
Revision Number	Date	Reason for Change
1.0	November 2009	Initial Revision

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	IMPAX RIS Interface Version 2005.3.4	5.7	May 2010¹
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Department System Scheduler / Order Filler	Departmental appointment Notification Instance Availability Notification	
	Performed Procedure Step Manager		
Patient Information Reconciliation (PIR)	Department System Scheduler / Order Filler		
	Performed Procedure Step Manager		
Basic Security (SEC)	Secure Node		
	Department System Scheduler / Order Filler		
Charge Posting (CHG)	Department System Scheduler / Order Filler		
	Performed Procedure Step Manager		
Patient Demographics Query (PDQ)	Patient Demographics Consumer		
Consistent Time (CT)	Time Client		
Audit Trail and Node Authentication (ATNA)	Secure Application		
	Secure Node		
Internet address for vendor's IHE information: www.agfa.com/healthcare/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	www.agfa.com/healthcare/hl7		
DICOM	www.agfa.com/healthcare/dicom		
Links to general information on IHE			
In North America: www.ihe.net		In Europe: www.ihe-europe.net	
		In Japan: www.jira-net.or.jp/ihe-j	

¹ This is the global release date. Please contact your local sales representative for the Regional release date.



Details as of PDF Creation Date

Document Metadata

Title:	001312 IMPAX RIS 5.7 IHE Integration Statement
Livelihood ID:	28864796
Version#:	3
Version Date:	2009/12/02 01:58 PM CET
Status:	Approved on 2009/12/10 12:55 PM CET
Owner:	Peter Buytaert (awabr)
Created By:	Peter Buytaert (awabr)
Created Date:	2009/11/30 10:21 AM CET
PDF Creation Date:	2009/12/10 12:55 PM CET

This document was approved by:

Signatures:

1. Sofie Neyrinck (amirm) on 2009/12/10 12:02 PM CET
2. Ann De Roose (amiqe) on 2009/12/03 05:17 PM CET
3. Bruno Laffin (nawyv) on 2009/12/03 04:14 PM CET

Detailed Approver History:

- Approval task originally assigned to and completed by Sofie Neyrinck (amirm)
- Approval task originally assigned to and completed by Bruno Laffin (nawyv)
- Approval task originally assigned to and completed by Ann De Roose (amiqe)

Version & Status History

Version#	Date Created	Status
3	2009/12/02 01:58 PM CET	Approved - 2009/12/10
2	2009/11/30 10:24 AM CET	Reviewed - 2009/12/01
1	2009/11/30 10:21 AM CET	