

AGFA HEALTHCARE IHE Integration Statement



IMPAX EE R20 V

Document No. 001314
Revision: 1.0
NodeID LiveLink : 29107314

When printed, this is NOT a controlled copy

Agfa HealthCare

Document Information

Service-related contact information worldwide	All service-related contact information is available on this URL→	http://www.agfa.com/en/he/support/support_service/index.jsp

Issued by:
Agfa HealthCare
SIV Connectivity
Septestraat 27
B-2640 Mortsel

tel: +32 3 444 7588
email: connectivity@agfa.com

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

Copyright © January, 10
Agfa HealthCare
All rights reserved

REVISION RECORD

IHE Integration Statement Template NodeID : 8818216		
Revision Number	Date	Reason for Change
1.0	28 Sep 2005	Initial revision
1.1	01 Feb 2007	Update on Agfa logo and Revision Record layout.
1.2	19 Mar 2008	Insert LiveLink NodeID in header and frontpage, adjust links.

IHE Integration Statement IMPAX EE R20 V		
Revision Number	Date	Reason for Change
1.0	07 Jan 2010	Initial version

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	IMPAX EE	R20 V	Nov. 3rd 2009
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Image Display Evidence Creator		
Post Processing Workflow	Image Display Evidence Creator		
Consistent Presentation of Images	Evidence Creator Image Display Print Composer		
Key Image Note	Evidence Creator Image Display		
Simple Image and Numeric Report Portable Data for Imaging	Report Reader Image Display Portable Media Creator Portable Media Importer Print Composer Report Reader	Web Content	
Access to Radiology Information	Image Display Report Reader	Multiple Sources Multiple Sources	
Import Reconciliation Workflow	Importer	Scheduled Importer	
Teaching File and Clinical Trial Export	Export Selector		
Cath Workflow	Image Display		
Echo Workflow	Image Display		
Consistent Time	Time Client	None	
Audit Trail and Node Authentication	Secure Node	Radiology Audit Trail	
Internet address for vendor's IHE information: www.agfa.com/healthcare/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	www.agfa.com/healthcare/hl7		
DICOM	www.agfa.com/healthcare/dicom		
Links to general information on IHE			
In North America: www.ihe.net		In Europe: www.ihe-europe.net	In Japan: www.jira-net.or.jp/ihe-j



Details as of PDF Creation Date

Document Metadata

Title:	001314_IMPAX EE R20 V IHE Integration Statement
Livelihood ID:	29107314
Version#:	3
Version Date:	2010/01/18 10:15 AM CET
Status:	Approved on 2010/01/18 05:54 PM CET
Owner:	Peter Buytaert (awabr)
Created By:	Peter Buytaert (awabr)
Created Date:	2010/01/12 02:15 PM CET
PDF Creation Date:	2010/01/18 06:24 PM CET

This document was approved by:

Signatures:

1. Bruno Laffin (nawyv) on 2010/01/18 01:51 PM CET
2. Peter Steiger (awpqz) on 2010/01/18 05:23 PM CET
3. Sabine Kolics (awppx) on 2010/01/18 02:56 PM CET

Detailed Approver History:

- Approval task originally assigned to and completed by Sabine Kolics (awppx)
- Approval task originally assigned to and completed by Bruno Laffin (nawyv)
- Approval task originally assigned to and completed by Peter Steiger (awpqz)

Version & Status History

Version#	Date Created	Status
3	2010/01/18 10:15 AM CET	Approved - 2010/01/18
2	2010/01/12 02:21 PM CET	
1	2010/01/12 02:15 PM CET	