

# AGFA HEALTHCARE IHE Integration Statement

## → IMPAX Cardiovascular 7.7 and 7.8

Status: Released  
Document No. 001315  
Revision: 1.1  
NodeID Livelink : 29200071

**When printed, this is NOT a controlled copy**

## Document Information

---

<b>Service-related contact information worldwide</b>	All service-related contact information is available on this URL→	<a href="http://www.agfa.com/en/he/support/support_service/index.jsp">http://www.agfa.com/en/he/support/support_service/index.jsp</a>

Issued by:  
Agfa HealthCare  
SIV Connectivity  
Septestraat 27  
B-2640 Mortsel

tel: +32 3 444 7588  
email: [connectivity@agfa.com](mailto:connectivity@agfa.com)

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

**Copyright © February, 10**  
**Agfa HealthCare**  
**All rights reserved**

## REVISION RECORD

---

IHE Integration Statement Template NodeID : 8818216		
Revision Number	Date	Reason for Change
1.0	28 Sep 2005	Initial revision
1.1	01 Feb 2007	Update on Agfa logo and Revision Record layout

IHE Integration Statement IMPAX Cardiovascular 7.7 and 7.8		
Revision Number	Date	Reason for Change
1.0	January 22 , 2010	Initial version
1.1	February 11, 2010	Update document numbering

# 1 IMPAX CARDIOVASCULAR 7.7

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	IMPAX Cardiovascular	7.7	September 2009
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Cardiac Cath Workflow (CATH)	Image Manager / Archive	Cardiac Cath, Intermittently Connected Modality	
	Image Display	None	
	PPS Manager	None	
Echo Workflow (ECHO)	Image Manager / Archive	Echocardiography, Intermittently Connected Modality	
	Image Display	Stress Echo	
	PPS Manager	None	
NM Image	Image Manager / Archive	None	
	Image Display	Cardiac NM Option	
Internet address for vendor's IHE information: <a href="http://www.agfa.com/healthcare/ihe">www.agfa.com/healthcare/ihe</a>			
Links to Standards Conformance Statements for the Implementation			
HL7	<a href="http://www.agfa.com/healthcare/hl7">www.agfa.com/healthcare/hl7</a>		
DICOM	<a href="http://www.agfa.com/healthcare/dicom">www.agfa.com/healthcare/dicom</a>		
Links to general information on IHE			
In North America: <a href="http://www.ihe.net">www.ihe.net</a>		In Europe: <a href="http://www.ihe-europe.net">www.ihe-europe.net</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>

## 2 IMPAX CARDIOVASCULAR 7.8

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
<b>AGFA HealthCare</b>	<b>IMPAX Cardiovascular</b>	<b>7.8</b>	<b>February 2010</b>
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
<b>Cardiac Cath Workflow (CATH)</b>	<b>Image Manager / Archive</b>	<b>Cardiac Cath, Intermittently Connected Modality</b>	
	<b>Image Display</b>	<b>None</b>	
	<b>PPS Manager</b>	<b>None</b>	
<b>Echo Workflow (ECHO)</b>	<b>Image Manager / Archive</b>	<b>Echocardiography, Intermittently Connected Modality</b>	
	<b>Image Display</b>	<b>Stress Echo</b>	
	<b>PPS Manager</b>	<b>None</b>	
<b>NM Image</b>	<b>Image Manager / Archive</b>	<b>None</b>	
	<b>Image Display</b>	<b>Cardiac NM Option</b>	
<b>Internet address for vendor's IHE information:</b> <a href="http://www.agfa.com/healthcare/ihe">www.agfa.com/healthcare/ihe</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	<a href="http://www.agfa.com/healthcare/hl7">www.agfa.com/healthcare/hl7</a>		
<b>DICOM</b>	<a href="http://www.agfa.com/healthcare/dicom">www.agfa.com/healthcare/dicom</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.ihe.net">www.ihe.net</a>		In Europe: <a href="http://www.ihe-europe.net">www.ihe-europe.net</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>



### Details as of PDF Creation Date

#### Document Metadata

<b>Title:</b>	001315_IMPAX Cardiovascular 7.7 and 7.8 IHE Integration Statement
<b>Livelink ID:</b>	29200071
<b>Version#:</b>	3
<b>Version Date:</b>	2010/02/11 04:24 PM CET
<b>Status:</b>	Approved on 2010/02/22 12:52 PM CET
<b>Owner:</b>	Paul Seifert (amica)
<b>Created By:</b>	Paul Seifert (amica)
<b>Created Date:</b>	2010/01/22 09:30 PM CET
<b>PDF Creation Date:</b>	2010/02/22 12:53 PM CET

**This document was approved by:**

#### Signatures:

1. Paul Seifert (amica) on 2010/02/11 05:04 PM CET
2. Bruno Laffin (nawyv) on 2010/02/17 01:43 PM CET

#### Detailed Approver History:

- Approval task originally assigned to and completed by Paul Seifert (amica)
- Approval task originally assigned to and completed by Bruno Laffin (nawyv)

#### Version & Status History

<b>Version#</b>	<b>Date Created</b>	<b>Status</b>
3	2010/02/11 04:24 PM CET	Approved - 2010/02/22
2	2010/02/11 04:22 PM CET	
1	2010/01/22 09:30 PM CET	Reviewed - 2010/02/10