

AGFA HEALTHCARE IHE Integration Statement



IMPAX SE Client/Standalone Workstation

Status: Released
Document No. 001168
Revision: 1.0

NodeID Livelink : 19225392

When printed, this is NOT a controlled copy

Document Information

Service-related contact information worldwide	All service-related contact information is available on this URL→	http://www.agfa.com/en/he/support/support_service/index.jsp

Issued by:
Agfa HealthCare
SI Connectivity
Septestraat 27
B-2640 Mortsel

tel: +32.3.444.7588
email: connectivity@agfa.com

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

Copyright © September, 07
Agfa HealthCare
All rights reserved

REVISION RECORD

IHE Integration Statement Template NodeID : 8818216		
Revision Number	Date	Reason for Change
1.0	28 Sep 2005	Initial revision
1.1	01 Feb 2007	Update on Agfa logo and Revision Record layout.

IHE Integration Statement IMPAX SE Client/Standalone		
Revision Number	Date	Reason for Change
1.0	12 Jul 2007	Initial version

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	IMPAX SE Client/Standalone	8.5	12 Jul 2007
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Consistent Presentation of Images	Image Display	None	
	Evidence Creator (GSPS Creator)	None	
	Print composer	None	
	Image Manager/Archive	None	
Portable Data for Imaging	Image Display	None	
	Portable Media Creator	None	
	Portable Media Importer	None	
	Report Reader	None	
Mammography Profile	Image Display	None	
	Print composer	None	
	Image Manager/Archive	None	
Internet address for vendor's IHE information: http://www.agfa.com/healthcare/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	http://www.agfa.com/healthcare/hl7		
DICOM	http://www.agfa.com/healthcare/dicom		
Links to general information on IHE			
In North America: http://www.ihe.net	In Europe: www.ihe-europe.org	In Japan: www.jira-net.or.jp/ihe-j	

**Document Metadata**

Title:	IMPAX SE Client/Standalone Workstation
Livelihood ID:	19225392
Version#:	2
Version Date:	2007/09/20 02:12 PM CET
Status:	Approved on 2007/09/20 04:55 PM CET
Owner:	Peter Buytaert (awabr)
Created By:	Peter Buytaert (awabr)
Created Date:	2007/09/20 02:09 PM CET

This document was approved by:

Signatures:

1. Bruno Laffin (NAWYV) on 2007/09/20 02:39 PM CET

Version & Status History

Version#	Date Created	Status
2	2007/09/20 02:12 PM CET	Approved - 2007/09/20
1	2007/09/20 02:09 PM CET	