


AGFA HEALTHCARE IHE Integration Statement



NX 1.0

Status: Released

Document no : 001086

Revision 1.0

NoteID Livelink : 10889599

When printed, this is NOT a controlled copy

Document Information

Service-related contact information worldwide	All service-related contact information is available on this URL→	http://www.agfa.com/en/he/support/support_service/index.jsp

Issued by:
Agfa HealthCare
<SIV Connectivity >
<Septestraat 27>
<B-2640 Mortsel>

tel: <+32.3.444.7338>
email: <connectivity@agfa.com>

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

Copyright © April, 06
Agfa HealthCare
All rights reserved

REVISION RECORD

Revision Number	Date	Reason for Change
1.0	23-March-2006	Initial version

Agfa HealthCare

IHE Integration Statement		Date
		23 March 2006
Vendor	Product Name	Version
AGFA Healthcare	NX	1.0
This product implements all transactions required in the IHE Technical Framework to support the IHE Profiles, Actors and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow	Acquisition Modality	Broad Worklist Query PPS Exception Management Billing and Material Management
Patient Info Reconciliation	Acquisition Modality	none
Consistent Presentation of Images	Acquisition Modality	none
	Print Composer	
Basic Security Integration	Acquisition Modality	None
	Secure Node	
Audit Trail and Node Authentication	Acquisition Modality	none
	Secure Node	
Portable Data for Imaging	Portable Media Creator	none
Internet address for vendor's IHE information: http://www.agfa.com/healthcare/ihe		
Links to Standards Conformance Statements for the Implementation		
HL7		
DICOM	http://www.agfa.com/healthcare/dicom	
Links to general information on IHE		
In North America: www.ihe.net	In Europe: www.ihe-europe.org	In Japan: http://www.ihe-j.org/en