

# AGFA HEALTHCARE IHE Integration Statement



## IMPAX 6.5.x Solution

Status: Released  
Document No. 001305  
Revision: 1.5  
NodeID Livelink: 28255601

**When printed, this is NOT a controlled copy**

Agfa HealthCare

## Document Information

---

<b>Service-related contact information worldwide</b>	All service-related contact information is available on this URL→	<a href="http://www.agfahealthcare.com/global/en/main/contact/index.jsp">http://www.agfahealthcare.com/global/en/main/contact/index.jsp</a>

Issued by:  
Agfa HealthCare  
SIV Connectivity  
Septestraat 27  
B-2640 Mortsel

tel: +32 3 444 7588  
email: [connectivity@agfa.com](mailto:connectivity@agfa.com)

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

**Copyright © January, 14**  
**Agfa HealthCare**  
**All rights reserved**

## REVISION RECORD

---

IHE Integration Statement Template NodeID : 8818216		
Revision Number	Date	Reason for Change
1.0	28 Sep 2005	Initial revision
1.1	01 Feb 2007	Update on Agfa logo and Revision Record layout

Revision Number	Date	Reason for Change
1.0	October 2009	Initial version
1.1	October 2010	Add support for ED, NM and REM profiles as Image Manager/Archive Actors
1.2	May 2011	- Add support for RWF profile as Image Manager/Archive and PPS Manager - Add IMPAX version 6.5.1
1.3	October 2011	Add IMPAX version 6.5.2
1.4	March 2012	Add IMPAX version 6.5.3
1.5	December 2013	Add IMPAX version 6.5.4 and 6.5.5

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
<b>AGFA HealthCare</b>	<b>IMPAX</b>	<b>6.5.0</b>	<b>October 2010</b>
		<b>6.5.1</b>	<b>June 2011</b>
		<b>6.5.2</b>	<b>November 2011</b>
		<b>6.5.3</b>	<b>March 2012</b>
		<b>6.5.4</b>	<b>April 2013</b>
		<b>6.5.5</b>	<b>December 2013</b>
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
<b>Scheduled Workflow (SWF)</b>	<b>Image Manager / Archive</b>	<b>None</b>	
	<b>Image Display</b>	<b>None</b>	
	<b>Evidence Creator</b>	<b>None</b>	
	<b>PPS Manager</b>	<b>None</b>	
<b>Patient information Reconciliation (PIR)</b>	<b>Image Manager / Archive</b>	<b>None</b>	
	<b>PPS Manager</b>	<b>None</b>	
<b>Consistent Presentation of Images (CPI)</b>	<b>Image Manager / Archive</b>	<b>None</b>	
	<b>Image Display</b>	<b>None</b>	
	<b>Evidence Creator</b>	<b>None</b>	
	<b>Print Composer</b>	<b>None</b>	
<b>Portable data for Imaging (PDI)</b>	<b>Portable Media Creator</b>	<b>None</b>	
	<b>Portable Media Importer</b>	<b>None</b>	
	<b>Image Display</b>	<b>None</b>	
	<b>Print Composer</b>	<b>None</b>	
<b>Access to Radiology Information (ARI)</b>	<b>Image Manager / Archive</b>	<b>None</b>	
	<b>Image Display</b>	<b>None</b>	
<b>Consistent Time (CT)</b>	<b>Time Client</b>	<b>None</b>	
<b>Mammography Image (MAMMO)</b>	<b>Image Manager / Archive</b>	<b>None</b>	
	<b>Image Display</b>	<b>None</b>	
	<b>Print Composer</b>	<b>None</b>	
<b>Cardiac Catheterization Workflow (CATH)</b>	<b>Image Display</b>	<b>None</b>	
<b>Simple Image &amp; Numeric Report (SINR)</b>	<b>Report Repository</b>	<b>None</b>	
<b>Reporting Workflow (RWF)</b>	<b>Image Manager / Archive</b>	<b>None</b>	
	<b>PPS Manager</b>	<b>None</b>	
<b>Key Image Notes (KIN)</b>	<b>Image Manager / Archive</b>	<b>None</b>	
	<b>Evidence Creator</b>	<b>None</b>	
	<b>Image Display</b>	<b>None</b>	
<b>Evidence Document (ED)</b>	<b>Image Manager / Archive</b>	<b>None</b>	

Agfa HealthCare

Integration Profiles Implemented	Actors Implemented	Options Implemented
NM Image (NM)	Image Manager / Archive	None
Radiation Exposure Monitoring (REM)	Image Manager / Archive	None
<b>Internet address for vendor's IHE information:</b> <a href="http://www.agfahealthcare.com/ihe">www.agfahealthcare.com/ihe</a>		
<b>Links to Standards Conformance Statements for the Implementation</b>		
HL7	<a href="http://www.agfahealthcare.com/hl7">www.agfahealthcare.com/hl7</a>	
DICOM	<a href="http://www.agfahealthcare.com/dicom">www.agfahealthcare.com/dicom</a>	
<b>Links to general information on IHE</b>		
IHE International: <a href="http://www.ihe.net">www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org/">www.ihe-j.org/</a>