

EASYLIFT

Unique wall-mounted, floor-standing lifting device especially designed for smooth and easy positioning in full body imaging.

- Designed for quick and easy lifting
- Ergonomic and user-friendly design
- Increased accuracy and ease of operation for improved workflow in full body imaging

Quick and easy positioning

The Easylift has been developed to work with the DR FLFS Detector Positioner or the CR Full Body Cassette Holder, allowing smooth and easy lifting at greatly reduced effort.

During installation, the Easylift can be modified for left or right hand operation. Its ergonomic design and user-friendly single handle adjustment system provide for fewer handling movements together with quick and easy height adjustment. The Easylift is a robust system that greatly improves handling and positioning.



Smooth and ergonomic workflow

The Easylift brings increased accuracy and ease of operation to Full Leg Full Spine exposures, enabling precise positioning of the DR FLFS Detector Positioner or the CR Full Body Cassette Holder with respect to the X-ray tube, regardless of patient height. This results in a smoother and more ergonomic workflow and enhanced productivity.

Key part of a unique complete solution

The Easylift brings the user all the benefits with the added value of easy and safe height manipulation. It is a key part of a complete range of DR and CR Full Leg Full Spine solutions, assuring ease of use and uniformly high quality images for all Full Leg Full Spine exposures in full body imaging.

Technical specifications

GENERAL

Dimensions

- W x D x H: 39 x 13.5 x 230 cm
(15 x 5.5 x 90.5")

Weight

- 73 kg (161 lbs)
- With DR FLFS Detector Positioner or CR Full Body Cassette Holder: approx. 100 kg (220 lbs)

Installation Requirements

- An Agfa technician will perform the installation of the Easylift, based on local guidance and instructions. Please refer to the installation manual for important installation requirements and preconditions.

Operation

- During installation, the Easylift can be modified for left or right hand operation

Transport details

- Packaging dimensions (W x D x H): 60 x 39 x 238 cm (24 x 15 x 94")
- Shipped in a wooden crate

SAFETY

Approvals

- CE, UL

For more information on Agfa, please visit our website on www.agfa.com ■

Agfa and the Agfa rhombus are trademarks of Agfa-Gevaert NV, Belgium, or its affiliates. All rights reserved. All information contained herein is intended for guidance purposes only, and characteristics of the products and services described in this publication can be changed at any time without notice. Products and services may not be available for your local area. Please contact your local sales representative for availability information. Agfa-Gevaert NV diligently strives to provide as accurate information as possible, but shall not be responsible for any typographical error.

© 2019 Agfa NV
All rights reserved
Published by Agfa NV
Septestraat 27 - 2640 Mortsel
Belgium

NFWD1 GB 00201908