

AGFA HEALTHCARE IHE Integration Statement

→ **IMPAX RIS 5.8**

Document No. 001364
Revision: 1.0
NodeID Livelink: 33326029

When printed, this is NOT a controlled copy

Document Information

Service-related contact information worldwide	All service-related contact information is available on this URL→	http://www.agfahealthcare.com/global/en/main/contact/index.jsp

Issued by:
Agfa HealthCare
SIV Connectivity
Septestraat 27
B-2640 Mortsel

tel: +32 3 444 7588
email: connectivity@agfa.com

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

Copyright © June, 11
Agfa HealthCare
All rights reserved

REVISION RECORD

IHE Integration Statement Template NodeID : 8818216		
Revision Number	Date	Reason for Change
1.0	28 Sep 2005	Initial revision
1.1	01 Feb 2007	Update on Agfa logo and Revision Record layout.
1.2	19 Mar 2008	Insert Livelink NodeID in header and front page, adjust links.

IHE Integration Statement IMPAX RIS 5.8		
Revision Number	Date	Reason for Change
1.0	May 2011	Initial Revision

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	IMPAX RIS Interface Version 2005.3.5	5.8	October 2011¹
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Department System Scheduler / Order Filler	Departmental appointment Notification Instance Availability Notification	
	Performed Procedure Step Manager		
Patient Information Reconciliation (PIR)	Department System Scheduler / Order Filler		
	Performed Procedure Step Manager		
Patient Administration Management (PAM)	Patient Demographics Consumer (PDC)		
	Patient Encounter Consumer (PEC)		
Charge Posting (CHG)	Department System Scheduler / Order Filler		
	Performed Procedure Step Manager		
Patient Demographics Query (PDQ)	Patient Demographics Consumer		
Consistent Time (CT)	Time Client		
Audit Trail and Node Authentication (ATNA)	Secure Application		
Internet address for vendor's IHE information: www.agfahealthcare.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	www.agfahealthcare.com/hl7		
DICOM	www.agfahealthcare.com/dicom		
Links to general information on IHE			
In North America: www.ihe.net	In Europe: http://www.ihe-europe.net/	In Japan: www.ihe-j.org/	

¹ This is the global release date. Please contact your local sales representative for the Regional release date.



Details as of PDF Creation Date

Document Metadata

Title:	001364_IMPAX RIS 5.8 IHE Integration Statement
Livelihood ID:	33326029
Version#:	3
Version Date:	2011/06/08 11:57 AM CET
Status:	Approved on 2011/06/14 11:54 AM CET
Owner:	Peter Buytaert (awabr)
Created By:	Peter Buytaert (awabr)
Created Date:	2011/05/19 05:16 PM CET
PDF Creation Date:	2011/06/14 11:55 AM CET

This document was approved by:

Signatures:

1. Bruno Laffin (nawyv) on 2011/06/08 05:01 PM CET
2. Peter Moreno (awlzk) on 2011/06/14 11:44 AM CET

Detailed Approver History:

- **Approval Workflow started on 2011/06/08 11:59 AM CET**
 - Approval task originally assigned to and completed by Peter Moreno (awlzk) on 2011/06/14 11:44 AM CET
 - Approval task originally assigned to and completed by Bruno Laffin (nawyv) on 2011/06/08 05:01 PM CET

Version & Status History

Version#	Date Created	Status
3	2011/06/08 11:57 AM CET	Approved - 2011/06/14
2	2011/05/19 05:17 PM CET	Reviewed - 2011/06/08
1	2011/05/19 05:16 PM CET	