

AGFA HEALTHCARE IHE Integration Statement

→ **IMPAX for Cardiology EMEA 1.0**

Document No.: 001440
Revision: 1.0
NodeID LiveLink: 40740498

When printed, this is NOT a controlled copy

Document Information

Service-related contact information worldwide	All service-related contact information is available on this URL→	http://www.agfahealthcare.com/global/en/main/contact/index.jsp
--	---	---

Issued by:
Agfa HealthCare
SIV Connectivity
Septestraat 27
B-2640 Mortsel

tel: +32 3 444 7588
email: connectivity@agfa.com

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

Copyright © June, 13
Agfa HealthCare
All rights reserved

REVISION RECORD

IHE Integration Statement Template NodeID : 8818216		
Revision Number	Date	Reason for Change
1.0	28 Sep 2005	Initial revision
1.1	01 Feb 2007	Update on Agfa logo and Revision Record layout.
1.2	19 Mar 2008	Insert Livelink NodeID in header and front page, adjust links.

IHE Integration Statement IMPAX for Cardiology EMEA 1.0		
Revision Number	Date	Reason for Change
1.0	2013-05-21	Initial version

Agfa HealthCare

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	IMPAX for Cardiology EMEA¹	1.0	May, 2013
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented²	Actors Implemented	Options Implemented	
Cardiac Catheterization Workflow (CATH)	Image Manager / Archive	none	
	PPS Manager	none	
	Order Placer	none	
Echocardiography Workflow (ECHO)	Image Manager / Archive	none	
	PPS Manager	none	
	Evidence Creator	none	
	Order Placer	none	
Internet address for vendor's IHE information: www.agfahealthcare.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	www.agfahealthcare.com/hl7		
DICOM	www.agfahealthcare.com/dicom		
Links to general information on IHE			
In North America: www.ihe.net	In Europe: http://www.ihe-europe.net/	In Japan: www.ihe-j.org/	

¹ IMPAX for Cardiology EMEA 1.0 is made of the components:

- CardioReport 2.9
- Image-Arena 4.6SP1
- Connectivity Manager 2.4.1
- IMPAX 6.5.3
- Mirth Connect 2.2.1

² The implemented actors and profiles are in addition to the ones supported by IMPAX 6.5.3
(<http://www.agfahealthcare.com/global/en/he/library/libraryopen?ID=28255601>)