

AGFA HEALTHCARE IHE Integration Statement

→ **IMPAX Cardiovascular 12.x**

Status: Released
Document No. 001463
Revision: 1.1
NodeID Livelink : 44199086

When printed, this is NOT a controlled copy

Agfa HealthCare

Document Information

Service-related contact information worldwide	All service-related contact information is available on this URL→	http://www.agfahealthcare.com/global/en/main/contact/index.jsp

Issued by:
Agfa HealthCare
SIV Connectivity
Septestraat 27
B-2640 Mortsel

tel: +32 3 444 7588
email: connectivity@agfa.com

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

Copyright © October, 14
Agfa HealthCare
All rights reserved

REVISION RECORD

IHE Integration Statement Template NodeID : 8818216		
Revision Number	Date	Reason for Change
1.0	28 Sep 2005	Initial revision
1.1	01 Feb 2007	Update on Agfa logo and Revision Record layout
1.2	19 Mar 2008	Insert Livelink NodeID in header and front page, adjust links.
1.3	25 Jan 2011	Update layout + links

IHE Integration Statement IMPAX Cardiovascular 12.x		
Revision Number	Date	Reason for Change
1.0	February 11, 2014	Initial version
1.1	October 2, 2014	Added CV 12.2.x version integration statement.

1 IMPAX CARDIOVASCULAR 12.0.1

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	IMPAX Cardiovascular	12.0.1	July 2013
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Cardiac Cath Workflow (CATH)	Image Manager / Archive	Cardiac Cath, Intermittently Connected Modality	
	Image Display	None	
	PPS Manager	None	
Echo Workflow (ECHO)	Image Manager / Archive	Echocardiography, Intermittently Connected Modality	
	Image Display	Stress Echo	
	PPS Manager	None	
NM Image	Image Manager / Archive	None	
	Image Display	Cardiac NM Option	
Internet address for vendor's IHE information: www.agfa.com/healthcare/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	www.agfa.com/healthcare/hl7		
DICOM	www.agfa.com/healthcare/dicom		
Links to general information on IHE			
In North America: www.ihe.net		In Europe: www.ihe-europe.net	In Japan: www.ihe-j.org/

2 IMPAX CARDIOVASCULAR 12.1.X

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	IMPAX Cardiovascular	12.1	December 2013
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Cardiac Cath Workflow (CATH)	Image Manager / Archive	Cardiac Cath, Intermittently Connected Modality	
	Image Display	None	
	PPS Manager	None	
Echo Workflow (ECHO)	Image Manager / Archive	Echocardiography, Intermittently Connected Modality	
	Image Display	Stress Echo	
	PPS Manager	None	
NM Image	Image Manager / Archive	None	
	Image Display	Cardiac NM Option	
Internet address for vendor's IHE information: www.agfa.com/healthcare/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	www.agfa.com/healthcare/hl7		
DICOM	www.agfa.com/healthcare/dicom		
Links to general information on IHE			
In North America: www.ihe.net		In Europe: www.ihe-europe.net	In Japan: www.ihe-j.org/

3 IMPAX CARDIOVASCULAR 12.2.X

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	IMPAX Cardiovascular	12.2	October 2014
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Cardiac Cath Workflow (CATH)	Image Manager / Archive	Cardiac Cath, Intermittently Connected Modality	
	Image Display	None	
	PPS Manager	None	
Echo Workflow (ECHO)	Image Manager / Archive	Echocardiography, Intermittently Connected Modality	
	Image Display	Stress Echo	
	PPS Manager	None	
NM Image	Image Manager / Archive	None	
	Image Display	Cardiac NM Option	
Internet address for vendor's IHE information: www.agfa.com/healthcare/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	www.agfa.com/healthcare/hl7		
DICOM	www.agfa.com/healthcare/dicom		
Links to general information on IHE			
In North America: www.ihe.net		In Europe: www.ihe-europe.net	In Japan: www.ihe-j.org/