

AGFA HEALTHCARE IHE Integration Statement



IMPAX 6.6.x Solution

Status: Released
Document No. 001471
Revision: 1.2
NodeID Livelink: 44908214

When printed, this is NOT a controlled copy

Document Information

Service-related contact information worldwide	All service-related contact information is available on this URL→	http://www.agfahealthcare.com/global/en/main/contact/index.jsp

Issued by:
Agfa HealthCare
SIV Connectivity
Septestraat 27
B-2640 Mortsel

tel: +32 3 444 7588
email: connectivity@agfa.com

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

Copyright © 2015
Agfa HealthCare
All rights reserved

REVISION RECORD

IHE Integration Statement Template NodeID : 8818216		
Revision Number	Date	Reason for Change
1.0	28 Sep 2005	Initial revision
1.1	01 Feb 2007	Update on Agfa logo and Revision Record layout

Revision Number	Date	Reason for Change
1.0	April 2 nd 2014	Initial version IMPAX 6.6
1.1	February 4 th 2015	Update document title to IMPAX 6.6.x Add version IMPAX 6.6.1
1.2	November 16 th 2015	Update layout Add support of DBT and IOCM (optional) as of IMPAX 6.6.1 SU3

Agfa HealthCare

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	IMPAX	6.6	July 2014
		6.6.1	February 2015
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented		Actors Implemented	Options Implemented
CARD – Cardiology Domain			
CATH - Cardiac Catheterization Workflow	Image Display	None	
ITI – IT Infrastructure Domain			
CT - Consistent Time	Time Client	None	
RAD – Radiology Domain			
ARI - Access to Radiology Information	Image Manager / Archive	None	
	Image Display	None	
CPI - Consistent Presentation of Images	Image Manager / Archive	None	
	Image Display	None	
	Evidence Creator	None	
	Print Composer	None	
DBT – Digital Breast Tomosynthesis ¹	Image Manager / Archive	None	
	Image Display	None	
ED - Evidence Document	Image Manager / Archive	None	
IOCM – Imaging Object Change Management ²	Change Requester	None	
KIN - Key Image Notes	Image Manager / Archive	None	
	Evidence Creator	None	
	Image Display	None	
MAMMO - Mammography Image	Image Manager / Archive	None	
	Image Display	None	
	Print Composer	None	
NM - NM Image	Image Manager / Archive	None	
PDI - Portable data for Imaging	Portable Media Creator	None	
	Portable Media Importer	None	
	Image Display	None	
	Print Composer	None	
PIR - Patient information Reconciliation	Image Manager / Archive	None	
	PPS Manager	None	
REM - Radiation Exposure Monitoring	Image Manager / Archive	None	
RWF - Reporting Workflow	Image Manager / Archive	None	
	PPS Manager	None	

¹ DBT is supported as of IMPAX 6.6.1 SU3² The IOCM module is an optional module provided by the Agfa Global Professional Services Team and is supported as of IMPAX 6.6.1 SU3

Agfa HealthCare

SINR - Simple Image & Numeric Report	Report Repository	None
SWF - Scheduled Workflow	Image Manager / Archive	None
	Image Display	None
	Evidence Creator	None
	PPS Manager	None
Internet address for vendor's IHE information: www.agfahealthcare.com/ihe		
Links to Standards Conformance Statements for the Implementation		
HL7	www.agfahealthcare.com/hl7	
DICOM	www.agfahealthcare.com/dicom	
Links to general information on IHE		
IHE International: www.ihe.net In Europe: www.ihe-europe.net In Japan: www.ihe-j.org/		