

AGFA HEALTHCARE IHE Integration Statement



IMPAX RIS 6.1.0

Document No. 001505
Revision: 1.1
NodeID Livelink: 48990605

When printed, this is NOT a controlled copy

Document Information

Service-related contact information worldwide	All service-related contact information is available on this URL→	http://www.agfahealthcare.com/global/en/main/contact/index.jsp

Issued by:
Agfa HealthCare
SIV Connectivity
Septestraat 27
B-2640 Mortsel

tel: +32 3 444 7588
email: connectivity@agfa.com

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

Copyright © 2015
Agfa HealthCare
All rights reserved

REVISION RECORD

IHE Integration Statement Template NodeID : 8818216		
Revision Number	Date	Reason for Change
1.0	28 Sep 2005	Initial revision
1.1	01 Feb 2007	Update on Agfa logo and Revision Record layout.
1.2	19 Mar 2008	Insert Livelink NodeID in header and front page, adjust links.

IHE Integration Statement IMPAX RIS 6.1.0		
Revision Number	Date	Reason for Change
1.0	March 2015	Initial Revision
1.1	March 20, 2015	Add PAM/PDC option Merge

Agfa HealthCare

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	IMPAX RIS Interface Version 6.1.0	6.1.0	February 2015¹
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Department System Scheduler / Order Filler	Departmental appointment Notification Availability of PPS-Referenced Instances	
	Performed Procedure Step Manager	None	
Patient Information Reconciliation (PIR)	Department System Scheduler / Order Filler	None	
	Performed Procedure Step Manager	None	
Patient Administration Management (PAM)	Patient Demographics Consumer (PDC)	None Merge	
	Patient Encounter Consumer (PEC)	None	
Charge Posting (CHG)	Department System Scheduler / Order Filler	None	
	Performed Procedure Step Manager	None	
Patient Demographics Query (PDQ)	Patient Demographics Consumer	None	
Consistent Time (CT)	Time Client	None	
Audit Trail and Node Authentication (ATNA)	Secure Application	None	
Internet address for vendor's IHE information: www.agfahealthcare.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	www.agfahealthcare.com/hl7		
DICOM	www.agfahealthcare.com/dicom		
Links to general information on IHE			
In North America: www.ihe.net	In Europe: http://www.ihe-europe.net/	In Japan: www.ihe-j.org/	

¹ This is the global release date. Please contact your local sales representative for the Regional release date.