

AGFA HEALTHCARE IHE Integration Statement



Enterprise Imaging CWP (Consolidated Web Platform) 8.0.x

Document No.: 001508
Revision: 1.1
NodeID Livelink: 49771363

When printed, this is NOT a controlled copy

REVISION RECORD

IHE Integration Statement Template NodeID : 8818216		
Revision Number	Date	Reason for Change
1.0	28 Sep 2005	Initial revision
1.1	01 Feb 2007	Update on Agfa logo and Revision Record layout.
1.2	19 Mar 2008	Insert Livelink NodeID in header and front page, adjust links.

IHE Integration Statement Agfa Healthcare Enterprise Imaging CWP 8.0.x		
Revision Number	Date	Reason for Change
1.0	2015-04-15	Initial version
1.1	2015-09-18	Update product name to CWP 8.0.x

Agfa HealthCare

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	CWP	8.0.x	April, 2015
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Image Display	None	
Consistent Presentation of Images (CPI)	Image Display	None	
Key Image Note (KIN)	Image Display	None	
Evidence Documents (ED)	Image Display	None	
Simple Image and Numeric Report (SINR)	Report Reader	None	
Access to Radiology Information (ARI)	Report Reader	None	
	Image Display	None	
Consistent Time (CT)	Time Client	None	
Audit Trail and Node Authentication (ATNA)	Secure Application	None	
Patient Identifier Cross-referencing (PIX,PIXv3)	Patient Identifier Cross-reference Consumer	None	
Patient Demographics Query (PDQ,PDQv3)	Patient Demographics Consumer	None	
Cross-Enterprise Document Sharing (XDS.b)	Document Consumer	None	
Cross-Enterprise Document Sharing for Imaging (XDS-I.b)	Imaging Document Consumer	None	
Cross-Enterprise Sharing of Scanned Documents (XDS-SD)	Document Consumer	None	
Internet address for vendor's IHE information: www.agfahealthcare.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	www.agfahealthcare.com/hl7		
DICOM	www.agfahealthcare.com/dicom		
Links to general information on IHE			
In North America: www.ihe.net	In Europe: http://www.ihe-europe.net/	In Japan: www.ihe-j.org/	