

AGFA HEALTHCARE IHE Integration Statement

Enterprise Imaging TRANSPORT Web 1.0

Status: Released
Document No. 001523
Revision: 1.0
Livelink NodeID: 51480964

When printed, this is NOT a controlled copy

Document Information

Service-related contact information worldwide	All service-related contact information is available on this URL →	http://www.agfahealthcare.com/global/en/main/contact/index.jsp

Issued by:
Agfa HealthCare
SIV Connectivity
Septestraat 27
B-2640 Mortsel

tel: +32 3 444 7588
email: connectivity@agfa.com

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

Copyright © 2015
Agfa HealthCare
All rights reserved

REVISION RECORD

IHE Integration Statement Template NodeID : 8818216		
Revision Number	Date	Reason for Change
1.0	28 Sep 2005	Initial revision
1.1	01 Feb 2007	Update on Agfa logo and Revision Record layout.
1.2	19 Mar 2008	Insert Livelink NodeID in header and front-page, adjust links.
1.3	25 Jan 2011	Update layout + links

IHE Integration Statement Enterprise Imaging TRANSPORT Web 1.0		
Revision Number	Date	Reason for Change
1.0	07 Sep 2015	Initial version

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	TRANSPORT Web	1.0	07 Sep 2015
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
ITI – IT Infrastructure Domain			
CT - Consistent Time	Time Client	None	
ATNA - Audit Trail and Node Authentication	Secure Node	None	
XDS.b – Cross Enterprise Document Sharing	Document Source	None	
Internet address for vendor's IHE information: www.agfahealthcare.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	www.agfahealthcare.com/hl7		
DICOM	www.agfahealthcare.com/dicom		
Links to general information on IHE			
In North America: www.ihe.net		In Europe: www.ihe-europe.net	In Japan: www.ihe-j.org/