

# AGFA HEALTHCARE IHE Integration Statement



## IMPAX 6.7.x Solution

Status: Released  
Document No. 001539  
Revision: 5  
NodeID Livelink: 54939764

**When printed, this is NOT a controlled copy**

## Document Information

---

<b>Service-related contact information worldwide</b>	All service-related contact information is available on this URL→	<a href="http://www.agfahealthcare.com/global/en/main/contact/index.jsp">http://www.agfahealthcare.com/global/en/main/contact/index.jsp</a>

Issued by:  
Agfa HealthCare  
SIV Connectivity  
Septestraat 27  
B-2640 Mortsel

tel: +32 3 444 9367  
email: [connectivity@agfa.com](mailto:connectivity@agfa.com)

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

**Copyright © 2020**  
**Agfa HealthCare**  
**All rights reserved**

## REVISION RECORD

---

IHE Integration Statement Template NodeID : 8818216		
Revision Number	Date	Reason for Change
1.0	28 Sep 2005	Initial revision
1.1	01 Feb 2007	Update on Agfa logo and Revision Record layout
1.2	19 Mar 2008	Insert Livelink NodeID in header and front page, adjust links.
1.3	25 Jan 2011	Update layout + links

IHE Integration Statement IMPAX 6.7.x		
Revision Number	Date	Reason for Change
1.0	20 May 2016	Initial version IMPAX 6.7
4	22 April 2020	Update document title to IMPAX 6.7.x Solution as the document is also valid for version 6.7.1
5	24 April 2020	All corrections made in review cycle

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
<b>AGFA HealthCare</b>	<b>IMPAX</b>	<b>6.7</b>	<b>May 2016</b>
		<b>6.7.1</b>	<b>July 2020</b>
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
<b>CARD – Cardiology Domain</b>			
<b>CATH</b> - Cardiac Catheterization Workflow	Image Display	None	
<b>ITI – IT Infrastructure Domain</b>			
<b>CT</b> - Consistent Time	Time Client	None	
<b>RAD – Radiology Domain</b>			
<b>ARI</b> - Access to Radiology Information	Image Manager / Archive	None	
	Image Display	None	
<b>CPI</b> - Consistent Presentation of Images	Image Manager / Archive	None	
	Image Display	None	
	Evidence Creator	None	
	Print Composer	None	
<b>DBT</b> – Digital Breast Tomosynthesis	Image Manager / Archive	None	
	Image Display	None	
<b>ED</b> - Evidence Document	Image Manager / Archive	None	
<b>IOCM</b> – Imaging Object Change Management <sup>1</sup>	Change Requester	None	
<b>KIN</b> - Key Image Notes	Image Manager / Archive	None	
	Evidence Creator	None	
	Image Display	None	
<b>MAMMO</b> - Mammography Image	Image Manager / Archive	None	
	Image Display	None	
	Print Composer	None	
<b>NM</b> – Nuclear Medicine Image	Image Manager / Archive	None	
<b>PDI</b> - Portable data for Imaging	Portable Media Creator	None	
	Portable Media Importer	None	
	Image Display	None	
	Print Composer	None	
<b>PIR</b> - Patient information Reconciliation	Image Manager / Archive	None	
	PPS Manager	None	
<b>REM</b> - Radiation Exposure Monitoring	Image Manager / Archive	None	
<b>RWF</b> - Reporting Workflow	Image Manager / Archive	None	
	PPS Manager	None	

<sup>1</sup> The IOCM module is an optional module provided by the Agfa Global Professional Services Team

<b>SINR</b> - Simple Image & Numeric Report	Report Repository	None	
<b>SWF</b> - Scheduled Workflow	Image Manager / Archive	None	
	Image Display	None	
	Evidence Creator	None	
	PPS Manager	None	
<b>Internet address for vendor's IHE information:</b> <a href="http://www.agfahealthcare.com/ihe">www.agfahealthcare.com/ihe</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	<a href="http://www.agfahealthcare.com/hi7">www.agfahealthcare.com/hi7</a>		
<b>DICOM</b>	<a href="http://www.agfahealthcare.com/dicom">www.agfahealthcare.com/dicom</a>		
<b>Links to general information on IHE</b>			
IHE International:	<a href="http://www.ihe.net">www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org/">www.ihe-j.org/</a>