

AGFA HEALTHCARE IHE Integration Statement

→ **NX 3.0.9000**

Document No. 001549
Revision: 2
NodeID Livelink: 57497642

When printed, this is NOT a controlled copy

Document Information

Service-related contact information worldwide	All service-related contact information is available on this URL→	http://www.agfahealthcare.com/global/en/main/contact/index.jsp

Issued by:
Agfa HealthCare
V&V Connectivity
Septestraat 27
B-2640 Mortsel

tel: +32 3 444 7588
email: connectivity@agfa.com

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

Copyright © 2016
Agfa HealthCare
All rights reserved

REVISION RECORD

For detailed version history and version numbers, see Livelink.

IHE Integration Statement, Global Template, NodeID : 8818216		
Revision Number	Date	Reason for Change
1.0	28 Sep 2005	Initial revision
1.1	01 Feb 2007	Update on Agfa logo and Revision Record layout.
1.2	19 Mar 2008	Insert Livelink NodeID in header and front-page, adjust links.
1.3	25 Jan 2011	Update layout + links

IHE Integration Statement NX 3.0.9000		
Revision Number	Date	Reason for Change
2	2016-11-30	Initial version – include NodeID

Agfa HealthCare

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	NX	3.0.9000	December, 2016
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
ITI – IT Infrastructure Domain			
ATNA - Audit Trail and Node Authentication	Secure Application	none	
CT - Consistent Time	Time Client	none	
RAD - Radiology Domain			
CPI - Consistent Presentation of Images	Acquisition Modality	none	
	Evidence Creator	none	
	Print Composer	none	
ED - Evidence Documents	Acquisition Modality	none	
MAMMO - Mammo Integration Profile	Acquisition Modality	none	
PDI - Portable Data for Imaging	Portable Media Creator	none	
PIR - Patient Info Reconciliation	Acquisition Modality	none	
REM - Radiation Exposure Monitoring	Acquisition Modality	none	
SWF - Scheduled Workflow	Acquisition Modality	Broad Worklist Query PPS Exception Management Billing and Material Management	
Internet address for vendor's IHE information: http://global.agfahealthcare.com/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	http://global.agfahealthcare.com/hl7		
DICOM	http://global.agfahealthcare.com/dicom		
Links to general information on IHE			
In North America: www.ihe.net		In Europe: www.ihe-europe.net	
		In Japan: www.ihe-j.org/	