

AGFA HEALTHCARE IHE Integration Statement

→ IMPAX RIS Connectivity 6.4.x

Document No.: 001631
Revision: 10
NodeID Livelink: 74497356

When printed, this is NOT a controlled copy

Document Information

Service-related contact information worldwide	All service-related contact information is available on this URL→	https://www.agfahealthcare.com/contact/
--	---	---

Issued by:
Agfa HealthCare
V&V Connectivity
Septestraat 27
B-2640 Mortsel

email: connectivity@agfa.com

Agfa HealthCare shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

Copyright © 2025
Agfa HealthCare
All rights reserved

REVISION RECORD

IHE Integration Statement Template NodeID: 8818216		
Revision Number	Date	Reason for Change
1.0	28 Sep 2005	Initial revision
1.1	01 Feb 2007	Update on Agfa logo and Revision Record layout.
1.2	19 Mar 2008	Insert Livelink NodeID in header and front page, adjust links.
1.3 (14)	25 Jan 2011	Update layout + links

IHE Integration Statement IMPAX RIS Connectivity 6.4.x		
Revision Number	Date	Reason for Change
2-4	November 19 to 30, 2020	Initial Revision and corrections made in review cycle
5-6	October 29 to November 10, 2021	Document is also valid for version 6.4.1
7	September 2022	Added version 6.4.2
8-9	January 2024	Added version 6.4.3
10	May 2025	Added version 6.4.4

IHE Integration Statement

Vendor	Product Name	Version	Date
AGFA HealthCare	IMPAX RIS Connectivity	6.4.0	December 2020
		6.4.1	November 2021
		6.4.2	November 2022
		6.4.3	March 2024
		6.4.4	June 2025

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Patient Administration Management (PAM)	Patient Demographics Consumer (PDC)	Merge
	Patient Encounter Consumer (PEC)	None
Charge Posting (CHG)	Department System Scheduler / Order Filler	None
	Performed Procedure Step Manager	None
Patient Demographics Query (PDQ)	Patient Demographics Consumer	None
Consistent Time (CT)	Time Client	None
Audit Trail and Node Authentication (ATNA)	Secure Application	None

Internet address for vendor's IHE information: <https://www.agfahealthcare.com/ihe>

Links to Standards Conformance Statements for the Implementation

HL7	https://www.agfahealthcare.com/hl7/
FHIR	https://www.agfahealthcare.com/fhir-conformance-statement/
DICOM	https://www.agfahealthcare.com/dicom-conformance/

Links to general information on IHE

International: www.ihe.net	In Europe: www.ihe-europe.net/	In Japan: www.ihe-j.org/
---	---	--