

AGFA HEALTHCARE IHE Integration Statement

→ Agfa HealthCare Engage Suite 3.6

Document No. 001636
Revision: 3
NodeID Livelihood: 74741533

When printed, this is NOT a controlled copy

Document Information

Service-related contact information worldwide	All service-related contact information is available on this URL→	http://global.agfahealthcare.com/main/contact/

Issued by:
Agfa HealthCare
V&V Connectivity
Septestraat 27
B-2640 Mortsel

tel: +32 3 444 9367
email: connectivity@agfa.com

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

Copyright © 2021
Agfa HealthCare
All rights reserved

REVISION RECORD

For detailed version history and version numbers, see Livelink.

IHE Integration Statement, Global Template, NodeID : 8818216		
Revision Number	Date	Reason for Change
1.0	28 Sep 2005	Initial revision
1.1	01 Feb 2007	Update on Agfa logo and Revision Record layout.
1.2	19 Mar 2008	Insert Livelink NodeID in header and front-page, adjust links.
1.3	25 Jan 2011	Update layout + links

IHE Integration Statement Engage Suite 3.6		
Revision Number	Date	Reason for Change
2-3	2021-01-14 to 01-29	Initial Revision and corrections made in review cycle

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	Engage Suite	3.6	February 2021
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
ATNA - Audit Trail and Node Authentication	Secure Application	None	
CT - Consistent Time	Time Client	None	
PIX / PIXV3 - Patient Identifier Cross-referencing	Patient Identifier Cross-reference Consumer	None	
PDQ / PDQV3 - Patient Demographics Query	Patient Demographics Consumer	None	
XDS - Cross-Enterprise Document Sharing	Document Consumer	None	
	Document Source	None	
Internet address for vendor's IHE information: https://global.agfahealthcare.com/iheintegration/			
Links to Standards Conformance Statements for the Implementation			
HL7	https://global.agfahealthcare.com/hl7/		
DICOM	https://global.agfahealthcare.com/dicomconformance/		
Links to general information on IHE			
In North America: www.ihe.net		In Europe: www.ihe-europe.net	In Japan: www.ihe-j.org/