

AGFA NV

IHE Integration Statement

→ **NX 3.0.25.00 / 4.0.25.00**

Document No. 005002
Revision: 4
NodeID Livelink: 82942843

When printed, this is NOT a controlled copy

Document Information

Service-related contact information worldwide	All service-related contact information is available on this URL→	https://medimg.agfa.com/main/contact/

Issued by:
Agfa NV
Septestraat 27
B-2640 Mortsel

email: connectivity@agfa.com

Agfa shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. Agfa reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of Agfa equipment.

Copyright © 2024
Agfa NV
All rights reserved

REVISION RECORD

IHE Integration Statement Template NodeID: 76934367		
Revision Number	Date	Reason for Change
2	June 2021	Initial revision
3	June 2021	All corrections made in review cycle

IHE Integration Statement NX 3.0.25.00 / 4.0.25.00		
Revision Number	Date	Reason for Change
4	2024-02-06	Version for NX 25 (LMatthieu)

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA NV	NX	3.0.25.00 / 4.0.25.00	January, 2024
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
ITI – IT Infrastructure Domain			
ATNA - Audit Trail and Node Authentication	Secure Application	none	
CT - Consistent Time	Time Client	none	
RAD - Radiology Domain			
CPI - Consistent Presentation of Images	Acquisition Modality	none	
	Evidence Creator	none	
	Print Composer	none	
ED - Evidence Documents	Acquisition Modality	none	
MAMMO - Mammo Integration Profile	Acquisition Modality	none	
PDI - Portable Data for Imaging	Portable Media Creator	none	
PIR - Patient Info Reconciliation	Acquisition Modality	none	
REM - Radiation Exposure Monitoring	Acquisition Modality	none	
SWF - Scheduled Workflow	Acquisition Modality	Broad Worklist Query PPS Exception Management Billing and Material Management	
Internet address for vendor’s IHE information: https://medimg.agfa.com/main/miscellaneous/interoperability/ihe_integration			
Links to Standards Conformance Statements for the Implementation			
HL7	https://medimg.agfa.com/main/miscellaneous/interoperability/hl7_conf_statements		
DICOM	https://medimg.agfa.com/main/miscellaneous/interoperability/dicom_connectivity		
Links to general information on IHE			
In North America: www.ihe.net		In Europe: www.ihe-europe.net	In Japan: www.ihe-j.org/



Details as of PDF Creation Date

Document Metadata	
Title:	005002_NX_3.0.25.00_4.0.25.00_IHE_Integration_Statement.docx
Livelink ID:	82942843
Version#:	4
Version Date:	2024-02-07 09:24 AM CET
Status:	Approved on 2024-02-08 04:54 PM CET
Owner:	Luc Matthieu (ameqm)
Created By:	Luc Matthieu (ameqm)
Created Date:	2024-02-05 09:12 AM CET
PDF Creation Date:	2024-02-08 04:54 PM CET

This document was approved by:

Signatures:

- 1. Luc Matthieu (ameqm) on 2024-02-07 09:48 AM CET
- 2. Eleonora Rossi (apypu) on 2024-02-07 02:44 PM CET

Detailed Approver History:

- Approval Workflow started on 2024-02-07 09:25 AM CET
 - Approval task originally assigned to and completed by Eleonora Rossi (apypu) on 2024-02-07 02:44 PM CET
 - Approval task originally assigned to and completed by Luc Matthieu (ameqm) on 2024-02-07 09:48 AM CET

Version & Status History

Version#	Date Created	Status
4	2024-02-07 09:24 AM CET	Approved - 2024-02-08
3	2024-02-06 05:09 PM CET	
2	2024-02-05 09:14 AM CET	Approval Cancelled - 2024-02-06
1	2024-02-05 09:12 AM CET	