

# **AGFA HEALTHCARE IHE Integration Statement**

## **AGFA HealthCare Enterprise Imaging 8.4.x**

Document No.: 001685  
Revision: 6  
NodeID Livelink: 83154606

**When printed, this is NOT a controlled copy**

## Document Information

---

<b>Service-related contact information worldwide</b>	All service-related contact information is available on this URL→	<a href="http://global.agfahealthcare.com/main/contact/">http://global.agfahealthcare.com/main/contact/</a>

Issued by:  
AGFA HealthCare  
V&V Connectivity  
Septestraat 27  
B-2640 Mortsel

email: [connectivity@agfa.com](mailto:connectivity@agfa.com)

AGFA HealthCare shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance or use of this publication. AGFA reserves the right to revise this publication and to make changes to its content at any time, without obligation to notify any person or entity of such revisions and changes. This publication may only be used in connection with the promotion, sales, installation and use of AGFA equipment.

**Copyright © 2024**  
**AGFA HealthCare**  
**All rights reserved**

## REVISION RECORD

---

IHE Integration Statement Template NodeID: 8818216		
Revision Number	Date	Reason for Change
1.0	28 Sep 2005	Initial revision
1.1	01 Feb 2007	Update on Agfa logo and Revision Record layout.
1.2	19 Mar 2008	Insert Livelihood NodeID in header and front page, adjust links.
1.3	25 Jan 2011	Update layout + links

IHE Integration Statement Enterprise Imaging 8.4.x		
Revision Number	Date	Reason for Change
2	March 2024	Initial revision for EI 8.4.0.000
3-4	May 2024	Changed version to EI 8.4.1.000
5-6	September 2024	Changed version to EI 8.4.2.000. Included WIA profile.

IHE Integration Statement			
Vendor	Product Name	Version	Date
AGFA HealthCare	Enterprise Imaging <sup>1</sup>	8.4.2.000	September 2024
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented		Actors Implemented	Options Implemented
<b>CARD – Cardiology Domain</b>			
CATH – Cardiac Catheterization Workflow	Image Display	none	
ECHO – Echocardiography Workflow	Image Display	none	
<b>ITI – IT Infrastructure Domain</b>			
ATNA – Audit Trail and Node Authentication	Audit Record Repository	none	
	Secure Application	Radiology Audit Trail	
CT – Consistent Time	Time Client	none	
PIX / PIXV3 - Patient Identifier Cross-referencing <sup>2</sup>	Patient Identifier Cross-reference Consumer	none	
PDQ / PDQV3 – Patient Demographics Query <sup>2</sup>	Patient Demographics Consumer	none	
XDS.b – Cross-Enterprise Document Sharing <sup>2</sup>	Document Consumer	none	
XDS-SD – Cross-Enterprise Sharing of Scanned Documents <sup>2</sup>	Document Consumer	none	
<b>RAD – Radiology Domain</b>			
ARI – Access to Radiology Information	Image Manager / Archive	none	
	Image Display <sup>3</sup>	Multiple Sources	
	Report Reader	none	
CPI – Consistent Presentation of Image	Image Manager / Archive	none	
	Image Display	none	
	Evidence Creator	none	
	Print Composer	none	
CXCAD – Chest X-Ray CAD Display	Image Manager / Archive	none	
DBT – Digital Breast Tomosynthesis	Image Manager / Archive	none	
	Image Display	none	
DIFF – MR Diffusion Imaging	Image Manager / Archive	none	
ED – Evidence Document	Image Manager / Archive	none	
	Image Display	none	
	Evidence Creator	none	
FUS – Image Fusion Integration Profile	Image Manager / Archive	none	
IOCM – Imaging Object Change Management	Change Requester	none	

<sup>1</sup> Enterprise Imaging 8.4.x comprises the following components:

- Enterprise Imaging Core Server (comprised in Enterprise Imaging / Radiology Suite / Cardiology suite / VNA)
- Enterprise Imaging Web Server (further referred to as Web Server for the server part and XERO Viewer for the display part)

<sup>2</sup> Supported by Web Server component only

<sup>3</sup> Some query keys not supported

	Image Manager/Archive	none
	Image Display <sup>4</sup>	none
<b>KIN</b> – Key Image Notes	Image Manager / Archive	none
	Image Display	none
	Evidence Creator	none
<b>MAMMO</b> – Mammography Image	Image Manager/Archive	none
	Image Display <sup>5</sup>	none
<b>NM</b> – Nuclear Medicine Image	Image Manager/Archive	none
	Image Display	none
	Evidence Creator	none
<b>PDI</b> – Portable Data for Imaging	Image Display	none
	Display	none
	Portable Media Creator	none
	Portable Media Importer	none
<b>PERF</b> – CT/MR Perfusion Imaging with Contrast	Image Manager / Archive	none
<b>PGP</b> – Presentation of Grouped Procedures	Image Manager / Archive	none
<b>PIR</b> – Patient Information Reconciliation	Image Manager / Archive	none
	PPS Manager	none
<b>REM</b> – Radiation Exposure Monitoring	Image Manager/Archive	none
<b>RWF</b> – Reporting Workflow	Image Manager/Archive	none
<b>SINR</b> – Simple Image and Numeric Report	Report Reader <sup>5</sup>	none
	Report Repository <sup>5</sup>	none
<b>SMI</b> – Stereotactic Mammography Image	Image Manager / Archive	none
<b>SWF</b> – Scheduled Workflow <sup>6</sup>	Image Manager / Archive	none
	PPS Manager	none
	Image Display <sup>7</sup>	none
	Evidence Creator	none
<b>SWF.b</b> – Scheduled Workflow.b	Image Display	none
<b>TCE</b> – Teaching File and Clinical Trial Export	Export Selector <sup>8</sup>	none
	Export Manager <sup>8</sup>	none
<b>WIA</b> – Web-based Image Access	Imaging Document Source <sup>9</sup>	none
<b>XCA-I</b> – Cross-Community Access for Imaging	Imaging Document Consumer <sup>2</sup>	none
<b>XDS-I.b</b> – Cross-Enterprise Document Sharing for Imaging	Imaging Document Consumer <sup>2</sup>	none
	Imaging Document Source	none
	Document Repository <sup>10</sup>	none
<b>RO – Radiation Oncology Domain</b>		
<b>BRTO-II</b> - Basic Radiation Therapy Objects II	Archive	none

<sup>4</sup> The 'configurable behavior' for handling images rejected for quality reasons is not supported

<sup>5</sup> Instance level query keys not supported

<sup>6</sup> SWF profile: 1 to 1 relationship between Accession Number and Requested procedure AND 1 to 1 relation between RP and SPS

<sup>7</sup> Some query keys not supported

<sup>8</sup> Both actors Export Selector and Export Manager are combined and cannot be dissociated

<sup>9</sup> Work in progress

<sup>10</sup> The Document Repository must be combined with the Imaging Document Source actor and cannot be dissociated.

**AGFA HealthCare**

---

<b>MMRO-III - Multimodality Image Registration for Radiation Oncology 2018</b>	Archive	none
<b>Internet address for vendor's IHE information:</b> <a href="https://www.agfahealthcare.com/ihe">https://www.agfahealthcare.com/ihe</a>		
<b>Links to Standards Conformance Statements for the Implementation</b>		
<b>HL7</b>	<a href="https://www.agfahealthcare.com/hl7">https://www.agfahealthcare.com/hl7</a>	
<b>DICOM</b>	<a href="https://www.agfahealthcare.com/dicom-conformance">https://www.agfahealthcare.com/dicom-conformance</a>	
<b>Links to general information on IHE</b>		
IHE USA: <a href="https://www.iheusa.org">https://www.iheusa.org</a>	In Europe: <a href="http://www.ihe-europe.net">www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">www.ihe-j.org</a>



### Details as of PDF Creation Date

#### Document Metadata

<b>Title:</b>	001685 Enterprise Imaging 8.4.x IHE Integration Statement
<b>Livelink ID:</b>	83154606
<b>Version#:</b>	6
<b>Version Date:</b>	2024-09-23 04:11 PM CET
<b>Status:</b>	Approved on 2024-09-25 03:03 PM CET
<b>Owner:</b>	Bruno Laffin (nawyv)
<b>Created By:</b>	Bruno Laffin (nawyv)
<b>Created Date:</b>	2024-03-19 08:14 AM CET
<b>PDF Creation Date:</b>	2024-09-25 03:03 PM CET

**This document was approved by:**

#### Signatures:

1. Bruno Laffin (nawyv) on 2024-09-24 12:23 PM CET
2. Laura Peters (axdhr) on 2024-09-23 08:05 PM CET

#### Detailed Approver History:

- **Approval Workflow started on 2024-09-23 04:12 PM CET**
  - Approval task originally assigned to and completed by Laura Peters (axdhr) on 2024-09-23 08:05 PM CET
  - Approval task originally assigned to and completed by Bruno Laffin (nawyv) on 2024-09-24 12:23 PM CET

#### Version & Status History

Version#	Date Created	Status
6	2024-09-23 04:11 PM CET	Approved - 2024-09-25
5	2024-09-13 02:01 PM CET	Reviewed - 2024-09-23
4	2024-05-24 09:41 PM CET	Published - 2024-06-04 Approved - 2024-06-04
3	2024-05-22 04:10 PM CET	Reviewed - 2024-05-24
2	2024-03-19 08:18 AM CET	Unpublished - 2024-06-04 Published - 2024-03-26 Approved - 2024-03-25 Reviewed - 2024-03-22
1	2024-03-19 08:14 AM CET	

**Applied Categories and Attributes:**

Document Status: Approved

Livelink ID: 83154606 Version: 6

<b>ITCo Library</b>	
Document Type:	Conformance Statement
Category > SubCategory > Item:	Connectivity > Connectivity & Application > To Agfa
Content Manager:	Bruno Laffin (nawyv)
Summary:	The document 001685 is the IHE Integration Statement for Enterprise Imaging 8.4.1.000.
Language:	English
Availability:	Internet
Language Master Document:	
Software Download URL:	