

CMYK DISTRIBUTORS ANNOUNCES

Filter Klear™

Fountain Recycling System



Eliminate Wastewater in the pressroom

Filter Klear utilizes advanced nanofiber technology featuring ceramic nanofibers on a micro-glass matrix. They provide outstanding reduction of virus, bacteria, cysts, endotoxin and many other submicron particles through both electro-adhesion and mechanical processes.

The result is ZERO discharge and Zero waste of incoming water



Keep your tank running clean

Benefits:

- ◆ Zero discharge
- ◆ Biofilter technology
 - cuts odors and bacteria
- ◆ Lightweight and portable
- ◆ Easy to install
 - connects directly to press via air supply hose
- ◆ Cost-effective
- ◆ Reduces press downtime and maintenance
- ◆ Enhances print quality
 - prevents dust and gunk in water trays
 - better hickey control
 - less plate wear
 - less scumming
- ◆ Decreases fountain solution consumption cost

The economical Litho Filter will reduce your costs and increase the productivity by extending the life of your fountain solution.



271 Essex Street • Weymouth, MA 02188
Exclusive distributors of FilterKlear™

Made in the USA

For information and sales contact:
CMYK DISTRIBUTORS, INC.
Sales: 973-459-0524
markwilliams@cmkydistributors.com

