



## Technical Data Sheet

### DYNAMIC DUO SYSTEMS CLEANER

#### READY TO USE, TWO PART FOUNTAIN SOLUTION SYSTEMS CLEANER

#### DESCRIPTION:

Dynamic Duo is an outstanding, two part ready to use fountain solution recirculation systems cleaner. The product helps to extend the life of fountain solution by cleaning and removing ink, coating, mold and paper from the press recirculation system. It even helps to eliminate calcium build up. Best of all, it is a safer to use technology than the competition since it does not have to be pre-mixed and it is a non-regulated product for shipping.

#### DIRECTIONS FOR USE:

Add **PART ONE** directly to the dirty fountain solution and allow it to circulate for 30 minutes. Do not drain the system after using **PART ONE**. Add **PART TWO** directly to the dirty fountain solution tank and let it circulate for 10 minutes. Drain and flush the system with clean water. Add fresh fountain solution for press start up.

One quart of each part will treat 5-8 gallons of dirty fountain solution. One gallon of each part will treat 20-32 gallons of dirty fountain solution.

**READ SDS BEFORE USE.**

#### FEATURES AND BENEFITS:

- \* Ready to Use. No pre-mixing or splashing.
- \* Non-regulated product for shipping improves cost savings.
- \* Outstanding cleaning power with excellent neutralization for fast start ups.
- \* Zero VOC formulation.
- \* Easy to use and takes only 40 minutes to clean the entire system.

#### SPECIFICATIONS:

Flashpoint: None

VOC: None

Odor: Detergent

#### HMIS RATING:

Health: 1 Flammability: 0 Reactivity: 0 Personal Protection: B

#### PACKAGING:

6 quarts per case (3 Part One and 3 Part Two) / 4 of gallons per case (2 Part One and 2 Part Two)